



sartorius stedim  
biotech

BACTair™



turning science into solutions

# Big Impact.

## A New Impact for Microbiological Air Monitoring

As regulations become more stringent, the microbial monitoring of ambient air is increasing importance in today's world. For example, microorganisms can have a major impact on product quality and the production process, and can even pose health hazards.

Air quality plays a key role in the pharmaceutical, biotechnological and food and beverage industries, hospitals and in the field of occupational and environmental protection.

The most frequently used method today for sampling airborne microorganisms is based on the Andersen principle, which traps particles on culture media plate by impaction. In this method, air is suctioned through a sieve, accelerated and directed against a culture medium plate. Due to their inertia, airborne organisms are prevented from

being swept away by the diverted stream of air and are impacted onto the culture medium plate. After sampling, the culture medium plate is incubated and the colonies grown are counted as colony-forming units/m<sup>3</sup> of air (cfu/m<sup>3</sup>).

Sartorius Stedim Biotech has developed a new system for sampling airborne microorganisms that allows impaction onto culture media plates, where the plates function directly as collection heads. This means that the collection properties are integrated right into the culture media plates. Metal sieve plates or metal collection heads with slots, which have to be sterilized for routine samplings on a regular basis are eliminated. Now, non-sterile sieves or slots have become a thing of the past.

The geometry of the culture medium plate and the 400 holes in the sieve plate yield exceptional sampling efficiency, which is generally higher than that of other impaction samplers.



This new method uses the AirPort MD8 air sampler to pump the air stream.



BACTair™ culture media plates are ready-to-connect to the AirPort MD8



BACTair™ culture media plate prior sampling

# Only 3 handling steps with BACTair™ provide you with an airborne microorganisms

1) Just connect your BACTair™ plate



2) Sample (Press START)



3) Remove your BACTair™ plate and incubate



## BACTair™ Features

- Gamma-sterile
- Integrated disposable sieve
- Pre-filled with agar medium
- Individually packaged

## BACTair™ Benefits

- No sterilization required
- No handling of re-usable sieves
- No preparation of media
- No desiccation effects

- Samples 1m<sup>3</sup> in just 8 min
- Optimized geometry
- Filled with sufficient amount of media

- Fast sampling
- High recovery efficiency (details on page 6)
- No effects due to evaporation (details on page 6)

- Protection with covers
- And after incubation:
- Optimized geometry
- No correlated sampling head

- The agar surface is protected
- No colony overlapping means no correction factor (details on page 6)
- No complicated correlation of sampling heads and devices

The detection of airborne microbial contamination has never been so easy and so reliable!



# exceptional recovery of

**BACTair™ saves your labor-time and guarantees reliable results:**

Preparation- and sterilization-free procedure reduce the risk of secondary contamination.

Culture media plate properties are maintained from purchase through to sampling.

Detects total viable airborne microorganisms in a very short sampling time.

Easy handling with no risk of secondary contamination.


Results are easy to access and evaluate.

Makes your calibration and sampler management fast and easy.

BACTair™  
culture media plate  
after sampling  
and incubation




# BACTair™ stands for optimized geometry



## Recovery Comparison Study

Aerosols of three different bacteria suspensions are released into a sampling channel under defined conditions. At the end of the channel the aerosols are sampled using four different air monitoring methods.


As a reference method two parallel impinger systems were used (100 % recovery). BACTair™ culture media plates show the highest recovery of bacteria due to optimized geometry and complete disposable design.



## Correction Factors

The recovery of colony-forming units depends on the number of holes in the sieve plate onto which the air is impacted. The lower the number of holes, the higher is the risk that more than one microorganism will pass into one single hole (colony overlapping).

A correction factor (K) can be calculated mathematically:  $K = x \cdot (\ln(x) - \ln(x-n))$   
 $x$  = number of holes in the sieve plate  
 $n$  = colony count.  
 BACTair™ culture media plates provide 400 air impaction holes. A correction factor is only relevant with very high CFU counts.



## Evaporation Effects

During the sampling period, agar media may dry out, thus inhibiting the growth of the collected microorganisms. The weight of several BACTair™ culture media plates was measured prior to and after sampling and compared to other impaction agar plates or strips.

Because of its optimal design BACTair™ shows the lowest desiccation effects and enables viable microorganisms to grow under optimal conditions.

The number of impaction-holes and their optimal position in relation to the sampling area are the key to reliable results.

# Ordering Information

Description	Order Number
<b>BACTair™ Culture Media Plates</b>	
BACTair™ – culture media plates Tryptic Soy Agar (TSA), 110 mm, individually, sterile packaged, 10 units	14320-110----ACD
BACTair™ – culture media plates Sabouraud agar (acc. USP), 110 mm, individually, sterile packaged, 10 units	14321-110----ACD
Other BACTair™ culture media types on request	
<b>Air Sampler</b>	
AirPort MD8 Air Sampler for BACTair™ culture media plates incl. charger	16757
<b>Accessories</b>	
Adapter for BACTair™ culture media plates on the AirPort MD8 air sampler	17803
Covers for BACTair™ culture media plates, 10 × 2 units individually, sterile packaged	1ZPX-D0002



# Sales and Service Contacts

For further contacts, visit [www.sartorius-stedim.com](http://www.sartorius-stedim.com)

## Europe

**Germany**  
Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen  
Phone +49.551.308.0  
Fax +49.551.308.3289  
[www.sartorius-stedim.com](http://www.sartorius-stedim.com)

Sartorius Stedim Systems GmbH  
Schwarzenberger Weg 73-79  
34212 Melsungen  
Phone +49.5661.71.3400  
Fax +49.5661.71.3702  
[www.sartorius-stedim.com](http://www.sartorius-stedim.com)

**France**  
Sartorius Stedim Biotech S.A.  
ZI Les Paluds  
Avenue de Jouques – BP 1051  
13781 Aubagne Cedex  
Phone +33.442.845600  
Fax +33.442.845619

Sartorius Stedim France SAS  
ZI Les Paluds  
Avenue de Jouques – CS 71058  
13781 Aubagne Cedex  
Phone +33.442.845600  
Fax +33.442.846545

**Austria**  
Sartorius Stedim Austria GmbH  
Franzosengraben 12  
A-1030 Vienna  
Phone +43.1.7965763.18  
Fax +43.1.796576344

**Belgium**  
Sartorius Stedim Belgium N.V.  
Leuvensesteenweg, 248/B  
1800 Vilvoorde  
Phone +32.2.756.06.80  
Fax +32.2.756.06.81

**Denmark**  
Sartorius Stedim Nordic A/S  
Hoerskaetten 6D, 1.  
DK-2630 Taastrup  
Phone +45.7023.4400  
Fax +45.4630.4030

**Italy**  
Sartorius Stedim Italy S.p.A.  
Via dell'Antella, 76/A  
50012 Antella-Bagno a Ripoli (FI)  
Phone +39.055.63.40.41  
Fax +39.055.63.40.526

**Netherlands**  
Sartorius Stedim Netherlands B.V.  
Edisonbaan 24  
3439 MN Nieuwegein  
Phone +31.30.6025080  
Fax +31.30.6025099

**Spain**  
Sartorius Stedim Spain SA  
C/Isabel Colbrand 10,  
Oficina 70  
Poligono Industrial de Fuencarral  
28050 Madrid  
Phone +34.90.2110935  
Fax +34.91.3589623

**Switzerland**  
Sartorius Stedim Switzerland GmbH  
Lerzenstrasse 21  
8953 Dietikon  
Phone +41.44.741.05.00  
Fax +41.44.741.05.09

**U.K.**  
Sartorius Stedim UK Limited  
Longmead Business Park  
Blenheim Road, Epsom  
Surrey KT19 9 QQ  
Phone +44.1372.737159  
Fax +44.1372.726171

## America

**USA**  
Sartorius Stedim North America Inc.  
5 Orville Drive  
Bohemia, NY 11716  
Toll-Free +1.800.368.7178  
Fax +1.631.254.4253

Sartorius Stedim SUS Inc.  
1910 Mark Court  
Concord, CA 94520  
Phone +1.925.689.6650  
Toll Free +1.800.914.6644  
Fax +1.925.689.6988

Sartorius Stedim Systems Inc.  
201 South Ingram Mill Road  
Springfield, MO 65802  
Phone +1.417.873.9636  
Fax +1.417.873.9275

**Argentina**  
Sartorius Argentina S.A.  
Int. A. Avalos 4251  
B1605ECS Munro  
Buenos Aires  
Phone +54.11.4721.0505  
Fax +54.11.4762.2333

**Brazil**  
Sartorius do Brasil Ltda  
Av. Dom Pedro I, 241  
Bairro Vila Pires  
Santo André  
São Paulo  
Cep 09110-001  
Phone +55.11.4451.6226  
Fax +55.11.4451.4369

**Mexico**  
Sartorius de México S.A. de C.V.  
Circuito Circunvalación Poniente No. 149  
Ciudad Satélite  
53100 Naucalpan, Estado de México  
Phone +52.5555.62.1102  
Fax +52.5555.62.2942

## Asia | Pacific

**Australia**  
Sartorius Stedim Australia Pty. Ltd.  
Unit 5, 7-11 Rodeo Drive  
Dandenong South Vic 3175  
Phone +61.3.8762.1800  
Fax +61.3.8762.1828

**China**  
Sartorius Stedim Beijing  
Representative Office  
No. 33, Yu'an Road,  
Airport Industrial Zone B, Shunyi District  
Beijing 101300  
Phone +86.10.80426516  
Fax +86.10.80426580

Sartorius Stedim Shanghai  
Representative Office  
Room 618, Tower 1, German Centre,  
Shanghai, PRC., 201203  
Phone +86.21.28986393  
Fax +86.21.28986392.11

Sartorius Stedim Guangzhou Office  
Room 704, Broadway Plaza,  
No. 233-234 Dong Feng West Road  
Guangzhou 510180  
Phone +86.20.8351.7921  
Fax +86.20.8351.7931

**India**  
Sartorius Stedim India Pvt. Ltd.  
#69/2-69/3, Jakkasandra  
Kunigal Road, Nelamangala Tq  
Bangalore – 562 123  
Phone +91.80.4350.5361  
Fax +91.80.4350.5253

**Japan**  
Sartorius Stedim Japan K.K.  
KY Building, 8-11  
Kita Shinagawa 1-chome  
Shinagawa-ku  
Tokyo 140-0001  
Phone +81.3.3740.5407  
Fax +81.3.3740.5406

**Malaysia**  
Sartorius Stedim Malaysia Sdn. Bhd.  
Lot L3-E-3B, Enterprise 4  
Technology Park Malaysia  
Bukit Jalil  
57000 Kuala Lumpur  
Phone +60.3.8996.0622  
Fax +60.3.8996.0755

**Singapore**  
Sartorius Stedim Singapore Pte. Ltd.  
10, Science Park Road, The Alpha  
#02-25, Singapore Science Park 2  
Singapore 117684  
Phone +65.6872.3966  
Fax +65.6778.2494