LABFLAME micro

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The smallest professional among the laboratory gas burners with button function, IR-Sensor, foot pedal and SCS safety package

Top Rate Safety

The Safety Control System SCS® demonstrates state-of-the-art safety technology which constantly analyzes potential hazards and, if necessary, initiates safety measures, such as an interruption of the gas supply. In addition to the ignition and flame control function and the overheating protection, the new safety concept also features continuous burner head control (BHC). BHC unfailingly identifies burner head clogging by liquids or solid substances and ensures the correct assembly of the burner head. Exceptional passive safety features: After prolonged use, the residual heat display will signal that the burner head is still hot in order to protect the user from burns.

Simple and Safe Operation

The *Labflame micro* can be operated with the button function, foot pedal or touch free IR-Sensor. For all options, different programs exist that have been developed from practical experience: flexible start-stop functions, continuous application for up to 60 minutes or programs for short flame sterilization – exact to the second.

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LABFLAME micro

Safety laboratory gas burner
with button function, IR-Sensor
and foot pedal

Art.-Nr./Art.-No.
LF6.004.000

Exclusive – The DoubleClick IR-Sensor

This safety function ensures that the burner can only be ignited by activating the DoubleClick IR-Sensor twice. As a result, unintentional ignition or ignition due to dropping or falling objects is virtually impossible, increasing safety considerably.

Robust

All external components are made of stainless steel; this guarantees the greatest possible working sterility. A tilt mechanism protects the burner tube against soiling. Spilled liquids simply escape through an integrated drain while solid substances can easily be removed by opening the burner shaft. It is even possible to dismantle the individual components of the removable burner head for in-depth cleaning.

Labflame micro - Ideal for use in cleanroom workbenches and laboratories.

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Labflame micro

Art.-No. LF6.004.000

with button function, stainless steel foot pedal and
DoubleClick IR-Sensor, adjustable IR-Sensor reaction distance
5 standard programs for button (function knob), IR-Sensor and foot pedal
SCS (Safety Control System) with BHC (Burner Head Control)
Removable and decomposable burner head,
Tilt mechanism, right / left
Holding device for 3 inoculation loop holders
Nozzles for natural gas, propane / butane gas
Turbo flame
Wrench 17 mm for gas connection
Screwdriver for burner head and cover of burner shaft
Tubing connector with swivel nut
Switching power supply (global)
Instruction manual and 2-year warranty

Technical Data:

Technology	Microprocessor
Programs	
Foot pedal:	Standard (flame during pressed foot pedal) Start-Stop with timer, 60 min
IR Sensor:	Start-Stop with timer, 60 min Auto-Off with learn function, 1 sec - 60 min
Button:	Start-Stop with timer, 60 min

Safety features

Safety Control System (SCS)
with gas safety cut off: ignition and flame control
temperature monitor
burner head clogging
and assembly monitor (BHC)
automatic unit switch off, 4 h
residual heat display

Gas supply and consumption

Gas connection: 1/4" left with gas filter
Gas types: I12ELL3B/P: natural gas E/LL, 18 - 25 mbar
liquid gas, 20 - 50 mbar
Connected load: 70 g/h liquid gas
Continuous cartridge operation: CV 360 - 40 min, Express 444 - 50 min,
CG 1750 - 150 min, C 206 - 170 min,
CP 250 - 210 min, CV 470 - 370 min

Temperatures

Flame temperature: 1350 °C on liquid gas
1300 °C on natural gas (E)
Temperature threshold level: 1.0 kW liquid gas, 1.0 kW natural gas

Electrical

Power consumption: 2 VA
Power connection: 100 - 240 V / 50/60 Hz / max. 0.3 A
9 V DC / 1 A

IR-Sensor coverage: 35 - 70 mm, adjustable
DoubleClick IR-Sensor: Time range 1 sec

Mechanical

Casing and operating controls: stainless steel / glass, UV and solvent resistant
Burner head: removable and decomposable, stainless steel
Cover of burner shaft: Ø 23 mm, with drains
Measurements (w x h x d): 85 x 49 x 86 mm
Weight: 700 g

Licenses

DIN-DVGW Reg.-No.: NG-2211AS0167
CE: EN 61326-1, EN61000-3-2, EN 61010
EEC guidelines: 89/336/EEC, 73/23/EEC

Turboflamme
im Handumdrehen eine extra
harte und stabile Flamme
Turbo flame
an extremely firm and
consistent flame with only a few
flicks of the wrist

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Fußpedal
aus Edelstahl
Foot pedal
made of
stainless steel

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Biotoool AG