Identification System Reagent

Reagents for the microbial identification systems.

DESCRIPTION

Identification System Reagent is a kit containing the reagents Alpha Naphtol, NaOH 40%, Hydrogen Peroxide, Methyl Red, Sulphanilic Acid, Alpha Naphthylamine, Kovac's Reagent, PYR Reagent, Ninhydrin and Vaseline Oil necessary for the proper use of systems ideated for the identification of bacteria such as Enterosystem 18R (ref. 71618), Streptosystem 12R (ref. 72560) and Listeria System 18R (ref. 71640).

KIT CONTENTS

- 1 x 10 ml vial of Alpha Naphtol.
- 1 x 10 ml vial of NaOH 40%.
- 1 x 5 ml vial of Hydrogen Peroxide.
- 1 x 10 ml vial of Methyl Red.
- 1 x 10 ml vial of Sulphanilic Acid.
- 1 x 10 ml vial of Alpha Naphthylamine.
- 1 x 10 ml vial of Kovac's Reagent.
- 1 x 10 ml vial of PYR Reagent.
- 1 x 10 ml vial of Ninhydrin.
- 1 x 10 ml vial of Vaseline Oil.
- 1 instruction sheet.

METHOD PRINCIPLE

Alpha Naphtol and NaOH 40% are used to carry out the Voges-Proskauer test.
 Hydrogen Peroxide is used to carry out the catalase test.
 Methyl Red is used to carry out the methyl red test.
 Sulphanilic Acid and Alpha Naphthylamine are used to carry out the nitrate test.
 Kovac's Reagent is used to carry out the indole test.
 PYR Reagent is used to carry out the peptidase test.
 Ninhydrin is used to carry out the hippurate test.
 Vaseline Oil is used for the maintenance of a sufficiently anaerobic environment.

REAGENTS

- Alpha Naphtol: 6% (w/v) alpha naphtol dissolved in ethyl alcohol.
- NaOH 40%: 40% (w/v) sodium hydroxide in aqueous solution.
- Hydrogen Peroxide: 7.5% (v/v) hydrogen peroxide solution.
- Methyl Red: 0.02% (w/v) methyl red dissolved in an aqueous solution containing 60% ethyl alcohol.
- Sulphanilic Acid: 0.8% (w/v) sulphanilic acid dissolved in 5N acetic acid.
- Alpha Naphthylamine: 0.5% (w/v) alpha naphthylamine dissolved in 5N acetic acid.
- Kovac's Reagent: 5% (w/v) *p*-dimethylaminobenzaldehyde dissolved in a solution of 25% hydrochloridric acid and 75% isobutyl alcohol.
- PYR Reagent: 0.2% (w/v) dimethylaminocinnamaldehyde dissolved in an aqueous solution containing 10% hydrochloric acid.
- Ninhydrin: 7% (w/v) ninhydrin dissolved in methoxyethanol.
- Vaseline Oil.

TEST PROCEDURE AND RESULTS INTERPRETATION

Refer to the package insert provided with the system being used.

QUALITY CONTROL FOR THE USER

For each reagent, check positive and negative results by using suitable organisms, according to bibliographics references.

PPRECAUTIONS

Identification System Reagent is classifiable as hazardous under current legislation; it is recommended that the Safety Data Sheet be consulted on its use. The product is intended for *in vitro* diagnostic use only and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE AND TRANSPORT CONDITIONS

Store at 2-8°C away from light in its original package, until the expiry date shown on the label. However, our stability studies have shown that the storage or transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, do not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

ELIMINATING USED MATERIAL

After use, used Identification System 18R Reagent and the material that has come into contact with the sample must be decontaminated and disposed of in accordance with the laboratory procedures for the decontamination and disposal of potentially infected material.

BIBLIOGRAPHY

1. Murray, Baron, Pfaller, Tenorev and Tenorev: Manual of Clinical Microbiology (1995).
2. Bayley and Scott's: Diagnostic Microbiology (1986).
3. Edwin H.Lenette: Manual of Clinical Microbiology (1995).

PRESENTATION

Product	Ref.	Content
Identification System Reagent	80260	100-200 tests

TABLE OF SYMBOLS

 LOT	Batch code	 IVD	<i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 REF	Catalogue number	 Temperature limitation	 Contains Sufficient for <n> tests	 Caution, consult accompanying documents	 Do not reuse	

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