

INSTRUCTIONS FOR USE



■ Inactivated Helix Elite™ Molecular Standards

INTENDED USE _____

Inactivated **Helix Elite™ Molecular Standards** are external controls intended for use as external positive controls for qualitative detection by molecular assays.

SUMMARY AND EXPLANATION _____

Molecular tests offer rapid and accurate results regarding the presence of an organism. Proper interpretation of a molecular test requires the use of a positive control. Inactivated **Helix Elite Molecular Standards** are easy-to-use process controls that can be used to monitor the extraction, amplification and detection of molecular assays or instruments. These independent controls may also be used in evaluation of laboratory proficiency and training, or determination of the lot-to-lot consistency of assay consumables as directed by various regulatory requirements and standards.

PRINCIPLES _____

Inactivated **Helix Elite Molecular Standards** are intact pathogens inactivated by chemical, radiological, or heat treatments. Target pathogens are stable at 2°C - 25°C through the expiration date. Each pellet is packaged in a single-use foil pouch. Users should follow assay manufacturers, or laboratory procedures for processing controls.

COMPOSITION _____

The preparation consists of inactivated pathogen(s) stabilized in a proprietary matrix of excipients.

WARNINGS AND PRECAUTIONS _____

- For In Vitro Diagnostic use only.
- For professional use only. To be used by personnel trained in the use of the assay.
- Do not open foil pouch until ready to use.
- Material has been inactivated, but should be handled and disposed of as infectious material.
- Wear proper personal protective equipment.
- Refer to the Safety Data Sheet (SDS) for more detailed information. The SDS can be located on the Microbiologics website at www.microbiologics.com or by contacting Technical Support at 320.229.7045 or U.S. Toll Free 1.866.286.6691.
- These products are not made with natural rubber latex.
- Customer is responsible for ensuring compatibility of the control with the assay in use.



MATERIALS REQUIRED BUT NOT PROVIDED

- Nucleic acid extraction kit and assay
- Instrumentation for detection
- Rehydration buffer such as PBS or viral transport medium
- Pipettors capable of delivering 0.5-1000µl volumes
- Nuclease-free aerosol barrier pipette tips
- Vortex
- Microcentrifuge

INSTRUCTIONS FOR USE

1. Open the pouch and ensure that the pellet is at the bottom of the vial before opening.
2. Rehydrate the lyophilized pellet with the appropriate buffer. A minimum volume of 100 µl is recommended.
3. Vortex the vial for 10 seconds at full speed to mix. Centrifuge to collect the rehydrated, inactivated target material at the bottom of the tube.
4. Use the appropriate volume for the assay being performed and follow laboratory protocols or manufacturer instructions for processing a sample.

Note: Each pellet is intended as a single use test. Dilutions may be performed and used immediately. Storage of the rehydrated or diluted material is not recommended.

STORAGE AND EXPIRATION

Store the Inactivated **Helix Elite Molecular Standards** at 2°C - 25°C in the original packaging up to the indicated expiration date. After opening the foil pouch, use the pellet immediately.

Inactivated **Helix Elite Molecular Standards** should not be used if:

- Stored improperly
- There is evidence of excessive exposure to heat or moisture
- The expiration date has passed

LIMITATIONS

This product may not be suitable for use with all kits and procedures. Only primers and probes that hybridize to sequences of the extracted nucleic acids of the organism should be expected to yield a positive reaction.

MICROBIOLOGICAL STATE

This product was prepared using suitable inactivation methods. While the product has been tested for innocuity, universal laboratory precautions are recommended and material should be treated as though it was a viable specimen.

KEY OF SYMBOLS



Batch Code (Lot)



Catalog Number



Caution, Consult Accompanying Documents



Contains sufficient for <n> tests



In Vitro Medical Device



Manufacturer



Temperature Limitation



Use By



Refer to Instructions for Use



Telephone Number

PRODUCT WARRANTY

- These products are warranted to meet the specifications and performance printed and illustrated in product inserts, instructions, and supportive literature.
- The warranty, expressed or implied, is limited when:
 - The procedures employed in the laboratory are contrary to printed and illustrated directions and instructions
 - The products are employed for applications other than the intended use cited in product inserts, instructions, and supportive literature

WEBSITE

Visit our website, www.microbiologics.com, for current technical information and product availability.

ACKNOWLEDGEMENTS



Microbiologics, Inc.
200 Cooper Avenue North
St. Cloud, MN 56303 USA
www.microbiologics.com

Customer Service

Tel. 320-253-1640
U.S. Toll Free 800-599-BUGS (2847)
Email: info@microbiologics.com

Technical Support

Tel. 320-229-7045
U.S. Toll Free 866-286-6691
Email: techsupport@microbiologics.com


ILLUSTRATED INSTRUCTIONS

1



Open the foil pouch and ensure that the pellet is at the bottom of the vial before opening.

2



Rehydrate the lyophilized pellet with the appropriate buffer.

A minimum volume of 100 μ l is recommended.

3



Vortex briefly to mix. Centrifuge to collect the rehydrated, inactivated target material at the bottom of the tube.

4

Use the appropriate volume for the assay being performed and follow laboratory protocols or manufacturer instructions for processing a sample.

Note: Each pellet is intended as a single use test. Dilutions may be performed and used immediately. Storage of the rehydrated or diluted material is not recommended.