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## Microsart® Funnel 100 | 250 & Microsart® Funnel Dispenser

A sterile, 100 ml and 250 ml respectively disposable funnel for vacuum filtration with Sartorius Stedim Biotech membrane filters. The system allows fast filtration and rapid sampling. Ideal for microbiological and analytical testing of beverages, pharmaceuticals, water and other liquids.

### Instruction for use:

1. Flame the Microsart® base and the inserted frit with the vacuum supply on.
2. Shut off the vacuum source.
3. Flame the blunt-tipped tweezers (16625) for membrane filters.
4. The Microsart® e.motion Dispenser fully automatically dispenses a membrane filter at the push of a button or via optical sensor. In case of single packed membrane filters, the transparent foil of the individually packed membrane filters need to be manually peeled back. Take the membrane filter with the tweezers and place it (yellow protective paper or grid should face upwards) onto the frit of the filtration apparatus.
5. If there is a yellow protective disc make sure to discard it before assembling the Microsart® Funnel 100 | 250 and the base.
6. Pass the lower part of the plastic-bag of the Microsart® Funnel 100 | 250 through the dispenser-opening (Microsart Funnel Dispenser 16A08) and tear open the bag at the perforated position.  
Hang the bag with the perforated hole on the upper part on the hook of the dispenser.
7. Remove a sterile funnel from the Funnel Dispenser.
8. Place the funnel on the Microsart® base. The leak tightness of the assembled unit is ensured by click-fit closure.
9. Add the sample to the funnel. For small sample volumes it is recommended to filter 10 – 20 ml of sterile buffer (e.g. 0.9 % NaCl solution) first.
10. Turn on the vacuum source to filter the sample.
11. Shut off the vacuum source. Rinse the inside of the funnel with 20 – 30 ml of sterile buffer. Turn on the vacuum source to filter.
12. Shut off the vacuum source. Flame the tweezers to remove the filter. Place the filter on a prepared Sartorius Stedim Biotech Nutrient Pad without entrapping any air bubbles.
13. Incubate the filtered sample under the appropriate conditions.

For further information, please refer to the Directions for Use of Nutrient Pad Sets.

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Z.A de Gesvrine – 4 rue Képler – B.P.4125  
44241 La Chapelle-sur-Erdre Cedex – France  
t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00  
commercial@humeau.com