sartorius stedim
biotech

Operating Instructions | Bedienungsanleitung | Mode d'emploi
Istruzioni per l'uso | Instrucciones de manejo

Laboratory Vacuum Pumps | Labor Vakuumpumpen | Pompes à vide de laboratoire | Pompe a vuoto da laboratorio | Bombas de vacío para laboratorio

16694-2-50-06

16694-1-60-06

85032-535-33

English	Page 3
Deutsch	Seite 13
Français	page 23
Italiano	pagina 33
Español	página 43

List of Contents

1.	Description, Operating Conditions	4	<p>You have selected a high-quality Sartorius Stedim Biotech product; the following tips will help you operate it safely, and reliably over a long period of time. Carefully study the operating instructions before using the pumps and observe at all times the relevant instructions to avoid dangerous situations.</p>
2.	Safety	4	
3.	Installation	4	
4.	Operation	5	
5.	Servicing	5	
6.	Cleaning	7	
7.	Trouble Shooting	7	
8.	Replacement Parts	8	
9.	Tables	8	
	Appendix	9	



Symbols

- ⓪ Position in the illustration
- Important point
- ❶ Task
- Advise to users
- ⚠ Warning

Sartorius Stedim Biotech reserve the right to make changes

Declaration of Conformity CE
 Following the EU Directive about electrical equipment for use within certain limits of voltage 73/23/EWG, appendix 3 and following the EC directive about electromagnetic compatibility 89/336/EWG, appendix 1

Company	Sartorius Stedim Biotech GmbH
Address	August-Spindler-Str. 11 D-37079 Göttingen Bundesrepublik Deutschland
	<p>We herewith declare that the device described below fulfills the relevant fundamental safety requirements and health regulations specified by the appropriate EU-Directive, with respect to its design and construction and to the version as commercialized</p> <p>This declaration becomes legally invalid if modifications are performed on the device which have not been certified by Sartorius Stedim Biotech GmbH.</p>
Designation of Device	Laboratory Vacuum Pump
Model, version	6 I
Cat.-No	16694-2-50-06 16694-1-60-06
Relevant directives of the EC	EU Directive about Electromagnetic Compatibility (89/336/EWG, appendix 3) EU Low Voltage Directive (73/23 EWG)
Applied harmonized standards	EN 61000-6-1; EN 61000-3-2/3 ; EN 61000-6-3; EN 55014; EN 61010-1
Applied national standards and technical specifications	
Date and Signature	 
Function of the Signatory	V. Niebel Senior Vice President Dr. Susanne Gerthausen Quality Management Representative

1. Description, Operating Conditions

The pumps transfer, evacuate and compress 100% oil-free.

1.1 Electrical Equipment

- See table 1 for full electrical data.
- The protection class of standard versions is IP 20.
- The pumps are fitted with a thermal switch to protect against overloading.

1.2 Operating Conditions

- Handling air, gases, and vapours at temperatures between +5°C ... +40°C.
- For maximum permissible operating pressure, ultimate vacuum, and flow capacity: see table 2.

The pumps must not be used in areas where there is a danger of explosion.

- The pumps must not be used for liquids.
- Before pumping a medium, the compatibility of materials of pump head, diaphragm and valves with the medium must be checked (for pump materials: see table 3).
- Pumps that may be employed as vacuum pump and compressor must not be used to produce vacuum and pressure at the same time.

If your potential application lies outside the above limits discuss it with our technical adviser (see last page for contact address).

1.3 Ambient Condition

When the pumps are operating the following ambient conditions must be maintained:

- Ambient temperature during operation: between +5°C ... +40°C.
- During operation an adequate supply of air for cooling must be provided.
- The pumps must not be used in areas where there is a danger of explosion.

1.4 Pump materials

See table 3.

2. Safety

- Note that the pumps may only be used for their intended purpose.

The pumps must not be used in areas where there is a danger of explosion.

- For vacuum pumps: The gas discharge at the pressure side must be drained off safely and reliably.

Components connected to the pumps must be designed to withstand the pneumatic performance of the pumps (see table 2).

- Plug the pump only into properly installed grounded outlets.

When the operation of the pump is interrupted by the thermal switch, the pump will re-start automatically after cooling down. Take all care necessary to prevent this leading to a dangerous situation.

Specific safety instructions and measures for the media being handled must be observed.

- Use only original Sartorius Stedim Biotech spare parts.
- The pumps conform to the safety regulations of the EC Low Voltage Directive 73/23 EEC, and of the EC Directive 89/336 EEC concerning Electromagnetic Compatibility. The requirements of the following harmonized standards are fulfilled: EN 61010 part 1, EN 61000-6-1, EN 61000-6-3.

3. Installation

- Choose a safe location (flat surface) for the pump.
- Install the pump so as to ensure adequate flow of air cooling.
- Fit the pump at the highest point in the system, so that condensate cannot collect in the head of the pump – that prolongs working life of structured diaphragm and pump.
- At the pump head, remove the protection plugs from the hose connectors.
- Mounting the accessory parts filter or silencer (if present): If the pump is used as a vacuum pump, mount a silencer at the pressure side if necessary. If the pump is used as a compressor, mount a filter at the suction side if necessary.
- Before mounting the filter or silencer, unscrew the corresponding hose connector from the threads in the pump head.
- Connect the suction and pressure lines (tube ID 4 mm). For flow direction see marking on the pump head.
- Arrange the suction and pressure lines so that condensate cannot run into the pump (sloping lines).
- For vacuum pumps: The gas discharge at the pressure side must be drained off safely and reliably.
- Plug the pump only into properly installed grounded outlets.
- Compare the supply data with the electrical data of the pump. The voltage must not vary by more than +10% and -10% from that shown on the type-plate.

4. Operation

- The pumps can be switched on and off using the power switch.

Specific safety instructions for the media being handled must be observed.

- Before pumping a medium, the compatibility of materials of pump head, diaphragm, and valves with the medium must be checked (for pump materials: see table 3).
- The pumps must not start against pressure or vacuum. This also applies when the pump restarts after the power has been cut off for a short period. If a pump starts under pressure or vacuum, then the thermal switch will be activated and switch the pump off.

When the operation of the pump is interrupted by the thermal switch, the pump will re-start automatically after cooling down. Take all care necessary to prevent this leading to a dangerous situation.

The maximum permissible operating pressure (see table 2) must not be exceeded.

- To prevent the maximum permissible operating pressure being exceeded, restriction or control of the air or gas flow should only be carried out in the suction line.
- If restriction or control of the air or gas flow is made on the pressure side ensure that the maximum permissible operating pressure is not exceeded. Excessive pressure (with all of the related hazards) can be prevented by placing a bypass line with a pressure-relief valve between the pressure and suction sides of the pump. For further information, contact our specialists (contact address: see last page).

- When the pump is at a standstill the inlet and exhaust must be at normal atmospheric pressure.
- Change the filter (accessory) if it is dirty.
- Pumps that may be employed as vacuum pump and compressor must not be used to produce vacuum and pressure at the same time.
- Diaphragm and valve plates are the only parts subject to wear. Wear is usually indicated by a drastic reduction in the pneumatic performance (vacuum, pressure, delivery). When replacing parts proceed as described in section 5.
- Ambient conditions: see chapter 1.3.

5. Servicing

Before working on the pump, disconnect it from the electricity supply by pulling the plug out.

- Diaphragm and valve plates are the only parts subject to wear. They are simple to change.
- Always change valve plates, diaphragm, and sealing rings at the same time. If the diaphragm and the valve plates are not changed at the same time the nominal performance of the pump is not guaranteed after the service.
- If aggressive, toxic or other types of gases hazardous to health have been pumped please observe:
 - 1.) Clean the pump and its components before servicing.
 - 2.) Ensure that the service personnel is not subject to a health hazard during diaphragm and valve plate changes. Apply the necessary safety measures (example: the use of protective gloves).
 - 3.) Ensure that the discarded parts and materials are safely and correctly disposed of.

Parts | tools required

- Service Set (see section 8)
- Philips-head screwdriver No. 1
- Small screwdriver
- Pencil

Changing the structured diaphragm, valve plates, and sealing rings in the following sequence:

- a) Remove pump head
- b) Change structured diaphragm
- c) Change valve plates and sealing rings
- d) Refit pump head

See figs. 4 and 5.

a) Removing pump head

- ❶ Mark the position of head plate ③, cover ⑤ and cover plate ⑥ relative to each other by a drawing line with a pencil. This helps avoid incorrect assembly later.
- ❷ Undo the 4 screws ④ in the head plate and lift the head plate with the cover ⑤ off the pump housing.
- ❸ Mark the position of intermediate plate ② and housing ① relative to each other by a drawing line with a pencil.
- ❹ Lift the intermediate plate ② off the housing ①.

b) Change structured diaphragm

- ❶ Using a small screwdriver, between the housing ① and the outer edge of the structured diaphragm ⑨, carefully lever the edge of the diaphragm lightly upwards.
- ❷ Grip the structured diaphragm ⑨ on opposite sides, unscrew it about two turns (anti-clockwise).
- ❸ Hold the pump with one hand, so that the head is pointing downwards. Turn the structured diaphragm ⑨ anti-clockwise to unscrew it.
- ❹ Take the diaphragm support ⑩ and diaphragm spacer(s) k off the threaded portion of the diaphragm and retain them.
- ❺ Check that all parts are free from dirt and clean them if necessary (see section 6. Cleaning).
- ❻ Put the diaphragm support ⑩ and diaphragm spacer(s) ⑪, in that order, on the threaded portion of the new structured diaphragm ⑨.
- ❼ Screw the new structured diaphragm ⑨, complete with diaphragm support ⑩ and diaphragm spacer(s) ⑪ into the connecting rod (clockwise) and tighten it by hand.

c) Changing the valve plates

- ❶ Remove the valve plates ⑦ and sealing rings ⑧ from the intermediate plate ② (sealing rings could cling to the head plate).
- ❷ Check that the valve seats in the head plate ③ and intermediate plate ② are clean. If scratches or distortion are evident on these parts they should be replaced.
- ❸ Lay the new valve plates ⑦ in the recesses in the intermediate plate ②. The valve plates for suction and pressure sides are identical, as are upper and lower sides of the plates.
- ❹ Check that the valve plates ⑦ are not deformed by moving them gently sideways in their recesses.
- ❺ Lay the new sealing rings ⑧ on the intermediate plate ②.

d) Refitting the pump head

- ❶ Place the intermediate plate ②, with valve plates ⑦ and sealing rings ⑧ on the housing ①, in the position indicated by the drawing line.
 - ❷ Place the head plate ③, with cover ⑤ on the housing ①, in the position indicated by the drawing line.
 - ❸ Check that the head plate ③ is centred by moving it gently sideways.
 - ❹ Tighten the screws ④, evenly and diagonally, first gently, then firmly.
- If you have any questions about servicing call our technical adviser (see last page for contact address).


Fig. 4: Pump parts


Fig. 5: Pump parts

Specification

- ❶ Housing
- ❷ Intermediate plate
- ❸ Head plate
- ❹ Screw
- ❺ Cover
- ❻ Cover plate
- ❼ Valve plate
- ❽ Sealing ring
- ❾ Structured diaphragm
- ❿ Diaphragm support
- ⓫ Diaphragm spacer(s)
- ⓬ Connecting rod
- ⓭ Eccentric
- ⓮ Counter weight

6. Cleaning

- When changing structured diaphragm and valve plates, inspect all parts for dirt before assembling the pump head and clean them if necessary.
- As far as possible, clean the parts with a dry cloth. Only use solvents for cleaning if the head materials cannot be attacked (check the resistance of the material!).

7. Trouble Shooting

Before working on the pump isolate the power supply securely, then check that the lines are not live.

- The following tips for fault-finding are best employed in the sequence shown.

Pump produces no flow

- Thermal switch has opened due to over-heating.
- Disconnect pump from mains and allow to cool. Trace cause of over-heating and eliminate it.
- Fuses in the pump are defective (only authorized | qualified personnel should investigate this problem)
- Remove the pump from the source of electrical power by pulling the power plug.
- At the pump housing, fuses are located at the IEC plug.
- Fuse ratings: see Table 1.
- Order numbers for fuses: see Chapter 8.
- Connections or lines are blocked.
- An external valve is closed, or a filter blocked.
- Liquid (condensate) has collected in the pump head.
 - Let the pump run for a few minutes pumping air.
 - Install the pump at the highest point in the system.
- Diaphragm or valve plates are worn.
 - Section 5 Servicing.

Flow, pressure, or vacuum too low

- Compare the actual performance with the figures in table 2.

- There is pressure on the pressure side, and at the same time vacuum, or a pressure above atmospheric, on the suction side.
 - Change pneumatic conditions.
- Liquid (condensate) has collected in the pump head.
 - Let the pump run for a few minutes pumping air.
 - Install the pump at the highest point in the system.
- The cross-section of pneumatic lines, or connected components is too small, or they are restricted.
 - To measure the performance, disconnect the pump from the system (small diameter tubing or a valve can significantly affect performance).
- There is a leak at a connector, in a line, or in the pump head.
- Diaphragm or valve plates are worn, or dirt is in the head.
 - Section 5 Servicing.
- If the pump does not operate properly and you cannot find any of the above faults, send it to the Sartorius Service Department.
- In order for Sartorius Stedim Biotech to repair the pump, the customer must provide a statement on the media which were pumped and on pump cleaning. Please fill out the corresponding Sartorius Stedim Biotech form, and submit it together with the pump. A sample statement for copying can be found in the Appendix of these operating instructions.

8. Replacement Parts

A Service Set contains all replacement parts needed for one complete service: 1 structured diaphragm, 2 valve plates and 2 sealing rings.

Order-No.
1ED----0054

Fuses
There are two fuses for each pump.

Electrical Version	Order-No. fuse (1 piece)
230 V	025250
115 V	029986
100 V	020085

9. Tables

Table 1: Electrical Data

	Electrical Type	
Voltage	230 V	115V
Frequency	50 Hz	60 Hz
Power consumption	65 W	60 W
Operating current	0.63 A	1.2 A
Fuse* (2×) T (A)		

* Order-No. see section 8

Table 2: Pneumatic Data

Delivery (l/min*)	Ultimate vacuum (mbar abs)	Maximum permissible operating pressure (bar g)
6	100	2.4

* Litre at SPT at atm. pressure

Table 3: Pump materials

Pump head	Structured diaphragm	Valves
PPS	EPDM	FPM

Material abbreviations according to DIN ISO 1629 and 1043.1

Appendix: Customer statement for repair order (sample statement for copying)

- In order for Sartorius Stedim Biotech to repair the pump, the customer must provide a statement on the media which were pumped and on pump cleaning. Please fill out the corresponding Sartorius Stedim Biotech form, and submit it together with the pump.

Statement | Certificate

We confirm that the pump model listed below (please specify)

Serial-No. (please specify)

was used to pump the following media:

and that the pump listed above was cleaned. There are no poisonous, aggressive, biological, radioactive or other dangerous media in the pump.

Company

Date | Signature

Information and Instructions on Disposal

Packaging that is no longer required must be disposed of at the local waste disposal facility. The packaging is made of environmentally friendly materials

that can be used as secondary raw materials.

The equipment, including accessories and empty non-rechargeable and rechargeable batteries, does not belong in your regular household waste; such equipment is manufactured from high-grade materials and can be recycled and reused. The European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) requires that electrical and electronic equipment be collected and disposed of separately from other unsorted municipal waste with the aim of recycling it. The crossed out waste bin symbol indicates that separate collection is required.

In Germany and many other countries, Sartorius AG takes care of the return and legally compliant disposal of its electrical and electronic equipment on its own. Such equipment may not be thrown out with household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators. Please contact our on-site service technicians or our Service Center in Goettingen:

Sartorius AG

Service Center
Weender Landstrasse 94-108
37075 Goettingen, Germany

Tel.: +49.551.308.3333 | 4444
Fax: +49.551.308.3730

E-Mail:
Jan.Naumann@sartorius.com

For equipment disposal in the following countries of the European Economic Area (EEA), please contact your local Sartorius subsidiary or distributor:

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In countries that are not members of the European Economic Area (EEA) or where no Sartorius affiliates, subsidiaries, dealers or distributors are located, please contact your local authorities or a commercial disposal operator.

Sartorius AG, its affiliates, subsidiaries, dealers and distributors will not take back equipment contaminated with hazardous materials (ABC contamination). Before returning any equipment to us, the sender must fill out the attached Non-hazardous Declaration completely and return it to us in advance. No repair or return is possible unless a completed Non-hazardous Declaration form has been submitted.

Prior to disposal or scrapping of the equipment, remove any rechargeable or non-rechargeable batteries and, if they are empty, dispose of them in local collection boxes. Only for customers in Germany: On request, Sartorius will have the Foundation of Manufacturers' Mutual Return System for Batteries (GRS, Stiftung Gemeinsames Rücknahmesystem Batterien) provide customers in Germany with battery collection boxes.

Note: In Germany, Sartorius AG is entered in the EAR Registry of the Foundation for Used Electrical and Electronic Equipment (Stiftung Elektro-Altgeräte Register). The WEEE registration number for Germany is DE 49923090.

Non-hazard Declaration

Declaration of Safety and Acceptability as Non-Hazardous to Health

The safety and health of our employees, the legal regulations governing the handling of materials that pose hazards to health and the environment and the regulations governing occupational safety mandate and require that this Non-hazardous Declaration be filled out completely for all equipment returned to the respective Sartorius Stedim Biotech affiliate, subsidiary, dealer or distributor listed for your country. No repair or return is possible unless a completed Non-hazardous Declaration form has been submitted.

Fax or mail a copy of a completed Non-hazardous Declaration form in advance to the respective Sartorius Stedim Biotech affiliate, subsidiary, dealer or distributor listed for your country (for the fax number or address, please refer to our list of countries). We require this information before the equipment | component arrives. An additional copy must accompany the equipment | component. If necessary, the shipping company must be notified.

Incomplete entries or non-compliance with this procedure will automatically lead to considerable delays in processing.

In the event of return, we do not differentiate between

- still utilizable and | or repairable used electrical and electronic equipment and no longer utilizable and | or repairable used electrical and electronic equipment or
- between waste that is subject to monitoring by the supervisory authorities and waste that is not subject to any special monitoring by the supervisory authorities (used electrical and electronic equipment).

Model: _____ Serial No.: _____

I | we guarantee that the following measures have been taken:

- The equipment | component has been freed of hazardous materials to ensure that persons involved with the handling | repair are not exposed to any hazard or danger whatsoever.
- The equipment has been packaged for safety and the complete labeling has been affixed to it.
- The shipping agent has been notified (as prescribed) about the hazardous nature of the shipment.

The person sending in the equipment | component is hereby notified that said person shall be liable to Sartorius Stedim Biotech, its affiliates, subsidiaries, dealers and distributors as well as to any third parties – in particular to any of employees of Sartorius Stedim Biotech, its affiliates, its subsidiaries, dealers or distributors involved with the handling | repair of the equipment | component – for any damage, caused by entries that are incomplete or incorrect on account of negligence, gross negligence or willful intent.

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

Declaration of Conformity CE
Following the EU Directive about electrical equipment for use within certain limits of voltage 73/23/EWG, appendix 3
and following the EC directive about electromagnetic compatibility 89/336/EWG, appendix 1

Company	Sartorius Stedim Biotech GmbH
Address	August-Spindler-Str. 11 D-37079 Göttingen Bundesrepublik Deutschland
<p>We herewith declare that the device described below fulfills the relevant fundamental safety requirements and health regulations specified by the appropriate EU-Directive, with respect to its design and construction and to the version as commercialized</p> <p>This declaration becomes legally invalid if modifications are performed on the device which have not been certified by Sartorius Stedim Biotech GmbH.</p>	
Designation of Device	Laboratory Vacuum Pump
Model, version	6 I
Cat.-No	16694-2-50-06 16694-1-60-06
Relevant directives of the EC	EU Directive about Electromagnetic Compatibility (89/336/EWG, appendix 3) EU Low Voltage Directive (73/23 EWG)
Applied harmonized standards	EN 61000-6-1; EN 61000-3-2/3 ; EN 61000-6-3; EN 55014; EN 61010-1
Applied national standards and technical specifications	
Date and Signature	 
Function of the Signatory	V. Niebel Senior Vice President Dr. Susanne Gerthausen Quality Management Representative

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

Declaration of Conformity CE
Following the EU Directive about electrical equipment for use within certain limits of voltage 73/23/EWG, appendix 3 and following the EC directive about electromagnetic compatibility 89/336/EWG, appendix 1

Company	Sartorius Stedim Biotech GmbH
Address	August-Spindler-Str. 11 D-37079 Göttingen Bundesrepublik Deutschland
<p>We herewith declare that the device described below fulfills the relevant fundamental safety requirements and health regulations specified by the appropriate EU-Directive, with respect to its design and construction and to the version as commercialized</p> <p>This declaration becomes legally invalid if modifications are performed on the device which have not been certified by Sartorius Stedim Biotech GmbH.</p>	
Designation of Device	Laboratory Vacuum Pump
Model, version	6 I
Cat.-No	16694-2-50-06 16694-1-60-06
Relevant directives of the EC	EU Directive about Electromagnetic Compatibility (89/336/EWG, appendix 3) EU Low Voltage Directive (73/23 EWG)
Applied harmonized standards	EN 61000-6-1; EN 61000-3-2/3 ; EN 61000-6-3; EN 55014; EN 61010-1
Applied national standards and technical specifications	
Date and Signature	 
Function of the Signatory	V. Niebel Dr. Susanne Gerthausen Senior Vice President Quality Management Representative

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
Company	Sartorius Stedim Biotech GmbH
Address	August-Spindler-Str. 11 D-37079 Göttingen Bundesrepublik Deutschland
<p>We herewith declare that the device described below fulfills the relevant fundamental safety requirements and health regulations specified by the appropriate EU-Directive, with respect to its design and construction and to the version as commercialized</p> <p>This declaration becomes legally invalid if modifications are performed on the device which have not been certified by Sartorius Stedim Biotech GmbH.</p>	
Designation of Device	Laboratory Vacuum Pump
Model, version	6 I
Cat.-No	16694-2-50-06 16694-1-60-06
Relevant directives of the EC	EU Directive about Electromagnetic Compatibility (89/336/EWG, appendix 3) EU Low Voltage Directive (73/23 EWG)
Applied harmonized standards	EN 61000-6-1; EN 61000-3-2/3 ; EN 61000-6-3; EN 55014; EN 61010-1
Applied national standards and technical specifications	
Date and Signature	 V. Niebel Dr. Susanne Gerthausen
Function of the Signatory	Senior Vice President Quality Management Representative

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Phone +49.551.308.0
Fax +49.551.308.3289

www.sartorius-stedim.com

Copyright by Sartorius Stedim Biotech GmbH, Goettingen, Germany. All rights reserved. No part of this publication may be reprinted or translated in any form or by any means without the prior written permission of Sartorius Stedim Biotech GmbH.

The status of the information, specifications and illustrations in this manual is indicated by the date given below. Sartorius Stedim Biotech GmbH reserves the right to make changes to the technology, features, specifications, and design of the equipment without notice.

Status:
June 2009,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany

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