
Reference Dryer

English

RD-8

Operating Instructions

The RD-8 was developed for the determining the remaining moisture of Milk powder.

Principle: An examination sample is dried under flowing, dried air. The mass loss of the examination sample is determined and listed in mass percentage.

This method is described in the international standard DIN/ISO 5537 or IDF 26.

The following parameters should be complied with precisely for drying milk powder in accordance with the standard above.

- | | | |
|-----------------|------------|------------|
| 1. Temperature: | 87°C | + 1°C |
| 2. Airflow: | 33 ml/min. | + 1ml/min. |
| 3. Drying time: | 5 h | + 0,2 h |

In addition, other powder-like samples can be examined or dried.

1. Set up and commissioning

1.1. Setup and electrical connection:

The device must be placed on a level, horizontal surface. Place the device on a level surface. It is connected to a grounding outlet (230V / 6A). The mains switch is located on the right front side.

1.2. Compressed air connection

For this, the compressed air hose that is included in the scope of delivery is plugged into the compressed air outlet located on the back. The connection is airtight immediately upon connection.

If the air hose is to be removed again, the small simultaneously pulled out. The other end of the compressed air hose is connected to the compressed air supply (either compressed air bottle or building supply). The air pressure can be between 2.5 bar and 7.5 bar.

2. Use of the RD-8

2.1 Temperature setting

The desired temperature (desired value) can be set by follows:

2.1.1. Press the "P" button
A desired temperature value appears on the display.

2.1.2. The displayed desired temperature can be reduced or increased with the "▼" or "▲" buttons.
The setting is automatically accepted after 2 seconds.

2.1.3. If no further adjustment takes place, the controller returns to normal display after approx. 30 seconds. This means that the respective actual temperature is displayed again. The last displayed desired temperature value is saved.

All settings are saved and are maintained even if the device is switched off or the network cable is unplugged.

You can find additional setting possibilities in the extensive operating instructions for the Jumo temperature controller.

2.2. Setting the flow amount

The flow amount can be set with rotary button (pressure regulator) located on the top of the device. For precise measurement of the flow amount, we recommend the use of a digital "flowmeter" such as the AD 1000.

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EG Declaration of Conformity in accordance with the EC-Directive for Machines 89/392 EWG, Appendix II A

We hereby declare that the design and construction of the machine described in the following and the model placed in circulation by us fulfills the basic safety and health requirements of the EC-Directive for Machines. The declaration becomes void if the machine is altered without consulting us.

Designation of the machine: **RD - 8**
Maschine typ: Drying unit for determining the remaining moisture of powdered substances, especially milk powder
Machine no.: 5700-
Applicable EC Directives: EG-Directive for Machines in the form 93/44/EWG
EC-Low Voltage Directive 73/23/EWG
EC-Directive for Elektromagnetic Compatibility 89/336/EWG

Applied harmonized standards, especially EN 292, Parts 1 and 2


Berlin, den
Funke-Dr.N.Gerber Labortechnik GmbH

Dipl.-Ing. K. Schäfer, Managing Director

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