**Instruction Manual**  
**Manuale di istruzioni**  
**Manuel d'instructions**  
**Manual de instrucciones**  
**Bedienungsanleitung**  
**操作指南**

**AREX 5 Heating Magnetic Stirrer**

F20500570, F20510570

**ARE 5 Heating Magnetic Stirrer**

F20500560, F20510560

**RC 5 Heating Plate**

F20700610, F20710610

**General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息**

Before using the unit, please read the following instruction manual carefully.

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Attention, surface chaude!

Prudencia, superficie caliente!

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.

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The product can be used with flammable liquids.

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**This unit must be used for indoor laboratory applications only.**  
The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a way not specified by the manufacturer or with non specified accessories, product's safety may be compromised.

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**IEC/EN 61010-1**  
**IEC/EN 61010-2-051**  
**IEC/EN 61010-2-010**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improve their quality.  
Nell'impegno di migliorare costantemente la qualità dei prodotti, VELP si riserva la facoltà di variarne le caratteristiche.  
Dans le but d'améliorer constamment la qualité de ses produits, VELP se réserve le droit d'apporter des modifications aux caractéristiques de ceux-ci.  
VELP se reserva el derecho de modificar las características de productos con el fin de mejorar constantemente su calidad.  
VELP behält sich zum Zwecke der ständigen Verbesserung der Produktqualität das Recht auf Änderung der Geräteeigenschaften vor.  
VELP 保留修改其产品特性的权利，以不断提高其质量。

## Safety Regulations / Norme di Sicurezza / Consignes de Sécurité / Advertencias de Seguridad / Sicherheitshinweise / 安全法规

The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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Use only the power cable provided with the instrument.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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The vessel must be made of a suitable material to withstand the foreseen temperature.

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Magnetic field doesn't interfere with the function of cardiac pacemakers or data media if they are farther than 20 cm to the instrument.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least) and between instruments of 10 cm (at least).

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
Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids or flammable liquids that have a fire point higher than 625°C and a flash point higher than 600°C.

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It is responsibility of the user to appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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The AREX / ARE 5 heating magnetic stirrer is used to heat and mix liquids inside a suitable container placed on the aluminium plate. The stirring is possible thanks to a magnetic drive stir bar placed into the vessel.

The RC 5 heating plate is used to heat liquids inside a suitable container placed on the aluminium plate.

The instrument is used in the laboratory for general use and for all those applications that require precise regulation of the heating plate temperature (and, for AREX / ARE 5, also of the stirring speed).

The Aluminum alloy heating plate ensures:

- Optimum heat distribution and a high specific power thanks to the circular configuration
- Temperature homogeneity
- High resistance to thermal stress and thermal shock

The AREX 5 heating plate has a special ceramic coating that ensures high resistance to chemical agents and surface abrasion. Using this heating plate at high temperatures may cause discoloring. This does not alter the thermal, mechanical and chemical resistance of the plate in any way.

For AREX / ARE 5, the powerful brushless motor is able to stir from 100 up to 1500 rpm. SpeedServo technology maintains constant speed as viscosity changes.

**Note:** *it's important to choose the most suitable magnetic stirring bar related to the type and the quantity of the liquid to be stirred as well as to the type of the beaker. Size and shape of the magnetic stir bar determines the stirring efficiency at any given speed. The stirring bar which satisfies most applications is code A00000356 (8 x 40 mm).*

The AREX 5 can work in combination with thermoregulator VTF for precise and accurate control of the sample temperature.

## 2. Installation

- Unpacking
  - Check the integrity of the unit after unpacking
- The box includes
  - AREX / ARE 5 heating magnetic stirrer or RC 5 heating plate
  - Instruction manual
  - Power cable
  - Magnetic stir bar Ø 8 x 40 mm (A00000356)
- Installation
  - Place the unit on a non-flammable surface
  - Make sure that the technical specification of the unit and of the power supply line are the same
  - Make sure that the electrical network is grounded
  - Make sure that the main switch is on position "0" (OFF)
  - Connect the instrument to an easily accessible socket (compliant with the current safety norms), using only the provided power cable

**Note:** *keep the power cable far away from the hot plate.*

- Place the container with the liquid and the stir bar (only for AREX / ARE 5) on the stirring plate

### 3. Operation

#### Power-up

- Turn the unit on using the main switch
- The right led between the two knobs lights up

#### Stirring (only for AREX / ARE 5)

- Rotate the right knob to set the stirring speed, the led over the speed knob lights up
- The real speed increases until setpoint achievement
- A microprocessor ensures constant speed even if the viscosity changes (SpeedServo)
- Turn the knob to the 0 positions to stop stirring

#### Heating

- Rotate the left knob to set the heating plate temperature
- If the set temperature is higher than the actual one, the led over the temperature knob lights up
- When the heating plate temperature exceeds 50 °C, the left led between the two knobs lights up until the temperature drops below 50 °C

**Note:** This warning is active even if the main switch is off but not if the instrument is not powered.

- Turn the knob to the 0 positions to stop heating

#### Operation with thermoregulator VTF (only for AREX 5)

- Turn off the instrument using the main switch
- Screw the threaded support rod into its seat on the back of the instrument and fasten the VTF thermoregulator onto the support rod. Place the temperature probe in the flask making sure that it is completely immersed into the sample.
- Plug the VTF into the dedicated socket on the back of the instrument
- Turn the unit on using the main switch. When the VTF is correctly connected, the display of the VTF lights up showing the temperature of the sample
- Select the temperature set point on the VTF
- Turn the temperature knob to the max value to start heating

**Note:** The AREX 5 always has primary control of the heating plate temperature. The temperature control function of the heating plate can also be used as a safety thermostat. In this case, the maximum temperature of the heating plate will not exceed the temperature setting on the magnetic stirrer, meaning that a longer heating time is required in order to reach the VTF thermoregulator temperature set, thus reducing temperature oscillation at the setpoint value.

### 4. Error and warning messages

When the control panel shows one of the following error messages, the instrument functions automatically stop.

Display indicator	Malfunction:
Flashing Heating led and Hot plate led ON	Heating plate overtemperature (T > 330°C)
Flashing Heating led and Hot plate led OFF	Heating plate safety probe overtemperature (T >370°C)
Flashing Stirring led (only for AREX / ARE 5)	The stirring system doesn't run correctly

To remove the error message, restart the instrument from the main switch. For the overtemperature of the safety probe is necessary to disconnect the instrument from the power supply.

If the alarm persists on the control panel, please contact VELP Scientifica's technical assistance service.

## 5. Maintenance

No routine or extraordinary maintenance is required except periodic cleaning.

### Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

### Repair

Repairs must be carried out by authorized Velp personnel only.

The transport of the instrument by freight forwarders, couriers or others must be carried out using the original shockproof packaging. Follow any instructions on the original packaging (e.g. palletizing).

It is the responsibility of the user to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

## 6. Technical data

		AREX 5		ARE 5		RC 5	
General features	Model	F20500570	F20510570	F20500560	F20510560	F20700610	F20710610
	Voltage	230 V 50/60 Hz	115 V 60 Hz	230 V 50/60 Hz	115 V 60 Hz	230 V 50/60 Hz	115 V 60 Hz
	Power input	630 W		630 W		600 W	
	Current consumption	2,8 A	5,5 A	2,8 A	5,5 A	2,6 A	5,2 A
	Dimensions (WxHxD)	160x85x270 mm (6,3x3,3x10,6 in)					
	Weight	1,7 kg (3,7 lb)				1,4kg (3,1 lb)	
	Construction material	Aluminum body –Technopolymer enclosure					
	Working in continuous	Admitted					
	Maximum load on the plate	25 kg					
	Noisiness	<< 80 dBa					
	Environmental temperature admitted	+5...+40 °C					
	Storage temperature admitted	-10...+60 °C					
	Level of electrical protection CEI EN60529	IP 42					
	Max humidity	80%					
	Overvoltage category	II					
	Pollution degree CEI EN61010-1	2					
	Max altitude	4000 m					
Heating	Heating plate power output	600 W					
	Heating plate dimensions	Ø 135 mm					
	Temperature range	50 ÷ 310 °C					
	Heat control type	Analogic					
	Heating plate construction material	CerAlTop™		Aluminum			
	Safety circuit	Separated with dedicated probe					
	Hot plate alarm	Over 50°C					
	Overtemperature alarm	Over 330°C					

<b>VTF</b>	Temperature range	-10 ÷ 300 °C (1 °C step)		
	Temperature resolution	0,2 °C		
	Accuracy	+/-0,5°C*		
	Electrical data	12Vdc – 1.2W (max)		
<b>Stir</b>	Stirring capacity	20 l H <sub>2</sub> O	15 l H <sub>2</sub> O	
	Speed range	100 ÷ 1500 rpm		
	Stir control type	Analogic		
	Motor type	BLDC		
	Motor rating output	10 W		

\* in the following conditions: 800 ml of water in 1 liter glass beaker (diameter 105 mm), stirring bar 8 x 40 mm, 600 rpm, 50°C.

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements
EN 61010-2-010	Particular requirements for laboratory equipment for the heating of material
EN 61010-2-051	Particular requirements for laboratory equipment for mixing and stirring
EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements

and satisfies the essential requirements of the following directives:

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2006/42/EC	Machinery directive
2014/30/EU	Electromagnetic compatibility directive
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (WEEE)	Waste of electric and electronic equipment

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung /  
加修改.

## 11. Declaration of conformity <sup>UK</sup>CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101      Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091      Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1      Safety requirements for electrical equipment for measurement, control, and laboratory use  
Part 1: General requirements

EN 61010-2-010    Particular requirements for laboratory equipment for the heating of material

EN 61010-2-051    Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1      Electrical equipment for measurement, control and laboratory use - EMC requirements -  
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597      Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032      Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic  
Equipment Regulations 2012

S.I. 2013/3113      Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

## Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO 14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortex mixers  
Dispersers  
COD Thermoreactors  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Overhead Shakers  
Turbidimeter  
Open Circulating Baths  
Pumps

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VELP Scientifica Srl