1.09204.0100  
1.09204.0103  
1.09204.0500  
1.09204.0503  
1.09204.1000  
1.09204.1022  
1.09204.2500  
1.09204.9025

## Microscopy

# Giemsa's azur eosin methylene blue solution

for microscopy

For professional use only

**IVD** In vitro Diagnostic Medical Device

CE 0123

### Intended purpose

This "Giemsa's azur eosin methylene blue solution - for microscopy" is used for human-medical cell diagnosis and serves the hematological, clinico-cytological, and histological investigation of sample material of human origin. It is a staining solution that when used together with other *in vitro* diagnostic products from our portfolio makes target structures evaluable for diagnostic purposes (by fixing, embedding, staining, counterstaining, mounting) in human-hematological, histological and clinico-cytological specimen materials, for example smears of whole blood and bone marrow, as well as paraffin sections.

Unstained structures are relatively low in contrast and are extremely difficult to distinguish under the light microscope. The images created using the staining solutions help the authorized and qualified investigator to better define the form and structure in such cases. Further tests must be carried out according to recognized, valid methods to reach a definitive diagnosis.

### Principle

When used in hematological applications, Giemsa's stain is frequently used in combination with other staining solutions, e.g. with May-Grünwald's solution for Pappenheim (MGG) overview staining. This staining solution generally stains the nuclei red, based on the molecular interaction between the eosin Y dye and an Azure B-DNA complex. Both dyes assemble to an Eosin Y - Azure B-DNA complex and the intensity of the resulting stain depends on the content of Azure B and the ratio of Azure B : Eosin Y. Furthermore, the resulting stain can vary depending on the influence of fixation, staining times, pH-value of the solutions or buffer substances. In histology and clinico-cytological applications, Giemsa's staining without additional stainings is used as an extended overview staining method. In this method, the color of the various cell components is influenced by pretreatment of the specimen material. Here, chromatin-containing structures (e.g. cell nuclei) appear in various blue shades, while the acidophilic components show up in a variety of red shades.

### Sample material

Sections of formalin fixed, paraffin embedded tissue (3 - 4 µm thick paraffin sections) or fresh, native whole blood and bone-marrow smears as well as clinical cytological material like urine sediment, sputum, smears from fine needle aspiration biopsies (FNAB), rinses, imprints are used as starting material.

### Reagents

Cat. No. 1.09204  
Giemsa's azur eosin methylene blue solution 100 ml, 500 ml,  
for microscopy 1 l, 2.5 l

### Also required:

Cat. No. 1.06009 Methanol for analysis EMSURE® 1 l, 2.5 l, 5 l  
ACS,ISO,Reag. Ph Eur  
Cat. No. 1.09468 Buffer tablets pH 7.2 100 tabs  
for preparing buffer solution acc. to WEISE  
for staining of blood smears

or

Cat. No. 1.11373 Buffer tablets pH 6.4 100 tabs  
for preparing buffer solution acc. to WEISE  
for the staining of blood smears

or

Cat. No. 1.11374 Buffer tablets pH 6.8 100 tabs  
for preparing buffer solution acc. to WEISE for  
the staining of blood smears

### for the staining of paraffin sections:

Cat. No. 1.00063 Acetic acid (glacial) 100% anhydrous 1 l, 2.5 l  
for analysis EMSURE®  
ACS,ISO,Reag. Ph Eur

Cat. No. 1.09634 2-Propanol for analysis EMSURE® 1 l, 2.5 l, 5 l  
ACS,ISO,Reag. Ph Eur

### for the Pappenheim staining:

Cat. No. 1.01424 May-Grünwald's eosine-methylene 100 ml, 500 ml,  
blue solution modified 1 l, 2.5 l  
for microscopy

### Sample preparation

The sampling must be performed by qualified personnel.

All samples must be treated using state-of-the-art technology.

All samples must be clearly labeled.

Suitable instruments must be used for taking samples and their preparation. Follow the manufacturer's instructions for application / use.

When using the corresponding auxiliary reagents, the corresponding instructions for use must be observed.

Deparaffinize and rehydrate paraffin sections in the conventional manner.

### Notes on Giemsa staining of paraffin sections

Always employ separate xylene or Neo-Clear™ (Cat. No. 1.09843) rinse baths when Giemsa staining paraffin sections as any ethanol traces in the solutions may result in the preparations being discolored.

### Pretreatment of bone marrow and iliac crest biopsy materials

Optimal results can be achieved using an OSTEOSOFT® mild decalcifier-solution (Cat. No. 1.01728).

To gently remove any calcification, the fixed biopsy materials are first placed in OSTEOSOFT® for 18 - 24 hours, after which they are transferred to histoprocessing. Blocks are carefully cut and, if required, are again treated with OSTEOSOFT® for an additional 20 minutes.

### Reagent preparation

#### Giemsa's azur eosin methylene blue solution

The solution is supplied as a concentrated staining solution and before use must be diluted with a buffer solution as described below. The diluted staining solution should be filtered before use.

#### Buffer solution

For preparation of approx. 1000 ml of solution, add and dissolve:

|                                                                                                                                           |          |
|-------------------------------------------------------------------------------------------------------------------------------------------|----------|
| Buffer tablet, Cat. No. 1.11373 (pH 6.4), Cat. No. 1.11374 (pH 6.8) or Cat. No. 1.09468 (pH 7.2) depending on the required reaction color | 1 tablet |
| Distilled water                                                                                                                           | 1000 ml  |

#### Dilute Giemsa's staining solution for manual staining

For preparation of approx. 200 ml solution mix:

|                                                              |        |
|--------------------------------------------------------------|--------|
| Giemsa's azur eosin methylene blue solution                  | 10 ml  |
| Buffer solution                                              | 190 ml |
| Mix well, leave to stand for 10 min, and filter if necessary |        |

#### Dilute Giemsa's staining solution for staining with automatic stainer

For preparation of approx. 300 ml solution mix:

|                                                              |        |
|--------------------------------------------------------------|--------|
| Giemsa's azur eosin methylene blue solution                  | 25 ml  |
| Buffer solution                                              | 275 ml |
| Mix well, leave to stand for 10 min, and filter if necessary |        |

In many cases precipitates of the dye form in the diluted staining solutions; these can be eliminated by repeating the filtration process.

#### Acetic acid 0.1 %, aqueous

For preparation of approx. 1000 ml solution mix:

|                   |         |
|-------------------|---------|
| Acetic acid 100 % | 1 ml    |
| Distilled water   | 1000 ml |

# Giemsa staining

## Procedure

### Air-dried smears

#### Staining in the staining cell / on the staining rack

The slides should be allowed to drip off well after the individual staining steps, as a measure to avoid any unnecessary cross-contamination of solutions.

The stated times should be adhered to in order to guarantee an optimal staining result.

|                                                            |        |
|------------------------------------------------------------|--------|
| Slide with air-dried smear                                 |        |
| Methanol                                                   | 3 min  |
| Dilute Giemsa's staining solution for manual staining      | 20 min |
| Buffer solution                                            | 1 min  |
| Buffer solution                                            | 1 min  |
| Air-dry (e.g. over night or at 50°C in the drying cabinet) |        |

#### Staining in the automatic stainer

The stated times should be adhered to in order to guarantee an optimal staining result.

|                                                                       | Time   | Station | Dip |
|-----------------------------------------------------------------------|--------|---------|-----|
| Slide with air-dried smear                                            |        |         |     |
| Methanol                                                              | 3 min  | 2       | on  |
| Dilute Giemsa's staining solution for staining with automatic stainer | 20 min | 3       | on  |
| Buffer solution                                                       | 1 min  | 4       | on  |
| Running tap water                                                     | 2 min  | 5       | on  |
| Dry                                                                   | 3 min  | 6       | -   |

All diluted solutions should be replaced after each working day.

Covering with non-aqueous mounting media (e.g. Neo-Mount™, DPX new, or Entellan™ new) and a cover glass is recommended for the storage of hematological specimens for several months.

After dehydration (ascending alcohol series) and clarification with xylene or Neo-Clear™, cytological samples can be mounted with water-free mounting agents (e.g. Entellan™ new, DPX new, or Neo-Mount™) and a cover glass and can then be stored.

The use of immersion oil is recommended for the analysis of stained slides with a microscopic magnification >40x.

## Result

|                          | Buffer solution pH 6.4 | Buffer solution pH 6.8 | Buffer solution pH 7.2 |
|--------------------------|------------------------|------------------------|------------------------|
| Nuclei resp. chromatin   | red to violet          | red to violet          | red to violet          |
| Cytoplasm of lymphocytes | blue                   | blue                   | blue                   |
| Cytoplasm of monocytes   | grey-blue              | grey-blue              | grey-blue              |
| Neutrophilic granules    | light violet           | light violet           | light violet           |
| Eosinophilic granules    | reddish to red-brown   | reddish to red-brown   | reddish to red-brown   |
| Basophilic granules      | dark violet            | dark violet            | dark violet            |
| Thrombocytes             | violet                 | violet                 | violet                 |
| Erythrocytes             | reddish                | reddish                | reddish-brownish       |

# Pappenheim staining

## with May-Grünwald's solution and Giemsa's solution

## Procedure

### Air-dried smears

#### Staining in the staining cell

The slides must be immersed and moved briefly in the solutions, simple immersion alone yields inadequate staining results.

The slides should be allowed to drip off well after the individual staining steps, as a measure to avoid any unnecessary cross-contamination of solutions.

The stated times should be adhered to in order to guarantee an optimal staining result.

|                                                            |        |
|------------------------------------------------------------|--------|
| Slide with air-dried smear                                 |        |
| May-Grünwald's eosine-methylene blue solution modified     | 3 min  |
| Dilute Giemsa's staining solution for manual staining      | 20 min |
| Buffer solution                                            | 1 min  |
| Buffer solution                                            | 1 min  |
| Air-dry (e.g. over night or at 50°C in the drying cabinet) |        |

#### Staining on the staining rack

The stated times should be adhered to in order to guarantee an optimal staining result.

|                                                            |      |                  |        |
|------------------------------------------------------------|------|------------------|--------|
| Slide with air-dried smear                                 |      |                  |        |
| May-Grünwald's eosine-methylene blue solution modified     |      | cover completely | 3 min  |
| Buffer solution                                            | 1 ml | mix              |        |
| Dilute Giemsa's staining solution for manual staining      |      | cover completely | 20 min |
| Buffer solution                                            |      | rinse            |        |
| Air-dry (e.g. over night or at 50°C in the drying cabinet) |      |                  |        |

All diluted solutions should be replaced after each working day. Only in the case of May-Grünwald's eosin methylene blue solution modified should the concentrated solution, when used daily, be renewed after one working week, at the latest, or otherwise when necessary. Concentrated May-Grünwald's eosin methylene blue solution modified (in the event of evaporation) may not be refilled, otherwise the concentration of the staining solution is no longer correct.

Covering with non-aqueous mounting media (e.g. Neo-Mount™, DPX new, or Entellan™ new) and a cover glass is recommended for the storage of hematological specimens for several months.

After dehydration (ascending alcohol series) and clarification with xylene or Neo-Clear™, cytological samples can be mounted with water-free mounting agents (e.g. Entellan™ new, DPX new, or Neo-Mount™) and a cover glass and can then be stored.

The use of immersion oil is recommended for the analysis of stained slides with a microscopic magnification >40x.

## Result

|                          | Buffer solution pH 6.4 | Buffer solution pH 6.8 | Buffer solution pH 7.2 |
|--------------------------|------------------------|------------------------|------------------------|
| Nuclei resp. chromatin   | red-violet             | purple to violet       | violet                 |
| Cytoplasm of lymphocytes | blue                   | blue                   | blue                   |
| Cytoplasm of monocytes   | grey-blue              | grey-blue              | grey-blue              |
| Neutrophilic granules    | light violet           | light violet           | violet                 |
| Eosinophilic granules    | brick-red              | brick-red              | red-brown              |
| Basophilic granules      | dark violet            | dark violet to black   | dark violet to black   |
| Thrombocytes             | violet                 | violet                 | violet                 |
| Erythrocytes             | reddish                | reddish                | reddish-grey           |

# Giemsa staining

## Procedure

### Paraffin sections of punched iliac crest specimens and detection of *Helicobacter pylori*

#### Staining in the staining cell

The slides should be allowed to drip off well after the individual staining steps, as a measure to avoid any unnecessary cross-contamination of solutions.

Deparaffinize histological slides in the conventional manner and rehydrate in a descending alcohol series.

The stated times should be adhered to in order to guarantee an optimal staining result.

Always employ separate xylene or Neo-Clear™ (Cat. No. 1.09843) rinse baths when Giemsa staining paraffin sections as any ethanol traces in the solutions may result in the preparations being discolored.

|                                                                                                                   |        |
|-------------------------------------------------------------------------------------------------------------------|--------|
| Slide with histological specimen                                                                                  |        |
| Distilled water                                                                                                   | 10 sec |
| Giemsa's azur eosin methylene blue solution (undiluted, filtered)                                                 | 15 min |
| Acetic acid 0.1 %                                                                                                 | 10 sec |
| Distilled water                                                                                                   | 10 sec |
| 2-Propanol                                                                                                        | 10 sec |
| 2-Propanol                                                                                                        | 10 sec |
| 2-Propanol                                                                                                        | 10 sec |
| Xylene or Neo-Clear™                                                                                              | 5 min  |
| Xylene or Neo-Clear™                                                                                              | 5 min  |
| Mount the Neo-Clear™-wet slides with Neo-Mount™ or the xylene-wet slides with e.g. Entellan™ new and cover glass. |        |

In the case of Giemsa's azur eosin methylene blue solution should the concentrated solution, when used daily, be renewed after one working week, at the latest, or otherwise when necessary. Concentrated Giemsa's azur eosin methylene blue solution (in the event of evaporation) may not be refilled, otherwise the concentration of the staining solution is no longer correct.

After dehydration (ascending alcohol series) and clarification with xylene or Neo-Clear™, histological slides can be covered with non-aqueous mounting agents (e.g. Entellan™ new, Neo-Mount™) and a cover glass and can then be stored.

The use of immersion oil is recommended for the analysis of stained slides with a microscopic magnification >40x.

## Result

|                                                                   |                   |
|-------------------------------------------------------------------|-------------------|
| Nuclei, cells                                                     | blue to dark blue |
| Collagen, osteoid                                                 | pale blue         |
| Eosinophilic granules                                             | red               |
| Acidic mucopolysaccharides, mastocytes granules, cartilage matrix | reddish-violet    |
| Acidophilic materials                                             | orange-red        |
| <i>Helicobacter pylori</i>                                        | blue to dark blue |

## Technical notes

The microscope used should meet the requirements of a medical diagnostic laboratory.

When using histoprocessors and automatic staining systems, please follow the instructions for use supplied by the supplier of the system and software. The diluted staining solution should be filtered before use.

Remove surplus immersion oil before filing.

## Analytical performance characteristics

"Giemsa's azur eosin methylene blue solution" stains and thereby visualizes biological structures, as described in the "Result" chapters of this IFU. The use of the product is only to be carried out by authorized and qualified persons, this includes, among other things, sample and reagent preparation, sample handling, histoprocessing, decisions regarding suitable controls and more.

The analytical performance of the product is confirmed by testing each production batch. The successful participation in international interlaboratory tests on a regular basis provide an additional and unaffiliated confirmation of analytical specificity and repeatability.

For the following stains, the analytical performance was confirmed in terms of specificity, sensitivity and repeatability of the product with a rate of 100 %:

|                             | Inter-assay Specificity | Inter-assay Sensitivity | Intra-assay Specificity | Intra-assay Sensitivity |
|-----------------------------|-------------------------|-------------------------|-------------------------|-------------------------|
| Hematological staining      |                         |                         |                         |                         |
| Erythrocytes                | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Nuclei                      | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Eosinophilic granules       | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Neutrophilic granules       | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Lymphocyte cytoplasm        | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Monocyte cytoplasm          | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Thrombocytes                | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Basophilic granules         | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Blood parasite nuclei       | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Histological staining       |                         |                         |                         |                         |
| Nuclei                      | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Cells                       | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Collagen                    | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Osteoid                     | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Eosinophilic granules       | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Acidic mucopoly-saccharides | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Mast cell granules          | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Cartilage matrix            | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| Acidophilic materials       | 20/20                   | 20/20                   | 14/14                   | 14/14                   |
| <i>Helicobacter pylori</i>  | 20/20                   | 20/20                   | 14/14                   | 14/14                   |

### Analytical performance results

Intra- (performed on the same batch) and inter-assay (performed on different batches) data list the number of correctly stained structures in relation to the number of performed assays.

## Clinical performance characteristics

Furthermore, the clinical performance of this product has been successfully proven in multiple scientific publications.

The diagnostic interpretation of the staining results, however, is to be carried out by qualified and authorized professionals, taking into account patient anamnesis, morphology, the use of adequate controls, and additional diagnostic tests, if appropriate. This method can be supplementarily used in human diagnostics.

The clinical performance of the use of the Giemsa's azur eosin methylene blue solution for the detection of *Helicobacter pylori*, for example, was compared with the immunohistochemical detection of *H. pylori*, since no "gold standard" for its diagnosis exists. The sensitivity and specificity of the Giemsa's azur eosin methylene blue solution for the detection of *H. pylori* in comparison to the immunohistochemical detection was determined as follows:

|             | Giemsa staining | Immunohistochemistry |
|-------------|-----------------|----------------------|
| Sensitivity | 13/15           | 15/15                |
| Specificity | 15/15           | 15/15                |

Sensitivity: 13 samples out of 15: 86.7 %

Specificity: 15 samples out of 15: 100 %

Positive predictive value (PPV): 100 %

Negative predictive value (NPV): 88.3 %

The results of this Performance Evaluation confirms that the product is suitable for the intended use and performs reliably.

## Diagnostics

Diagnoses are to be made only by authorized and qualified personnel. Valid nomenclatures must be used. This method can be supplementarily used in human diagnostics. Further tests must be selected and implemented according to recognized methods. Suitable controls should be conducted with each application in order to avoid an incorrect result.

## Storage

Store the Giemsa's azur eosin methylene blue solution - for microscopy at +15 °C to +25 °C.

## Shelf-life

The Giemsa's azur eosin methylene blue solution - for microscopy can be used until the stated expiry date.

After first opening of the bottle, the contents can be used up to the stated expiry date when stored at +15 °C to +25 °C.

The bottles must be kept tightly closed at all times.

## Capacity

3500 - 5000 stainings / 500 ml

## Additional instructions

### For professional use only.

In order to avoid errors, the application must be carried out by qualified personnel only. National guidelines for work safety and quality assurance must be followed. Microscopes equipped according to the standard must be used.

## Protection against infection

Effective measures must be taken to protect against infection in line with laboratory guidelines.

## Instructions for disposal

The package must be disposed of in accordance with the current disposal guidelines.

Used solutions and solutions that are past their shelf-life must be disposed of as special waste in accordance with local guidelines. Information on disposal can be obtained under the Quick Link "Hints for Disposal of Microscopy Products" at [www.microscopy-products.com](http://www.microscopy-products.com). Within the EU the currently applicable REGULATION (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006 applies.

## Auxiliary reagents

|                  |                                                                                          |                                                                         |
|------------------|------------------------------------------------------------------------------------------|-------------------------------------------------------------------------|
| Cat. No. 1.00063 | Acetic acid (glacial) 100% anhydrous for analysis EMSURE® ACS,ISO,Reag. Ph Eur           | 1 l, 2.5 l                                                              |
| Cat. No. 1.00496 | Formaldehyde solution 4%, buffered, pH 6.9 (approx. 10% Formalin solution) for histology | 350 ml and 700 ml (in bottle with wide neck), 5 l, 10 l, 10 l Titripac® |
| Cat. No. 1.00579 | DPX new non-aqueous mounting medium for microscopy                                       | 500 ml                                                                  |
| Cat. No. 1.00974 | Ethanol denatured with about 1 % methyl ethyl ketone for analysis EMSURE®                | 1 l, 2.5 l                                                              |
| Cat. No. 1.01424 | May-Grünwald's eosine-methylene blue solution modified for microscopy                    | 100 ml, 500 ml, 1 l, 2.5 l                                              |
| Cat. No. 1.01728 | OSTEOSOFT® mild decalcifier-solution for histology                                       | 1 l, 10 l Titripac®                                                     |
| Cat. No. 1.03699 | Immersion oil Type N acc. to ISO 8036 for microscopy                                     | 100-ml dropping bottle                                                  |
| Cat. No. 1.04699 | Immersion oil for microscopy                                                             | 100-ml dropping bottle, 100 ml, 500 ml                                  |
| Cat. No. 1.06009 | Methanol for analysis EMSURE® ACS,ISO,Reag. Ph Eur                                       | 1 l, 2.5 l, 5 l                                                         |
| Cat. No. 1.07961 | Entellan™ new rapid mounting medium for microscopy                                       | 100 ml, 500 ml, 1 l                                                     |
| Cat. No. 1.08298 | Xylene (isomeric mixture) for histology                                                  | 4 l                                                                     |
| Cat. No. 1.09016 | Neo-Mount™ anhydrous mounting medium for microscopy                                      | 100-ml dropping bottle, 500 ml                                          |

|                  |                                                                                                    |                                |
|------------------|----------------------------------------------------------------------------------------------------|--------------------------------|
| Cat. No. 1.09203 | Giemsa's azur-eosin-methylene blue for microscopy                                                  | 25 g, 100 g                    |
| Cat. No. 1.09468 | Buffer tablets pH 7.2 for preparing buffer solution acc. to WEISE for staining of blood smears     | 100 tabs                       |
| Cat. No. 1.09634 | 2-Propanol for analysis EMSURE® ACS, ISO, Reag. Ph Eur                                             | 1 l, 2.5 l, 5 l                |
| Cat. No. 1.09843 | Neo-Clear™ (xylene substitute) for microscopy                                                      | 5 l                            |
| Cat. No. 1.11373 | Buffer tablets pH 6.4 for preparing buffer solution acc. to WEISE for the staining of blood smears | 100 tabs                       |
| Cat. No. 1.11374 | Buffer tablets pH 6.8 for preparing buffer solution acc. to WEISE for the staining of blood smears | 100 tabs                       |
| Cat. No. 1.11609 | Histosec™ pastilles solidification point 56-58°C embedding agent for histology                     | 1 kg, 10 kg (4x 2.5 kg), 25 kg |
| Cat. No. 1.15161 | Histosec™ pastilles (without DMSO) solidification point 56-58°C embedding agent for histology      | 10 kg (4x 2.5 kg), 25 kg       |

## Hazard classification

Cat. No. 1.09204

Please observe the hazard classification printed on the label and the information given in the safety data sheet. The safety data sheet is available on the website and on request.

## Main components of the product

|                             |         |
|-----------------------------|---------|
| Cat. No. 1.09204            |         |
| C.I.52015 + Azure           | 4.1 g/l |
| C.I.45380                   | 2.4 g/l |
| contains CH <sub>3</sub> OH |         |
| 1 l = 0.99 kg               |         |

## Other IVD products

|                  |                                                                                                                      |                       |
|------------------|----------------------------------------------------------------------------------------------------------------------|-----------------------|
| Cat. No. 1.00869 | Entellan™ new for cover slipper for microscopy                                                                       | 500 ml                |
| Cat. No. 1.01383 | Wright's eosin methylene blue solution for microscopy                                                                | 100 ml, 500 ml, 2.5 l |
| Cat. No. 1.02439 | Eosin Y-solution 0.5%, alcoholic for microscopy                                                                      | 500 ml, 2.5 l         |
| Cat. No. 1.03999 | Formaldehyde solution min. 37% free from acid stabilized with about 10% methanol and calcium carbonate for histology | 1 l, 2.5 l, 25 l      |
| Cat. No. 1.05174 | Hematoxylin solution modified acc. to Gill III for microscopy                                                        | 500 ml, 1 l, 2.5 l    |
| Cat. No. 1.05175 | Hematoxylin solution modified acc. to Gill II for microscopy                                                         | 500 ml, 2.5 l         |
| Cat. No. 1.05387 | Leishman's eosin methylene blue solution modified for microscopy                                                     | 500 ml                |
| Cat. No. 1.09844 | Eosin Y-solution 0.5% aqueous for microscopy                                                                         | 1 l, 2.5 l            |
| Cat. No. 1.11661 | Hemacolor® Rapid staining of blood smear staining kit for microscopy                                                 | 1 set                 |
| Cat. No. 1.17081 | Eosin Y solution 1%, alcoholic for microscopy                                                                        | 1 l                   |

## General remark

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

## Literature

- Atlas der klinischen Hämatologie, Löffler, Rastetter, Haferlach, 2004, Springer Verlag 6. Auflage
- Theory and Practice of Histological Techniques, John d. Bancroft, Marilyn Gamble, 2008, Churchill Livingstone Elsevier, Sixth Edition
- Histological & Histochemical Methods, J. A. Kiernan, 1990, Pergamon Press, Second Edition
- Romeis - Mikroskopische Technik, Editors: Maria Mulisch, Ulrich Welsch, 2015, Springer Spektrum, 19. Auflage
- Sobotta, Lehrbuch Histologie, Welsch, 2006, Urban&Fischer, 2. Auflage
- Atlas of clinical cytology, Paul Lopez Cardozo, EM editon medizin
- Klinische Hämatologie, Herbert Begemann, 1975, Georg Thieme Verlag, 2. Auflage
- Conn's Biological Stains: A Handbook of Dyes, Stains and Fluorochromes for Use in Biology and Medicine, 10th Edition, (ed. Horobin, R.W. and Kiernan, J.A). Bios, 2002

H225: Highly flammable liquid and vapor.

H301 + H311 + H331: Toxic if swallowed, in contact with skin or if inhaled.

H317: May cause an allergic skin reaction.

H370: Causes damage to organs (Eyes, Central nervous system).

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P233: Keep container tightly closed.

P280: Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301 + P310: IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P304 + P340 + P311: IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.

Consult instructions  
for use

Manufacturer

Catalog number

Batch code

Caution, consult  
accompanying documents

Use by  
YYYY-MM-DD

Temperature  
limitation

Status: 2022-Oct-04

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly available resources.

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