E-Gel™ Sizing DNA Ladder

PRODUCT INFORMATION SHEET

Pub. No. MAN0017359

Rev. A.0



Contents

Catalog No. 10488100

Amount 100 applications



(1) Kit contents



Storage


Product is shipped at ambient temperature.

 Store at room temperature or at 4°C for up to 6 months, or at -20°C for long term storage.



Product description

- The Invitrogen[™] E-Gel[™] Sizing DNA Ladder is designed for sizing and quantification of double stranded DNA on 2% E-Gel[™] agarose gels.
- The E-Gel[™] Sizing DNA Ladder consists of 19 individual chromatography-purified DNA fragments ranging in size from 50 bp to 1,500 bp.
- Reference bands at 200 bp, 500 bp, and 1,000 bp are included for easy orientation.
- The ladder is supplied with 1X E-Gel[™] Sample Loading Buffer for sample DNA.





- Visit our product pages for additional information and protocols.
- Go online to view related DNA ladders and markers.
- For support, visit thermofisher.com/support.



Required materials

- E-Gel[™] SizeSelect[™] II or E-Gel[™] EX Agarose Gel (See Choosing the right DNA ladder for your E-Gel[™] agarose gel)
- TE Buffer (Cat. No. AM9858)
- Ultrapure[™] DNase/RNase-Free Distillated Water (Cat. No. 10977023)



Important guidelines

- Do not heat the E-Gel[™] Sizing DNA Ladder before loading.
- Load the same volume of DNA sample and DNA ladder.
- For quantification, adjust the concentration of the sample to equalize it approximately with the amount of DNA in the nearest band of the ladder.
- Dilute sample DNA in TE buffer to avoid degradation of DNA sample.
- Choosing the right DNA ladder for your E-Gel™ agarose gel
- Troubleshooting
- Limited product warranty and disclaimer details

Prepare DNA ladders and samples for electrophoresis

This protocol provides a brief description of how to use the DNA ladder with E-GelTM agarose gels. For detailed instructions on using specific types of E-GelTM agarose gels, go to thermofisher.com or contact Technical Support.

Step			Action
1		Prepare DNA ladder	 a. Thaw, mix and briefly centrifuge DNA ladder before use. b. Prepare DNA ladder. For E-Gel™ SizeSelect™ II Agarose Gels, mix and use 25 μL of the ladder without dilution. For E-Gel™ EX Agarose Gels, mix and use 20 μL of the ladder without dilution.
2		Prepare samples	 a. Dilute your sample 2- to 10-fold with TE Buffer (Cat. No. AM9858), 1X E-Gel™ Sample Loading Buffer (Cat No. 10482055), or water. b. Mix gently.
3		Load samples and DNA ladders	 a. Load DNA ladders and DNA samples into the appropriate wells of the E-Gel™ agarose gel. Add 25 µL for E-Gel™ SizeSelect™ II Agarose Gels. Note: All wells of E-Gel™ SizeSelect™ II Agarose Gels must be pre-filled with 50 µL of water. Add 20 µL for E-Gel™ EX Agarose Gels. b. Add water to any empty wells, so that all wells contain an equal volume of liquid.
4		Perform electrophoresis	a. Choose the appropriate E-Gel™ run protocol for your gel type based on the electrophoresis device being used. Gel type Program Recommended run time E-Gel™ Power Snap Electrophoresis Device (Cat. No. G8100) E-Gel™ SizeSelect™ II Agarose Gel (2%) E-Gel EX 1-2% 10 min (20 min max) E-Gel™ EX Agarose Gel (2%) b. Run the program to start electrophoresis.
5		Visualize agarose gel	Visualize DNA ladder and samples. • Use the E-Gel™ Power Snap Camera (Cat. No. G8200), E-Gel™ Imager (Cat. No. 466612), or other blue light imager to detect DNA bands stained with SYBR™ stains. • UV transilluminator to detect DNA bands stained with ethidium bromide.