**LANGE** 

DOC012.98.90429

# LOC 100

09/2012, Edition 1A

Ръководство за потребителя	bg
Uživatelská příručka	cs
Brugervejledning	da
Bedienungsanleitung	de
Εγχειρίδιο Χρήστη	el
User Manual	en
Manual del usuario	es
Manuel d'utilisation	fr
Korisnički priručnik	hr
Felhasználói kézikönyv	hu
Manuale utente	it
Basisgebruikershandleiding	nl
Instrukcja obsługi	pl
Manual do utilizador	pt
Manual al utilizatorului	ro
Руководство пользова	ru
Návod na použitie	sk
Uporabniški priročnik	sl
Bruksanvisning	sv
Kullanım Kılavuzu	tr

Български .....	3
Čeština .....	15
Dansk .....	27
Deutsch .....	39
Ελληνικά .....	51
English .....	63
Español .....	75
Français .....	87
Hrvatski .....	99
Magyar .....	111
Italiano .....	123
Nederlands .....	135
Polski .....	147
Português .....	159
Română .....	171
Русский .....	183
Slovenčina .....	195
Slovenščina .....	207
Svenska .....	219
Türkçe .....	231

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# Specifications

Subject to change without notice.

Performance specifications		
<b>Display</b>		128 × 64 pixels
<b>Input</b>		3 touch keys
<b>Power supply</b>		2 rechargeable AA Ni-MH batteries (included within the scope of delivery)
	<b>Operating hours</b>	1000 work cycles, approx. 2 months
	<b>Charging time required for full charge</b>	18 hours
<b>Interface</b>		Mini USB
<b>Real-time clock</b>	<b>Accuracy</b>	1 second/tag
<b>Size</b>		74 × 30 × 168 mm
<b>Weight</b>		150 g
<b>Housing rating</b>		IP 54
<b>Protection class</b>		Class III

## Safety information

Please read this entire manual before unpacking, setting up, or operating this equipment. Pay attention to all danger and warning statements. Failure to do so could result in serious injury to the operator or damage to the instrument.

To make sure that the protection provided by this instrument is not impaired, do not use or install this instrument in any manner other than that specified in these operating instructions.

### **DANGER**

Indicates a potentially or imminently hazardous situation that, if not avoided, results in death or serious injury.

### **WARNING**

Indicates a potentially or imminently hazardous situation that, if not avoided, may result in death or serious injury.

### **CAUTION**

Indicates a potentially hazardous situation that may result in minor or moderate injury.

### **NOTICE**

Indicates a situation that, if not avoided, could result in damage to the instrument. Information that requires special emphasis.

**Note:** Information that supplements points in the main text.

## Warning signs

Read all labels and tags attached to the instrument. Failure to do so may result in personal injury or damage to the instrument. If noted on the instrument, corresponding symbols will be included in the relevant safety information in the user manual.

	This symbol may be noted on the instrument and refers to operational and/or safety information in the user manual.
	Electrical equipment marked with this symbol may not be disposed of via European domestic or public disposal systems after 12 August 2005. In accordance with current regulations (EU Directive 2002/96/EC), users within the EU must now return old electrical equipment to the manufacturer for disposal. The user bears no costs for this process. <b>Note:</b> Please contact the manufacturer or supplier for instructions on how to return end-of-life equipment, manufacturer-supplied electrical accessories and all auxiliary items for proper disposal or recycling.

## RFID module

RFID technology is a radio application. Radio applications are subject to national conditions of authorization. The use of the LOC 100 is currently permitted in the following countries: EU, Norway, Switzerland, Iceland, Liechtenstein, Croatia, Turkey, USA, Canada, UAE, Australia and New Zealand. The manufacturer advises that the use of the LOC 100 outside of the above-mentioned regions may contravene national laws. The manufacturer reserves the right also to obtain authorization in other countries. In case of doubt, please contact the distributor.

The LOC 100 contains an RFID module to receive and transmit information and data. The RFID module operates with a frequency of 13.56 MHz.

## Certification (Integrated RFID module)

FCC ID: YCB-ZBA987

IC: 5879A-ZBA987

This device complies with Part 15 of the FCC rules and Industry Canada license exempt RSS standard(s).

Operation is subject to the following two conditions:

1. this device may not cause interference, and
2. this device must accept any interference, including interference that may cause undesired operation of the device."

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Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.

3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above

### **WARNING**

The instrument may not be used in dangerous environments.

The manufacturer and its suppliers reject any express or indirect guarantee for use with high-risk activities.

Follow the following safety information, in addition to any local guidelines in force.

Safety information for the correct use of the instrument:

- Do not operate the instrument in hospitals or comparable establishments in the vicinity of medical equipment, such as pace makers or hearing aids.
- Do not operate the instrument in the vicinity of highly flammable substances, such as fuels, highly flammable chemicals and explosives.
- Do not operate the instrument near combustible gases, vapors or dust.
- Do not vibrate or jolt the instrument.
- The instrument can cause interference in immediate proximity to televisions, radios and computers.
- Do not open the instrument.

Improper use of the instrument voids the warrantee.

## Batteries

### DANGER

Risk of fire and explosion.

Only use the NiMH rechargeable batteries (Type AA, 1.2 V / min. 2450 mAh) provided in the measuring instrument and make sure that the batteries are correctly inserted into the battery compartment. Incorrect insertion of the batteries could cause damage to the instrument, fire or explosions.

Do not use nonrechargeable batteries under any circumstances. The use of nonrechargeable batteries could cause damage to the instrument, fire or explosions.

## Scope of delivery

The LOC 100 scope of delivery comprises the following items:

- LOC 100 RFID read/write device
- USB cable
- Two AA NiMH batteries
- LOC 100 user manual

**Note:** *If any of these items are missing or damaged, please contact the manufacturer or the responsible sales representative immediately.*

## Product description

The LOC 100 RFID read/write device can be used to write the operator ID and location ID onto a sample bottle label together with the time and date.

Every sample bottle then clearly has the following information assigned to it:

- Location — where was the sample taken?
- Time — when was the sample taken?
- Operator — who took the sample?

This information (sample ID) can then be read by an RFID-compatible photometer so that it can be saved together with the test results. This allows the sample location, sample time and the operator to be documented consistently.


The RFID starter set for sample identification comprises:

- The LOC 100 RFID read/write device
- Two operator RFID tags
- Five location RFID tags
- Five sets of sample bottle labels with sample RFID tags, in red (3×), blue (3×), yellow (3×), green (3×) and black (3×).

All RFID tags can be written and read out as many times as required using the LOC 100 or the RFID module in the photometer.


The following describes the writing and reading process of the individual components in detail.

**Figure 1 Front view of LOC 100**

1. Left key (Read)	4. Mini-USB port
2. On/Off key	5. Display
3. Right key (Write)	

Figure 2 Rear view of LOC 100

1. Screw

2. Battery compartment cover

3. RFID module

4. AA rechargeable battery (2x)

# Commissioning

## Prepare the LOC 100 RFID read/write device

### DANGER


Risk of fire and explosion.

Only use the NiMH rechargeable batteries (type AA, 1.2 V / min. 2450 mAh) provided in the measuring instrument and make sure that the batteries are correctly inserted into the battery compartment. Incorrect insertion of the batteries could cause damage to the instrument, fire or explosions.

Do not use nonrechargeable batteries under any circumstances. The use of nonrechargeable batteries could cause damage to the instrument, fire or explosions.

1. Open the battery compartment by releasing the four screws on the back side.
2. Insert two type AA batteries into the battery compartment; note the poles.
3. Close the battery compartment and secure it with the four screws.
4. Switch on the LOC 100 using the On/Off key.
  - a. To charge the batteries, connect the LOC 100 to the photometer via a USB cable. Alternatively, the batteries can be charged via a USB cable using a standard PC or laptop.
  - b. It is also possible to charge the batteries in a standard external charger for NiMH rechargeable batteries.
  - c. The battery charge is shown via a status display on the LOC 100 screen.

The nominal capacity of the batteries is not reached until repeated charge/discharge cycles have been completed.
5. Press and hold down the On/Off key to open the Configuration menu. Use the three keys to adjust the brightness, contrast, date, time, tone and reader of the operator ID.

**Note:** If the LOC 100 is connected to the photometer via a USB cable, the rechargeable batteries are charged and the date and time are transferred from the photometer to the LOC 100.

**Note:** All settings can also be made from the PC. To do so, connect the LOC 100 to the PC using a USB cable. The LOC 100 is displayed as an RFID mass storage device; the settings for brightness, contrast, date, time and tone can be changed in the *RFID.exe* file.

## Prepare sample bottles

Figure 3 Sample ID label on a sample bottle

Glue the sample ID label onto the sample bottle so that the lower edge of the label is approximately 2 cm above the bottom of the bottle.

## Initialize operator and location RFID tag on the photometer

1. Specify the operator ID (refer to photometer operating instructions).

2. Select **Initialize RFID Tag**.


3. Follow the onscreen instructions and hold the operator RFID tag in front of the photometer RFID module.

A sound confirms that the writing process has been successful.


4. The operator RFID tag can be rewritten at any time. To do so, confirm the "Overwrite Operator ID" message with **OK** and enter the password, if the old operator ID is protected with a password.
5. Specify a sample ID (refer to photometer operating instructions).

6. Press **Initialize RFID Tag**.

7. Follow the onscreen instructions and hold the location RFID tag in front of the photometer RFID module.

A sound confirms that the writing process has been successful.

8. The location RFID tag can be rewritten at any time. To do so, confirm the "Overwrite location RFID tag" message with **OK**.

## Read data and write on label

1. Switch on the LOC 100.
2. Read in the operator ID initialized on the photometer by holding the LOC 100 above the operator RFID tag and briefly pressing the left key. The LOC 100 attempts to read the tag for five seconds. Once the tag has been read successfully, the operator ID is shown on the screen. (This step is optional. The option to read out the operator ID can be deactivated in the Configuration menu).

3. Read in the sample ID initialized on the photometer by holding the LOC 100 above the location RFID tag and briefly pressing the left key. The LOC 100 attempts to read the tag for five seconds. Once the tag has been read successfully, the location ID is shown in the display.

4. Write on the empty sample RFID ID tag by holding the LOC 100 over the label and pressing the right key. A successful transfer is shown in the display.


**Note:** This method can be used to write sample RFID tags as many times as required. Sample RFID tags are dishwasher-proof. However, repeated cleaning, particularly at high temperatures, can impair the function of the tag. Read out the RFID tag using the photometer (*Read out the sample RFID tag on the photometer*) to check the sample RFID tag.

**Note:** One long tone indicates that the tag has been read or written successfully. An error, e. g. because the distance from the tag is too great, is indicated by five short tones and an error message.

**Figure 4** Read and write data using the LOC 100

- |    |  |
|----|--|
| 1. | Left key (Read)                        |
| 2. | Right key (Write)                      |
| 3. | Operator RFID tag or location RFID tag |

## Read out the sample RFID tag on the photometer

1. Hold the sample bottle label in front of the photometer RFID module.  
The sample ID is read in and displayed as a button on the toolbar in the photometer
2. To determine the sample parameters, prepare the sample in accordance with the process specification for the relevant tests. Insert the prepared cell into the photometer cell compartment.
3. This action will start the measurement process and the results will be displayed.
4. Press the Sample ID key in the toolbar.  
The read-in sample ID is assigned to the measurement.
5. Remove the cell.  
The measurement result is saved with the sample ID in the photometer.

**Note:** *The sample ID will be shown as a key in the toolbar until the photometer is switched off.*

## Error messages

Displayed errors	Error description	Troubleshooting
E-01	RFID communication error: Tag sends no answer. Distance too great or tag defective.	Position the RFID tag closer to the instrument or replace the tag.
E-02	RFID collision error: Several tags within range.	Remove all tags from the reading area, except the one that needs to be read.
E-03	Unknown or incorrect tag: Attempts to write on a tag intended for reading or vice versa. Reads in the operator ID, even though this function was deactivated.	Check whether the correct tag is being used for the procedure.
E-04	Data format error: The data structure on the tag is incorrect.	Tag is damaged and must be replaced.

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80539-0389 U.S.A.
(U.S.A. only)