
LANGE 

DOC012.98.90328

SIP 10

03/2012, Edition 3A

Ръководство за потребителя	bg
Uživatelská příručka	cs
Brugervejledning	da
Bedienungsanleitung	de
Εγχειρίδιο Χρήστη	el
User Manual	en
Manual del usuario	es
Manuel d'utilisation	fr
Korisnički priručnik	hr
Felhasználói kézikönyv	hu
Manuale utente	it
Basisgebruikershandleiding	nl
Instrukcja obsługi	pl
Manual do utilizador	pt
Manual al utilizatorului	ro
Návod na použitie	sk
Uporabniški priročnik	sl
Bruksanvisning	sv
Kullanım Kılavuzu	tr

Български	3
Čeština	19
Dansk	33
Deutsch	47
Ελληνικά	61
English	77
Español	91
Français	105
Hrvatski	119
Magyar	133
Italiano	147
Nederlands	161
Polski	175
Português	189
Română	203
Slovenský jazyk	217
Slovenščina	231
Svenska	245
Türkçe	259

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Technical data

Subject to changes.

Performance specifications		
Cell		Dual Path Length Pour-Thru Cell, 1 inch/1 cm, plastic
		Pour-Thru Cell, 1 cm, quartz glass
Wavelength range	Plastic	340–900 nm
	Quartz glass	190–900 nm
Rinsing volume		Minimum of 25 mL
Flow speed		1 mL/second
Interface		USB type A
Power supply		Via USB cable, 530 mA, 5 V
Ambient storage requirements		–10–60 °C (14–140 °F), 85 % relative humidity (without condensate formation)
Ambient operating requirements		10–40 °C (50–104 °F), 80 % relative humidity (without condensate formation)
Dimensions		120 × 85 × 200 mm
Weight		0.5 kg
Enclosure rating		IP30

General information

Safety notes

Please read the entire manual before the equipment is unpacked, set up or put into operation. Observe all hazard and warning notes. Failure to do so could lead to serious injury of the operator or to damage to the device.

To make sure that the device safety features are not impaired, this device must not be used or installed in any manner other than that specified in this user manual.

Hazard notes in this manual

DANGER

Indicates a potentially or imminently dangerous situation that, if it is not avoided, leads to death or to serious injuries.

WARNING

Indicates a potentially or imminently dangerous situation that, if it is not avoided, can lead to death or to serious injuries.

CAUTION

Indicates a possible dangerous situation that can have minor or moderate injuries as a result.

NOTICE

Indicates a situation that, if not avoided, could result in damage to the device. Information that particularly should be emphasized.

Note: Information that supplements aspects from the main text.

Warning labels

Observe all marks and labels that are attached to the device. Failure to do so can result in personal injury or damage to the device.

	This symbol is a warning triangle. Follow all safety notes that follow this symbol to prevent possible injuries. If this symbol is attached to the device, it refers to information in the operation and/or safety sections of the user manual.
	This symbol can be attached to a housing or a barrier in the product and shows that electric shock risk and/or the risk of a death through electric shock exists.
	Electrical equipment marked with this symbol may as of August 12, 2005 Europe-wide no longer be disposed of in unsorted house or industrial waste. According to valid provisions (EU Directive 2002/96/EC), from this point consumers in the EU must return old electrical devices to the manufacturer for disposal. This is free for the consumer. Note: <i>Instructions on the correct disposal of all (marked and unmarked) electrical products supplied or manufactured by HACH LANGE can be obtained from your local HACH LANGE sales office.</i>

Chemical and biological safety

WARNING

Chemical hazards. Chemicals can cause serious injuries if handled without the relevant protection.

Always wear personal protective equipment to detach and replace tubes in accordance with the safety data sheets applicable to these chemicals.

CAUTION

Chemical hazards. Do not use the Pour-Thru Cell in tests that require the use of organic solvents such as alcohols, toluene, chloroform, trichloroethane or cyclohexanone. It is possible that these solvents may attack the plastic components of the Pour-Thru Cell and therefore cause damage to the device and chemical exposure for the user.

With normal operation of the device, the use of chemicals that pose a risk health or biologically harmful samples can be required.

- Before handling these substances, read all danger notes and safety information printed on the containers of the original solution and in the safety data sheet.
- Dispose of all consumed solutions in accordance with the national regulations and laws.
- Select the type of protective equipment according to the concentration and quantity of the dangerous substance at the respective work place.
- Always purge the cell and tubes thoroughly after use with deionized water.

Product overview

The SIP 10 Sipper module is an optional accessory for the DR 3900 and DR 6000 spectrophotometer. The measurement accuracy is improved by working with the SIP 10. The same optical characteristics apply to both the zero compensation and the measurements. All test solutions flow through a cell and errors due to optical differences between different cells are eliminated. A peristaltic pump pumps a constant amount of fluid through the Pour-Thru Cell.

DR 3900 product contents

Check that the order is complete. If something is missing or damaged, contact the manufacturer or distribution partner immediately.

- SIP 10 Sipper module, pump tubing pre-assembled
- Dual Path Length Pour-Thru Cell, 1 inch/1 cm
- USB connection cable
- Tygon® tubing, 1.70 m (5.6 ft)
- Stainless steel tube adapter
- User manual, multilingual

DR 6000 product contents

Check that the order is complete. If something is missing or damaged, contact the manufacturer or distribution partner immediately.

- SIP 10 Sipper module, pump tubing pre-assembled
- Pour-Thru Cell according to order variant:
 - Dual Path Length Pour-Thru Cell, 1 inch/1 cm, plastic
 - Pour-Thru Cell, 1 cm, quartz glass
- USB connection cable
- Tube set consisting of
 - 2× Tygon tubes with fitting
 - Tygon tube
 - 2× Pharmed® tubes with fitting
 - Rubber tube guide
- Stainless steel tube adapter
- Cell compartment cover
- Tray for storage
- User manual, multilingual

Installation

NOTICE

Make sure that tubes are not bent.
Check the fill level of the waste vessel regularly.

NOTICE

There may be residual amounts of fluid in the tubes. Always wear suitable protective clothing and goggles when changing the tubes or when performing maintenance work in order to prevent contamination of the skin or eyes.


NOTICE

Connect the device directly to the photometer via the USB cable.
Do not use a USB hub connection.

DR 3900 installation

1. Position the SIP 10 to the left of the photometer.
2. Clean the Pour-Thru Cell with a lint-free cloth.
3. Place the Pour-Thru Cell into the cell compartment of the photometer according to the required path length.
4. Cut a piece of tubing no longer than 40 cm (1.3 ft) for the feed tube.
5. Connect one end of the feed tube to the IN connection of the Pour-Thru Cell.
6. Connect the other end of the feed tube to the stainless steel tube adapter and place this in the zero-solution/sample vessel.
7. Cut a piece of tubing approx. 30 cm (1 ft) long for the connection tube.
8. Connect the connection tube to the OUT connection of the Pour-Thru Cell and the input connection of the peristaltic pump.
9. Use the remaining tubing as a discharge tube, and connect one end to the output connection of the peristaltic pump.
10. Place the other end of the discharge tube in a suitable waste vessel.
11. Connect the SIP 10 USB cable to a USB port on the photometer.

Figure 1 SIP 10 and DR 3900 fully installed

1	SIP 10 Sipper module	6	Pour-Thru Cell
2	Discharge tube	7	Feed tube
3	Peristaltic pump	8	Photometer
4	Connection tube	9	Stainless steel tube adapter
5	USB cable		

DR 6000 installation

1. Position the SIP 10 on the storage tray to the left of the photometer.
2. Open the cell compartment cover up to the stop.
3. Push the rubber tube guide into the groove in the top left of the cell compartment.
4. Push the two Pharmed tubes with the free ends onto both connections of the Pour-Thru Cell.
5. Clean the Pour-Thru Cell with a lint-free cloth.
6. Place the Pour-Thru Cell into the cell compartment of the photometer according to the required path length.
7. Push the Pharmed tube with the connection to the IN-connection of the Pour-Thru Cell into the bottom guide in the tube guide.
8. Push the Pharmed tube with the connection to the OUT-connection of the Pour-Thru Cell into the top guide in the tube guide.

Note: Make sure that the tubes are sitting securely in the guide so that they do not become squashed.

The fitting ends of two Pharmed tubes are now protruding from the side of the photometer.

9. Install the cell compartment cover, starting at the bottom right corner. Close the cover.
10. Connect the fitting end of the bottom Pharmed tube to the fitting side of a Tygon tube.
This is the feed tube.
11. Connect the other end of the feed tube to the stainless steel tube adapter and place this in the zero-solution/sample vessel.
12. Connect the fitting end of the top Pharmed tube to the fitting side of the second Tygon tube.
This is the connection tube.
13. Connect the other end of this connection tube to the input connection of the peristaltic pump.
14. Connect one end of the Tygon tube (without fitting) to the output connection of the peristaltic pump.


This is the discharge tube.

15. Place the other end of the discharge tube in a suitable waste vessel.

16. Connect the SIP 10 USB cable to a USB port on the photometer.

Note: Shorten the feed tube and connection tube as much as possible in order to optimize the pump volume of the system.

Figure 2 SIP 10 and DR 6000 fully installed


1	SIP 10 Sipper module on the tray for storage	8	Rubber tube guide
2	Discharge tube: Tygon tube (without fitting)	9	Pour-Thru Cell
3	Peristaltic pump	10	Feed tube: Pharmed tube with fitting
4	USB cable	11	Feed tube: Tygon tube with fitting
5	Connection tube: Tygon tube with fitting	12	Photometer
6	Connection tube: Pharmed tube with fitting	13	Stainless steel tube adapter
7	Cell compartment cover		

Operation

SIP 10 setup

When a SIP 10 module is connected to the photometer, an additional **SIPPER** button is shown in the toolbar.

1. Press **SIPPER** on the toolbar.
A Sipper options menu opens.
2. Press **ON** to modify the settings for a Sipper cycle.

- **SIP TIME:** The sip time determines the sample volume that is pumped through the cell.
In order to make sure that all of the liquid is replaced in the Dual Path Length Pour-Thru Cell, 1 inch/1 cm, the minimum time is 25 seconds.
Default setting: 25 seconds The flow rate is 1 mL/sec.

Note: If you use other Pour-Thru Cells, the SIP time must be determined individually. One method to determine this is to use two solutions with different concentrations and measure the absorbance continuously.

- **SETTLE TIME:** The settle time defines the duration of the standstill phase between the pump process and measurement process. During this time, air bubbles that have formed during the pump process can escape and sample turbulence can cease.
Default setting: 5 seconds
- **PURGE TIME:** The purge time determines the volume of rinse agent that is pumped through the cell after each measurement. This additional purging can be deactivated.
Default setting: off
- **PURGE START:** The start of the purge time can be triggered automatically after the measurement process, or manually. To start the cycle manually, press **PURGE** to trigger the purge cycle.
Default setting: auto

Note: The **PURGE START** button is only active when the purge time setting is ON.

3. Configure the required settings and press **OK** to confirm.

Execution of a stored program


1. Access the required test via **STORED PROGRAMS**.
2. Place the feed tube in the zero solution and press **ZERO**.

The zero solution is pumped in according to the selected SIP time. Once the settle time has passed, the zero measurement is taken.

Note: The remaining time is shown on the display. Press **CANCEL** to stop the operation.

Note: If the **PURGE START: AUTO** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

The zero measurement is shown on the display.

- Place the feed tube in the sample solution and press **READ**.
The sample solution is pumped in over the selected SIP time. Once the settle time has passed, the measurement is taken.

Note: The remaining time is shown on the display. Press **CANCEL** to stop the operation.

Note: If the **PURGE START: AUTO** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

Note: Either deionized water or the next sample can be used as a rinse agent.

The measurement result is shown on the display.

For the tests listed in [Table 1](#), an additional purge cycle of the Pour-Thru Cell with deionized water is required between samples.

Table 1 Tests for which purging with deionized water is required

Aluminum, Aluminon	Chlorine dioxide, LR	Cobalt, PAN
Copper, Porphyrin	Hardness, Calmagite	Manganese, LR, PAN
Nickel, PAN	Nitrate, MR	Nitrate, HR

The Dual Path Length Pour-Thru Cell, 1 inch/1 cm can also be used for the Nessler Method for nitrogen and TKN if cleaned correctly. For this, pour a few sodium thiosulfate crystals into the cell to clean it. Rinse out the crystals with deionized water.

The use of the Pour-Thru Cell is not possible for the tests listed in [Table 2](#), since chemical problems or other complications can occur. Use the sample cell specified in the procedural instructions for these tests.

Table 2 Tests for which the Pour-Thru Cell is not suitable

Aluminum ECR	Arsenic	Barium	Boron, Carmine
Cyanuric acid	Fluoride	Formaldehyde	Lead, LeadTrak
Mercury	Nickel, Heptoxime	Nitrite, HR	PCB
Phenols	Potassium	Selenium	Silver
Suspended solids	Sulfate	TPH	Volatile acids
Zinc	Surfactants, Anionic (Detergents)		

The Pour-Thru Cell can also be used for other tests. However, please note that the Pour-Thru Cell requires a minimum sample volume of 25 mL to make sure that the previous solution is fully rinsed out of the cell. Increase the reagent and sample volumes accordingly.

Execution of a single or multi-wavelength measurement


1. Press **SINGLE WAVELENGTH** or **MULTI-WAVELENGTH**.
2. Press **OPTIONS** to adjust the settings within single/multi-wavelength mode.
3. Place the feed tube in the zero solution and press **ZERO**.

The zero solution is pumped in according to the selected SIP time. Once the settle time has passed, the zero measurement is taken.

Note: The remaining time is shown on the display. Press **CANCEL** to stop the operation.

Note: If the **PURGE START: AUTO** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

The zero measurement is shown on the display.

4. Place the feed tube in the sample solution and press **READ**.
The sample solution is pumped in over the selected SIP time. Once the settle time has passed, the measurement is taken.
Note: The remaining time is shown on the display. Press **CANCEL** to stop the operation.
Note: If the **PURGE START: AUTO** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.
Note: Either deionized water or the next sample can be used as a rinse agent.

The measurement result is shown on the display.


Execution of a time course measurement

1. Press **TIME COURSE**.
2. Press **OPTIONS** to adjust the settings within time course mode.
3. Place the feed tube in the zero solution and press **ZERO**.

The zero solution is pumped in according to the selected SIP time. Once the settle time has passed, the zero measurement is taken.

Note: If the **PURGE START: AUTO** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

The zero measurement is shown on the display.

4. Place the feed tube in the sample solution and press **READ**.
The sample solution is pumped in over the selected SIP time. Once the settle time has passed, the measurement is taken.
The measurement is taken in accordance with the time course settings. The **ZERO/READ** buttons change to **MARK/STOP** during the time course. To abort the measurement process, press **STOP**.

Once the entire time has passed, the measurement is shown as a curve on the display.

Note: If the **PURGE START: AUTO** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

Maintenance

Cleaning

⚠ CAUTION

Personal injury hazard. Only qualified experts may perform the tasks described in this section of the manual.

NOTICE

There may be residual fluid in the tubes. Always wear suitable protective clothing and goggles when changing the tubes or when performing maintenance work, in order to prevent contamination of the skin or eyes.

SIP 10

Clean the SIP 10 with a damp cloth. Do not use any type of solvent (for example, acetone).

Tray for storage (for DR 6000 only)

If you no longer require the SIP 10 for current measurements, you can store all of the accessories in the tray:

- Cell compartment cover and USB connection cable to the right of the SIP 10
- Fitted tubes and rubber tube guide behind the SIP 10
- Pour-Thru Cell in the SIP 10

Pour-Thru Cell

Purge the Pour-Thru Cell before and after every series of tests with deionized water.

If the Pour-Thru Cell is heavily contaminated, repeat the purge cycle several times.

Before the Pour-Thru Cell is put into storage, perform the purge cycle several times with air to remove any residual fluid. Hold the Pour-Thru Cell with the OUT connection slanted downwards during this procedure.

Check the Pour-Thru Cell viewing windows. If they are dirty or fogged up, place the Pour-Thru Cell in a soap solution or in dilute acid. Then rinse the Pour-Thru Cell thoroughly with deionized water.

Tube cleaning

⚠ CAUTION

Chemical exposure hazard. Familiarise yourself with the necessary safety procedures and the correct handling of the chemicals before the work and read and follow all relevant safety data sheets.

Clean the tubes with deionized water after each series of measurements.

The tubes are exposed to chemicals and must be replaced at least every 12 months.


Replace the pump tubing

The pump tubing is exposed to mechanical loads and chemicals and must be replaced at least every 12 months.

1. Perform the purge cycle several times with air to remove any residual fluid from the system.
2. Disconnect the USB cable from the photometer.
3. Remove the connection tube and discharge tube from the input and output connections on the peristaltic pump.
4. Tilt back the lever behind the peristaltic pump (Figure 3 step 1)
The peristaltic pump cover will move back.
5. Lift up the peristaltic pump cover (Figure 3 step 2) and dispose of the cover, pump tubing and connections.
6. Position the new cover with pre-assembled pump tubing and connections on the peristaltic pump.
7. Push down the cover on the peristaltic pump and tilt up the lever behind the peristaltic pump (Figure 3 steps 3 and 4).
The peristaltic pump cover will move forward.
8. Connect the connection tube and the discharge tube to the input and output connections on the peristaltic pump (refer to Figure 1 and Figure 2).

9. Connect the SIP 10 USB cable to a USB port on the photometer.

Figure 3 Change pump tubing

Troubleshooting

Error messages

Potential SIP 10 errors are displayed by the photometer.

Table 3 Error messages

Error displayed	Cause	Resolution
Sipper module not connected. Please check cable.	Connection between SIP 10 and photometer is disrupted.	Check the USB cable. The cable length must not exceed 1 m (3.3 ft) and the cable must be directly connected to the photometer. Remove other devices that are connected where necessary.
Please check sipper and tube.	Pump tubing is not correctly inserted.	Loosen the cover and position again. Where required, run the pump briefly before the lever is moved to make sure that the tube is correctly positioned around the rollers.

Replacement parts

Designation	Order no.
SIP 10 Sipper module set for DR 3900 complete with tube set and 1 inch/1 cm Pour-Thru Cell, EU	LQV157.99.10001
SIP 10 Sipper module set for DR 6000 complete with tray, tube set and 1 inch/1 cm Pour-Thru Cell, EU	LQV157.99.20001
SIP 10 Sipper module set DR 6000 complete with tray, tube set and 1 cm quartz glass Pour-Thru Cell, EU	LQV157.99.30001
SIP 10 Sipper module; including pump tubing, EU	LQV157.99.00001
USB cable, type AB 1 m (3.3 ft)	LZQ104
Dual Path Length Pour-Thru Cell, plastic, 1 inch/1 cm	LZV876
UV Pour-Thru Cell, quartz glass, 1 cm	LZV510
Pump tubing (Lagoprene®) pre-assembled with cover and connections	LZV877
Complete tube set for DR 3900, including Tygon tubing, 1.70 m (5.6 ft), pump tubing (Lagoprene pre-assembled with cover and connections), stainless steel tube adapter	LZV875
Complete tube set for DR 6000 and Dual Path Length Pour-Thru Cell, including Tygon tubing and Pharmed tubing, pump tube (Lagoprene pre-installed with cover and connections), stainless steel tube adapter	LZQ102
Complete tube set for DR 6000 "drinking water application," including Pharmed tubing, pump tube (Lagoprene pre-installed with cover and connections), stainless steel tube adapter	LZQ100
Tray for storage	HTT057
Cell compartment cover	LZQ105
User manual, multilingual SIP 10, EU	DOC012.98.90328
Information on Tygon tubes, EC no. 1907/2006 (REACH)	DOC063.98.90334

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80539-0389 U.S.A.
(U.S.A. only)