RS232 Interface

Instruction Manual

Interfase RS232

Manual de instrucciones

Interface RS232

Manuel d'instructions

RS232 Schnittstelle

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Interfaccia RS232

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
INTRODUCTION

This Interface Kit is for use with the following Ohaus products: Scout Pro, Traveler and Navigator.


INTERFACE INSTALLATION

Install the Interface module on the underside of the balance as shown.


Note: The appearance of your model may be different.


1. Slide Blank Cover to rear.

2. Remove the Blank Cover.

3. Install Interface Module.

4. Slide Module forward.

SETUP

- A) The Ohaus Interface is preset to communicate using the following settings: 2400 baud, 7 bit, no parity, stop bits 2, no handshake. Determine the RS232 parameters required for the printer or computer that is to be connected. See computer or printer documentation if assistance is required. If the parameters do not match, it will be necessary to change either the balance settings, or the computer / printer settings.
- B) Upon installation of the interface, the balance will recognize the RS232 Interface and add relevant items to the menu. Configure the balance to the desired RS232 and printing parameters; refer to the balance User Manual for assistance in using the menus.

PRINT

Stable

On, Off

A-Print

Cont, On.Stbl, On.Acc*, 5sec, 15sec,
30sec, 60sec, off

End

Yes, No

*Note: Print and RS232 menu
selections may vary between
the different balance series.

RS232

RS232

On / off

Baud

600, 1200, 2400, 4800, 9600, 19200

Parity

7-even, 7-odd, 7-none, 8-none

Handshake

None, Xon-Xoff, RTS-CTS

END

Yes, No

SETUP DEFINITIONS

| | |
|---------------------------------|--|
| PRINT / Stable - On | Balance transmits stable data only. |
| PRINT / Auto Print – Continuous | Balance repeatedly sends data as fast as possible. |
| PRINT / Auto Print – On Stable | Balance automatically outputs data upon stability. |
| PRINT / Auto Print – (xx) sec | Balance sends data every (xx) seconds. |
| PRINT / Auto Print – Off | Only sends data when PRINT is pressed. |
| RS232 / RS232 –Off | Turns Interface off, saving battery power. |
| RS232 / Baud, Parity, Handshake | Set parameters to match printer or computer. |

CONNECTION

The RS232 Interface Kit cable terminates with a 9 pin Sub-D female connector.

Active pins: Pin 2 = TXD, Pin 3 = RXD, Pin 5 = Ground, Pins 4 and 6 are connected for hardware handshake, Pin 7 = CTS, Pin 8 = RTS.

The Interface connector will connect directly to a computer, an Ohaus printer and many other printers.

OPERATIONRS232 OUTPUT

With the Adapter installed the balance will operate in one of the following three ways:

- RS232 = on, Auto Print = off, Stable = on or off
Pressing PRINT will send the display data.
If Stable is set to On, the balance will wait for a stable reading before sending data.
- RS232 = on, Auto Print = on, Stable = on or off
The balance will automatically send data based on the settings in the Menu.
If Stable is set to On, only stable values will be sent.
- RS232 = off
The interface is turned off. During battery operation, this saves significant power.

Data sent from the Interface is in standard ASCII format terminated with a [CRLF] (carriage return-line feed). The output format is as follows:

Scout Pro (SPx) Output:**Output Format 1 (for the models 303/123/202/402/602/2001/6001/401FZH/601FZH/6000FZH):**Output String:

| | |
|-----------------------|---|
| [weight] | 12 characters (right justified) |
| [space] | 1 character |
| [unit] | 5 characters |
| [space] | 1 character |
| [stability indicator] | 1 character; "?" when unstable, blank when stable |
| [space] | 1 character |
| [Legend] | 10 character |

Note: All of the fields have fixed length.

Output examples (space indicated by '*'):

```
*****0.00*g*****
*****12.73*g*****?
*****12.72*g*****?
*****0.4500*oz****
*****0.5950*oz****?
*****3*PCS**?
*****3*PCS**
*****0.85*oz*****WET*WT
*****0.85*oz*****AIR*WT
```

Output Format 2 (for the models 401/601/6000):

Output String:

| | |
|-----------------------|---|
| [weight] | 11~12 characters (right justified) |
| [space] | 1 character |
| [unit] | 1~5 characters |
| [space] | 1 character |
| [stability indicator] | 1 character; "?" when unstable, blank when stable |
| [space] | 1 character |
| [Legend] | 10 character |

Note: The unit field length varies with different units. The weight field could be 11 or 12, depending on if the weight string has an additional decimal point.

Output example (space indicated by '*'):

```
*****200*g**
*****311*g*?
*****0.211*kg*?
```

Navigator (NV/NVL/NVT) Output:

Output String:

| | |
|-----------------------|---|
| [weight] | 10 characters (right justified) |
| [space] | 1 character |
| [unit] | 1~5 characters (left justified) |
| [space] | 1 character |
| [stability indicator] | 1 character; "?" when unstable, blank when stable |
| [space] | 1 character |
| [NET] | 0 or 3 character |
| [Legend] | 0~16 character |

Note: The unit field length varies with different units, the Legend field is only printed in certain modes, and the time string only printed at internal print.

Output example (space indicated by '*'):

```
*****200*g***
*****15.06*g***
*****15*g***NET
*****124*g*?*NET
*****15*g***NET**00:00:02
*****122*g*?*NET**00:00:03
*****15*g***NET*UNDER
***5:10.75*1b:oz*?*NET*ACCEPT*00:00:05
***0:*0.50*1b:oz*?*NET*UNDER*00:00:06
```

Traveler (Tax/SE) Output:

Output String:

| | |
|-----------------------|---|
| [weight] | 11 characters (right justified) |
| [space] | 1 character |
| [unit] | 1~5 characters |
| [space] | 1 character |
| [stability indicator] | 1 character; "?" when unstable, blank when stable |
| [space] | 1 character |
| [Legend] | 10 character |

Note: The unit field length varies with different units.

Output example (space indicated by '*'):

```
*****-0.01*g*?
*****-0.04*g*?
*****4.20*g*?
*****34.91*g*?
```

RS232 INPUT

The balance will respond to various commands sent via the interface adapter. Terminate the following commands when with a [CR] or [CRLF].

Scout Pro and Traveler Commands

| | |
|---------------|--|
| ? | print current mode |
| 0A | turn Auto-print off |
| SA | Auto-print on, prints on stability |
| CA | Continuous Auto-print |
| (n)A | Auto-print on 1 to 3600 second intervals (n = 1 to 3600) |
| C | perform span calibration |
| L | perform linearity calibration |
| 0M | gram mode |
| 1M | ounce mode |
| 2M | troy ounce mode |
| 3M | pennyweight mode |
| 4M | parts counting mode |
| 5M | pound mode |
| T | tare balance, same as pressing ON-ZERO |
| V | print software version |
| (Esc)R | reset balance to factory defaults |
| P | same as pressing PRINT |
| LE | print last error code, i.e. [Err 0] |
| 0S | print unstable data |
| 1S | print stable data only |

Navigator Commands

| | |
|-------------|--|
| P | same as pressing Print |
| SP | print stable weight only |
| IP | immediate print of displayed weight (stable or unstable) |
| CP | Continuous print of weights |
| SLP | Auto-print stable non-zero weight only |
| SLZP | Auto-print stable non-zero weight and zero reading |
| xP | Auto-print on 1 to 3600 second intervals (x = 1 to 3600) |
| 0P | turns auto-print off |
| PM | print current mode |
| M | advance to the next enabled mode |
| PU | print current unit |
| U | advance to the next enabled unit |
| T | same as pressing Tare |
| Z | same as pressing Zero |
| PV | print software version |

AUTO-PRINT OPERATION

Once Auto-Print is activated in the menu, the balance will send data as required. To temporarily stop Auto-Printing, press the PRINT key. If there is data in the print buffer the printer will finish printing this data. A second press will resume Auto-Printing.

ACCESSORIES

For a complete listing of Ohaus printers and other accessories, contact Ohaus Corporation or visit www.Ohaus.com

COMPLIANCE

This accessory has been tested and complies with the approvals listed in the applicable Instruction Manual.

Disposal

In conformance with the European Directive 2002/96 EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

For disposal instructions in Europe, refer to www.ohaus.com/weee.

Thank you for your contribution to environmental protection.

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



On, Off

A-Print

Cont, On.Stbl, On.Acc*, 5sec, 15sec, 30sec,
60sec, off

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P/N 83032207 C © 2020 Ohaus Corporation, all rights reserved / todos los derechos reservados / tous droits réservés / Alle Rechte vorbehalten / tutti i diritti riservati / все права защищены.