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When a Peltier Module is installed, **Temperature Control** are available in the "Instrument Setup" menu. Temperature Modes can be modified in this menu.

1. Touch **Temperature Control** in the "Instrument Setup" menu.

2. Enter the desired Parameter for

- **Temperature:** Turn temperature control off or turn temperature control on and input the set temperature.
- **Warning:** The warning alerts the user if the actual temperature deviates from the set temperature. Use this soft key to turn the warning off or to turn the warning on by inputting the deviation from set temperature that will cause the instrument to warn the user.
- **Unit:** Select temperature units of Celsius or Fahrenheit.
- **Idle:** The idle causes the instrument to automatically turn temperature control off after a user input amount of time. Use this key to turn the idle off or to turn the idle on by inputting the amount of time in hours to turn temperature control off.

Note: Alternatively, if a Peltier Module is installed, you can touch Temperature Control option icon directly on the Measurement screen. The "Temperature Control" option screen will appear where you can change the temperature.

Determination of the duration of the change of temperature

Carrying out the analysis in a 1 cm quartz cell

Duration of the temperature stabilization to the new value:

Approx. 10 min with an accuracy of $\pm 1^\circ\text{C}$

Approx. 10 min with an accuracy of $\pm 0.5^\circ\text{C}$

The following procedure can be used to determine the exact duration of the temperature stabilization:

Pipette 3 mL sample into the cell. Insert the cell in the module and place a temperature probe in the cell.

Note the initial temperature of the sample and set the temperature that has to be reached. With the help of a timer, you can determine the exact time needed to reach the desired temperature.

Note: *It takes longer for the temperature to stabilize in plastic cells, due to the lower thermal conductivity of plastic.*

Contact Information

HACH Company World Headquarters

P.O. Box 389
Loveland, Colorado
80539-0389 U.S.A.
Tel (800) 227-HACH
(800) -227-4224
(U.S.A. only)
Fax (970) 669-2932
orders@hach.com
www.hach.com

Repair Service in the United States:

HACH Company
Ames Service
100 Dayton Avenue
Ames, Iowa 50010
Tel (800) 227-4224
(U.S.A. only)
Fax (515) 232-3835

Repair Service in Canada:

Hach Sales & Service
Canada Ltd.
1313 Border Street, Unit 34
Winnipeg, Manitoba
R3H 0X4
Tel (800) 665-7635
(Canada only)
Tel (204) 632-5598
Fax (204) 694-5134
canada@hach.com

Repair Service in Latin America, the Caribbean, the Far East, Indian Subcontinent, Africa, Europe, or the Middle East:

Hach Company World
Headquarters,
P.O. Box 389
Loveland, Colorado,
80539-0389 U.S.A.
Tel +001 (970) 669-3050
Fax +001 (970) 669-2932
intl@hach.com

HACH LANGE GMBH

Willstätterstraße 11
D-40549 Düsseldorf
Tel. +49 (0)2 11 52 88-320
Fax +49 (0)2 11 52 88-210
info@hach-lange.de
www.hach-lange.de

HACH LANGE LTD

Pacific Way
Salford
GB-Manchester, M50 1DL
Tel. +44 (0)161 872 14 87
Fax +44 (0)161 848 73 24
info@hach-lange.co.uk
www.hach-lange.co.uk

HACH LANGE LTD

Unit 1, Chestnut Road
Western Industrial Estate
IRL-Dublin 12
Tel. +353(0)1 46 02 5 22
Fax +353(0)1 4 50 93 37
info@hach-lange.ie
www.hach-lange.ie

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HACH LANGE SU ANALİZ SİSTEMLERİ LTD.ŞTİ.