

300

Spectroquant® NOVA 30

M

Contents

1	Photometers	5
1.1	Photometry	5
1.2	The Photometers	6
2	Photometric Test Kits	6
2.1	Basic Principle	6
2.1.1	Spectroquant® Cell Tests	7
2.1.2	Spectroquant® Reagent Tests	7
2.2	Notes for Practical Use	8
2.2.1	Measuring Range	8
2.2.2	Influence of pH	10
2.2.3	Influence of Temperature	10
2.2.4	Time Stability	10
2.2.5	Influence of Foreign Substances	11
2.2.6	Dosing of Reagents	11
2.2.7	Shelf-life of the Reagents	12
3	Sample Preparation	12
3.1	Taking Samples	12
3.2	Preliminary Tests	13
3.3	Dilution	13
3.4	Filtration	14
3.5	Homogenization	15
3.6	Decomposition	15
4	Pipetting System	17
5	Analytical Quality Assurance (AQA)	18
5.1	Quality Control at the Manufacturer	18
5.2	Quality Control for the User	19
5.2.1	Checking the Photometer	20
5.2.2	Checking the Overall System	20
5.2.3	Checking the Pipettes	21
5.2.4	Checking Thermoreactors	21
5.2.5	Testing for Handling Errors	22
5.3	Determination of Sample Influences	22
5.4	Definition of Errors	23


1 Photometers

1.1 Photometry

When a beam of light is transmitted through a colored solution, then this beam loses its intensity, in other words a part of the light is absorbed by the solution. Depending on the substance in question, this absorption occurs at specific wavelengths.

Monochromators (e.g. narrow-band interference filters, lattices) are used to select the wavelength from the total spectrum of a tungsten-halogen lamp (VIS spectrum), a deuterium lamp (UV spectrum) or, respectively, a xenon lamp.

The intensity of the absorption can be characterized using the transmittance T (or, respectively, T in percent).



$$T = I/I_0$$

I_0 = Initial intensity of the light

I = Intensity of the transmitted light

If the light is not absorbed at all by a solution, then this solution has a transmittance of 100 %; a complete absorption of the light in the solution means 0 % transmittance.

The measure generally used for the absorption of light is the absorbance (A), since this correlates directly with the concentration of the absorbing substance. The following connection exists between absorbance and transmittance:

$$A = -\log T$$

Experiments by BOUGUER (1698–1758) and LAMBERT (1728–1777) showed that the absorbance is dependent on the thickness of the absorbing layer of the cell used. The relationship between the absorbance and the concentration of the analyte in question was discovered by BEER (1825–1863). The combination of these two natural laws led to the derivation of *Lambert-Beer's law*, which can be described in the form of the following equation:

$$A = \varepsilon_{\lambda} \cdot c \cdot d$$

ε_{λ} = Molar absorptivity, in $l/mol \times cm$

d = Path length of the cell, in cm

c = Concentration of the analyte, in mol/l

1 Photometers

1.2 The Photometers

The photometers that belong to the Spectroquant® Analysis System differ from conventional photometers in the following important aspects:

- The calibration functions of all test kits are electronically stored.
- The measurement value can be immediately read off from the display in the desired form.
- The method for the test kits (Cell Tests **and** reagent tests) belonging to the Spectroquant® analysis system is automatically selected via the scanning of the bar code.
- All cells formats used are automatically identified and the correct measuring range is selected automatically.
- Instrument-supported AQA ensures that measurement results can be used as secure, reproducible, and recognized analytical results.
- New methods can be downloaded from the internet site www.service-test-kits.com and permanently stored in the instrument.

For technical data and instructions for use please refer to the section “Function description” or can also be found on the internet.

2 Photometric Test Kits


2.1 Basic Principle

By means of reagents, the component of a sample to be analyzed is converted into a colored compound in a specific reaction. The reagents or reagent mixtures contain – in addition to the reagent selective for a parameter to be determined – a number of auxiliary substances that are essential for the course of the reaction. These include, for example, buffers for adjusting the pH to the optimal value for the reaction, and masking agents that suppress or minimize the influence of interfering ions.

The color reactions are in most cases based on standardized analytical methods specifically optimized in terms of ease of use, a low working effort, and shorter reaction times. Furthermore, methods cited in the literature or developed by ourselves are also used. Details on the respective reference procedures are stated in the package insert or else in the parameter overview.

2 Photometric Test Kits

2.1.1 Spectroquant® Cell Tests





Additional reagent(s)

Certain cell tests, e. g. COD or nitrite, already contain all necessary reagents in the cells, and the sample must merely be added with a pipette. In other tests, however for reasons of chemical compatibility it is necessary to separate the test into two or three different reagent mixtures. In such cases, besides the sample a metered reagent must also be added.

2.1.2 Spectroquant® Reagent Tests

The principle behind the reagent tests is that the reagents necessary for the color reaction are combined in the form of liquid concentrates or solid-substance mixtures. A few drops of the reagent concentrate are added to the sample. This means that there is no need to dilute the sample, which in turn enhances the sensitivity of the detection. The procedure generally used in classical photometry by which the sample is made up to a defined volume in a volumetric flask is dispensed with.



The method is selected automatically by means of the scanning of the bar code by the AutoSelector.

All cells formats used are automatically identified and the correct measuring range is selected automatically.

Subsequently the result is automatically shown on the display.

2 Photometric Test Kits

2.2 Notes for Practicle Use

2.2.1 Measuring range


The intensity of the color of a solution, measured as the absorbance, is proportional to the concentration of the respective analyte only within a specific range. This measuring range (effective range) is electronically stored in the photometers for each individual test kit .

Below the specified measuring range, either a different cell or else another procedure must be used. The **lower limit of the measuring range** either takes the form of nonlinearity of the calibration curve, as shown in the figure, or else is given by the method detection limit. The **method detection limit** of an analytical method is the lowest concentration of the analyte in question that can be measured quantitatively with a defined degree of probability (e.g. 99 %).

The **upper limit of the measuring range** is the point at which the linear correlation between the concentration and the absorbance ends. In such a case the sample must be diluted accordingly so that it lies ideally in the middle of the effective range (least-error measurement).

In photometry it is conventional practice to measure against the reagent blank value. Here the analysis is carried out "blind", i.e. without any analyte added. Instead of the sample volume, the corresponding quantity of distilled or DI water is used. This **reagent blank value is prestored** in the photometers belonging to the Spectroquant® analysis system, which means that - due to the high batch reproducibility - it is possible to dispense with a separate measurement of the reagent blank. At the lower limit of the measuring range, the accuracy of the determination can be enhanced by performing the measurement against a separately prepared reagent blank.

In some cases the intensity of the color of the solution and thus the absorbance can drop again when **very high concentrations of the analyte** are present (see package insert).



2 Photometric Test Kits

2.2.2 Influence of pH

Chemical reactions follow an optimal course only within a certain pH range. The reagents contained in the test kits produce an adequate buffering of the sample solutions and ensure that the pH optimal for the reaction in question is obtained.

Strongly acidic ($\text{pH} < 2$) and strongly alkaline ($\text{pH} > 12$) sample solutions can prevent the pH from being adjusted to an optimal range, since under certain circumstances the buffering capacity of the test-kit reagents may not be sufficient. Any necessary correction is made by the dropwise addition of diluted acid (reduces the pH) or diluted lye (raises the pH), testing the pH with suitable indicator strips after each drop is added. The addition of the acid or lye results in a dilution of the test solution. When up to five drops are added to 10 ml of sample, the change in the volume can be neglected, since the resultant error is lower than 2 %. The addition of larger quantities should be duly considered by adjusting the sample volume accordingly.


The specified pH values for the sample solution and, wherever applicable, for the measurement solution are defined in the respective package inserts and in the analysis instructions in chapter 3 of the manual.

2.2.3 Influence of Temperature

The temperature of the sample solution and the reagents may have an effect on the color reaction and thus on the measurement result. The typical temperature course is illustrated in the figure.

If the sample temperature is lower than 15 °C, false-low results must be reckoned with. Temperatures exceeding 30 °C generally influence the stability of the compound that is formed in the reaction. The optimal temperature for the color reaction is stated in the package inserts of the respective Spectroquant® test kits.


Attention! After thermic decomposition procedures, the determination of COD or total contents of nitrogen, phosphorus, or metal, a sufficient waiting time must be allowed for to permit the solution cool to room temperature.



2.2.4 Time Stability

Most of the color reactions require a certain time to reach the maximum color intensity. The solid curve in the figure at the right gives a schematic impression of a typical time course. The behavior of relatively instable color reactions with time is shown by the dotted curve.

The reaction time specified in the working instructions refers to the period of time from the addition of the last reagent until the actual measurement. In addition, the package inserts for the individual test kits also state the time interval in which the measurement value does not change. The maximum time interval is 60 minutes; this time should not be exceeded, even in the case of stable color reactions.



2 Photometric Test Kits

2.2.5 Influence of Foreign Substances

Foreign substances in the sample solution can

- raise the measurement value as a result of an amplification of the reaction
- lower the measurement value as a result of a prevention of the reaction.

A quantification of this effects is stated in tabular form in the respective package inserts for the most important foreign ions. The tolerance limits have been determined for the individual ions; they may not be evaluated cumulatively.

Suitability for use in seawater

A tabular survey (see appendix 1) provides information on the suitability of the tests in connection with seawater and also on the tolerances for salt concentrations.

2.2.6 Dosing the Reagents




Small amounts of liquids are dosed by counting the number of drops from a leakproof bottle.



When using dropper bottles it is extremely important that the bottle be held vertically and that the drops be added slowly (approx. 1 drop per second). If this is not observed, the correct drop size and thus the correct amount of reagent are not achieved.




A positive-displacement pipette should be used for larger quantities of liquid or for the exact dosage of smaller reagent quantities. In these cases the reagent bottles are not fitted with a dropper insert.



Solid substances are dosed either with the dose-metering cap or with microspoons that are integrated into the screw cap of the respective reagent bottle. The dose-metering cap is used for solid reagents or reagent mixtures that are free-flowing.

In all other cases the substances are dosed with the microspoon.

In this case it is necessary to add only level microspoonfuls. To this end the spoon must be drawn over the brim of the reagent bottle.



At the first use replace the black screw cap of the reagent bottle by the dose-metering cap.

Hold the reagent bottle vertically and, at each dosage, press the slide all the way into the dose-metering cap. Before each dosage ensure that the slide is completely retracted.



Reclose the reagent bottle with the black screw cap at the end of the measurement series, since the function of the reagent is impaired by the absorption of atmospheric moisture.

2 Photometric Test Kits

2.2.7 Shelf-life of the Reagents

The Spectroquant® test kits are in most cases stable for 3 years when stored in a cool, dry place. A few test kits have a lower shelf-life of 18 or 24 months or must else be stored in a refrigerator.

COD Cell Tests must be stored protected from light.

The expiry date of the package unit is printed on the outer label. The shelf-life may become reduced when the reagent bottles are not reclosed tightly after use or when the test kit is stored at temperatures higher than those specified.

3 Sample Preparation

Sample preparation covers all the steps necessary before the actual analysis can be performed.

3.1 Taking Samples

The taking of samples is the first and **most important step** on the way to obtaining the correct analysis result. Not even the most exact method of analysis can correct any mistakes made in the taking of the sample. The objective of the sampling procedure is to gain a sample with a representative composition. The most important precondition for **gaining a representative sample** is the identification of the suitable sampling site. Here it must be borne in mind that the solution to be investigated can display varying concentrations in different places at different times.

In sampling, a distinction is made between manual and automatic methods. In many cases a true picture of the average composition of the sample can be obtained only once several individual samples have been collected; this can be done manually or with an automatic sampler.

Clean plastic containers with a volume of 500 or 1000 ml are suitable for collecting samples. They should be rinsed several times, under vigorously shaken, with the water to be investigated, and then filled free of air bubbles and immediately closed tightly. The containers must be protected against the effects of air and heat and then be forwarded for the further analytical steps as soon as possible. In exceptional cases, preservation measures in the form of short-term refrigeration at +2 to +5 °C and chemical conservation can be taken.

Parameter	Preservation
COD	+2 to +5 °C max. 24 h or -18 °C max. 14 days
N compounds: NH ₄ -N, NO ₃ -N, NO ₂ -N	analyze immediately, only in exceptional case +2 to +5 °C max. 6 h
P compounds: PO ₄ -P, P total	short-term storage, no preservation; with nitric acid to pH 1, max. 4 weeks
Heavy metals	short-term storage, no preservation; with nitric acid to pH 1, max. 4 weeks

3 Sample Preparation

3.2 Preliminary Tests

Correct measurement results can be obtained only within the measuring range specified for each individual parameter. When dealing with sample solutions of an unknown concentration, it is advisable to establish whether the sample concentration is indeed within the specified measuring range, ideally roughly in the middle of the range.

Preliminary tests enhance the analytical reliability and make the determination of the necessary dilution ratios in the case of high concentrations easier. **MQuant™ Test Strips** are very well suited for preliminary tests.

3.3 Dilution

Dilution of samples is necessary for two reasons:

- The concentration of the parameter under investigation is too high, i. e. it lies outside the measuring range.
- Other substances contained in the sample interfere with the determination (matrix interference); false-high or false-low results may ensue.

The following auxiliaries are absolute prerequisites for the dilution of the sample:

- Volumetric flasks of varying sizes (e. g. 50, 100 and 200 ml)
- Positive-displacement pipette
- Distilled or DI water.

Only dilutions carried out with these auxiliary products are of sufficient reliability in the area of trace analysis, to which photometry belongs (for the simplified procedure see page 14).

An important aspect here is that once the volumetric flask has been filled up to the mark with distilled water the flask is closed and the contents are thoroughly mixed.

The **dilution factor (D_F)** resulting from the dilution procedure is calculated as follows:

$$D_F = \frac{\text{Final volume (total volume)}}{\text{Initial volume (sample volume)}}$$

The analytical result is subsequently multiplied by the dilution factor.

A calculation can be dispensed with when the dilution is programmed into the photometer. The **dilution number** (see the table on page 14) is entered and the measurement value is subsequently calculated correctly and immediately displayed.

3 Sample Preparation

All dilutions should be made in such a way that the measurement value lies in the middle of the measuring range. As a rule, the dilution factor should never be higher than 100. In the event that yet higher dilutions become necessary all the same, then this must be done in two separate steps.

Example

Step 1: Make up 2 ml of sample to 200 ml with distilled water;
 $D_F = 100$, dilution number 1+99

Step 2: Take 5 ml of the above solution and make up to 100 ml;
 $D_F = 20$, dilution number 1+19

The dilution factor for the total dilution is calculated by multiplying the individual dilutions:

$$D_{F_{\text{total}}} = D_{F_1} \times D_{F_2} = 100 \times 20 = 2000, \text{ dilution number } 1+1999$$

Simplified procedure

Dilutions up to 1:10 can also be prepared without volumetric flasks in a glass beaker, measuring the volumes of the sample and the dilution water using a previously calibrated positive-displacement pipette (see table for instructions).

Desired dilution	Volume of sample [ml]	Volume of distilled water [ml]	Dilution factor	Dilution number
1:2	5	5	2	1+1
1:3	5	10	3	1+2
1:4	2	6	4	1+3
1:5	2	8	5	1+4
1:10	1	9	10	1+9

3.4 Filtration

Strongly turbid samples require pretreatment before they can be determined in a photometer, since the effect of turbidity can result in considerable variations in the measurement values and in false-high readings. Care must be taken here to ensure that the substance to be determined is not contained in the suspended material, in which case a sample decomposition must be carried out.


Compounds that always occur in dissolved form (for example ammonium, nitrate, nitrite, chlorine, chloride, cyanide, fluoride, orthophosphate, and sulfate) permit a previous filtration, even when the sample solution is strongly turbid.

Weak turbidity is eliminated by the **automatic turbidity-correction** feature built into the photometer (see Function description, "Device set-up/Correction function"); in such cases it is not necessary to filter the sample before analysis.


As a measure to distinguish between dissolved and undissolved water-borne substances, the water sample can be filtered through a simple paper filter. Following the recommendations stated in the reference methods, membrane filters with a pore size of 0.45 µm are required for fine filtration.

3 Sample Preparation


Procedure for microfiltration




Draw out the liquid to be filtered with the syringe.



Screw the syringe tightly into the front side of the membrane-filter attachment.



Hold the syringe upright and slowly depress the piston upwards until the membrane-filter is fully wetted free of air bubbles.




Filter the contents of the syringe into the intended glass vessel.

3.5 Homogenization

As a measure to ensure that a representative sample can be taken in the presence of suspended matter in the water sample in question, for certain parameters - e.g. COD and the total content of heavy metals - the sample must be homogenized. This must be carried out using a high-speed blender (2 minutes at 5000 – 20 000 rpm and taking the sample while stirring).


3.6 Decomposition



Water-borne substances can be present in the sample for investigation in a variety of forms: as the ion, bound more or less solidly in a complex, or as a solid substance.

3 Sample Preparation

The manner in which the sample is pretreated enables the three proportions to be distinguished from each other. This can be illustrated using a copper-containing wastewater sample as an example.



Decomposition converts the substance to be determined into an analyzable form. In most cases, decomposition agents take the form of acids in combination with oxidizing agents; in exceptional cases (e. g. in the determination of total nitrogen) an alkaline decomposition is more effective. The type of decomposition procedure used depends on the analyte to be determined and the sample matrix.

The ready-to-use sample-decomposition products **Spectroquant® Crack Set 10** and **20** are suited for the preparation of the sample materials for the determinations stated in the table.

The decomposition processes are carried out in the **Spectroquant® thermoreactor** (capacity: 12 or 24 decomposition cells) at 120 °C or, respectively, 100 °C. Details regarding the heating times and further treatment can be found in the package inserts contained in the **Spectroquant® Crack Set** packs.

Determination of	Sample preparation with
Total phosphorus*	Crack Set 10 / 10C**
Total chromium* [= sum of chromate and chromium(III)]	Crack Set 10 / 10C
Total metal [= sum of free and complex-bound metal]	Crack Set 10 / 10C
Total nitrogen*	Crack Set 20


* The decomposition reagents are already contained in the packs of the respective cell tests.

** Decomposition cells are included in the pack; empty cells are required for the decomposition for Crack Sets 10 and 20.

3 Sample Preparation

In the event that the sample to be analyzed is a highly contaminated material (high proportion of organic substances) or water-insoluble samples, decomposition using concentrated acids and other agents is indispensable. Corresponding examples from the **collection of applications** for real samples are available on request.

The necessity for decomposition can be checked according to the following diagram:



For wastewater with a consistent composition, this check as a rule need be carried out only once. It is, however, advisable to check the result periodically.

4 Pipetting System

Positive-displacement pipettes permit

- an exact dosage of the sample volume
- a precise measurement of sample and reagent volumes and of the volumes of water for dilution purposes.

Pipettes of varying volumes and also ones with a fixed volume are available.

Sources of error and hints on how to avoid them:

- Closely follow the instructions for use contained with the pipette in question.
- Check the pipetted volumes
 - a) by weighing using analytical scales (weighing accuracy ± 1 mg),
1 ml of water at 20 °C = 1.000 g ± 1 mg
 - b) using Spectroquant® PipeCheck;
this is a photometric check of the pipette, and scales are not necessary (see section "AQA").
- Avoidance of spread effects by rinsing the pipette several times with the solution to be pipetted.
- Always exchange the pipette tip.
- Draw up the liquid slowly and depress piston completely to discharge the liquid.

5 Analytical Quality Assurance (AQA)

The objective of analysis must always be to determine the true content of the analyte in question as accurately and precisely as possible.

Analytical Quality Assurance represents a suitable and indispensable method by which the quality of the user's own work can be assessed, errors in the measurement system diagnosed, and the comparability with the results obtained using the respective reference methods demonstrated.

Details regarding the necessity of AQA can be found in the in Memorandum A 704 of the German Association for the Water Sector, Wastewater, and Waste Materials (Deutsche Vereinigung für Wasserwirtschaft, Abwasser und Abfall e.V., DWA) and in the corresponding self-control/self-monitoring regulations of the German federal states (available in english).

Causes for errors can include:

- the working materials used
- the handling
- the sample under investigation.

These errors have effects on both the accuracy and precision of the results obtained.

5.1 Quality Control at the Manufacturer

Photometers and photometric test kits possess specifications that are adhered to and above all else also documented by the manufacturer.

The **certificate for the photometer** enclosed with each device documents the quality of the measuring device.

Prüfprotokoll Test record		Seite 1 von 1 page 1 of 1
Photometer / photometer	Modell / Model Screen-Nr. / Serial no.	SO NOVA 60 Akku 07340913
Transponderfunktion / Funktion of transponder	Korrekt Einlesen eines Testdatenatzes / correct reading of test data set	ok
Selbsttest / Self Check	Signalabgleich ohne Küvette / signal adjustment without cell	ok
LS-Check / LS-Check	Korrekte Erkennung von Test-Barcodes / correct identification of test barcodes	ok
Nullabgleich Rechteckküvette 10 mm / Zero adjustment rectangular cell 10 mm		ok
Nullabgleich Rechteckküvette 20 mm / Zero adjustment rectangular cell 20 mm		ok
Nullabgleich Rundküvette / Zero adjustment round cell		ok
Bei allen Filterpositionen Abgleich auf Extinktion E = 0 mit entionisiertem Wasser als Messlösung / at all filterpositions adjustment to absorbance A = 0 using deionised water as measurement solution		
Photometrische Richtigkeit / photometric accuracy	Extinktion E einer Testlösung in Rundküvette / absorbance A of test solution in round cell	
Wellenlänge / wavelength (nm): 605	Schwert E / nominal value A: 0.395	
Toleranz E / toleranz A: ±0.020	Messwert E / measured value A: 0.379	ok
Linearität / linearity	Extinktionsergebnisse von 2-Phasfiltern in separater (E1, E2) und kombinierter Anordnung (E12) / absorbance data of 2 plane filters in separate (E1, E2) and combined configuration (E12)	
Messwerte E / Measured values A (445nm):	E1 = 0.986 E2 = 1.015 E12 = 2.001	
Anforderung / requirement:	±0.02 ±(E1 + E2 - E12) ≤ 0.02	ok
Elektrische Sicherheit nach IEC 1010 / electrical safety according to IEC 1010	Keine visuellen Mängel, keine Gerüche, keine losen Teile und Befestigungen / no visual flaws, no burrs, no loose parts and fastenings	ok
Datum / Date:	30.01.2007	Prüfer / Tester: Michael Dobry
Dieses Dokument wurde mittels EDV erstellt und ist ohne Unterschrift gültig. This document has been generated using electronic data processing and is valid without signature.		

5 Analytical Quality Assurance (AQA)

Chargenzertifikat
Lot Certificate / Certificado del lote

M

Spectroquant® CSB-Küvettestest
Spectroquant® COD Cell Test / Spectroquant® Test en cubetas DQO

Act.No. / Cat.No. / Art. No.	14888-0201	n =	10
Measuring Range / Intervalo de medida	1.0 - 40.0 mg/l CSB/COD/DQO	Measured Result / Resultado (Standard / Percent) / CSB/COD/DQO	32.0
Change No. / Lot no. / Lote no.	001/0001	Target value / Valor nominal	32.0
Expiry date / Fecha de caducidad	31.10.2012	Relative Error / Error relativo	0.0 %
Standard / Standard / Patrón	Potassium hydrogen citrate 1.0000	Standard Deviation / Desviación estándar	0.5
Wavelength / Wavelength / Longitud de onda	440 nm	Confidence Interval (95%) / Intervalo de confianza (95 % de probabilidad)	± 0.3 mg/l
Cell / Cuveta	10 mm quartz (quartz) reference	Standard Deviation of the Method / Desviación estándar del procedimiento	± 0.3 mg/l
Supplier / Fabricante	Dr. Spectro	Variation Coefficient of the Method / Coeficiente de variación del procedimiento	± 0.9 %
Batch / Lote / Referencia	148880201 / 1481/0001 / 01		

Calibration Function / Función de calibración	Y = 0.8000 X - 0.0000	Calibration Error / Error de calibración	0.0 %
Designing / Diseño / Plantilla	1	Tolerance / Tolerancia	± 0.5 %
Optimization / Optimización / Determinación según el método en ordenada	1	Tolerance / Tolerancia	± 0.5 %
Linearity / Linealidad / Verificación de la linealidad	1	Tolerance / Tolerancia	± 0.5 %
Verification / Verificación / Verificación de la precisión	1	Tolerance / Tolerancia	± 0.5 %
Standard Deviation of the Method / Desviación estándar del procedimiento	± 0.3 mg/l		
Variation Coefficient of the Method / Coeficiente de variación del procedimiento	± 0.9 %		

Merck KGaA
Laborleiter / Head of Lab.
Jefe de laboratorio

Qualitätskontrolle
Quality control / Control de calidad

The **certificate for the test kit**, available for each lot produced, documents the quality of the reagents contained in the test kit.

Calibration function:

The calculated function must agree, within specified tolerances, with the function electronically stored in the photometer.

Confidence interval:

Maximum deviation from the desired value over the entire measuring range; every measurement value can be affected by this deviation; this parameter is a measure for the accuracy.

Standard deviation for the procedure:

Measurement for the dispersion of the measurement values over the entire measuring range, expressed in \pm mg/l.

Coefficient of variation for the procedure:


Measurement for the dispersion of the measurement values over the entire measuring range, expressed in %. The smaller the standard deviation/ coefficient of variation for the procedure, the more pronounced the linearity of the calibration curve.

5.2 Quality Control for the User

A complete check comprises the entire system, i.e. the working equipment and the mode of operation. The photometer offers an optimum degree of support in this regard, in the form of the different quality mode. The instrument, or the whole system (including reagents and all accessories) will be checked, depending on which quality mode selected. All of checking operations can thus be supported by the photometer and the check values accordingly documented as per GLP (Good Laboratory Practice) recommendations (see Function description, "Analytical Quality Assurance").


The following diagram provides an overview regarding internal quality-assurance aspects:

Checking the working equipment



= Test for the overall system

Checking the handling operations



Influence of the sample

Test for recovery

5 Analytical Quality Assurance (AQA)

5.2.1 Checking the Photometer



As soon as the photometer is activated it is running a Self-Check. This means the hardware and the software of the photometer is checked and compared with internal standards.

As soon as the photometer is activated it is running a Self-Check. This means the hardware and the software of the photometer is checked and compared with internal standards.

The photometer itself is checked in the **AQA 1 mode** with the **Spectroquant® PhotoCheck**: the pack includes round cells containing stable test solutions (**secondary standards**) for checking the photometer at the 445, 525, and 690 nm wavelengths. The test solutions are measured in a **reference photometer** monitored with **primary standards**, and the certificate stating the absorbance values is enclosed with the package unit. These desired values with the permissible tolerances are entered into the photometer or else handwritten into the control chart. For the measurement the cell is placed in the compartment for the round cell and identified by the photometer via the bar code, and the measured absorbance is compared with the desired value. The absorbance is shown on the display and can be entered into the corresponding control chart.

The measurement of four cells for a given wavelength tests – in addition to the wavelength accuracy – also the linearity of the absorbance over the effective range.

The verification of the instrument, as it is required by DIN/ISO 9000 or GLP, can be easily performed by using the Spectroquant® PhotoCheck. The PhotoCheck hence offering the possibility to check the instrument. All of the corresponding documentation, required by these certification guidelines, is done by the photometer automatically.



5.2.2 Checking the Overall System

Test for the overall system includes checking the working equipment and checking the handling operations.

The **overall system** can be checked using standard solutions of a known content, preferably with the Spectroquant® CombiCheck; this corresponds with the **AQA 2 mode** in the photometer.

Spectroquant® CombiCheck are ready-to-use standard solutions that in terms of the analyte concentration are finely adjusted to the individual test kits. They contain a mixture of several analytes that do not interfere with each other. The standard solution (R-1) is used in the same way as a sample. A double determination is recommended as a measure to diagnose any random errors.

Standard solutions for photometric applications (CRM) are ready-to-use standard solutions that in terms of the analyte concentration are finely adjusted to the individual test kits. The standard solution is used in the same way as a sample. A double determination is recommended as a measure to diagnose any random errors.

In addition to the CombiCheck and the standard solutions for photometric applications, it is also possible to use **CertiPUR® standard solutions** for this checking procedure. These contain 1000 mg of the respective analyte per liter of solution.

They can be diluted to different final concentrations, which should preferably lie approximately in the middle of the measuring range of the respective test kit. The table presented in Appendix 2 provides an overview of the available CombiCheck and ready-to-use standard solutions.

5 Analytical Quality Assurance (AQA)

Due to limited shelf-life characteristics, there are no CombiCheck or ready-to-use standard solutions for certain parameters. Appendix 3 is a compilation of **standard working procedures** necessary to make your own solutions of a defined concentration. This allows the control of parameters where there are no simple to prepare solutions available.


If the test for the overall system shows that all requirements are fulfilled, the individual results are flagged as AQA2. If not, an error message is given and the individual components of the instrument have to be checked in detail.

5.2.3 Checking the Pipettes



The **Spectroquant® PipeCheck** is used to check the pipettes. The pack contains cells filled with color-dye concentrates. After the addition of a predefined volume of water using the pipette in question, the cell is measured against a corresponding reference cell also contained in the pack. The difference in the absorbance values of the measurement cell and reference cell may not exceed the tolerances given in the package insert. If the tolerances are exceeded, the instructions given in the section “Pipetting system” must be followed accordingly.

5.2.4 Checking Thermoreactors



This is checked by means of the thermosensor. The thermoreactor is pre-heated as described in the Instructions for use. When the control lamp goes out, the temperature is measured in any one of the bores of the thermoreactor. The following desired temperatures must be achieved:

Block temperature 100 °C = desired temperature 100 ±3 °C
Block temperature 120 °C = desired temperature 120 ±3 °C
Block temperature 148 °C = desired temperature 148 ±3 °C

The even distribution of the temperature over all bores can also be documented using the thermosensor.

5 Analytical Quality Assurance (AQA)

5.2.5 Testing for Handling Errors

The user's own mode of operation must also be subjected to an exact analysis.

The following questions may serve as a guide in this regard:

- Is the test kit optimal for the measurement assignment in question?
- Is the test kit's measuring range suitable?
- Were the operating instructions for the test followed?
- Was the sample volume correct?
- Was the pipette handled properly?
- Was a new pipette tip used?
- Is the pH of the sample and measurement solution correct?
- Was the reaction time adhered to?
- Does the sample and reagent temperature lie within the correct range?
- Is the cell clean and free from scratches?
- Has the expiry date for the test kit been exceeded?

5.3 Determination of Sample Influences (matrix effects)

The influence of other substances contained in the sample may, under certain circumstances, be so great that their recovery rates lie in the region of several percent. It is recommended to check for any influence by using the addition solution contained in the Spectroquant® CombiCheck pack.

A defined quantity of the **addition solution** (R-2), which contains a known concentration of the respective analyte, is added to the sample and the recovery rate is determined. The following difference is then calculated:

Result (sample + addition solution) – Result (sample)

If the calculated difference is equal to the concentration of analyte of addition solution that was added, the recovery rate is 100 %. If the difference is less than 90 %, then a matrix interference is present.

5 Analytical Quality Assurance (AQA)

5.4 Definition of Errors

It is obvious that measurement results as a rule may be associated with errors. This applies equally to standardized methods of analysis (reference methods) and to routine analysis. The discovery and the minimization of errors must be the objective here.

A distinction is made between systematic errors and random errors.


Systematic errors are present when all the results of an analysis deviate from the true value with the same algebraic sign. Examples here include: a wrong sample volume, a wrong pH, a wrong reaction time, a sample-matrix influence, etc. Systematic errors thus affect the **accuracy** of the method of analysis.

Accuracy = Deviation of the measured concentration from the true concentration


Random errors manifest themselves in the form of a wide range of deviation of the results of a given sample. These can be kept to a minimum by ensuring good operating techniques and multiple determination with calculation of the mean values. Random errors make the result of the analysis unreliable; they influence the **precision**.

Precision = Dispersion of the results among each other


The following diagram illustrates the aspects of accuracy and precision:




Accuracy: poor
Precision: poor
Major errors have been made!



Accuracy: good
Precision: poor
Calculation of the mean values from at least three – or better even more – parallel determinations yields an approximation of the true value.




Accuracy: poor
Precision: good
The high degree of precision mistakenly indicates a correct value!



Accuracy: good
Precision: good
The ideal objective!

enno



Spectroquant® **NOVA 30**

Description of Function

General instructions

Notes on this operating manual

To ensure that you become rapidly acquainted with your photometer, the first chapter contains an overview and a short manual of the meter. The second chapter contains notes for the safe operation of the photometer.

Chapter 3 describes the commissioning of the photometer. The remaining chapters provide a comprehensive description of the functions and technical data of the photometer.

Symbols used



indicates notes that you must read – for your own safety, the safety of others and to protect your meter from being damaged.



indicates notes that draw your attention to special features.

Scope of delivery

- Photometer
- Power pack
- Product documentation

Warranty

The designated meter is covered by a warranty of 2 years from the date of purchase. The meter warranty extends to manufacturing faults that are determined within the period of warranty. The warranty excludes components that are replaced during maintenance, such as batteries, accumulators, lamps etc.

The warranty claim extends to restoring the meter to readiness for use but not, however, to any further claim for damages. Improper handling or unauthorized opening of the instrument invalidates any warranty claim.

To ascertain the warranty liability, return the meter and proof of purchase together with the date of purchase freight paid or prepaid.

Accuracy when going to press

The use of advanced technology and the high quality standard of our instruments are the result of continuous development. This may result in differences between this operating manual and your meter. We cannot guarantee that there are absolutely no errors in this manual. We are sure you will understand that we cannot accept any legal claims resulting from the data, figures or descriptions. The information in this manual is subject to change without notice.


1. Overview	29
1.1 Description of the operating elements	29
1.2 Identifying the connectors	29
1.3 Short manual	30
1.3.1 Measuring the concentration	30
1.4 Selecting and calling up the menu items	31
2. Safety	32
2.1 Authorized use	32
2.2 General instructions	32
2.2.1 Labeling of notes	32
2.2.2 Dangers of disregarding the safety instructions	32
2.2.3 Qualification of the personnel	32
2.2.4 Technical state of the meter	32
3. Commissioning	34
3.1 Preparing the photometer	34
3.2 Switching on the photometer	34
4. Measuring the Concentration	35
4.1 Display of concentration and absorption	35
4.2 Measuring using cell tests	36
4.3 Measuring using tests without barcode (manual method selection)	37
5. Measuring the Absorbance/Transmission	38
5.1 Switching to the Absorbance/ Transmission measuring mode	38
5.2 Measuring the absorbance or transmission	38
5.3 Measuring using cell tests	39
5.4 Measuring using tests without barcode	39
6. Documentation	40
6.1 Resetting the number of the measured value	41
6.2 Download memory	41
6.3 Download of the methods list	43
7. Method Parameters	44
7.1 Citation form	45
7.1.1 Changing the citation form	45
7.1.2 Performing a difference measurement	46
7.2 Selecting the unit	48
8. Analytical Quality Assurance (AQA)	49
8.1 Activating AQA	49
8.1.1 Activating AQA via the menu guide	49
8.1.2 Changing AQA intervals	51
8.1.3 Locking the system	53
8.1.4 Changing the password	54
8.1.5 Performing an AQA reset	54

Contents

8.2	Photometer monitoring (AQA1)	55
8.2.1	Entering PhotoCheck standards	55
8.2.2	Download of PhotoCheck standards	57
8.2.3	Erasing PhotoCheck standards	57
8.2.4	Performing Photometer monitoring	58
8.3	Total system monitoring with standard solutions (AQA2)	60
8.3.1	Entering standards	60
8.3.2	Output of standards	62
8.3.3	Erasing standards	63
8.3.4	Monitoring of the total system using standard solutions	63
9.	Correction functions	65
9.1	Blank value	65
9.1.1	Activating the blank value measurement	66
9.1.2	Measuring the blank value	67
9.1.3	Erasing blank values	67
9.1.4	Recalling blank values	68
9.2	Turbidity correction	69
10.	Zero adjustment	70
11.	Meter Setup	71
11.1	Selecting the language	71
11.2	Setting the date/time	72
11.3	Reset	72
11.4	System info	73
12.	Updating method data	74
13.	RS 232 C interface	76
13.1	Principle course of the remote control	76
13.2	Command list	76
13.3	Output format of measured values	77
13.4	Data transmission	77
13.5	Pin assignment	77
14.	Maintenance, cleaning, disposal	78
14.1	Maintenance - Changing the lamp	78
14.2	Cleaning - Actions to take if a cell is broken	78
14.3	Disposal	79
15.	Technical Data	80
16.	What to do if...	82


1. Overview

1.1 Description of the operating elements



1.2 Identifying the connectors

- ① Power pack
- ② Connection for power pack
- ③ RS 232 interface



1. Overview


1.3 Short manual

The short manual lists all of the steps necessary to determine the concentration of a sample and to activate AQA2 at a glance.


1.3.1 Measuring the concentration

– To switch on the photometer, open the cover.

The photometer performs a check (*Self-Check*) of the entire system and then switches automatically to the *concentration measuring mode*.




Measuring mode, *concentration*



– Insert the round cell with barcode in the round cell shaft until it clicks into place.

Align the line mark to the notch of the photometer. The message *measuring... appears*.

i If the *select method* menu is displayed, align the line mark of the round cell to the notch of the photometer.




The measured value appears on the display. Measured values outside the specified measuring range are output in small numerals.

Repeat the measurement:

– Press .

1. Overview

1.4 Selecting and calling up the menu items

- To switch on the photometer, open the cover.
- Press .

The following display appears:

```
Setup
└─ documentation
  method parameter
  Meter Setup
```

The following display appears:

```
Setup
  documentation
  method parameter
  └─ Meter Setup
```

```
meter setup
  return
  └─ AQA functions
    correction funct.
    adjust zero
    set date/time
```


Example:

The *documentation* menu item is preselected in the *setup* menu (▶).



Select a menu item, e. g. *meter setup*:

- Press .




The *meter setup* menu item is preselected (▶).

- Call up the *meter setup* submenu by pressing .



The required menu item is

- selected using .
- called up using .

Selection lists:

- Changes to the settings are accepted after confirmation by pressing .
- Current settings are marked by "◆".
- Change to other configuration levels by
 - Selecting the menu item, *return*
 - Pressing .
- Scroll with .

Character input:

- by using ,
- the character to be input is shown in reverse video.
- Confirm each input with .

2. Safety

This operating manual contains basic instructions to be followed in the commissioning, operation and maintenance of the meter. Consequently, all responsible personnel must read this operating manual before

working with the meter. The operating manual must always be available in the vicinity of the meter.

2.1 Authorized use

The photometer was developed for use in the laboratory for water analysis. Observe the technical specifications of the cells according to chapter 15 TECHNICAL DATA.

Any other use is considered **unauthorized**.

2.2 General instructions

The photometer is constructed and tested according to the EN 61010-1 safety regulations for electronic measuring instruments. It left the factory in a safe and secure technical condition.

The smooth functioning and operational safety of the photometer can only be guaranteed under the climatic conditions specified in chapter 15 TECHNICAL DATA of this operating manual.

Opening the photometer or adjustment, maintenance and repair work must only be performed by personnel authorized by the manufacturer.

The only exceptions to this are the activities described in chapter 14 MAINTENANCE, CLEANING, DISPOSAL. Non-

compliance results in the loss of warranty claims. Follow the points listed below when operating the photometer:

- Follow local safety and accident prevention regulations.
- Observe the enclosed instructions concerning reagents and accessories.
- Observe the regulations when dealing with dangerous substances.
- Follow the operating instructions at the workplace.
- Use only original spare parts.

2.2.1 Labeling of notes



indicates notes that you must read – for your own safety, the safety of others and to protect your meter from being damaged.



indicates notes that draw your attention to special features.

2.2.2 Dangers of disregarding the safety instructions

Disregarding the safety instructions can adversely affect the safety of both the user and the environment as well as the equipment.

Non-compliance with the safety instructions will result in the loss of any warranty claims.

2.2.3 Qualification of the personnel

The personnel responsible for the commissioning, operation and maintenance must have the necessary qualifications for this work. If the personnel do not have the required skills they have to be instructed.

Furthermore, it must be ensured that the personnel read and completely understand the present operating manual.

2.2.4 Technical state of the meter

2. Safety

It is the responsibility of the operator to continuously observe the overall technical condition (externally recognizable deficits and damage as well as alterations to the operational behavior) of the meter. If safe operation is no longer possible, the equipment must be taken out of service and secured against inadvertent operation.

Safe operation is no longer possible if

- the equipment has been damaged in transport
- the equipment has been stored under adverse conditions for a lengthy period of time
- the equipment is visibly damaged
- the equipment no longer operates as prescribed.

If you are in any doubt, please contact the supplier of the photometer.

3. Commissioning

The photometer operates at an environmental temperature of +5 °C to +40 °C. During transport from cold to warm surroundings, condensation can form resulting in the malfunction of the meter.

Before putting the photometer into service, wait until it has adapted to the new environmental conditions (see also chapter 15 TECHNICAL DATA).

3.1 Preparing the photometer

- Place the photometer on a hard, flat surface and protect it against intensive light and heat.

Line operation


- Plug the original power pack into the socket on the photometer
- Plug the power pack into the line socket
- Switch on the photometer (open the cover).

Battery operation

- Charge the battery for approx. 5 hours before the initial commissioning. To do this:
 - Plug the original power pack into the socket on the photometer
 - Plug the power pack into the line socket and then the battery will be charged.


During battery operation or when the meter is at a standstill for longer periods of time, the battery runs down. This can result in your photometer no longer being ready for operation.

When the following symbol is displayed, charge the

battery: 


3.2 Switching on the photometer

- To switch on the photometer, open the cover. The photometer performs a check (*Self-Check*) of the entire system and then switches automatically to the *concentration* measuring mode.




Self-check of the photometer

After approx.
5 s:



Automatic change to the measuring mode,
concentration


4. Measuring the Concentration

- Call up the *concentration* measuring mode by actuating .

```
Concentration
-----
insert cell
```

Measuring mode, *concentration*

4.1 Display of concentration and absorption

- Press  to call up the *setup* menu
- Call up the *meter setup* submenu.

```
Meter Setup
-----
return
AQA Functions
Correction Funct.
adjust zero
▶ conc. / abs.
```

- In the *meter setup* menu, call up the *conc. / ext.* submenu.

```
conc. / abs.
-----
conc.
▶ conc. + abs.
return
```

Selection of the measuring mode:



- *conc.*
- *conc. + abs.*

Example

```
14729          P04-P
-----
          0.998 E
          12.3 mg/l
```

4. Measuring the Concentration


4.2 Measuring using cell tests



- Insert the round cell with barcode into the round cell shaft until it clicks into place.
Align the line mark to the notch of the photometer.

The photometer reads the barcode of the round cell and automatically selects the relevant method.

After approx.
2 s:



The measured value appears on the display.




If the *select method* menu is displayed, align the line mark of the round cell to the notch of the photometer.

4. Measuring the Concentration



4.3 Measuring using tests without barcode (manual method selection)

When measuring using cell tests without barcode, the method must be selected manually.



```
select method
-----
method: 186
                14729
                PO4-P
#          0.5-25.0 mg/l
```

The last method set up manually appears on the display.

- Select the required method with 
- Confirm with .

```
14729          PO4-P
-----
measuring...
```


After approx.
2 s:

```
14729          PO4-P
-----
          12.3 mg/l
```

The measured value appears on the display.

5. Measuring the Absorbance/Transmission

5.1 Switching to the Absorbance/ Transmission measuring mode

- Call up the *setup* measuring mode by actuating .

```
Setup
-----
documentation
method Parameter
▶abs./trm. %
meter setup
```


- In the *setup* menu, call up the *abs./trm. %* submenu.

```
absorbance +
transmission
return
```

Selection of the measuring mode:

- *absorbance*
- *transmission*

5.2 Measuring the absorbance or transmission

- Call up the *absorbance* or *transmission* measuring mode (depending on the selection in the *abs./trm. %* menu) by actuating .

```
Absorbance
-----
insert cell
```

Measuring mode, *absorbance*

```
transmission
-----
insert cell
```



Measuring mode, *transmission*

i The transmission measurement is not described separately in the following example as it operates in exactly the same way as the absorbance measurement. However, the result of the measurement is displayed as % Transmission instead of A for Absorbance.


i A measured reference absorbance is also effective in the measuring mode, *transmission*. It is displayed as reference absorbance.

5. Measuring the Absorbance/Transmission

5.3 Measuring using cell tests





- Insert the round cell with barcode into the round cell shaft until it clicks into place. Align the line mark to the notch of the photometer.



- The measured value for the wavelength displayed at the top right appears. This measured value is automatically stored. If necessary, call up further wavelengths:

– with .





- The measured value for the selected wavelength appears and can be stored and output to the interface with .



- Sample display for transmission measurement

5.4 Measuring using tests without barcode



- The last wavelength measured appears on the display.

- Select the wavelength:


with .

- Start the measurement:

Press .

6. Documentation

The measured values can be documented as follows:

- Storage in the measured value memory
 - Output to a connected printer via the serial interface (automatic when a printer is connected)
 - Transmission to a PC for further processing (by using the relevant software, e.g. Multi/ACHATII or – less conveniently – by means of a terminal program).
- To switch on the photometer, open the cover.
- Press .

The following display appears:

```
Setup
-----
▶documentation
  method parameter
  Meter Setup
```

```
documentation
-----
▶no. of meas. value
  download memory
  output methods
  return
```

- Call up the *documentation* menu with .

The following functions can be selected:

- *no. of meas. value*
 - reset the number
- *download memory*
 - total
 - from date
- *output methods*
 - all

The current settings are marked by "◆" in the selection lists of the respective submenus.

6. Documentation



6.1 Resetting the number of the measured value

```
documentation
├──no. of meas. value
│   ├──download memory
│   ├──output methods
│   └──return
```

```
no. of meas. value
reset number:
├──Yes   +
│   └──No
└──return
```

- Call up the *no. of meas. value* submenu.

- *yes*
The numbering of the measured values starts again with 001 (default)
- *no*
Consecutive numbering of the measured values (from 001 to 999)

- Select the menu item with 
- Confirm with 


6.2 Download memory

The measured value storage can be selectively downloaded to either the display or serial interface. The selection of the output medium is made after the specification of the sorting criteria.



```
documentation
├──no. of meas. value
│   ├──download memory
│   ├──output methods
│   └──return
```

```
download memory
├──total
│   ├──from date
│   └──return
```

- Call up the *download memory* submenu.

 The *download memory* menu item only appears after at least one measurement has been performed.

The following sorting criteria can be set:



- *total* – all stored measured values
 - *from date* – all measured values from a special date
- Select the menu item with 
 - Confirm with 

6. Documentation

Selecting "total"




```
download memory
└to display
  to printer/PC
  return
```

Select the output medium:

- *to display*
- *to printer/PC* (serial interface).
- Select the menu item with 
- Confirm with  to start the memory download.



Selecting "from date"

```
download memory
from date: 07.02.98
↓
```

- Input the date using 
- Erase the input using 
- Confirm with 

```
download memory
└to display
  to printer/PC
  return
```

Select the output medium:

- *to display*
- *to printer/PC* (serial interface).
- Select the menu item with 
- Confirm with  to start the memory download.


Memory download to display

```
download memory
009 07.05.97 17:24
Feed 14554 Ni
      3.66 mg/l
↓ return AQA2
```

Each data record appears individually on the display beginning with the data record just measured.

The display shows:

- no. of meas. value
- date/time
- I. D. number
- method designation
- citation
- meas. value
- unit
- Where necessary, AQA ID, e.g. AQA2.


- Scroll with 

6. Documentation

Memory download to printer/PC

```
download memory
data transmission
runs:
           121
↵ cancel
```

Memory download to the serial interface:

- Display of the transmitted no. of measured value (continuation display) beginning with the last measured value.
- Cancel with .

Sample printout:

003	14541	10.02.98	11:56:33	t	80	mg/l	COD
002	14541	10.02.98	11:54:21	t	70	mg/l	COD
001	14729	03.02.98	18:30:53	*	0.3	mg/l	PO4-P

6.3 Download of the methods list



The stored methods are downloaded to the printer/PC via the serial interface.

```
documentation
no. of meas.value
download memory
▶ output methods
return
```

- Call up the *output methods* submenu.

```
output methods
▶ all
return
```

The following parameters can be set:


- *all* – Download of all stored methods
- Select the menu item with .
- Start the download with .

7. Method Parameters

The following parameters can be set in the *method parameters* menu:

- citation
- unit

– To switch on the photometer, open the cover.

– Press .

The following display appears:


```
Setup
-----
documentation
▶method parameter
AQA-Check
Meter Setup
```

```
method parameter
-----
method: 286
                14729
                P04-P
↵ 0.5-25.0 mg/l
```

```
method parameter
-----
▶Citation
Dilution
return
```

– Call up the *method parameters* submenu.

– Input the method number

– Confirm with .

– Select the menu item with .

– Call up the parameter by pressing .

7. Method Parameters

7.1 Citation form

7.1.1 Changing the citation form

Example:

Change the citation form from "NH₄-N" to "NH₄".

```
method parameter
└─ Citation
  Dilution
  return
```


– Call up the *citation* submenu.

```
Citation 14739
└─ NH4-N ✦
  NH4
  return
```

The current setting: NH₄-N (✦).

```
Citation 14739
  NH4-N ✦
└─ NH4
  return
```

– Using , scroll to NH₄

– Confirm with .


```
Citation 14739
  NH4-N
└─ NH4 ✦
  return
```

– Citation form NH₄ is set (✦).

7. Method Parameters

7.1.2 Performing a difference measurement


Difference measuring is possible for some methods (e.g. Iron II/III, Ca-/Mg Hardness).

 For more information on this, see part, "Analysis specifications".

Example:

Determination of iron (II) and iron (III).

```
method parameter
-----
method: 106
                        14896
                        Fe
↓      1.0-50.0 mg/l
```



- Enter method 106
- Confirm with .

```
method parameter
-----
▶Citation
Dilution
return
```

- Call up the *citation* menu item.


```
Citation                14896
-----
▶Fe  +
FeII,FeIII  Δ
return
```


The current setting: Fe

- Using  scroll to Fe II, Fe III Δ
- Confirm with .

```
Citation                14896
-----
Fe
▶FeII,FeIII  Δ +
return
```

- Citation form Fe II, Fe III Δ (+) is set.

- Change to measuring by pressing .



```
14896  FeII,FeIII Δ
-----
Σ Fe
measuring...
```

- Start the 1st measurement by inserting cell 1.


7. Method Parameters


After approx.
2 s:

14896	FeII,FeIII Δ
3.2 mg/l	
Σ Fe	
\leftarrow FeII	

The 1st measured value appears on the display: Σ Fe.

– Remove cell 1

– Press .



14896	FeII,FeIII Δ
FeII	
measuring...	

– Start the 2nd measurement by inserting cell 2.

After approx.
2 s:

14896	FeII,FeIII Δ
2.1 mg/l	
FeII	
\leftarrow FeII,FeIII	

The 2nd measured value appears on the display:
Iron II.

– Continue to the display of both measured values

using .

14896	FeII,FeIII Δ
FeII	2.1 mg/l
FeIII	1.1 mg/l

Display of both measured values as a summary.

7. Method Parameters

7.2 Selecting the unit



The preset unit is "mg/l".
It can be changed to "mmol/l".

```
method parameter
-----
Citation
unit
return
```

- Call up the *unit* submenu.

```
unit 14729
-----
mg/l  +
mmol/l
return
```

The current setting: *mg/l* (⊕)

- Using  scroll to *mmol/l*
- Confirm with .

```
unit 14729
-----
mg/l  +
mmol/l
return
```

- Unit *mmol/l* (⊕) is set.

8. Analytical Quality Assurance (AQA)

Analytical quality assurance (AQA) can be performed in two steps:

- **AQA1** – Photometer monitoring
- **AQA2** – Total system monitoring with standard solutions



The total system monitoring (AQA2) is a method-specific check using standard solutions.

If this is performed successfully, it also includes photometer monitoring (AQA1).

See also part "General information" for further information on Analytical Quality Assurance (AQA).

The AQA mode must be activated in the photometer. In the delivery state it is switched off.

The AQA mode is activated by using a menu to select

- monitoring of the photometer (AQA1)
- monitoring of the total system using standard solutions (AQA2)

8.1 Activating AQA

- To switch on the photometer, open the cover.

8.1.1 Activating AQA via the menu guide

- Press

```
Setup
-----
documentation
method parameter
▶Meter Setup
```

- Call up the *meter setup* submenu.

```
meter setup
-----
return
▶AQA functions
correction funct.
adjust zero
set date/time
```

The *meter setup* submenu appears with the *AQA functions* menu item preselected.

- Confirm with

A password request appears:

```
AQA-Configuration
-----
input password:

  0 0 0 0
```

A separate password protects settings of the AQA-configuration against unauthorized access (Changing the password see section 8.1.5).

- Input the password with :
Only **numeric** characters are allowed.
Default: 0000
- Confirm with

If the input was incorrect:

8. Analytical Quality Assurance (AQA)

```
AQA-Configuration
-----
wrong Password
```

After the password has been successfully input, the *AQA configuration* submenu appears:

```
AQA-Configuration
-----
return
▶AQA-Mode
  AQA-Standards
  AQA-Intervals
  System locked
```

```
AQA-Mode
-----
▶off *
  n weeks
  n measurements
  return
```

```
Setup
-----
documentation
method Parameter
▶AQA-Check
  Meter Setup
```

```
AQA-Check
-----
  Meter
▶system
  return
```


- Repeat the input.



If you have forgotten the password, contact the service department.

- Call up the *AQA mode* function.

Default: *off* (no monitoring)

- Select *AQA mode*:
 - *off*
 - *n weeks*
 - *n measurements*
- Confirm with .
- In the *setup* menu, call up the *AQA check* submenu.

Selection of the *AQA mode*:

- *meter*
- *system*



The menu item, *meter*, only appears after the corresponding PhotoCheck standards have been input (see section 8.2.1).

8. Analytical Quality Assurance (AQA)

8.1.2 Changing AQA intervals

AQA intervals specify the interval between two AQA checks. A fixed time interval (*n weeks*) or a number of measurements (*n measurements*) can be specified as the interval.

The respective values that were input remain stored even if they are not activated.

Additionally, two separate intervals can be set up for both photometer monitoring (AQA1) and system monitoring (AQA2).



For the total system monitoring (AQA2), a change of the time interval (*n weeks*) even retroactively applies to monitoring processes that are already running.

Changing the number of measurements (*n measurements*) does not affect monitoring processes already running.

Thus, individual numbers of measurements can be set for different methods.

```
AQA-Configuration
-----
return
AQA-Mode
AQA-Standards
▶AQA-Intervals
System locked
```

After an interval has expired, the following consequences become effective:

- Warning and loss of AQA identification
- Locking of the method for concentration measurements (as long as the locking is active).

Setting ranges:

- Photometer monitoring (AQA1):
 - 1 to 52 weeks (default: 12 weeks) or
 - 1 to 9999 measurements (default: 1500)
- Monitoring of the total system using standard solutions (AQA2):
 - 1 to 52 weeks (default: 4 weeks) or
 - 1 to 9999 measurements (default: 100)



With the *n measurements* setting, a difference measurement (see section 7.1.2) is counted as one measurement only.

- In the *AQA configuration* menu, call up the *AQA intervals* submenu. According to the selection in the *AQA mode* menu, a fixed time interval (*n weeks*) or a number of measurements (*n measurements*) is set in the *AQA intervals* menu.



If the *AQA mode* function is switched off, the *AQA intervals* submenu is not available.

AQA interval, "n weeks"

The AQA interval, *n weeks*, is only effective if the *n weeks* setting is active for the *AQA mode* function.

The specified number of *n weeks* applies to:

- the photometer with AQA1
- all methods with AQA2.

- In the *AQA intervals* menu, call up the *n weeks* submenu.

```
AQA-Intervals
-----
AQA-Meter: 12 w
AQA-System: 04 w
← confirm
```

- To return without change, press three times
- Enter the time interval for *AQA meter* with , confirm with .

8. Analytical Quality Assurance (AQA)

AQA interval, "n measurements"

The AQA interval, *n measurements*, is only effective if the *n measurements* setting is active for the *AQA mode* function.

The AQA2 check starts the monitoring for one method at a time.

The specified number of *measurements* applies to:

- the instrument with AQA1 (total number of measurements performed, independent of whether AQA2 is active for some parameters)
- each method an AQA check will then be performed for with AQA2.
Thus, it is possible to define individual numbers of measurements for different methods.
The measurements are counted separately for each monitored method.

The monitoring intervals of AQA2 monitoring processes already started for other methods are not affected by changing the number of *measurements*. Thus the number of *measurements* can be set for further methods no matter which monitoring processes were started before.








When an AQA2 check is performed, the number of *measurements* last set in the *AQA intervals* menu is automatically taken over.

Therefore, you should check and, if necessary, change the currently set number of *measurements* before each AQA2 check.

The currently set number of *measurements* for the AQA2 check is saved for the active method and output in the report individually (section 8.3.4).

- In the *AQA intervals* menu, call up the *n measurements* submenu.

```
AQA-Intervals
-----
AQA meter:
  1500 measurements
AQA system:
  0100 measurements
← confirm
```

- To return without change, press  three times
- Enter the number of measurements for *AQA meter* with , confirm with 
- Enter the number of measurements for *AQA system* with , confirm with 

8. Analytical Quality Assurance (AQA)

8.1.3 Locking the system

The function *system locked* is effective if, for a monitored method,



- no AQA check was performed,
- the AQA check “system” has expired.

As a result, a concentration measurement is not possible for this method.

```
AQA-Configuration
return
AQA-Mode
AQA-Standards
AQA-Intervals
▶System locked
```

- Call up the *system locked* submenu.

```
System locked
off
▶on      +
return
```


- Select the menu item with 
- Confirm with .

8. Analytical Quality Assurance (AQA)



8.1.4 Changing the password

When delivered, the default password is 0000. This password can be changed as follows:


```
AQA-Configuration
-----
AQA-Standards
AQA-Intervals
System locked
change Password
reset
```

- Call up the *change password* submenu.
- Confirm with .

```
AQA-Password
-----
input password:
(0000)
  0 0 0 0
```

- Input the required password, e.g. 0100, with .
- Confirm with .


```
AQA-Password
-----
confirm Password:
(0100)
  0 0 0 0
```

- Input the password once again:
- Confirm with .

8.1.5 Performing an AQA reset


If the Analytical Quality Assurance is to be switched off completely or reset to the delivery state, this can be made via the *reset* function in the *AQA configuration* submenu.

```
AQA-Configuration
-----
AQA-Intervals
System locked
change Password
reset
return
```

- Call up the *reset* submenu
- Confirm with .

```
AQA-Configuration
-----

reset
cancel
```

- Select the *reset* menu item
 - Confirm with .
- An AQA reset is performed.

8. Analytical Quality Assurance (AQA)

8.2 Photometer monitoring (AQA1)

8.2.1 Entering PhotoCheck standards



A Spectroquant® PhotoCheck is required to perform the photometer monitoring (AQA1). **At least 1 standard** must be input. We recommend, however, to input all available standards.

- Press to call up the *setup* menu
- Call up the *meter setup* submenu.
- Call up the *AQA functions* submenu.
- Input the password
- Call up the *AQA standards* submenu and the following display appears:

```
AQA-Standards
├─PhotoCheck
│  standard solution
│  return
```

- Call up the *PhotoCheck* submenu.

```
PhotoCheck-Standards
├─input
│  output
│  erase
│  return
```

Select between

- *input*
Input the theoretical value (absorbance) from the lot certificate of Spectroquant® PhotoCheck
- *output*
Print/display theoretical values
- *erase*
Erase theoretical values.





The *erase* and *output* menu items only appear after at least one standard has been input.

8. Analytical Quality Assurance (AQA)


Example:

445-1 nm, theoretical value (absorbance) 0.200,
admissible tolerance ± 0.020

```
PhotoCheck-Standards
return
▶445-1
  445-2
  445-3
  445-4
```



- Select with 
- Quit via the menu item, *return*
- Confirm with .

```
PhotoCheck      445-1
theor.val.: 0.200 A
↵confirm
```

- Input the theoretical value, 445-1
- Confirm with .


If the standard is already stored, this value appears on the display.

```
PhotoCheck      445-1
theor.val.: 0.200 A
Tolerance: ±0.020 A
↵confirm
```

- Input the tolerance with 
- Confirm with .

```
PhotoCheck-Standards
return
▶445-1 ✓
  445-2
  445-3
  445-4
```

PhotoCheck standard 445-1 is input.

- Select the next one with 
- Input all PhotoCheck standards in this way.

8. Analytical Quality Assurance (AQA)



8.2.2 Download of PhotoCheck standards

```
PhotoCheck-Standards
input
output
▶output
erase
return
```

- In the *PhotoCheck standards* submenu, call up the *output* menu item.

```
download PhotoCheck
▶to display
to printer/PC
return
```

Select the output medium:

- to display
- to printer/PC (serial interface).
- Select with 
- Confirm with  to start the download.

Example: Report output

AQA check meter
26.08.97

AQA1
13:19

AQA interval

12 weeks

test sol.	unit	theor. val.	tolerance	AQA date
445-1	A	0.200	0.020	26.08.97

8.2.3 Erasing PhotoCheck standards



At least 1 standard must still be stored to be able to perform the AQA check function (meter monitoring).

```
PhotoCheck-Standards
input
output
▶erase
return
```

- In the *PhotoCheck standards* submenu, call up the *erase* menu item.

```
erase PhotoCheck
445-2
445-3
▶445-4
return
```

Displays the stored PhotoCheck standards:

- Select with 
- Quit via *return*
- Erase with .


8. Analytical Quality Assurance (AQA)

8.2.4 Performing Photometer monitoring

Photometer monitoring (AQA1) includes a check of the


- Light barriers using the L1/L2 cells (contained within the scope of delivery of the Spectroquant® PhotoCheck)
- Absorbance measurement using PhotoCheck

standards.

- Press  to call up the *setup* menu
- Call up the *AQA check* submenu
- Call up the *meter* submenu.

The following display appears:

```
L-Check
-----
use L1
      ↵ cancel
```




After approx.
1 s:

```
L-Check
-----
L1 ok
```

After approx.
4 s:

```
L-Check
-----
use L2
      ↵ cancel
```



After approx.
1 s:

```
L-Check
-----
L2 ok
```

- Insert the L1 cell.



If the *error* message appears, clean the cell shaft with a damp, lint-free cloth and repeat the check.


If the message reoccurs, inform the service department.

- Insert the L2 cell.

8. Analytical Quality Assurance (AQA)

After successful light barrier testing, the PhotoCheck standards (test solutions) are measured.

Example:



```
PhotoCheck      445-1
-----
use test
solution      445-1
↓ cancel
```

- Insert a cell with the test solution, 445-1. The photometer measures the absorbance of the test solution and compares the result with the value entered.

Absorbance test OK...

After approx.
3 s:

```
PhotoCheck      445-1
-----
0.211 A
ok
```

- Insert the next test solution
- Cancel:
To cancel the check means no release for the next "meter" AQA interval!

...or error message

```
PhotoCheck      445-1
-----
---- A
Error
```

Error elimination:

1. Repeat the measurement (insert the cell again)
2. If necessary, perform a zero adjustment and repeat the check
3. Exchange the test solution (each packet contains two identical test solutions)
4. Use a new Spectroquant® Photo-Check packet
5. Quit and have the photometer checked in the factory

The absorbance test is terminated if an error message occurs and the meter is **not released**. On switching on, the warning message "AQA interval expired" appears until the AQA was successfully performed or the AQA mode was switched off.

Example: Report output

AQA check meter	26.08.97	operator:	AQA1	10:23	
AQA interval	AQA check AQA1	L check	12 weeks	ok	ok
test sol.	meas. value	unit	theor. val.	tolerance	result
445-1	0.211	A	0.200	0.020	ok

8. Analytical Quality Assurance (AQA)

8.3 Total system monitoring with standard solutions (AQA2)

8.3.1 Entering standards



The standards compiled in the table "Spectroquant® CombiCheck and standard solutions" (see part "General information") are already stored method-specifically in the photometer. These values can be overwritten.

For **total system monitoring** (AQA2), only one standard per test can be stored at a time. The input of a standard is only complete with the input of the tolerances for finding it again, i.e. it is then first stored (no premature quitting).

- Press to call up the *setup* menu
- Call up the *meter setup* submenu.
- Call up the *AQA functions* submenu
- Input the password
- Call up the *AQA standards* submenu and the following display appears:

```
AQA-Standards
PhotoCheck
▶standard solution
return
```

```
standard solution
▶input
output
erase
return
```

```
input standard
method: 386
14729
PO4-P
↵ 0.5-25.0 mg/l
```

- Call up the *standard solutions* submenu.

Select between

- *input*
Enter standards
- *output*
Print/display standards
- *erase*
Erase standards.

Displays the last selected method.

- Select the method with
- Confirm with
- Input the standards.


8. Analytical Quality Assurance (AQA)

Example:


Method 14729 with a preset theoretical value of 15.0 mg/l and tolerance of 1.0 mg/l (CombiCheck 80).

Change to: theoretical value = 8 mg/l, tolerance = 0.7 mg/l (CombiCheck 20).

```
input standard
method: 86
14729
PO4-P
0.5-25.0 mg/l
```

– Confirm with .


```
standard 14729
theor.val.: 15.0 mg/l
(06.3-18.8 mg/l)
confirm
```


– Enter the new theoretical value, e.g. 8.0 mg/l, with .

Values in parentheses indicate the range in which the theoretical value should move.

– Confirm with .


```
standard 14729
theor.val.: 08.0 mg/l
Tolerance: ±1.0 mg/l
confirm
```

– Input the tolerance (0.7 mg/l) with .

– Confirm with .

```
standard 14729
theor.val.: 08.0 mg/l
Tolerance: ±0.7 mg/l
confirm
```

Both standard and tolerance values have been overwritten.


– Confirm with .

8. Analytical Quality Assurance (AQA)

8.3.2 Output of standards

The current list of stored standards is output via the RS 232 interface (PC/printer) or via the display.

```
standard solution
input
output
erase
return
```

- Select the *output* submenu
- Confirm with .

```
download Standard
to display
to printer/PC
return
```

Select the output medium:

- *to display*
- *to printer/PC* (serial interface).

- Select with 
- Confirm with  to start the download.

Example: Report output

AQA check system		AQA2		
26.08.97		13:57		
system locked		on		
method	unit	theor. val.	tolerance	AQA date
14554	mg/l	2.00	0.20	24.08.97
14555	mg/l	5000	400	26.08.97

8. Analytical Quality Assurance (AQA)



8.3.3 Erasing standards

Erasing the method-specific standard solutions leads to the change of the measured value identification from AQA2 to AQA1 (with activated AQA mode).



```
AQA-Standards
PhotoCheck
▶standard solution
return
```

- Call up the *standard solutions* submenu.

```
standard solution
input
output
▶erase
return
```

- Select the *erase* menu item, *erase* with 
- Confirm with 

```
erase standard
▶14560
14729
return
```


- Select the standard to be erased with 
- Erase with 

8.3.4 Monitoring of the total system using standard solutions

The AQA2 check can be performed after it has been activated (see section 8.1).

The following display appears:

```
AQA-Check
insert cell
```



- i** For AQA2 with the setting, *n measurements*, we recommend to check and, if necessary, change the currently set number of *measurements* before each AQA check (8.1.2 CHANGING AQA INTERVALS).

- Insert cell with prepared solution ready to be measured (e.g. using Spectroquant® Combi-Check). The photometer reads the barcode, identifies the method and performs the AQA2 check.

8. Analytical Quality Assurance (AQA)

After approx.
2 s:

AQA check OK ...

AQA-Check	14554
1.93 mg/l ok	

...or error message

AQA-Check	14554
3.45 mg/l Error	

- Repeat the check
If the error is repeated, perform troubleshooting of the error. See "Analytical Quality Assurance" in part "General information".



The *system* AQA2 check must be performed **separately** for each method monitored. The release is stored with the date and the specified interval. The AQA2 interval *system* set up for the respective method begins again.


Example: Report output (AQA mode: *n weeks*)

AQA check system			AQA2		
26.08.97			11:02		
operator:					
AQA interval			4 weeks		
method	meas. value	unit	theor. val.	tolerance	result
14554	1.95	mg/l	2.00	0.20	ok

Example: Report output (AQA mode: *n measurements*)

AQA check system			AQA2		
26.08.97			11:02		
operator:					
AQA interval			100 measurements		
method	meas. value	unit	theor. val.	tolerance	result
14554	1.95	mg/l	2.00	0.20	ok

9. Correction functions

- To switch on the photometer, open the cover.
- Press .
- In the *setup* menu, call up the *meter setup* submenu.
The following display appears:


```
Meter Setup
-----
return
AQA Functions
▶Correction Funct.
adjust zero
set date/time
```

- Call up the *correction funct.* submenu.

The following display appears:

```
Correction Funct.
-----
▶Blank Value
Turbidity Correct.
return
```

Select the correction function:

- blank value
- turbidity correct.
- Confirm with .

9.1 Blank value

The blank value (= reagent blank value) for each method is stored in the photometer. When the *blank value* function is active, the stored value is ignored and the measured value of a self-prepared reagent blank solution is used instead.

This procedure increases the measuring accuracy for some tests (for more information, see part "Analytical procedures").

A blank value is always stored for the method that was just called up.

A blank value remains stored until it is erased (menu item, *erase blank value*) or overwritten.

The *reset setup* function sets the *blank value* to *off*. The stored blank values, however, remain stored.


The *reset total* function resets all settings and blank values at once.


If a measured blank value is stored and the *blank value* function is active for a method, this blank value is used for determining the measured value and the measured value is documented accordingly.

The *blank value* function is not active when delivered.

9. Correction functions

Measuring the concentration with a blank value


- Press  to call up the *concentration* measuring mode.





The value measured against the prepared blank solution is displayed.

9.1.1 Activating the blank value measurement

- In the *correction funct.* menu, call up the *blank value* submenu. The following display appears:




The *blank value meas.* function appears:

- Select the *on* menu item with 
- Confirm with .

i The stored blank values determined from blank solutions prepared by the user can be deactivated by switching off the blank value measurement. When doing so, the blank values remain stored in the memory and can be reactivated later.

Activating or deactivating the blank value function applies to all measurements using methods a blank value was stored for in the memory.

The *blank value* function is active and appears in the *setup* menu:



- To measure the blank value, call up the *blank value* submenu in the *setup* menu.

9. Correction functions

9.1.2 Measuring the blank value

```
Blank Value
┌───────────┐
│ meas. blank value  
│ erase blank value  
│ recall blank values  
│ return
```

- Call up the *meas. blank value* menu item.

i The menu items, *erase blank value* and *recall blank values* first appear after at least one blank value has been measured.

```
meas. blank value
┌───────────┐
│ insert cell
```

- Insert a cell with blank solution to start a measurement.
The message, *measuring...*, appears on the display.


After approx.
2 s:

```
Blank Value 14729
┌───────────┐
│ 0.033 A  
└───────────┘
│ ↵ return
```

9.1.3 Erasing blank values

A measured blank value is erased via the menu item, *erase blank value*.

```
Blank Value
┌───────────┐
│ meas. blank value  
│ erase blank value  
│ recall blank values  
│ return
```

- Select the *erase blank value* menu item
- After confirming with  the *erase blank value* menu opens.

i The *erase blank value* menu item first appears after a blank value has been measured.

```
erase blank value
┌───────────┐
│ all  
│ single  
│ return
```

Select between





- *all*
Erase all stored blank values
- *single*
Erase individual stored blank value

9. Correction functions

```
erase blank value
12.01.2004 14:57
14758
0.100 A
←
```

```
erase blank value
12.01.2004 14:57
14758
erase
cancel
```





i Each stored blank value is displayed with the date of the blank value measurement and the relevant method designation.

- Select the blank value with 
- Erase the displayed blank value with 
- Select the *erase* menu item with 
- Confirm with 

9.1.4 Recalling blank values

```
Blank Value
meas. blank value
erase blank value
recall blank values
return
```

```
recall blank values
12.01.2004 14:57
14758
0.100 A
←return
```

- Select the *recall blank values* menu item with 
- Confirm with 
- Select the blank value with 
- Return with 

9. Correction functions

9.2 Turbidity correction

Turbidity correction is used in sample solutions that contain finely distributed suspended particles. The suspended particles cause a light absorption. This leads to incorrect (too high) measured values. The function remains permanently switched on after it has been activated. Values that were measured using turbidity correction are given an identifier in the **display** and in the **documentation** (printout and storage).


The *turbidity correct.* function is not active when delivered.



This function is not necessary, or useful, in all methods. If the turbidity correction is active, the photometer automatically decides whether to perform the function or not depending on the method.

- In the *correction funct.* menu, call up the *turbidity correct.* submenu.


The following display appears:



The *turbidity correct.* function appears:


- Select the *on* menu item with
- Confirm with

- Press to call up the *concentration* measuring mode.



- Insert the measuring cell.


After approx.
2 s:



Display of the measured value with turbidity correction switched on: Identified by *TC*.

Warning of excessive turbidity:

If the turbidity absorbance of $0.100 A$ is exceeded, the meter displays the measured value together with a warning.



10. Zero adjustment

Zero adjustment is necessary

- after changing the lamp
- after the error message, *PhotoCheck* (AQA1) occurs
- on initial commissioning
- if the photometer was mechanically stressed, e.g. percussion, transport
- if the ambient temperature changed by more than 5 °C since the last zero adjustment
- at least every six months.



Only perform the zero adjustment against distilled water in an optically perfect cell.

- Press
 - In the *setup* menu, call up the *meter setup* submenu.
- The following display appears:

```
meter setup
-----
return
AQA functions
correction funct.
▶adjust zero
set date/time
```

```
adjust zero
-----
insert cell
```

After approx.
2 s:

```
adjust zero
-----
round ok
```

When performing the zero adjustment observe the following points:


- Only use a clean, scratch-free round cell with distilled water. A prepared zero cell is provided with your photometer. In addition, a prepared zero cell is contained in the scope of delivery of the *PhotoCheck* (article 14693).
- If the round cell is visibly contaminated, or at least every 24 months, clean and refill it (minimum filling level 20 mm). Then check the cell for scratches.

- Call up the *zero adjustment* submenu with

- Insert a cell with distilled water.
- The message, *measuring...*, appears on the display.

Successful zero adjustment

11. Meter Setup

- To switch on the photometer, open the cover.
- Press 
- In the *setup* menu, call up the *meter setup* submenu.
The following display appears:

```
meter setup
-----
return
▶AQA functions
correction funct.
adjust zero
set date/time
```

This chapter describes four functions of the *meter setup* menu:

- *select language*
- *set date/time*
- *Performing a meter reset*
- *system info*

11.1 Selecting the language

The following languages are stored in the photometer:

- Deutsch (German)
- English
- Français (French)
- Italiano (Italian)
- Português (Portuguese)
- Polski (Polish)
- Dansk (Danish)
- Svenska (Swedish)
- Español (Spanish)
- Nederlands (Dutch)
- Indonesia (Indonesian)
- Čeština (Czech)
- Magyar (Hungarian)
- Russkij (Russian)
- Türkçe (Turkish)
- Brasil (Brazilian)



This is the order in which the available languages appear in the *select language* menu.

The available languages are listed in the language of the respective country in the photometer.



When *Russkij* is selected as the language, the Cyrillic alphabet is used for the user guidance. Method designation and ID numbers are always displayed in Latin script.

For output to the RS 232 C interface, Cyrillic characters are converted to Latin characters according to GOST.

```
meter setup
-----
correction funct.
adjust zero
set date/time
▶select language
system info
```

- Call up the *select language* menu item.

```
select language
-----
return
Deutsch
▶English
Français
Italiano
```

- Select a language, e.g. English
- Confirm with 
- Press the  key again:
Return to the *meter setup* submenu.
The displays appear in English.


11. Meter Setup

11.2 Setting the date/time


```
Meter Setup
-----
AQA Functions
Correction Funct.
adjust zero
▶set date/time
select language
```

```
Date/Time
-----
Date      01.01.98
          (dd.mm.yy)
Time      16:45
          (hh:mm)
◀ confirm
```

– Call up the *set date/time* menu item.

– Input the date using 


– Confirm with 

– Input the time with 

– Confirm with .

11.3 Reset

It is possible to reset the photometer to its factory settings (delivery state) in single steps. The *reset total* function resets all settings and blank values at once.

 All AQA functions are retained when *meter setup* is used. See section 8.1.5 for AQA reset.

```
Meter Setup
-----
set date/time
select language
system info
▶reset
return
```

```
reset
-----
▶total
meas.storage
Setup
return
```

– Call up the *reset* menu item.

Select between

- *total*
Erase the measured value storage and reset the settings to the delivery state
- *meas. storage*
Erase the measured value storage
- *setup*
Reset all settings to the delivery state.


Example: Performing a total reset

```
reset
-----

reset total

▶reset
cancel
```

– Select the *reset* menu item

– Confirm with .

A meter reset is performed (measured value memory and setup).

11. Meter Setup

11.4 System info

```
Meter Setup
-----
adjust zero
set date/time
select language
▶system info
reset
```

– Call up the *system info* menu item.

```
Meter Setup
-----
Software: 3.15
methods: 33.00

◀ return
```

Sample display

12. Updating method data

You will always find the latest method data for your photometer on the Internet. A method update contains all new test sets and methods respectively. Additionally, minor modifications of already existing methods are transferred with it. With a method update, you receive all new methods and, at the same time, can easily and conveniently update all method data.

The software provided for downloading contains the program file and method data. It can be downloaded from our homepage with a mouse click.

The files are packed in a self-decompressing archive file (*.exe) or in a zip file (*.zip) and can be decompressed after the download.

Carry out the update as follows:

To download and update the photometer method data via the built-in RS232 interface, you need the following:

- PC (Win 95 or higher) with Internet connection
- PC cable (available as an accessory)
- An *.exe or *.zip file from the Internet; contains the "UpdateMethodData.exe" program file and 6 method data files (pls6md.xxx, pls12md.xxx, plspekmd.xxx, nova30md.xxx, nova60md.xxx, nova400md.xxx; xxx = version).

- Switch on the photometer (open the cover).
- Switch on the PC.
- Download from the Internet the software including the method data (*.exe or *.zip) and copy it into a separate directory or on a floppy disk.
- Decompress the *.exe file with a double-click or decompress the *.zip file with Winzip.
- Connect the serial interfaces of the PC and photometer with the cable.
- Start the "UpdateMethodData.exe" program file by double-clicking. The "Update Method Data" window appears. In the upper half of the window there is the name of your photometer (among other things), behind it there is the method version in brackets (e.g. 8.00).



All method data are reloaded into the photometer with the update. The old method data are overwritten by this.

- Click on the "Search meter" button. The program automatically recognizes the connected photometer. Another "Update Method Data" window appears.
- Click on the "Start" button to start the method download. The process takes approx. 3 minutes. You can terminate it at any time by clicking on the "Cancel" button. In this case, however, the download has to be carried out once again completely so that the photometer can save the method data and is operative.

12. Updating method data

During the download, the following display appears on the photometer screen:

```
remote
```

- After the download, confirm the "Data successfully downloaded" message. The download is finished. The photometer returns to the *concentration* measuring mode.



You can check whether the new method data are stored in the photometer. To do so, proceed as follows:

```
Meter Setup
-----
adjust zero
set date/time
select language
system info
reset
```

- In the setup menu, call up the meter setup submenu.
- Call up the *system info* menu item.

```
Meter Setup
-----
Software: 3.15
methods: 33.00
return
```

Sample display (the software version is irrelevant here). The method version (here: 33.00) has to agree with the method version for your photometer in the "Update Method Data" window during the download.

Error messages

Message	Meaning	Remedy
No meter found	Connection PC - photometer out of order or not available	<ul style="list-style-type: none"> – Tightly connect the cable to the serial interfaces of the PC and photometer. – Use the correct cable
	Photometer not recognized	<ul style="list-style-type: none"> – Select the photometer manually

13. RS 232 C interface

Via the interface, data can be

- output to a printer and
- exchanged with a personal computer (PC)

For this, the following items are available as accessories:

- Printer cable
- Printer
- Interface cable
- Communication software.

13.1 Principle course of the remote control

String to meter	Reply from meter	Operating mode
S <CR>	> <CR>	Remote (remote control)
Command xx (see 15.2 command list)	Reply string command xx <CR>	Remote (remote control)
.		
.		
.		
CLOC <CR>		Concentration measurement



The keyboard of the photometer is locked in the *remote* operating mode.

13.2 Command list

Command	Function
S	Begin communication
CLOC	Switchover to normal operation (concentration measurement)
CDAT [anz]	Reads out stored measured values; [anz] = number of the measured values to be output
CMES [MMM]	Measurement and transmission of the concentration value with date/time; [MMM] = method number (e.g. 086 for method 14729)
CEXT [LLL]	Measurement and transmission of the absorbance value for the wavelength; [LLL] = wavelength
CBLA [MMM]	Measurement and transmission of the sample blank value; [MMM] = method number
CCLB [MMM]	Erase measured sample blank values; [MMM] = method number



The error message, *Invalid command*, appears if commands are unknown or cannot be carried out (e. g. if optional parameters do not agree with the cell coding). Optional parameters [MMM] and [LLL] need only be input for uncoded cells.

13. RS 232 C interface

13.3 Output format of measured values

Character	Meaning
3	consecutive number (not required for interface commands CMES, CEXT and CBLA)
5	method designation
6	I. D. number
17	date and time
4	special characters
9	meas. value
10	unit
12	citation
4	AQA ID (AQA2/AQA1)

Notes:

Data fields are separated by spaces.
Character set: IBM, code page 437


Meaning of the special characters:

! = Measuring with blank value (concentration) or reference absorbance (absorbance)
t/T = Measurement with turbidity correction/with high turbidity
* = Measured value outside the measuring range
Q = AQA measurement

13.4 Data transmission

Baud rate	4800
Data bits:	8
Stop bits:	1
Parity:	none
Handshake:	Hardware
Max. cable length	15 m


13.5 Pin assignment



Photometer 9-pin socket	Computer 9-pin socket	25 pin plug	Printer with RS 232 C interface
1	4	20	-
2	3	2	TXD
3	2	3	RXD
4	1 and 6	6	-
5	5	7	SG
6	4	20	-
7	8	5	-
8	7	4	DTR (if not available: short-circuit CTS and RTS)
9		-	-

14. Maintenance, cleaning, disposal

14.1 Maintenance - Changing the lamp



- Switch off the photometer and disconnect it from the power line
- Carefully turn up the photometer and park it safely
- Screw off the lamp cover on the underside of the photometer



Let the lamp of the photometer cool down.


- Pull out the plug ①
- Unscrew the screw ②
- Remove the lamp with its holder ③ by pulling it gently upwards



Do not touch the new light bulb of the photometer.

- Insert a new preset lamp and screw it tight using the screw ②

- Connect the plug ① of the new lamp
- Screw the lamp cover on again
- Set up the photometer again and connect it to the power line
- Press and hold
- Switch on the meter (open the cover) and after the following display appears, release



- Press

14.2 Cleaning - Actions to take if a cell is broken



Do not rotate the photometer to pour out the liquid!

The photometer has a draining mechanism under the cell shaft that, when operated correctly, prevents any liquid coming into contact with electronic components.

- Switch off the photometer (close the cover) and disconnect it from the line power
- Let the liquid drain off
- Carefully remove any pieces of glass, e.g. using tweezers
- Carefully clean the cell shaft with a damp, lint-free cloth

- Let the cell shaft dry
After it is dry, check the photometer:
- Perform a photometer monitoring (see section 8.2).

▶ 14. Maintenance, cleaning, disposal

14.3 Disposal

Packing

The measuring instrument is sent out in a protective transport packing.

We recommend: Keep the packing material in case you have to send the measuring instrument back for service. The original packing prevents the measuring instrument from being damaged during transport.

Measuring instrument

Dispose of the measuring instrument as electronic waste at an appropriate collection point. It is illegal to dispose of it in household refuse.

Within the European Union, the batteries are removed at a specialized treatment center at the instrument's end of life. The instruments are taken to one of those specialized treatment centers via the recycling system set up for this purpose.

15. Technical Data

Optical measuring principle	Filter photometer with reference beam absorption measurement; simultaneous recording of all wavelengths
Light source	Tungsten halogen lamp, preset
Receiver	6 x photo diode array
Optical filters	340 nm, 445 nm, 525 nm, 550 nm, 605 nm, 690 nm, Accuracy: ± 2 nm; Half width: 340 nm = 30 nm ± 2 nm; all others = 10 nm ± 2 nm
Photometric reproducibility	0.001 A at 1.000 A
Photometric resolution	0.001 A
Warm-up time	none
Measuring time	approx. 2 s
Types of measurement	Concentration (method dependent, selectable display form), absorbance, transmission
Measuring range absorbance	-0.300 A to 3.200 A
Measuring range transmission	0.1 % to 1000 %
Balancing	Permanently stored
Drift correction	Automatic on each Self-Check
Retrofitting of new methods	via the Internet
Bar code recognition	automatic selection of the method; automatic recognition of the reagents lot
Cell recognition	automatic
Self-Check	<i>Test:</i> Memory, optics, electronic measured value recording, barcode recognition, cell recognition <i>Automatic calibration:</i> Optics, electronic measured value recording, barcode recognition
Time/Date	Real-time clock in the photometer
Dimensions	H: 140 mm, D: 270 mm, W: 260 mm

Weight	approx. 2.3 kg (battery version: 2.8 kg)
Guidelines and norms used	are defined in a separate document: Declaration of Conformity
Power pack	FRIWO FW 7555O/15 Friwo Part. No. 1822367 Input: 100 ... 240 V ~ / 50 ... 60 Hz / 400 mA Output: 15 V DC / 1 A
Power consumption in line operation	max. 1300 mA
Batteries	
● Backup battery	1 x 3,0 V Lithium battery, soldered in the instrument
● Battery (optional)	Built-in battery: NiMH rechargeable battery 7.2 V/2500 mAh, operating time with new, fully charged battery: typical 40 hours with 10 measurements per hour, trickle charging in line operation, approx. 5 h charging time for a discharged battery, total discharge protection
Climatic class	2, VDI/VDE 3540
Ambient temperature	Storage: -25 °C to +65 °C Operation: +5 °C to +40 °C
Allowable relative humidity	Annual mean: 75 % 30 days/year: 95 % other days: 85 %
Test certificate	CE
Operating elements	On/off switch actuated by opening/closing the lid of the cell shaft cover Silicon keyboard with 4 function keys Cell shaft – for round cells (flat cell floor, external/internal diameter 16 mm / 13.8 mm) –
Display	Graphical display 128 x 64 pixels
Connections	
● Digital interface	RS 232 C 9-pin socket to connect to PC or printer

15. Technical Data

- **Power supply** 2-pin socket to connect the plug-in power supply unit

Data storage Cyclical memory to record 500 measured values

Equipment Statement

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Software settings when delivered

Measured value number:	1
blank value is:	Off
turbidity correct.:	Off
language:	English
Date of the last valid AQA1 check:	invalid (not yet measured)
AQA1 interval:	12 weeks
AQA2 interval:	4 weeks
AQA password:	0000
AQA mode:	Off
Lock measurement if AQA2 expired:	Off
Checks to be measured with AQA1:	none
AQA2 values:	none

Settings after reset - total

Measured value storage and setup reset

Settings after reset - meas. storage

Meas. value number:	1
Measured values:	none


Settings after reset - setup

Measured value number:	1
blank value:	Off
reference absorbance:	Off
turbidity correct.:	Off
Language:	unchanged


Settings after reset - AQA

Date of the last valid AQA1 check:	invalid (not yet measured)
AQA1 interval:	12 weeks
AQA2 interval:	4 weeks
AQA password:	0000
AQA mode:	Off
Lock measurement if AQA2 expired:	Off
Checks to be measured with AQA1:	none (Input theoretical values and tolerances are not erased and are offered again with the next input).
AQA2 values:	none (theoretical values and tolerances of all methods are set to default values according to the "Spectroquant [®] CombiCheck and standard solutions" table in the part "General information".)

16. What to do if...

The display remains blank when switched on	Connect the photometer to the line power via the power pack. In the case of battery operation: Battery empty, charging required (approx. 5h); line operation is possible without restrictions during charging time.
 appears	Battery nearly empty. Charging required (see chapter 3 COMMISSIONING).
Date/time is lost when switched off	The backup battery of the real time clock is empty and has to be replaced. Send the photometer to the service department for this.
Password forgotten	Inform the service department.
Photometer does not react	The connected printer is off line. Switch on the printer or pull out the interface cable.
Error messages:	
<i>remove cell</i>	The message remove cell appears on the display although no cell is inserted. Clean the cell shaft with a damp, lint-free cloth. If the error message still appears, return the photometer to the service department.
<i>lamp defective</i>	Replace the lamp (see chapter chapter 14 MAINTENANCE, CLEANING, DISPOSAL).
<i>no zero adjustment</i>	No zero adjustment is stored in the meter for the cell. Perform zero adjustment (see chapter chapter 11 ZERO ADJUSTMENT).
<i>method invalid</i>	No data is stored in the photometer for the selected method. Update method data (see chapter chapter 12 UPDATING METHOD DATA).
<i>wrong method</i>	During a difference measurement, the method was changed between the first and second measurement. During a difference measurement, the method must remain identical.
<i>E_0</i>	Hardware error: Send the photometer to the service department.
<i>E_1, E_2 or E_3</i>	Replace the lamp (see chapter chapter 14 MAINTENANCE, CLEANING, DISPOSAL). If the error message remains, send the meter to the service department.

enno



Spectroquant® **NOVA 30**

Analytical Procedures
Appendices

Contents

Table – **Available photometric test kits**

Analytical Procedures

Appendix 1 – **Suitability of Test Kits for Testing Seawater**

Appendix 2 – **Spectroquant® CombiCheck and Standard Solutions**

Appendix 3 – **Instructions for the Preparation of Standard Solutions**

Available photometric test kits

The following methods are programmed into the photometer and measurements can be made without any further adjustments. Method selection is achieved through a barcode on the cell (for cell tests) or through a barcode on the AutoSelector (for reagent tests). The method number listed in column 1 is for manual selection. The total range relates to the cited test in column 2.

Method No.	Determination		Total Range	Method
208	Acid Capacity Cell Test to pH 4.3 (total alkalinity)	101758	0.40 – 8.00 mmol/l	Indicator reaction
196	Aluminium Cell Test*	100594	0.02 – 0.50 mg/l Al	Chromazurole S
104	Ammonium Cell Test	114739	0.010 – 2.000 mg/l NH ₄ -N	Indophenol blue
051	Ammonium Cell Test	114558	0.20 – 8.00 mg/l NH ₄ -N	Indophenol blue
052	Ammonium Cell Test	114544	0.5 – 16.0 mg/l NH ₄ -N	Indophenol blue
053	Ammonium Cell Test	114559	4.0 – 80.0 mg/l NH ₄ -N	Indophenol blue
156	AOX Cell Test*	100675	0.05 – 2.50 mg/l AOX	Oxidation to chloride
157	BOD Cell Test*	100687	0.5 – 3000 mg/l O ₂	Modification of Winkler method
067	Cadmium Cell Test	114834	0.025 – 1.000 mg/l Cd	Cation derivative
165	Calcium Cell Test*	100858	10 – 250 mg/l Ca	Phthalein purple
095	Chloride Cell Test*	114730	5 – 125 mg/l Cl	Iron(III)-thiocyanat
218	Chloride Cell Test*	101804	0.5 – 15.0 mg/l Cl	Iron(III)-thiocyanat
141	Chlorine Cell Test* (free chlorine)	100595	0.03 – 6.00 mg/l Cl ₂	S-DPD
142	Chlorine Cell Test* (free and total chlorine)	100597	0.03 – 6.00 mg/l Cl ₂	S-DPD
194	Chlorine Cell Test* (free and total chlorine)	100086/100087/ 100088	0.03 – 6.00 mg/l Cl ₂	DPD
039	Chromate Cell Test*	114552	0.05 – 2.00 mg/l Cr	Diphenylcarbazide
039	Chromate Cell Test* (total chromium)	114552	0.05 – 2.00 mg/l Cr	Peroxodisulfate oxidation, diphenylcarbazide
020	Chromium Baths		20 – 400 g/l CrO ₃	Inherent color
031	COD Cell Test*	114560	4.0 – 40.0 mg/l COD	Chromosulfuric acid oxidation, chromate determination
211	COD Cell Test*	101796	5.0 – 80.0 mg/l COD	Chromosulfuric acid oxidation, chromate determination
014	COD Cell Test*	114540	10 – 150 mg/l COD	Chromosulfuric acid oxidation, chromate determination
105	COD Cell Test*	114895	15 – 300 mg/l COD	Chromosulfuric acid oxidation, chromate determination
093	COD Cell Test*	114690	50 – 500 mg/l COD	Chromosulfuric acid oxidation, chromate determination
023	COD Cell Test*	114541	25 – 1500 mg/l COD	Chromosulfuric acid oxidation, chromium(III) determination
094	COD Cell Test*	114691	300 – 3500 mg/l COD	Chromosulfuric acid oxidation, chromium(III) determination
024	COD Cell Test*	114555	500 – 10000 mg/l COD	Chromosulfuric acid oxidation, chromium(III) determination
209	COD Cell Test*	101797	5000 – 90000 mg/l COD	Chromosulfuric acid oxidation, chromium(III) determination
137	COD Cell Test (Hg free)*	109772	10 – 150 mg/l COD	Chromosulfuric acid oxidation, chromate determination
138	COD Cell Test (Hg free)*	109773	100 – 1500 mg/l COD	Chromosulfuric acid oxidation, chromium(III) determination
220	COD Cell Test for seawater*	117058	5.0 – 60.0 mg/l COD	Chloride depletion, chromosulfuric acid oxidation, chromate determination
221	COD Cell Test for seawater*	117059	50 – 3000 mg/l COD	Chloride depletion, chromosulfuric acid oxidation, chromium(III) chromate determination
026	Copper Cell Test*	114553	0.05 – 8.00 mg/l Cu	Cuprizone
083	Copper Baths		10.0 – 50.0 g/l Cu	Inherent color
228	Cyanide Cell Test* (free cyanide)	102531	0.010 – 0.500 mg/l CN	Barbituric acid and pyridinecarboxylic acid
075	Cyanide Cell Test* (free cyanide)	114561	0.010 – 0.500 mg/l CN	Barbituric acid and pyridinecarboxylic acid
075	Cyanide Cell Test* (readily liberated cyanide)	114561	0.010 – 0.500 mg/l CN	Citric acid, barbituric acid, and pyridinecarboxylic acid
028	Formaldehyde Cell Test*	114500	0.10 – 8.00 mg/l HCHO	Chromotropic acid
	Hardness see Total Hardness or Residual Hardness			

* turbidity correction possible

Available photometric test kits

Method No.	Determination		Total Range	Method
037	Iron Cell Test	114549	0.05 – 4.00 mg/l Fe	Triazine
106	Iron Cell Test*	114896	1.0 – 50.0 mg/l Fe (Fe(II) and Fe(III))	2,2'-Dipyridyl
066	Lead Cell Test*	114833	0.10 – 5.00 mg/l Pb	PAR
158	Magnesium Cell Test*	100815	5.0 – 75.0 mg/l Mg	Phthalein purple
159	Manganese Cell Test*	100816	0.10 – 5.00 mg/l Mn	Formaloxime
017	Nickel Cell Test*	114554	0.10 – 6.00 mg/l Ni	Dimethylglyoxime
057	Nickel Baths		10 – 120 g/l Ni	Inherent color
059	Nitrate Cell Test*	114542	0.5 – 18.0 mg/l NO ₃ -N	Nitrospectral
030	Nitrate Cell Test*	114563	0.5 – 25.0 mg/l NO ₃ -N	2,6-Dimethylphenol
107	Nitrate Cell Test*	114764	1.0 – 50.0 mg/l NO ₃ -N	2,6-Dimethylphenol
151	Nitrate Cell Test*	100614	23 – 225 mg/l NO ₃ -N	2,6-Dimethylphenol
035	Nitrite Cell Test*	114547	0.010 – 0.700 mg/l NO ₂ -N	Griess reaction
197	Nitrite Cell Test*	100609	1.0 – 90.0 mg/l NO ₂ -N	Iron(II) ethylenediammonium sulfate
068	Nitrogen (total) Cell Test	114537	0.5 – 15.0 mg/l N	Peroxodisulfate oxidation, nitrospectral
153	Nitrogen (total) Cell Test*	100613	0.5 – 15.0 mg/l N	Peroxodisulfate oxidation, 2,6-dimethylphenol
108	Nitrogen (total) Cell Test	114763	10 – 150 mg/l N	Peroxodisulfate oxidation, 2,6-dimethylphenol
092	Oxygen Cell Test*	114694	0.5 – 12.0 mg/l O ₂	Modification of Winkler method
186	pH Cell Test	101744	6.4 – 8.8	Phenol red
212	Phosphate Cell Test	100474	0.05 – 5.00 mg/l PO ₄ -P	Phosphomolybdenum blue
055	Phosphate Cell Test	114543	0.05 – 5.00 mg/l PO ₄ -P	Phosphomolybdenum blue
055	Phosphate Cell Test (total phosphorus)	114543	0.05 – 5.00 mg/l P	Peroxodisulfate oxidation, phosphomolybdenum blue
213	Phosphate Cell Test	100475	0.5 – 25.0 mg/l PO ₄ -P	Phosphomolybdenum blue
086	Phosphate Cell Test	114729	0.5 – 25.0 mg/l PO ₄ -P	Phosphomolybdenum blue
086	Phosphate Cell Test (total phosphorus)	114729	0.5 – 25.0 mg/l P	Peroxodisulfate oxidation, phosphomolybdenum blue
152	Phosphate Cell Test	100616	3.0 – 100.0 mg/l PO ₄ -P	Phosphomolybdenum blue
214	Phosphate Cell Test	100673	3.0 – 100.0 mg/l PO ₄ -P	Phosphomolybdenum blue
214	Phosphate Cell Test (total phosphorus)	100673	3.0 – 100.0 mg/l P	Peroxodisulfate oxidation, phosphomolybdenum blue
069	Phosphate Cell Test*	114546	0.5 – 25.0 mg/l PO ₄ -P	Vanadatomoxybdate
103	Potassium Cell Test	114562	5.0 – 50.0 mg/l K	Kalignost, turbidimetric
150	Potassium Cell Test	100615	30 – 300 mg/l K	Kalignost, turbidimetric
098	Residual Hardness Cell Test*	114683	0.50 – 5.00 mg/l Ca	Phthalein purple
168	Sodium Cell Test in nutrient solutions*	100885	10 – 300 mg/l Na	indirectly as chloride
229	Sulfate Cell Test	102532	1.0 – 50.0 mg/l SO ₄	Bariumsulfate, turbidimetric
064	Sulfate Cell Test	114548	5 – 250 mg/l SO ₄	Bariumsulfate, turbidimetric
154	Sulfate Cell Test	100617	50 – 500 mg/l SO ₄	Bariumsulfate, turbidimetric
082	Sulfate Cell Test	114564	100 – 1000 mg/l SO ₄	Bariumsulfate, turbidimetric
193	Surfactants (nonionic) Cell Test*	101787	0.10 – 7.50 mg/l n-Ten	TBPE
182	Suspended Solids		50 – 750 mg/l SusS	
172	TOC Cell Test	114878	5.0 – 80.0 mg/l TOC	Peroxodisulfate oxidation, indicator
173	TOC Cell Test	114879	50 – 800 mg/l TOC	Peroxodisulfate oxidation, indicator
178	Total Hardness Cell Test*	100961	5 – 215 mg/l Ca	Phthalein purple
	Water hardness see Total Hardness or Residual Hardness			
191	Volatile Organic Acids Cell Test*	101763	50 – 3000 mg/l HOAc	Esterification
222	Volatile Organic Acids Cell Test*	101749	50 – 3000 mg/l CH ₃ COOH	Esterification
223	Volatile Organic Acids Test*	101809	50 – 3000 mg/l CH ₃ COOH	Esterification
174	Zinc Cell Test	100861	0.025 – 1.000 mg/l Zn	PAR
074	Zinc Cell Test	114566	0.20 – 5.00 mg/l Zn	PAR

* turbidity correction possible


Acid Capacity to pH 4.3 (Total Alkalinity)

101758

Cell Test

Measuring range: 0.40 – 8.00 mmol/l

20 – 400 mg/l CaCO₃




Pipette 4.0 ml of **AC-1** into a round cell.



Add 1.0 ml of the sample with pipette, close the cell with the screw cap, and mix.



Add 0.50 ml of **AC-2** with pipette, close the cell with the screw cap, and mix.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a sodium hydroxide solution 0.1 mol/l, Cat.No. 109141, can be used after diluting accordingly (see section “Standard solutions”).

Aluminium

100594

Cell Test

Measuring 0.02 – 0.50 mg/l Al

range: Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 3 – 10. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Pipette 6.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 1 level blue microspoon of **Al-1K**, close with the screw cap.




Shake the cell vigorously to dissolve the solid substance.



Add 0.25 ml of **Al-2K** with pipette, close with the screw cap, and mix.



Reaction time:
5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use aluminium standard solution Certipur®, Cat.No. 119770, concentration 1000 mg/l Al can be used after diluting accordingly.

Ammonium

114739

Cell Test


Measuring	0.010 – 2.000 mg/l NH ₄ -N
range:	0.01 – 2.58 mg/l NH ₄
	0.010 – 2.000 mg/l NH ₃ -N
	0.01 – 2.43 mg/l NH ₃
	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 4 – 13. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Pipette 5.0 ml of the sample into a reaction cell close with the screw cap, and mix.




Add 1 dose of **NH₄-1K** using the blue dose-metering cap, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 15 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

Very high ammonium concentrations in the sample produce turquoise-colored solutions (measurement solution should be yellow-green to green) and false-low readings are yielded. In such cases the sample must be diluted (plausibility check).

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 50, Cat.No. 114695, or the Standard solution for photometric applications, CRM, Cat.No. 125022 and 125023.

Ready-for-use ammonium standard solution Certipur®, Cat.No. 119812, concentration 1000 mg/l NH₄⁺, can also be used after diluting accordingly.


To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 50) is highly recommended.

Ammonium

114558

Cell Test


Measuring	0.20 – 8.00 mg/l NH ₄ -N
range:	0.26 – 10.30 mg/l NH ₄
	0.20 – 8.00 mg/l NH ₃ -N
	0.24 – 9.73 mg/l NH ₃
	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 4 – 13
If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Pipette 1.0 ml of the sample into a reaction cell close with the screw cap, and mix.




Add 1 dose of **NH₄-1K** using the blue dose-metering cap, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 15 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

Very high ammonium concentrations in the sample produce turquoise-colored solutions (measurement solution should be yellow-green to green) and false-low readings are yielded. In such cases the sample must be diluted (plausibility check).

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 10, Cat.No. 114676, or the Standard solution for photometric applications, CRM, Cat.No. 125022, 125023, 125024, and 125025.

Ready-for-use ammonium standard solution Certipur®, Cat.No. 119812, concentration 1000 mg/l NH₄⁺, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 10) is highly recommended.

Ammonium

114544

Cell Test


Measuring	0.5 – 16.0 mg/l NH ₄ -N
range:	0.6 – 20.6 mg/l NH ₄
	0.5 – 16.0 mg/l NH ₃ -N
	0.6 – 19.5 mg/l NH ₃
	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 4 – 13.
If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Pipette 0.50 ml of the sample into a reaction cell close with the screw cap, and mix.




Add 1 dose of **NH₄-1K** using the blue dose-metering cap, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 15 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

Very high ammonium concentrations in the sample produce turquoise-colored solutions (measurement solution should be yellow-green to green) and false-low readings are yielded. In such cases the sample must be diluted (plausibility check).

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 20, Cat.No. 114675, or the Standard solution for photometric applications, CRM, Cat.No. 125023, 125024, 125025, and 125026.

Ready-for-use ammonium standard solution Certipur®, Cat.No. 119812, concentration 1000 mg/l NH₄⁺, can also be used after diluting accordingly.


To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 20) is highly recommended.

Ammonium

114559

Cell Test


Measuring	4.0 – 80.0 mg/l NH ₄ -N
range:	5.2 – 103.0 mg/l NH ₄
	4.0 – 80.0 mg/l NH ₃ -N
	4.9 – 97.3 mg/l NH ₃
	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 4 – 13. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Pipette 0.10 ml of the sample into a reaction cell close with the screw cap, and mix.




Add 1 dose of **NH₄-1K** using the blue dose-metering cap, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 15 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

Very high ammonium concentrations in the sample produce turquoise-colored solutions (measurement solution should be yellow-green to green) and false-low readings are yielded. In such cases the sample must be diluted (plausibility check).

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 70, Cat.No. 114689, or the Standard solution for photometric applications, CRM, Cat.No. 125025, 125026, and 125027.

Ready-for-use ammonium standard solution Certipur®, Cat.No. 119812, concentration 1000 mg/l NH₄⁺, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 70) is highly recommended.

AOX


Adsorbable Organic Halogens (x)

100675


Cell Test

Measuring range: 0.05–2.50 mg/l AOX


Preparation of the adsorption column:




Place the column in an empty cell. Fill 1 level blue microspoon of **AOX-1** into the column using the glass funnel.



Run 3 separate 1-ml portions of **AOX-2** through the column. Discard the wash solution.



Run 3 separate 1-ml portions of **AOX-3** through the column. Discard the wash solution.




Close the bottom end of the column with the stopper. Apply to the column 1 ml of **AOX-3**. Close the top end of the column with the stopper and swirl to eliminate air bubbles. Remove the stopper on the top end and fill the column to the brim with **AOX-3**.


Sample enrichment:




Check the pH of the sample, specified range: pH 6 – 7. If required, add dilute sodium hydroxide solution or nitric acid drop by drop to adjust the pH.




Attach the glass reservoir to the prepared column (closed at the bottom end).



Fill 100 ml of the sample and 6 drops of **AOX-4** into the reservoir.



Remove the stopper from the column outlet and run the sample through completely.



Detach the column from the reservoir. Apply 3 separate 1-ml portions of **AOX-3**. Discard the wash solution.


AOX

Adsorbable Organic Halogens (x)

100675

Cell Test


Digestion:




Fill the 10-ml syringe with 10 ml of reagent **AOX-5** and attach the syringe with the column outlet using the connector. Place the top end of the column on an empty cell and rinse the charcoal filling of the column into an empty 16-mm cell.




Add 2 level green microspoons of **AOX-6**, close the cell with the screw cap, and mix.



Heat the cell at 120 °C in the thermoreactor for 30 minutes.



Remove the cell from the thermoreactor and place in a test-tube rack to cool to room temperature.




Add 5 drops of **AOX-4**, close the cell and mix; clear supernatant: **pretreated sample**.

Determination:




Pipette 0.20 ml of **AOX-1K** into a reaction cell, and mix.



Add 7.0 ml of **pretreated sample** with glass pipette, close the cell with the screw cap, and mix.



Reaction time: 15 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) Spectroquant® AOX Standard, Cat.No. 100680, concentration 0.2 – 2.0 mg/l can be used.

BOD

Biochemical Oxygen Demand


100687

Cell Test

Measuring 0.5 – 3000 mg/l O₂

range: Expression of results also possible in mmol/l.

Preparation and incubation:



Check the pH of the sample, specified range: pH 6 – 8. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Fill 2 oxygen reaction bottles each with **pretreated sample** and 2 glass beads to overflowing. Close bubble-free with the slanted ground-glass stoppers.



Fill 2 oxygen reaction bottles each with **inoculated nutrient-salt solution** and 2 glass beads to overflowing. Close bubble-free with the slanted ground-glass stoppers.

Measurement of initial oxygen concentration

= **Result 1**
(measurement sample)
= **Result 1**
(blank)



Use one bottle of **pretreated sample** and one of **inoculated nutrient-salt solution** for the measurement of the initial oxygen concentration.

Incubate one bottle of **pretreated sample** and one of **inoculated nutrient-salt solution** closed in a thermostatic incubation cabinet at $20 \pm 1^\circ\text{C}$ for 5 days.

Determination:

Measurement of final oxygen concentration

= **Result 2**
(measurement sample)
= **Result 2**
(blank)



Add 5 drops of **BOD-1K** and then 10 drops of **BOD-2K**, close bubble-free, and mix for approx. 10 seconds.



Reaction time:
1 minute




Add 10 drops of **BOD-3K**, reclose, and mix.



Fill the solution into a round cell.

After incubation, use one bottle of **pretreated sample** and one of **inoculated nutrient-salt solution** for the measurement of the final oxygen concentration.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Calculation:

BOD of measurement sample:
Result 1 – Result 2 (measurement sample) = A in mg/l

BOD of blank:
Result 1 – Result 2 (blank) = B in mg/l

BOD of original sample in mg/l = A • dilution factor – B

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) Spectroquant BOD Standard (acc. to EN 1899), Cat.No. 100718, can be used.


Cadmium

114834

Cell Test

Measuring 0.025 – 1.000 mg/l Cd

range: Expression of results also possible in mmol/l.



Check the pH of the sample, specified range: pH 3 – 11. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.




Pipette 5.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 0.20 ml of **Cd-1K** with pipette, close the cell with the screw cap, and mix.




Add 1 level green microspoon of **Cd-2K**, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time:
2 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

For the determination of **total cadmium** a pretreatment with Crack Set 10C, Cat.No. 114688 or Crack Set 10, Cat.No. 114687, and thermoreactor is necessary.

Result can be expressed as sum of cadmium (Σ Cd).

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 30, Cat.No. 114677.

Ready-for-use cadmium standard solution Certipur®, Cat.No. 119777, concentration 1000 mg/l Cd, can also be used after diluting accordingly.


To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 30) is highly recommended.

Calcium

100858

Cell Test

Measuring	10–250 mg/l Ca
range:	14–350 mg/l CaO
	25–624 mg/l CaCO ₃
	Expression of results also possible in mmol/l.



Check the pH of the sample, specified range: pH 3 – 9.
If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.



Pipette 1.0 ml of the sample into a reaction cell, close with the screw cap, and mix.




Add 1.0 ml of **Ca-1K** with pipette, close the cell with the screw cap, and mix.



Reaction time:
exactly 3 minutes



Add 0.50 ml of **Ca-2K** with pipette, close the cell with the screw cap, and mix.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a freshly prepared standard solution can be used (see section “Standard solutions”).

Chloride

114730

Cell Test

Measuring 5–125 mg/l Cl

range: Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 1 – 12. If required, add dilute ammonia solution or nitric acid drop by drop to adjust the pH.



Pipette 0.50 ml of **CI-1K** into a reaction cell, close with the screw cap, and mix.



Add 1.0 ml of the sample with pipette, close with the screw cap, and mix.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 10 and 20, Cat.Nos. 114676 and 114675.

Ready-for-use chloride standard solution Certipur®, Cat.No. 119897, concentration 1000 mg/l Cl⁻, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck) is highly recommended.

Chloride

101804

Cell Test

Measuring 0.5–15.0 mg/l Cl


range: Expression of results also possible in mmol/l.



Check the pH of the sample, specified range: pH 3 – 11. If required, add dilute ammonia solution or nitric acid drop by drop to adjust the pH.




Pipette 10 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 0.25 ml of **Cl-1K** with pipette, close with the screw cap, and mix.



Reaction time:
10 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use chloride standard solution Certipur®, Cat.No. 119897, concentration 1000 mg/l Cl⁻, can be used after diluting accordingly.

Chlorine

100595


Determination of free chlorine

Cell Test


Measuring	0.03–6.00 mg/l Cl ₂
range:	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 4 – 8. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Pipette 5.0 ml of the sample into a round cell.




Add 1 level blue micro-spoon of Cl₂-1, close with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 1 minute



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

Very high chlorine concentrations in the sample produce yellow-colored solutions (measurement solution should be red) and false-low readings are yielded. In such cases the sample must be diluted (plausibility check).

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a freshly prepared standard solution can be used (see section “Standard solutions”).

Chlorine

100597

Determination of free chlorine and total chlorine

Cell Test

Measuring 0.03–6.00 mg/l Cl₂

range: Expression of results also possible in mmol/l and also in free Cl₂ [Cl₂(f)], combined Cl₂ [Cl₂(b)], and total Cl₂ [Cl₂(t)].


Determination of free chlorine




Check the pH of the sample, specified range: pH 4 – 8.
If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Pipette 5.0 ml of the sample into a round cell.




Add 1 level blue micro-spoon of Cl₂-1, close with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 1 minute



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Determination of total chlorine

Same preparation as described above, add 2 drops of Cl₂-2, close the cell with the screw cap, and mix after dissolving solid.

A differentiation between free and combined chlorine [Cl₂(f) and Cl₂(b)] can be performed on the photometer. Prior to measuring, select the differentiation measurement and choose the corresponding citation form. Then measure the free chlorine, press enter, remove the cell, add 2 drops of Cl₂-2, close with the screw cap, mix, and measure the total chlorine. After pressing enter, the individual measuring values for free and combined chlorine are shown on the display.

Important:

Very high chlorine concentrations in the sample produce yellow-colored solutions (measurement solution should be red) and false-low readings are yielded. In such cases the sample must be diluted (plausibility check).
After each determination of total chlorine rinse the cell with sulfuric acid 25 % and subsequently several times with distilled water.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a freshly prepared standard solution can be used (see section "Standard solutions").

Chlorine (with liquid reagents)

100086/100087/
100088

Determination of free chlorine and total chlorine

Cell Test


Measuring 0.03–6.00 mg/l Cl₂

range: Expression of results also possible in mmol/l and also in free Cl₂ [Cl₂(f)], combined Cl₂ [Cl₂(b)], and total Cl₂ [Cl₂(t)].

Determination of free chlorine



Check the pH of the sample, specified range: pH 4 – 8.
If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Place 6 drops of **Cl₂-1** into a round cell.




Add 3 drops of **Cl₂-2**, close with the screw cap, and mix.



Add 10 ml of the sample with pipette, close with the screw cap, and mix.



Reaction time: 1 minute



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Determination of total chlorine

Same preparation as described above, add 2 drops of **Cl₂-3**, close with the screw cap, and mix after the end of the reaction time.

A differentiation between free and combined chlorine [Cl₂(f) and Cl₂(b)] can be performed on the photometer. Prior to measuring, select the differentiation measurement and choose the corresponding citation form. Then measure the free chlorine, press enter, remove the cell, add 2 drops of Cl₂-3, close with the screw cap, mix, and measure the total chlorine. After pressing enter, the individual measuring values for free and combined chlorine are shown on the display.

Important:

Very high chlorine concentrations in the sample produce yellow-colored solutions (measurement solution should be red) and false-low readings are yielded. In such cases the sample must be diluted (plausibility check).
After each determination of total chlorine rinse the cell with sulfuric acid 25 % and subsequently several times with distilled water.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a freshly prepared standard solution can be used (see section "Standard solutions").


Chromate

114552


Determination of chromium(VI)

Cell Test


Measuring	0.05 – 2.00 mg/l Cr
range:	0.11 – 4.46 mg/l CrO ₄
	Expression of results also possible in mmol/l.



Check the pH of the sample, specified range: pH 1 – 9. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Add 6 drops of **Cr-3K** into a reaction cell, close with the screw cap.




Shake the cell vigorously to dissolve the solid substance and leave to stand for **1 minute**.



Add 5.0 ml of the sample with pipette, close the cell with the screw cap, and mix.



Reaction time: 1 minute



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use chromate standard solution Certipur®, Cat.No. 119780, concentration 1000 mg/l CrO₄²⁻, can be used after diluting accordingly.

Chromate

Determination of total chromium
= sum of chromium(VI) and chromium(III)

114552

Cell Test

Measuring 0.05–2.00 mg/l Cr

range: 0.11–4.46 mg/l CrO₄


Expression of results also possible in mmol/l and also in Cr total (Σ Cr), Cr(III), and Cr(VI).




Check the pH of the sample, specified range: pH 1 – 9. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.




Pipette 10 ml of the sample into an empty round cell (Empty cells, Cat.No. 114724).



Add 1 drop of **Cr-1K**, close with the screw cap, and mix.




Add 1 dose of **Cr-2K** using the blue dose-metering cap, close the reaction cell with the screw cap.




Heat the cell in the thermoreactor at 120 °C (100 °C) for 1 hour.



Remove the cell from the thermoreactor and place in a test-tube rack to cool to room temperature: **pretreated sample**.



Add 6 drops of **Cr-3K** into a reaction cell, close the cell with the screw cap.




Shake the cell vigorously to dissolve the solid substance and leave to stand for **1 minute**.



Add 5.0 ml of the **pretreated sample** with pipette, close with the screw cap, and mix.



Reaction time:
1 minute



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

A differentiation between chromium(VI) and chromium(III) can be performed on the photometer. Prior to measuring, select the differentiation measurement and choose the corresponding citation form. Then measure the total chromium, press enter and measure the chromium(VI) (see analytical procedure for chromium(VI)). After pressing enter, the individual measuring values for Cr VI and Cr III are shown on the display.


Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use chromate standard solution Certipur[®], Cat.No. 119780, concentration 1000 mg/l CrO₄²⁻, can be used after diluting accordingly.


Chromium in electroplating baths

Inherent color


Measuring range: 20–400 g/l CrO₃




Pipette 5.0 ml of the sample into a 100-ml volumetric flask, fill to the mark with distilled water and mix thoroughly.




Pipette 4.0 ml of the dilute sample into a 100-ml volumetric flask, fill to the mark with distilled water and mix thoroughly.



Pipette 5.0 ml of the 1:500 dilute sample into an empty round cell (Empty cells, Cat. No. 114724).



Add 5.0 ml of **sulfuric acid 40%**, close the cell with the screw cap, and mix.



Place the cell into the cell compartment. Select method no. **20**.


COD

Chemical Oxygen Demand

114560

Cell Test


Measuring	4.0–40.0 mg/l COD or O ₂
range:	Expression of results also possible in mmol/l.




Suspend the bottom sediment in the cell by swirling.




Carefully pipette 3.0 ml of the sample into a reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**




Heat the reaction cell in the thermoreactor at 148 °C for 2 hours.




Remove the cell from the thermoreactor and place in a test-tube rack to cool.



Swirl the cell after 10 minutes.



Replace the cell in the rack for complete cooling to room temperature. **Very important!**



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 50, Cat.No. 114695, or the Standard solution for photometric applications, CRM, Cat.No. 125028.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 50) is highly recommended.


COD

Chemical Oxygen Demand


101796

Cell Test


Measuring	5.0–80.0 mg/l COD or O ₂
range:	Expression of results also possible in mmol/l.




Suspend the bottom sediment in the cell by swirling.




Carefully pipette 2.0 ml of the sample into a reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**




Heat the reaction cell in the thermoreactor at 148 °C for 2 hours.




Remove the cell from the thermoreactor and place in a test-tube rack to cool.



Swirl the cell after 10 minutes.



Replace the cell in the rack for complete cooling to room temperature. **Very important!**



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 50, Cat.No. 114695, or the Standard solution for photometric applications, CRM, Cat.No. 125028.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 50) is highly recommended.


COD

Chemical Oxygen Demand

114540

Cell Test


Measuring 10–150 mg/l COD or O₂
range: Expression of results also possible in mmol/l.




Suspend the bottom sediment in the cell by swirling.




Carefully pipette 3.0 ml of the sample into a reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**




Heat the reaction cell in the thermoreactor at 148 °C for 2 hours.




Remove the cell from the thermoreactor and place in a test-tube rack to cool.



Swirl the cell after 10 minutes.



Replace the cell in the rack for complete cooling to room temperature. **Very important!**



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 10, Cat.No. 114676, or the Standard solution for photometric applications, CRM, Cat.No. 125029.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 10) is highly recommended.


COD

Chemical Oxygen Demand

114895

Cell Test


Measuring 15–300 mg/l COD or O₂
range: Expression of results also possible in mmol/l.




Suspend the bottom sediment in the cell by swirling.




Carefully pipette 2.0 ml of the sample into a reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**




Heat the reaction cell in the thermoreactor at 148 °C for 2 hours.




Remove the cell from the thermoreactor and place in a test-tube rack to cool.



Swirl the cell after 10 minutes.



Replace the cell in the rack for complete cooling to room temperature. **Very important!**



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 60, Cat.No. 114696, or the Standard solution for photometric applications, CRM, Cat.No. 125029 and 125030.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 60) is highly recommended.


COD

Chemical Oxygen Demand

114690

Cell Test


Measuring 50–500 mg/l COD or O₂
range: Expression of results also possible in mmol/l.




Suspend the bottom sediment in the cell by swirling.




Carefully pipette 2.0 ml of the sample into a reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**




Heat the reaction cell in the thermoreactor at 148 °C for 2 hours.




Remove the cell from the thermoreactor and place in a test-tube rack to cool.



Swirl the cell after 10 minutes.



Replace the cell in the rack for complete cooling to room temperature. **Very important!**



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 60, Cat.No. 114696, or the Standard solution for photometric applications, CRM, Cat.No. 125029, 125030, and 125031.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 60) is highly recommended.


COD

Chemical Oxygen Demand

114541

Cell Test


Measuring	25–1500 mg/l COD or O ₂
range:	Expression of results also possible in mmol/l.




Suspend the bottom sediment in the cell by swirling.




Carefully pipette 3.0 ml of the sample into a reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**




Heat the reaction cell in the thermoreactor at 148 °C for 2 hours.




Remove the cell from the thermoreactor and place in a test-tube rack to cool.



Swirl the cell after 10 minutes.



Replace the cell in the rack for complete cooling to room temperature. **Very important!**



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 20, Cat.No. 114675, or the Standard solution for photometric applications, CRM, Cat.No. 125029, 125030, 125031, and 125032.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 20) is highly recommended.


COD

Chemical Oxygen Demand

114691

Cell Test


Measuring	300–3500 mg/l COD or O ₂
range:	Expression of results also possible in mmol/l.




Suspend the bottom sediment in the cell by swirling.




Carefully pipette 2.0 ml of the sample into a reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**




Heat the reaction cell in the thermoreactor at 148 °C for 2 hours.




Remove the cell from the thermoreactor and place in a test-tube rack to cool.



Swirl the cell after 10 minutes.



Replace the cell in the rack for complete cooling to room temperature. **Very important!**



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 80, Cat.No. 114738, or the Standard solution for photometric applications, CRM, Cat.No. 125031, 125032, and 125033.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 80) is highly recommended.


COD

Chemical Oxygen Demand

114555

Cell Test


Measuring 500–10000 mg/l COD or O₂
range: Expression of results also possible in mmol/l.




Suspend the bottom sediment in the cell by swirling.




Carefully pipette 1.0 ml of the sample into a reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**




Heat the reaction cell in the thermoreactor at 148 °C for 2 hours.




Remove the cell from the thermoreactor and place in a test-tube rack to cool.



Swirl the cell after 10 minutes.



Replace the cell in the rack for complete cooling to room temperature. **Very important!**



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 70, Cat.No. 114689, or the Standard solution for photometric applications, CRM, Cat.No. 125032, 125033, and 125034.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 70) is highly recommended.


COD

Chemical Oxygen Demand

101797

Cell Test


Measuring 5000–90000 mg/l COD or O₂
range: Expression of results also possible in mmol/l.




Suspend the bottom sediment in the cell by swirling.




Carefully pipette 0.10 ml of the sample into a reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**




Heat the reaction cell in the thermoreactor at 148 °C for 2 hours.




Remove the cell from the thermoreactor and place in a test-tube rack to cool.



Swirl the cell after 10 minutes.



Replace the cell in the rack for complete cooling to room temperature. **Very important!**



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use the Standard solution for photometric applications, CRM, Cat.No. 125034 and 125035.

COD (Hg-free)

Chemical Oxygen Demand


109772

Cell Test


Measuring 10–150 mg/l COD or O₂
range: Expression of results also possible in mmol/l.




Carefully pipette 2.0 ml of the sample into a reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**




Heat the reaction cell in the thermoreactor at 148 °C for 2 hours.




Remove the cell from the thermoreactor and place in a test-tube rack to cool.



Swirl the cell after 10 minutes.



Replace the cell in the rack for complete cooling to room temperature. **Very important!**



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use the Standard solution for photometric applications, CRM, Cat.No. 125028 and 125029.

COD (Hg-free)

Chemical Oxygen Demand

109773


Cell Test

Measuring 100–1500 mg/l COD or O₂


range: Expression of results also possible in mmol/l.




Carefully pipette 2.0 ml of the sample into a reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**




Heat the reaction cell in the thermoreactor at 148 °C for 2 hours.




Remove the cell from the thermoreactor and place in a test-tube rack to cool.



Swirl the cell after 10 minutes.



Replace the cell in the rack for complete cooling to room temperature. **Very important!**



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use the Standard solution for photometric applications, CRM, Cat.No. 125029, 125030, 125031, and 125032.

COD


Chemical Oxygen Demand
for seawater / high chloride contents

117058

Cell Test

Measuring range: 5.0–60.0 mg/l COD or O₂ 16-mm cell

Chloride depletion:




Pipette with glass pipette 20 ml of the sample into a 300-ml Erlenmeyer flask with NS 29/32.



Pipette with glass pipette 20 ml of distilled water (Water for chromatography LiChrosolv®, Cat.No. 115333, is recommended) into a second 300-ml Erlenmeyer flask with NS 29/32.




Add to each a magnetic stirring rod, and cool in the ice bath.




Add **slowly** to each Erlenmeyer flask 25 ml of **Sulfuric acid for the determination of COD** (Cat. No. 117048) with glass pipette **under cooling and stirring**.




Cool both Erlenmeyer flasks to room temperature in the ice bath.




Fill 6 - 7 g each of **Sodalime with indicator** (Cat. No. 106733) into two absorption tubes (Cat. No. 115955).



Close the absorption tubes with the glass stoppers, and attach to the top of the Erlenmeyer flasks.



Stir at 250 rpm for 2 h at room temperature: depleted sample / depleted blank



Check the chloride content of the depleted sample using MColorTest™ Chloride Test (Cat. No. 111132) according to the application (see the website):
Specified value
<2000 mg/l Cl⁻.

Chloride determination (acc. to application - brief version):

Fill 5.0 ml of sodium hydroxide solution 2 mol/l, Cat. No. 109136, into the test vessel of the MColorTest™ Chloride Test, Cat. No. 111132.

Carefully allow to run from the pipette 0.5 ml of depleted sample down the inside of the tilted test vessel into the sodium hydroxide solution and mix (**Wear eye protection! The test vessel becomes hot!**).

Add 2 drops of reagent Cl-1 and swirl. The sample directly turns yellow in color. (Reagent Cl-2 is not required.)

Holding the reagent bottle vertically, slowly add reagent Cl-3 dropwise to the sample while swirling until its color changes from yellow to blue-violet. Shortly before the color changes, wait a few seconds after adding each drop.

Result in mg/l chloride = number of drops x 250


COD

Chemical Oxygen Demand
for seawater / high chloride contents


117058

Cell Test

Determination:




Suspend the bottom sediment in two cells by swirling.




Carefully pipette 5.0 ml of the **depleted sample** into a reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**




Carefully pipette 5.0 ml of the **depleted blank** into a second reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**
(Blank cell)




Heat both cells in the thermoreactor at 148 °C for 2 hours.



Remove both cells from the thermoreactor and place in a test-tube rack to cool.




Swirl both cells after 10 minutes.




Replace both cells in the rack for complete cooling to room temperature. **(Very important!)**



Configure the photometer for blank-measurement.



Place the blank cell into the cell compartment. Align the mark on the cell with that on the photometer.



Place the cell containing the sample into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a COD/chloride standard solution must be prepared from Potassium hydrogen phthalate, Cat.No. 102400 and Sodium chloride, Cat.No. 106404 (see section "Standard solutions").

COD


Chemical Oxygen Demand
for seawater / high chloride contents

117059

Cell Test

Measuring range: 50–3000 mg/l COD or O₂ 16-mm cell

Chloride depletion:




Pipette with glass pipette 20 ml of the sample into a 300-ml Erlenmeyer flask with NS 29/32.



Pipette with glass pipette 20 ml of distilled water (Water for chromatography LiChrosolv®, Cat.No. 115333, is recommended) into a second 300-ml Erlenmeyer flask with NS 29/32.




Add to each a magnetic stirring rod, and cool in the ice bath.




Add **slowly** to each Erlenmeyer flask 25 ml of **Sulfuric acid for the determination of COD** (Cat. No. 117048) with glass pipette **under cooling and stirring**.




Cool both Erlenmeyer flasks to room temperature in the ice bath.




Fill 6 - 7 g each of **Sodalime with indicator** (Cat. No. 106733) into two absorption tubes (Cat. No. 115955).



Close the absorption tubes with the glass stoppers, and attach to the top of the Erlenmeyer flasks.



Stir at 250 rpm for 2 h at room temperature: depleted sample / depleted blank



Check the chloride content of the depleted sample using the MColorTest™ Chloride Test (Cat. No. 111132) as per the application instructions (see the website): specified value <250 mg/l Cl⁻.

Chloride determination (acc. the application instructions - abridged version):

Fill 5.0 ml of sodium hydroxide solution 2 mol/l, Cat. No. 109136, into the test vessel of the MColorTest™ Chloride Test, Cat. No. 111132.

Carefully allow to run from the pipette 0.5 ml of depleted sample down the inside of the tilted test vessel onto the sodium hydroxide solution and mix (**Wear eye protection! The cell becomes hot!**).

Add 2 drops of reagent Cl-1 and swirl. The sample directly turns yellow in color. (Reagent Cl-2 is not required.)

Holding the reagent bottle vertically, slowly add reagent Cl-3 dropwise to the sample while swirling until its color changes from yellow to blue-violet. Shortly before the color changes, wait a few seconds after adding each drop.

Result in mg/l chloride = number of drops x 250


COD

Chemical Oxygen Demand
for seawater / high chloride contents


117059

Cell Test

Determination:




Suspend the bottom sediment in two cells by swirling.




Carefully pipette 3.0 ml of the **depleted sample** into a reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**




Carefully pipette 3.0 ml of the **depleted blank** into a second reaction cell, close tightly with the screw cap, and mix vigorously. **Caution, the cell becomes hot!**
(Blank cell)




Heat both cells in the thermoreactor at 148 °C for 2 hours.



Remove both cells from the thermoreactor and place in a test-tube rack to cool.




Swirl both cells after 10 minutes.




Replace both cells in the rack for complete cooling to room temperature. **(Very important!)**



Configure the photometer for blank-measurement.



Place the blank cell into the cell compartment. Align the mark on the cell with that on the photometer.



Place the cell containing the sample into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a COD/chloride standard solution must be prepared from Potassium hydrogen phthalate, Cat.No. 102400 and Sodium chloride, Cat.No. 106404 (see section "Standard solutions").

Copper

114553

Cell Test


Measuring	0.05–8.00 mg/l Cu
range:	Expression of results also possible in mmol/l.



Check the pH of the sample, specified range: pH 4 – 10.
If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.




Pipette 5.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 5 drops of **Cu-1K**, close the cell with the screw cap, and mix.



Reaction time:
5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

Very high copper concentrations in the sample produce turquoise-colored solutions (measurement solution should be blue) and false-low readings are yielded. In such cases the sample must be diluted (plausibility check).

For the determination of **total copper** a pretreatment with Crack Set 10C, Cat.No. 114688, or Crack Set 10, Cat.No. 114687 and thermoreactor is necessary.

Result can be expressed as sum of copper (Σ Cu).

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 30, Cat.No. 114677.


Ready-for-use copper standard solution Certipur®, Cat.No. 119786, concentration 1000 mg/l Cu, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 30) is highly recommended.


Copper in electroplating baths

Inherent color

Measuring range: 10.0–50.0 g/l Cu




Pipette 25 ml of the sample into a 100-ml volumetric flask, fill to the mark with distilled water and mix thoroughly.



Pipette 5.0 ml of the 1:4 dilute sample into an empty round cell (Empty cells, Cat.No. 114724).



Add 5.0 ml of **sulfuric acid 40%**, close the cell with the screw cap, and mix.



Place the cell into the cell compartment. Select method no. **83**.

Cyanide

102531

Determination of free cyanide

Cell Test

Measuring 0.010–0.500 mg/l CN


range: Expression of results also possible in mmol/l and cyanide free [CN(f)].




Check the pH of the sample, specified range: pH 4.5 – 8.0. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Pipette 5.0 ml of the sample into a reaction cell, close with the screw cap, and dissolve the solid substance.




Add 1 level blue microspoon of **CN-1K**, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 10 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use cyanide standard solution Certipur[®], Cat.No. 119533, concentration 1000 mg/l CN⁻, can be used after diluting accordingly.

Cyanide

114561

Determination of free cyanide

Cell Test

Measuring 0.010–0.500 mg/l CN

range: Expression of results also possible in mmol/l and cyanide free [CN(f)].




Check the pH of the sample, specified range: pH 4.5 – 8.0. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Pipette 5.0 ml of the sample into a reaction cell, close with the screw cap, and dissolve the solid substance.




Add 1 level blue microspoon of **CN-3K**, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 10 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use cyanide standard solution Certipur[®], Cat.No. 119533, concentration 1000 mg/l CN⁻, can be used after diluting accordingly.

Cyanide


114561

Determination of readily liberated cyanide


Cell Test

Measuring 0.010–0.500 mg/l CN


range: Expression of results also possible in mmol/l and cyanide readily liberated [CN(v)].




Check the pH of the sample, specified range: pH 4.5 – 8.0. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.




Pipette 10 ml of the sample into an empty round cell (Empty cells, Cat.No. 114724).




Add 1 dose of **CN-1K** using the green dose-metering cap, close the cell with the screw cap.




Heat the cell in the thermoreactor at 120 °C (100 °C) for 30 minutes.




Remove the cell from the thermoreactor and place in a test-tube rack to cool to room temperature.



Swirl the cell before opening.




Add 3 drops of **CN-2K**, close with the screw cap, and mix: **pretreated sample**.



Pipette 5.0 ml of the **pretreated sample** into a reaction cell, close with the screw cap, and dissolve the solid substance.




Add 1 level blue micro-spoon of **CN-3K**, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 10 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use cyanide standard solution Certipur[®], Cat.No. 119533, concentration 1000 mg/l CN⁻, can be used after diluting accordingly.


Formaldehyde

114500

Cell Test

Measuring 0.10–8.00 mg/l HCHO


range: Expression of results also possible in mmol/l.



Check the pH of the sample, specified range: pH 0 – 13.



Add 1 level green micro-spoon of **HCHO-1K** into a reaction cell, close with the screw cap.




Shake the cell vigorously to dissolve the solid substance.



Add 2.0 ml of the sample with pipette, close the cell with the screw cap, and mix. **Caution, cell becomes hot!**



Reaction time: 5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:


To check the measurement system (test reagents, measurement device, and handling) a formaldehyde standard solution must be prepared from Formaldehyde solution 37%, Cat.No. 104003 (see section "Standard solutions").

Iron

114549

Cell Test

Measuring	0.05 – 4.00 mg/l Fe
range:	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 1 – 10. If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.



Pipette 5.0 ml of the sample into a reaction cell, close with the screw cap, and mix.




Add 1 level blue microspoon of **Fe-1K**, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 3 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

For the determination of **total iron** a pretreatment with Crack Set 10C, Cat.No. 114688, or Crack Set 10, Cat.No. 114687 and thermoreactor is necessary.

Result can be expressed as sum of iron (Σ Fe).

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 30, Cat.No. 114677.

Ready-for-use iron standard solution Certipur®, Cat.No. 119781, concentration 1000 mg/l Fe, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 30) is highly recommended.

Iron

114896


Determination of iron(II) and iron(III)

Cell Test

Measuring 1.0–50.0 mg/l Fe

range: Expression of results also possible in mmol/l and also in Fe(II), Fe(III).

Determination of iron (II)




Check the pH of the sample, specified range: pH 3 – 8.
If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.



Pipette 1.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Reaction time:
5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.


Determination of iron (II + III)




Check the pH of the sample, specified range: pH 3 – 8.
If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.



Pipette 1.0 ml of the sample into a reaction cell, close with the screw cap, and mix.




Add 1 dose of **Fe-1K** using the blue dose-metering cap, close the reaction cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time:
5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

A differentiation between iron(II) and iron(III) can be performed on the photometer. Prior to measuring, select the differentiation measurement and choose the corresponding citation form.

Then measure the iron(II + III), press enter and measure the iron(II). After pressing enter, the individual measuring values for Fe II and Fe III are shown on the display.

Important:

For the determination of **total iron** a pretreatment with Crack Set 10C, Cat.No. 114688, or Crack Set 10, Cat.No. 114687, and thermoreactor is necessary.

Result can be expressed as sum of iron (Σ Fe).

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use iron standard solution Certipur®, Cat.No. 119687, concentration 1000 mg/l Fe(III), can be used after diluting accordingly.

Lead

114833

Cell Test


Measuring 0.10–5.00 mg/l Pb

range: Expression of results also possible in mmol/l.


Samples of total hardness 0–10 °d




Check the total hardness of the sample.




Check the pH of the sample, specified range: pH 3–6.
If required, add dilute ammonia solution or nitric acid drop by drop to adjust the pH.



Add 5 drops of **Pb-1K** into a reaction cell and mix.




Add 5.0 ml of the sample with pipette, close the cell with the screw cap, and mix.




Place the cell into the cell compartment. Align the mark on the cell with that on the photometer = **Result A**


Samples of total hardness > 10 °d




Add 1 level grey micro-spoon of **Pb-2K** to the already measured cell, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer = **Result B**



Important:

For the determination of **total lead** a pretreatment with Crack Set 10C, Cat.No. 114688, or Crack Set 10, Cat.No. 114687, and thermoreactor is necessary.

Result can be expressed as sum of lead (Σ Pb).

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 40, Cat.No. 114692.

Ready-for-use lead standard solution Certipur®, Cat.No. 119776, concentration 1000 mg/l Pb, can also be used after diluting accordingly.


To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 40) is highly recommended.

Magnesium


100815

Cell Test


Measuring	5.0 – 75.0 mg/l Mg
range:	Expression of results also possible in mmol/l.



Check the pH of the sample, specified range: pH 3 – 9. If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.




Pipette 1.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 1.0 ml of **Mg-1K** with pipette, close the cell with the screw cap, and mix.



Reaction time: **exactly 3 minutes**



Add 3 drops of **Mg-2K**, close the cell with the screw cap and mix.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a freshly prepared standard solution can be used (see section “Standard solutions”).


Manganese

100816

Cell Test

Measuring 0.10–5.00 mg/l Mn


range: Expression of results also possible in mmol/l.



Check the pH of the sample, specified range: pH 2 – 7. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.




Pipette 7.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 2 drops of **Mn-1K**, close the cell with the screw cap, and mix.




Reaction time: 2 minutes



Add 3 drops of **Mn-2K**, close the cell with the screw cap, and mix.



Reaction time: 5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 30, Cat.No. 114677.

Ready-for-use manganese standard solution Certipur®, Cat.No. 119789, concentration 1000 mg/l Mn, can also be used after diluting accordingly.


To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 30) is highly recommended.

Nickel

114554

Cell Test

Measuring	0.10–6.00 mg/l Ni
range:	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 3–8. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.




Pipette 5.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Reaction time:
1 minute




Add 2 drops of **Ni-1K**, close with the screw cap, and mix.



Add 2 drops of **Ni-2K**, close the cell with the screw cap, and mix.



Reaction time:
2 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

For the determination of **total nickel** a pretreatment with Crack Set 10C, Cat.No. 114688, or Crack Set 10, Cat.No. 114687 and thermoreactor is necessary.

Result can be expressed as sum of nickel (Σ Ni).

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 40, Cat.No. 114692.

A nickel standard solution Titrisol®, Cat.No. 109989, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 40) is highly recommended.

Nickel in electroplating baths

Inherent color


Measuring range: 10–120 g/l Ni



Pipette 5.0 ml of the sample into an empty round cell (Empty cells, Cat.No. 114724).



Add 5.0 ml of **sulfuric acid 40%**, close the cell with the screw cap, and mix.



Place the cell into the cell compartment. Select method no. **57**.

Nitrate


114542

Cell Test

Measuring	0.5 – 18.0 mg/l NO ₃ -N
range:	2.2 – 79.7 mg/l NO ₃
	Expression of results also possible in mmol/l.



Add 1 level yellow micro-spoon of **NO₃-1K** into a reaction cell and close with the screw cap.




Shake the cell **vigorously for 1 minute** to dissolve the solid substance.



Add very slowly 1.5 ml of the sample with pipette, close with the screw cap, and mix **briefly**.
Caution, cell becomes hot!



Reaction time:
10 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 20, Cat.No. 114675, or the Standard solution for photometric applications, CRM, Cat.No. 125037 and 125038.

Ready-for-use nitrate standard solution Certipur®, Cat.No. 119811, concentration 1000 mg/l NO₃⁻, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 20) is highly recommended.

Nitrate


114563

Cell Test

Measuring	0.5 – 25.0 mg/l NO ₃ -N
range:	2.2 – 110.7 mg/l NO ₃
	Expression of results also possible in mmol/l.




Pipette 1.0 ml of the sample into a reaction cell, **do not mix**.



Add 1.0 ml of **NO₃-1K** with pipette, close the cell with the screw cap, and mix. **Caution, cell becomes hot!**



Reaction time:
10 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 20, Cat.No. 114675, or the Standard solution for photometric applications, CRM, Cat.No. 125037 and 125038.

Ready-for-use nitrate standard solution Certipur®, Cat.No. 119811, concentration 1000 mg/l NO₃⁻, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 20) is highly recommended.

Nitrate


114764

Cell Test

Measuring	1.0 – 50.0 mg/l NO ₃ -N
range:	4 – 221 mg/l NO ₃
	Expression of results also possible in mmol/l.




Pipette 0.50 ml of the sample into a reaction cell, **do not mix**.



Add 1.0 ml of **NO₃-1K** with pipette, close the cell with the screw cap, and mix. **Caution, cell becomes hot!**



Reaction time:
10 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 80, Cat.No. 114738, or the Standard solution for photometric applications, CRM, Cat.No. 125037, 125038, and 125039.

Ready-for-use nitrate standard solution Certipur®, Cat.No. 119811, concentration 1000 mg/l NO₃⁻, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 80) is highly recommended.

Nitrate

100614

Cell Test

Measuring	23 – 225 mg/l NO ₃ -N
range:	102 – 996 mg/l NO ₃
Expression of results also possible in mmol/l.	




Pipette 1.0 ml of **NO₃-1K** into a reaction cell, **do not mix**.



Add 0.10 ml of the sample with pipette, close the cell with the screw cap, and mix. **Caution, cell becomes hot!**



Reaction time: 5 minutes, **measure immediately**.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:


To check the measurement system (test reagents, measurement device, and handling) ready-for-use nitrate standard solution Certipur[®], Cat.No. 119811, concentration 1000 mg/l NO₃⁻, can be used after diluting accordingly as well as the Standard solution for photometric applications, CRM, Cat.No. 125039 and 125040.

Nitrite


114547

Cell Test


Measuring	0.010–0.700 mg/l NO ₂ -N
range:	0.03 –2.30 mg/l NO ₂
Expression of results also possible in mmol/l.	



Check the pH of the sample, specified range: pH 2 – 10. If required, add dilute sulfuric acid drop by drop to adjust the pH.




Pipette 5.0 ml of the sample into a reaction cell, close with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time:
10 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use nitrite standard solution Certipur[®], Cat.No. 119899, concentration 1000 mg/l NO₂⁻, can be used after diluting accordingly as well as the Standard solution for photometric applications, CRM, Cat.No. 125041.

Nitrite

100609

Cell Test

Measuring	1.0 – 90.0 mg/l NO ₂ -N
range:	3 – 296 mg/l NO ₂
	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 1 – 12. If required, add dilute sulfuric acid drop by drop to adjust the pH.



Add 2 level blue microspoons of **NO₂-1K** into a reaction cell.




Add 8.0 ml of the sample with pipette and close with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 20 minutes, **measure immediately**. **Do not shake or swirl** the cell before the measurement.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use nitrite standard solution Certipur[®], Cat.No. 119899, concentration 1000 mg/l NO₂⁻, can be used after diluting accordingly as well as the Standard solution for photometric applications, CRM, Cat.No. 125042.


Nitrogen (total)

114537


Cell Test

Measuring 0.5 – 15.0 mg/l N


range: Expression of results also possible in mmol/l.




Pipette 10 ml of the sample into an empty round cell (Empty cells, Cat.No. 114724).




Add 1 level blue micro-spoon of **N-1K**.




Add 6 drops of **N-2K**, close the cell with the screw cap, and mix.



Heat the cell in the thermoreactor at 120 °C (100 °C) for 1 hour.




Remove the cell from the thermoreactor and place in a test-tube rack to cool to room temperature: **pretreated sample**.



Swirl the cell after 10 minutes.



Add 1 level yellow micro-spoon of **N-3K** into a **reaction cell**, close the cell with the screw cap.




Shake the cell **vigorously for 1 minute** to dissolve the solid substance.



Add very slowly 1.5 ml of the **pretreated sample** with pipette, close the cell with the screw cap, and mix **briefly**. **Caution, cell becomes hot!**



Reaction time: 10 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 50, Cat.No. 114695, or the Standard solution for photometric applications, CRM, Cat.No. 125043 and 125044.


To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 50) is highly recommended.

Nitrogen (total)


100613

Cell Test


Measuring 0.5 – 15.0 mg/l N
range: Expression of results also possible in mmol/l.




Pipette 10 ml of the sample into an empty round cell (Empty cells, Cat.No. 114724).




Add 1 level blue micro-spoon of **N-1K**.




Add 6 drops of **N-2K**, close the cell with the screw cap, and mix.



Heat the cell in the thermoreactor at 120 °C (100 °C) for 1 hour.




Remove the cell from the thermoreactor and place in a test-tube rack to cool to room temperature: **pretreated sample**.




Swirl the cell after 10 minutes.




Pipette 1.0 ml of the **pretreated sample** into a reaction cell, **do not mix!**



Add 1.0 ml of **N-3K** with pipette, close the cell with the screw cap, and mix. **Caution, cell becomes hot!**



Reaction time: 10 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 50, Cat.No. 114695, or the Standard solution for photometric applications, CRM, Cat.No. 125043 and 125044.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 50) is highly recommended.


Nitrogen (total)

114763


Cell Test

Measuring 10–150 mg/l N


range: Expression of results also possible in mmol/l.




Pipette 1.0 ml of the sample into an empty round cell (Empty cells, Cat.No. 114724).




Add 9.0 ml of distilled water (Water for analysis EMSURE®, Cat.No. 116754, is recommended) with pipette.




Add 1 level blue micro-spoon of **N-1K**.




Add 6 drops of **N-2K**, close the cell with the screw cap, and mix.




Heat the cell in the thermoreactor at 120 °C (100 °C) for 1 hour.




Remove the cell from the thermoreactor and place in a test-tube rack to cool to room temperature: **pretreated sample**.



Swirl the cell after 10 minutes.




Pipette 1.0 ml of the **pretreated sample** into a reaction cell, **do not mix!**



Add 1.0 ml of **N-3K** with pipette, close the cell with the screw cap, and mix. **Caution, cell becomes hot!**



Reaction time: 10 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 70, Cat.No. 114689, or the Standard solution for photometric applications, CRM, Cat.No. 125044 and 125045.


To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 70) is highly recommended.

Oxygen

114694

Cell Test

Measuring 0.5–12.0 mg/l O₂
range: Expression of results also possible in mmol/l.



Check the pH of the sample, specified range: pH 6 – 8. If required, add dilute sodium hydroxide solution or nitric acid drop by drop to adjust the pH.



Fill watersample into a reaction cell to overflowing and make sure, that no air bubbles are present.



Place the filled cell in a test-tube rack.



Add with microspoon 1 glass bead.



Add 5 drops of O₂-1K.




Add 5 drops of O₂-2K, close the cell with the screw cap, and shake for 10 seconds.



Reaction time:
1 minute



Add 10 drops of O₂-3K, close the cell with the screw cap, mix, and clean from outside.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a oxygen standard solution must be prepared (application see the website).


Measuring range: pH 6.4 – 8.8



Pipette 10 ml of the sample into a round cell.



Add 4 drops of **pH-1**, close the cell with the screw cap, and mix.
Attention!
The reagent bottle must be held **vertically by all means!**



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) buffer solution pH 7.00 Certipur[®], Cat.No. 109407, can be used.


Phosphate

100474

Determination of orthophosphate

Cell Test


Measuring	0.05 – 5.00 mg/l PO ₄ -P
range:	0.2 – 15.3 mg/l PO ₄
	0.11 – 11.46 mg/l P ₂ O ₅
	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 0 – 10. If required, add dilute sulfuric acid drop by drop to adjust the pH.




Pipette 5.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 5 drops of **P-1K**, close the cell with the screw cap, and mix.




Add 1 dose of **P-2K** using the blue dose-metering cap, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time:
5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

For the determination of **total phosphorus = sum of orthophosphate, polyphosphate and organophosphate** either Phosphate Cell Test, Cat. No. 114543, 114729, and 100673 or Phosphate Test, Cat. No. 114848 in conjunction with Crack Set 10/10C, Cat. No. 114687 resp. 114688 can be used.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 10, Cat.No. 114676.

Ready-for-use phosphate standard solution Certipur®, Cat.No. 119898, concentration 1000 mg/l PO₄³⁻, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 10) is highly recommended.

Phosphate

114543

Determination of orthophosphate

Cell Test


Measuring	0.05 – 5.00 mg/l PO ₄ -P
range:	0.2 – 15.3 mg/l PO ₄
	0.11 – 11.46 mg/l P ₂ O ₅
	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 0 – 10. If required, add dilute sulfuric acid drop by drop to adjust the pH.




Pipette 5.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 5 drops of **P-2K**, close the cell with the screw cap, and mix.




Add 1 dose of **P-3K** using the blue dose-metering cap, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time:
5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 10, Cat.No. 114676.

Ready-for-use phosphate standard solution Certipur®, Cat.No. 119898, concentration 1000 mg/l PO₄³⁻, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 10) is highly recommended.

Phosphate

Determination of total phosphorus
= sum of orthophosphate, polyphosphate, and organophosphate

114543


Cell Test

Measuring 0.05 – 5.00 mg/l P

range: 0.2 – 15.3 mg/l PO₄

0.11 – 11.46 mg/l P₂O₅


Expression of results also possible in mmol/l and also in P total (Σ P), and P org* [P(o)].




Check the pH of the sample, specified range: pH 0 – 10.
If required, add dilute sulfuric acid drop by drop to adjust the pH.




Pipette 5.0 ml of the sample into a reaction cell, close with the screw cap, and mix.




Add 1 dose of **P-1K** using the green dose-metering cap, close the cell with the screw cap.




Heat the cell in the thermoreactor at 120 °C (100 °C) for 30 minutes.




Remove the cell from the thermoreactor and place in a test-tube rack to cool to room temperature.



Add 5 drops of **P-2K**, close the cell with the screw cap, and mix.




Add 1 dose of **P-3K** using the blue dose-metering cap, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time:
5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

A differentiation between orthophosphate (PO₄-P) and P org* (P(o)) can be performed on the photometer. Prior to measuring, select the differentiation measurement and choose the corresponding citation form. Then measure the P total, press enter and measure the orthophosphate (see analytical procedure for orthophosphate). After pressing enter, the individual measuring values for PO₄-P and P(o) are shown on the display.

* P org is the sum of polyphosphate and organophosphate.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 10, Cat.No. 114676, or the Standard solution for photometric applications, CRM, Cat.No. 125046 and 125047.

Ready-for-use phosphate standard solution Certipur®, Cat.No. 119898, concentration 1000 mg/l PO₄³⁻, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 10) is highly recommended.

Phosphate

100475

Determination of orthophosphate

Cell Test


Measuring	0.5–25.0 mg/l PO ₄ -P
range:	1.5–76.7 mg/l PO ₄
	1.1–57.3 mg/l P ₂ O ₅
	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 0–10. If required, add dilute sulfuric acid drop by drop to adjust the pH.




Pipette 1.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 5 drops of **P-1K**, close the cell with the screw cap, and mix.




Add 1 dose of **P-2K** using the blue dose-metering cap, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time:
5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

For the determination of **total phosphorus = sum of orthophosphate, polyphosphate and organophosphate** either Phosphate Cell Test, Cat. No. 114543, 114729, and 100673 or Phosphate Test, Cat. No. 114848 in conjunction with Crack Set 10/10C, Cat. No. 114687 resp. 114688 can be used.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 20 and 80, Cat.Nos. 114675 and 114738.

Ready-for-use phosphate standard solution Certipur®, Cat.No. 119898, concentration 1000 mg/l PO₄³⁻, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck) is highly recommended.

Phosphate

114729

Determination of orthophosphate

Cell Test


Measuring	0.5–25.0 mg/l PO ₄ -P
range:	1.5–76.7 mg/l PO ₄
	1.1–57.3 mg/l P ₂ O ₅
	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 0–10. If required, add dilute sulfuric acid drop by drop to adjust the pH.




Pipette 1.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 5 drops of **P-2K**, close the cell with the screw cap, and mix.




Add 1 dose of **P-3K** using the blue dose-metering cap, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time:
5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 20 and 80, Cat.Nos. 114675 and 114738.

Ready-for-use phosphate standard solution Certipur®, Cat.No. 119898, concentration 1000 mg/l PO₄³⁻, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck) is highly recommended.

Phosphate

Determination of total phosphorus
= sum of orthophosphate, polyphosphate, and organophosphate

114729


Cell Test

Measuring 0.5–25.0 mg/l P

range: 1.5–76.7 mg/l PO₄

1.1–57.3 mg/l P₂O₅


Expression of results also possible in mmol/l and also in P total (ΣP), and P org* [P(o)].




Check the pH of the sample, specified range: pH 0–10. If required, add dilute sulfuric acid drop by drop to adjust the pH.




Pipette 1.0 ml of the sample into a reaction cell, close with the screw cap, and mix.




Add 1 dose of **P-1K** using the green dose-metering cap, close the cell with the screw cap.




Heat the cell in the thermoreactor at 120 °C (100 °C) for 30 minutes.




Remove the cell from the thermoreactor and place in a test-tube rack to cool to room temperature.



Add 5 drops of **P-2K**, close the cell with the screw cap, and mix.




Add 1 dose of **P-3K** using the blue dose-metering cap, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

A differentiation between orthophosphate (PO₄-P) and P org* (P(o)) can be performed on the photometer. Prior to measuring, select the differentiation measurement and choose the corresponding citation form. Then measure the P total, press enter and measure the orthophosphate (see analytical procedure for orthophosphate). After pressing enter, the individual measuring values for PO₄-P and P(o) are shown on the display.

* P org is the sum of polyphosphate and organophosphate.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 20 and 80, Cat.No. 114675 and 114738, or as well as the Standard solution for photometric applications, CRM, Cat.No. 125047 and 125048.

Ready-for-use phosphate standard solution Certipur®, Cat.No. 119898, concentration 1000 mg/l PO₄³⁻, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck) is highly recommended.


Phosphate

100616

Determination of orthophosphate

Cell Test

Measuring	3.0 – 100.0 mg/l PO ₄ -P
range:	9 – 307 mg/l PO ₄
	7 – 229 mg/l P ₂ O ₅
	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 0 – 10. If required, add dilute sulfuric acid drop by drop to adjust the pH.




Pipette 0.20 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 5 drops of **PO₄-1K**, close the cell with the screw cap, and mix.




Add 1 dose of **PO₄-2K** using the blue dose-metering cap, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time:
5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

For the determination of **total phosphorus = sum of orthophosphate, polyphosphate and organophosphate** either Phosphate Cell Test, Cat. No. 114543, 114729, and 100673 or Phosphate Test, Cat. No. 114848 in conjunction with Crack Set 10/10C, Cat. No. 114687 resp. 114688 can be used.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use phosphate standard solution Certipur®, Cat.No. 119898, concentration 1000 mg/l PO₄³⁻, can be used after diluting accordingly.

Phosphate

100673

Determination of orthophosphate

Cell Test


Measuring	3.0 – 100.0 mg/l PO ₄ -P
range:	9 – 307 mg/l PO ₄
	7 – 229 mg/l P ₂ O ₅
	Expression of results also possible in mmol/l.




Check the pH of the sample, specified range: pH 0 – 10. If required, add dilute sulfuric acid drop by drop to adjust the pH.




Pipette 0.20 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 5 drops of **P-2K**, close the cell with the screw cap, and mix.




Add 1 dose of **P-3K** using the blue dosing cap, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time:
5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use phosphate standard solution Certipur®, Cat.No. 119898, concentration 1000 mg/l PO₄³⁻, can be used after diluting accordingly.

Phosphate

Determination of total phosphorus
= sum of orthophosphate, polyphosphate, and organophosphate

100673


Cell Test

Measuring 3.0 – 100.0 mg/l PO₄-P

range: 9 – 307 mg/l PO₄

7 – 229 mg/l P₂O₅


Expression of results also possible in mmol/l and also in P total (ΣP), and P org* [P(o)].




Check the pH of the sample, specified range: pH 0 – 10. If required, add dilute sulfuric acid drop by drop to adjust the pH.




Pipette 0.20 ml of the sample into a reaction cell, close with the screw cap, and mix.




Add 1 dose of **P-1K** using the green dose-metering cap, close the cell with the screw cap.




Heat the cell in the thermoreactor at 120 °C (100 °C) for 30 minutes.




Remove the cell from the thermoreactor and place in a test-tube rack to cool to room temperature.



Add 5 drops of **P-2K**, close the cell with the screw cap, and mix.




Add 1 dose of **P-3K** using the blue dose-metering cap, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

A differentiation between orthophosphate (PO₄-P) and P org* (P(o)) can be performed on the photometer. Prior to measuring, select the differentiation measurement and choose the corresponding citation form. Then measure the P total, press enter and measure the orthophosphate (see analytical procedure for orthophosphate). After pressing enter, the individual measuring values for PO₄-P and P(o) are shown on the display.

* P org is the sum of polyphosphate and organophosphate.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use phosphate standard solution Certipur[®], Cat.No. 119898, concentration 1000 mg/l PO₄³⁻, can be used after diluting accordingly as well as the Standard solution for photometric applications, CRM, Cat.No. 125047, 125048, and 125049.

Phosphate

114546

Determination of orthophosphate

Cell Test


Measuring	0.5 – 25.0 mg/l PO ₄ -P
range:	1.5 – 76.7 mg/l PO ₄
	1.1 – 57.3 mg/l P ₂ O ₅
	Expression of results also possible in mmol/l.



Check the pH of the sample, specified range: pH 0 – 10.
If required, add dilute sulfuric acid drop by drop to adjust the pH.



Pipette 5.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

For the determination of **total phosphorus = sum of orthophosphate, polyphosphate and organophosphate** either Phosphate Cell Test, Cat. No. 114543, 114729, and 100673 or Phosphate Test, Cat. No. 114848 in conjunction with Crack Set 10/10C, Cat. No. 114687 resp. 114688 can be used.

Quality assurance:


To check the measurement system (test reagents, measurement device, and handling) ready-for-use phosphate standard solution Certipur[®], Cat.No. 119898, concentration 1000 mg/l PO₄³⁻, can be used after diluting accordingly.

Potassium


114562

Cell Test

Measuring	5.0 – 50.0 mg/l K
range:	Expression of results also possible in mmol/l.



Filter turbid samples.




Check the pH of the sample, specified range: pH 3 – 12.
If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Pipette 2.0 ml of the sample into a reaction cell, close with the screw cap, and mix.




Check the pH, specified range: pH 10.0 – 11.5.



Add 6 drops of **K-1K**, close the cell with the screw cap, and mix.




Add 1 level blue micro-spoon of **K-2K**, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time:
5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use potassium standard solution Certipur®, Cat.No. 170230, concentration 1000 mg/l K, can be used after diluting accordingly.


Potassium

100615


Cell Test

Measuring 30–300 mg/l K

range: Expression of results also possible in mmol/l.



Filter turbid samples.




Check the pH of the sample, specified range: pH 3 – 12.
If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.



Pipette 0.50 ml of the sample into a reaction cell, close with the screw cap, and mix.




Check the pH, specified range: pH 10.0 – 11.5.



Add 6 drops of **K-1K**, close the cell with the screw cap, and mix.




Add 1 level blue micro-spoon of **K-2K**, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time:
5 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use potassium standard solution Certipur®, Cat.No. 170230, concentration 1000 mg/l K, can be used after diluting accordingly.


Residual Hardness

114683

Cell Test

Measuring	0.50 – 5.00 mg/l Ca
range:	0.070 – 0.700 °d
	0.087 – 0.874 °e
	0.12 – 1.25 °f

Measuring	0.70 – 7.00 mg/l CaO
range:	1.2 – 12.5 mg/l CaCO ₃
Expression of results also possible in mmol/l.	



Check the pH of the sample, specified range: pH 5–8.
If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.




Pipette 4.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 0.20 ml of **RH-1K**, close the cell with the screw cap, and mix.



Reaction time: 10 minutes, **measure immediately**.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use calcium standard solution Certipur[®], Cat.No. 119778, concentration 1000 mg/l Ca, can be used after diluting accordingly. (Pay attention to pH value!)

Sodium

in nutrient solutions

100885

Cell Test

Measuring 10–300 mg/l Na

range: Expression of results also possible in mmol/l.




Pipette 0.50 ml of **Na-1K** into a reaction cell and mix.



Add 0.50 ml of the sample with pipette, close the cell with the screw cap, and mix.



Reaction time:
1 minute



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use chloride standard solution Certipur[®], Cat.No. 119897, concentration 1000 mg/l Cl⁻ (corresponds to 649 mg/l Na), can be used after diluting accordingly (see section "Standard solutions").


Sulfate

102532


Cell Test

Measuring 1.0–50.0 mg/l SO₄


range: Expression of results also possible in mmol/l.



Filter turbid samples.




Check the pH of the sample, specified range: pH 2–10. If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.



Pipette 10 ml of the sample into a reaction cell, close with the screw cap, and mix.




Add 1 level green micro-spoon of SO₄-1K, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 2 minutes, **measure immediately**.




Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.


Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use sulfate standard solution Certipur[®], Cat.No. 119813, concentration 1000 mg/l SO₄²⁻, can be used after diluting accordingly.


Measuring	5–250 mg/l SO ₄
range:	Expression of results also possible in mmol/l.



Filter turbid samples.




Check the pH of the sample, specified range: pH 2–10. If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.



Pipette 5.0 ml of the sample into a reaction cell, close with the screw cap, and mix.




Add 1 level green microspoon of SO₄-1K, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 2 minutes, **measure immediately**.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 10, Cat.No. 114676, or the Standard solution for photometric applications, CRM, Cat.No. 125050 and 125051.

Ready-for-use sulfate standard solution Certipur®, Cat.No. 119813, concentration 1000 mg/l SO₄²⁻, can also be used after diluting accordingly.


To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 10) is highly recommended.

Sulfate


100617

Cell Test

Measuring 50 – 500 mg/l SO_4
range: Expression of results also possible in mmol/l.



Filter turbid samples.




Check the pH of the sample, specified range: pH 2–10. If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.



Pipette 2.0 ml of the sample into a reaction cell, close with the screw cap, and mix.




Add 1 level green microspoon of $\text{SO}_4\text{-1K}$, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 2 minutes, **measure immediately**.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 10, Cat.No. 114676, or the Standard solution for photometric applications, CRM, Cat.No. 125051 and 125052.

Ready-for-use sulfate standard solution Certipur®, Cat.No. 119813, concentration 1000 mg/l SO_4^{2-} , can also be used after diluting accordingly.


To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 10) is highly recommended.

Sulfate


114564

Cell Test


Measuring 100–1000 mg/l SO₄
range: Expression of results also possible in mmol/l.



Filter turbid samples.




Check the pH of the sample, specified range: pH 2–10. If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.



Pipette 1.0 ml of the sample into a reaction cell, close with the screw cap, and mix.




Add 1 level green micro-spoon of SO₄-1K, close the cell with the screw cap.



Shake the cell vigorously to dissolve the solid substance.



Reaction time: 2 minutes, **measure immediately**.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 20, Cat.No. 114675, or the Standard solution for photometric applications, CRM, Cat.No. 125051, 125052 and 125053.

Ready-for-use sulfate standard solution Certipur®, Cat.No. 119813, concentration 1000 mg/l SO₄²⁻, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 20) is highly recommended.

Surfactants (nonionic)

101787

Cell Test

Measuring 0.010–7.50 mg/l surfactants (nonionic)


range: (calculated as Triton® X-100)



Check the pH of the sample, specified range: pH 3–9. If required, add dilute sodium hydroxide solution or sulfuric acid drop by drop to adjust the pH.




Pipette 4.0 ml of the sample into a reaction cell. Close with the screw cap.




Shake the cell for **1 minute vigorously**.



Reaction time: 2 minutes



Swirl the cell before measurement.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:


To check the measurement system (test reagents, measurement device, and handling) a surfactants standard solution must be prepared from Triton® X-100, Cat.No. 112298 (see section “Standard solutions”).

Suspended Solids


Measuring range: 50 – 750 mg/l of suspended solid



Homogenize 500 ml of sample for 2 minutes in a mixer running at high speed.



Transfer the solution into a cell.



Place the cell into the cell compartment, select method no. **182**.

TOC

Total Organic Carbon

114878


Cell Test

Measuring range: 5.0 – 80.0 mg/l TOC


Removal of inorganic bound carbon (TIC):




Check the pH of the sample, specified range: pH 2– 12.
If required, add dilute sulfuric acid drop by drop to adjust the pH.



Place 25 ml of the sample into a suitable glass vessel.



Add 3 drops of **TOC-1K** and mix.



Check the pH, specified range pH < 2.5.



Stir for 10 minutes.


Preparation of measurement sample :




Pipette 3.0 ml of stirred sample into a reaction cell.




Add 1 level grey micro-spoon of **TOC-2K**. **Immediately** close the cell tightly with an **aluminium cap** (Cat.No. 173500).



Heat the cell, standing on its head, at 120 °C in the thermoreactor for 2 hours.



Remove the cell from the thermoreactor and let it, **standing on its head**, to cool for 1 hour.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a TOC standard solution Certipur®, Cat.No. 109017, concentration 1000 mg/l TOC, can be used after diluting accordingly.

TOC

Total Organic Carbon

114879


Cell Test

Measuring range: 50 – 800 mg/l TOC


Removal of inorganic bound carbon (TIC):




Check the pH of the sample, specified range: pH 2– 12. If required, add dilute sulfuric acid drop by drop to adjust the pH.



Pipette 1.0 ml of the sample and 9.0 ml of distilled water (Water for chromatography LiChrosolv®, Cat.No. 115333, is recommended) into a suitable glass vessel.



Add 2 drops of **TOC-1K** and mix.



Check the pH, specified range pH < 2.5



Stir for 10 minutes.


Preparation of measurement sample :




Pipette 3.0 ml of stirred sample into a reaction cell.




Add 1 level grey micro-spoon of **TOC-2K**. **Immediately** close the cell tightly with an **aluminium cap** (Cat.No. 173500).



Heat the cell, standing on its head, at 120 °C in the thermoreactor for 2 hours.



Remove the cell from the thermoreactor and let it, **standing on its head**, to cool for 1 hour.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a TOC standard solution Certipur®, Cat.No. 109017, concentration 1000 mg/l TOC, can be used after diluting accordingly.

Total Hardness

100961

Determination of total hardness

Cell Test

Measuring 5 –215 mg/l Ca

range: 0.7 – 30.1 °d

0.9 – 37.6 °e


1.2 – 53.7 °f

Measuring 7 –301 mg/l CaO

range: 12 –537 mg/l CaCO₃

Expression of results also possible in mmol/l

and also in mg/l Mg .



Check the pH of the sample, specified range: pH 3 – 9.
If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.




Pipette 1.0 ml of the sample into a reaction cell, close with the screw cap, and mix.



Add 1.0 ml of **H-1K** with pipette, close the cell with the screw cap, and mix.



Reaction time:
3 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a freshly prepared standard solution can be used (see section “Standard solutions”).

Total Hardness

100961


Differentiation between Ca- and Mg-hardness

Cell Test

Measuring	0.12 – 5.36 mmol/l
range:	0.7 – 30.1 °d
	0.9 – 37.6 °e
	1.2 – 53.7 °f

Differentiation possible only in mmol/l.

A differentiation between calcium- and magnesium-hardness can be performed on the photometer. Prior to measuring, select the differentiation measurement and choose the corresponding citation form.



Check the pH of the sample, specified range: pH 3 – 9.
If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.




Pipette 1.0 ml of the sample into a reaction cell, close with the screw cap, and mix.




Add 1.0 ml of **H-1K** with pipette, close the cell with the screw cap, and mix.



Reaction time:
3 minutes




Place the cell into the cell compartment. Align the mark on the cell with that on the photometer = **Result total hardness**




Press enter,
remove the cell.



Add 3 drops of **H-2K** to the already measured cell, close the cell with the screw cap, and mix.



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer = **Result magnesium**



After pressing enter, the individual measuring values for Ca- and Mg-hardness are shown on the display.


Volatile Organic Acids

101763


Cell Test

Measuring 50 – 3000 mg/l volatile organic acid


range: (calculated as acetic acid)



Check the pH of the sample, specified range: pH 2– 12.




Pipette 0.75 ml of **OA-1** into a round cell.




Add 2 drops of **OA-2**.




Add 0.50 ml of the sample with pipette, close with the screw cap, and mix.



Heat the cell in the thermoreactor at 100 °C for 10 minutes. Then cool to room temperature under running water.



Add 5 drops of **OA-3**.



Add 0.50 ml of **OA-4** with pipette, close the cell with the screw cap, and mix.




Reaction time: 3 minutes



Add 5.0 ml of **OA-5** with pipette, close the cell with the screw cap, and shake vigorously.



Reaction time: 10 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:


To check the measurement system (test reagents, measurement device, and handling) a standard solution must be prepared from sodium acetate anhydrous, Cat.No. 106268 (see section “Standard solutions”).

Volatile Organic Acids


101749

Cell Test

Measuring	50 – 3000 mg/l volatile organic acid	(calculated as acetic acid)
range:	71 – 4401 mg/l volatile organic acid	(calculated as butyric acid)




Check the pH of the sample, specified range: pH 2– 12.



Pipette 0.50 ml of **OA-1** into a round cell.




Add 0.50 ml of the sample with pipette, close with the screw cap, and mix.



Heat the cell in the thermoreactor at 100 °C for 15 minutes. Then cool to room temperature under running water.



Add 1.0 ml of **OA-2** with pipette.




Add 1.0 ml of **OA-3** with pipette, close the cell with the screw cap, and mix.



Add 1.0 ml of **OA-4** with pipette, close the cell with the screw cap, and shake vigorously.



Reaction time:
1 minute



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a standard solution must be prepared from sodium acetate anhydrous, Cat.No. 106268 (see section “Standard solutions”).

Volatile Organic Acids


101809

Test


Measuring	50 – 3000 mg/l volatile organic acid	(calculated as acetic acid)
range:	71 – 4401 mg/l volatile organic acid	(calculated as butyric acid)



Check the pH of the sample, specified range: pH 2– 12.




Pipette 0.75 ml of **OA-1** into a round cell.



Add 0.50 ml of **OA-2** with pipette.



Add 0.50 ml of the sample with pipette, close with the screw cap, and mix.



Heat the cell in the thermoreactor at 100 °C for 15 minutes. Then cool to room temperature under running water.



Add 1.0 ml of **OA-3** with pipette.




Add 1.0 ml of **OA-4** with pipette, close the cell with the screw cap, and mix.



Add 1.0 ml of **OA-5** with pipette, close the cell with the screw cap, and shake vigorously.



Reaction time:
1 minute



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) a standard solution must be prepared from sodium acetate anhydrous, Cat.No. 106268 (see section “Standard solutions”).

Zinc

100861


Cell Test

Measuring 0.025 – 1.000 mg/l Zn


range: Expression of results also possible in mmol/l.



Check the pH of the sample, specified range: pH 1–7.
If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.



Pipette 10 ml of sample into a glass vessel.




Add 1 level green microspoon of **Zn-1K** and shake to dissolve the solid substance: **sample-reagent mixture**.



Pipette 0.50 ml of **Zn-2K** into a reaction cell, close with the screw cap, and mix.




Add 2.0 ml of the **sample-reagent mixture** with pipette, close the cell with the screw cap, and mix.



Add 5 drops of **Zn-3K**, close the cell with the screw cap, and mix.



Reaction time: 15 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

For the determination of **total zinc** a pretreatment with Crack Set 10C, Cat.No. 114688, or Crack Set 10, Cat.No. 114687, and thermoreactor is necessary.

Result can be expressed as sum of zinc (Σ Zn).

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) ready-for-use zinc standard solution Certipur[®], Cat.No. 119806, concentration 1000 mg/l Zn, can be used after diluting accordingly.

Zinc

114566


Cell Test

Measuring 0.20–5.00 mg/l Zn

range: Expression of results also possible in mmol/l.



Check the pH of the sample, specified range: pH 3 – 10. If required, add dilute sodium hydroxide solution or hydrochloric acid drop by drop to adjust the pH.



Add 5 drops of **Zn-1K** into a reaction cell, close with the screw cap, and mix.




Add 0.50 ml of the sample with pipette, close the cell with the screw cap, and mix.



Add 5 drops of **Zn-2K**, close the cell with the screw cap, and mix.



Reaction time: 15 minutes



Place the cell into the cell compartment. Align the mark on the cell with that on the photometer.

Important:

For the determination of **total zinc** a pretreatment with Crack Set 10C, Cat.No. 114688, or Crack Set 10, Cat.No. 114687, and thermoreactor is necessary.

Result can be expressed as sum of zinc (Σ Zn).

Quality assurance:

To check the measurement system (test reagents, measurement device, and handling) we recommended to use Spectroquant® CombiCheck 40, Cat.No. 114692.

Ready-for-use zinc standard solution Certipur®, Cat.No. 119806, concentration 1000 mg/l Zn, can also be used after diluting accordingly.

To check for sample-dependent effects the use of addition solutions (e.g. in CombiCheck 40) is highly recommended.

Suitability of Test Kits for Testing Seawater and Tolerance Limits of Neutral Salts

Test kit	Cat. No.	Seawater	Limit of tolerance, salts in %		
			NaCl	NaNO ₃	Na ₂ SO ₄
Acid Capacity Cell Test	101758	no	–	–	–
Aluminium Cell Test	100594	yes	20	20	20
Ammonium Cell Test	114739	no	5	5	5
Ammonium Cell Test	114558	yes	20	10	15
Ammonium Cell Test	114544	yes	20	15	20
Ammonium Cell Test	114559	yes	20	20	20
AOX Cell Test	100675	no	0.4	20	20
BOD Cell Test	100687	yes	20	20	20
Cadmium Cell Test	114834	no	1	10	1
Calcium Cell Test	100858	no	2	2	1
Chloride Cell Test	114730	yes	–	20	1
Chloride Cell Test	101804	no	–	0.5	0.05
Chlorine Cell Test	100595	no	10	10	10
Chlorine Cell Test	100597	no	10	10	10
Chlorine reagents (liquid) (free and total)	100086/100087/ 100088	no	10	10	10
Chromate Cell Test (chromium(VI))	114552	yes	10	10	10
Chromate Cell Test (chromium total)	114552	no	1	10	10
COD Cell Test	114560	no	0.4	10	10
COD Cell Test	101796	no	0.4	10	10
COD Cell Test	114540	no	0.4	10	10
COD Cell Test	114895	no	0.4	10	10
COD Cell Test	114690	no	0.4	20	20
COD Cell Test	114541	no	0.4	10	10
COD Cell Test	114691	no	0.4	20	20
COD Cell Test	114555	no	1.0	10	10
COD Cell Test	101797	no	10	20	20
COD Cell Test (Hg free)	109772	no	0	10	10
COD Cell Test (Hg free)	109773	no	0	10	10
COD Cell Test (seawater)	117058	yes	35	10	10
COD Cell Test (seawater)	117059	yes	35	10	10
Copper Cell Test	114553	yes	15	15	15
Cyanide Cell Test	102531	no	10	10	10
Cyanide Cell Test	114561	no	10	10	10
Formaldehyde Cell Test	114500	no	5	0	10
Hardness, see Total Hardness Cell Test					
Iron Cell Test	114549	yes	20	20	20
Iron Cell Test	114896	no	5	5	5
Lead Cell Test	114833	no	20	20	1
Magnesium Cell Test	100815	yes	2	2	1
Manganese Cell Test	100816	no	20	20	20
Nickel Cell Test	114554	no	20	20	20
Nitrate Cell Test	114542	no	0.4	–	20
Nitrate Cell Test	114563	no	0.2	–	20
Nitrate Cell Test	114764	no	0.5	–	20
Nitrate Cell Test	100614	no	2	–	20
Nitrite Cell Test	114547	yes	20	20	15
Nitrite Cell Test	100609	yes	20	20	15
Nitrogen (total) Cell Test	114537	no	0.5	–	10
Nitrogen (total) Cell Test	100613	no	0.2	–	10
Nitrogen (total) Cell Test	114763	no	2	–	20
Oxygen Cell Test	114694	no	10	5	1
pH Cell Test	101744	yes	–	–	–
Phosphate Cell Test	100474	yes	5	10	10
Phosphate Cell Test (orthophosphates)	114543	yes	5	10	10
Phosphate Cell Test (phosphorus total)	114543	no	1	10	10
Phosphate Cell Test	100475	yes	20	20	20
Phosphate Cell Test (orthophosphates)	114729	yes	20	20	20
Phosphate Cell Test (phosphorus total)	114729	yes	5	20	20
Phosphate Cell Test	100616	yes	20	20	20
Phosphate Cell Test (orthophosphates)	100673	yes	20	20	20
Phosphate Cell Test (phosphorus total)	100673	yes	20	20	20
Phosphate Cell Test	114546	yes	20	20	20

Suitability of Test Kits for Testing Seawater and Tolerance Limits of Neutral Salts

Test kit	Cat. No.	Seawater	Limit of tolerance, salts in %		
			NaCl	NaNO ₃	Na ₂ SO ₄
Potassium Cell Test	114562	yes	20	20	20
Potassium Cell Test	100615	yes	20	20	20
Residual Hardness Cell Test	114683	no	0.01	0.01	0.01
Sodium Cell Test	100885	no	–	10	1
Sulfate Cell Test	102532	no	2	0.007	–
Sulfate Cell Test	114548	yes	10	0.1	–
Sulfate Cell Test	100617	yes	10	0.1	–
Sulfate Cell Test	114564	yes	10	0.5	–
Surfactants (nonionic) Cell Test	101787	no	2	5	2
TOC Cell Test	114878	no	0.5	10	10
TOC Cell Test	114879	no	5	20	20
Total Hardness Cell Test	100961	no	2	2	1
Volatile Organic Acids Cell Test	101763	no	20	20	10
Volatile Organic Acids Cell Test	101749	no	20	20	10
Volatile Organic Acids Test	101809	no	20	20	10
Zinc Cell Test	100861	no	20	20	1
Zinc Cell Test	114566	no	10	10	10

Spectroquant® CombiCheck and Standard Solutions

Test kit, Cat. No. or method	Evalu- ation as	CombiCheck, Cat. No.	Confidence interval		Diluted and ready-to-use standard solutions, CRM			Ready-to-use standard solution, Cat. No.
			Spec. value for the standard	max. working tolerance	Cat. No.	concen- tration	expanded measurement uncertainty	
Acid Capacity Cell Test, 101758	OH	–	5.00 mmol/l*	± 0.50 mmol/l	–	–	–	see prep. instr.
Aluminium Cell Test, 100594	Al	–	0.25 mg/l*	± 0.03 mg/l	–	–	–	119770
Ammonium Cell Test, 114739	NH ₄ -N	CombiCheck 50, 114695	1.00 mg/l	± 0.10 mg/l	125022	0.400 mg/l	± 0.012 mg/l	–
					125023	1.00 mg/l	± 0.04 mg/l	119812
Ammonium Cell Test, 114558	NH ₄ -N	CombiCheck 10, 114676	4.00 mg/l	± 0.30 mg/l	125022	0.400 mg/l	± 0.012 mg/l	–
					125023	1.00 mg/l	± 0.04 mg/l	–
					125024	2.00 mg/l	± 0.07 mg/l	–
					125025	6.00 mg/l	± 0.13 mg/l	119812
Ammonium Cell Test, 114544	NH ₄ -N	CombiCheck 20, 114675	12.0 mg/l	± 1.0 mg/l	125023	1.00 mg/l	± 0.04 mg/l	–
					125024	2.00 mg/l	± 0.07 mg/l	–
					125025	6.00 mg/l	± 0.13 mg/l	–
					125026	12.0 mg/l	± 0.4 mg/l	119812
Ammonium Cell Test, 114559	NH ₄ -N	CombiCheck 70, 114689	50.0 mg/l	± 5.0 mg/l	125025	6.00 mg/l	± 0.13 mg/l	–
					125026	12.0 mg/l	± 0.4 mg/l	–
					125027	50.0 mg/l	± 1.2 mg/l	119812
AOX Cell Test, 100675	AOX	–	1.00 mg/l*	± 0.10 mg/l	–	–	–	100680
BOD Cell Test, 100687	O ₂	–	210 mg/l	± 20 mg/l	–	–	–	100718
Cadmium Cell Test, 114834	Cd	CombiCheck 30, 114677	0.500 mg/l	± 0.060 mg/l	–	–	–	119777
Calcium Cell Test, 100858	Ca	–	75 mg/l*	± 7 mg/l	–	–	–	see prep. instr.
Chloride Cell Test, 114730	Cl	CombiCheck 20, 114675	60 mg/l	± 10 mg/l	–	–	–	–
		CombiCheck 10, 114676	25 mg/l	± 6 mg/l	–	–	–	119897
Chloride Cell Test, 101804	Cl	–	7.5 mg/l*	± 0.8 mg/l	–	–	–	119897
Chlorine Cell Test, 100595	Cl ₂	–	3.00 mg/l*	± 0.30 mg/l	–	–	–	see prep. instr.
Chlorine Cell Test, 100597	Cl ₂	–	3.00 mg/l*	± 0.30 mg/l	–	–	–	see prep. instr.
Chlorine Cell Test (liquid reagent), 100086/100087	Cl ₂	–	3.00 mg/l*	± 0.30 mg/l	–	–	–	see prep. instr.
Chlorine Cell Test (liquid reagent), 100086/100087/100088	Cl ₂	–	3.00 mg/l*	± 0.30 mg/l	–	–	–	see prep. instr.
Chromate Cell Test, 114552	Cr	–	1.00 mg/l*	± 0.10 mg/l	–	–	–	119780
COD Cell Test, 114560	COD	CombiCheck 50, 114695	20.0 mg/l	± 4.0 mg/l	125028	20.0 mg/l	± 0.7 mg/l	see prep. instr.
COD Cell Test, 101796	COD	CombiCheck 50, 114695	20.0 mg/l	± 2.0 mg/l	125028	20.0 mg/l	± 0.7 mg/l	see prep. instr.
COD Cell Test, 114540	COD	CombiCheck 10, 114676	80 mg/l	± 12 mg/l	125029	100 mg/l	± 3 mg/l	see prep. instr.
COD Cell Test, 114895	COD	CombiCheck 60, 114696	250 mg/l	± 20 mg/l	125029	100 mg/l	± 3 mg/l	–
					125030	200 mg/l	± 4 mg/l	see prep. instr.
COD Cell Test, 114690	COD	CombiCheck 60, 114696	250 mg/l	± 25 mg/l	125029	100 mg/l	± 3 mg/l	–
					125030	200 mg/l	± 4 mg/l	–
					125031	400 mg/l	± 5 mg/l	see prep. instr.
COD Cell Test, 114541	COD	CombiCheck 20, 114675	750 mg/l	± 75 mg/l	125029	100 mg/l	± 3 mg/l	–
					125030	200 mg/l	± 4 mg/l	–
					125031	400 mg/l	± 5 mg/l	–
					125032	1000 mg/l	± 11 mg/l	see prep. instr.
COD Cell Test, 114691	COD	CombiCheck 80, 114738	1500 mg/l	± 150 mg/l	125031	400 mg/l	± 5 mg/l	–
					125032	1000 mg/l	± 11 mg/l	–
					125033	2000 mg/l	± 32 mg/l	see prep. instr.
COD Cell Test, 114555	COD	CombiCheck 70, 114689	5000 mg/l	± 400 mg/l	125032	1000 mg/l	± 11 mg/l	–
					125033	2000 mg/l	± 32 mg/l	–
					125034	8000 mg/l	± 68 mg/l	see prep. instr.
COD Cell Test, 101797	COD	–	50000 mg/l*	± 5000 mg/l	125034	8000 mg/l	± 68 mg/l	–
					125035	50000 mg/l	± 894 mg/l	see prep. instr.
COD Cell Test, 109772	COD	–	80 mg/l*	± 12 mg/l	125028	20.0 mg/l	± 0.7 mg/l	–
					125029	100 mg/l	± 3 mg/l	see prep. instr.
COD Cell Test, 109773	COD	–	750 mg/l*	± 75 mg/l	125029	100 mg/l	± 3 mg/l	–
					125030	200 mg/l	± 4 mg/l	–
					125031	400 mg/l	± 5 mg/l	–
					125032	1000 mg/l	± 11 mg/l	see prep. instr.
COD Cell Test, 117058	COD	–	30.0 mg/l*	± 3.0 mg/l	–	–	–	see prep. instr.
COD Cell Test, 117059	COD	–	1500 mg/l*	± 150 mg/l	–	–	–	see prep. instr.
Copper Cell Test, 114553	Cu	CombiCheck 30, 114677	2.00 mg/l	± 0.20 mg/l	–	–	–	119786
Cyanide Cell Test, 102531	CN	–	0.250 mg/l*	± 0.030 mg/l	–	–	–	119533
Cyanide Cell Test, 114561	CN	–	0.250 mg/l*	± 0.030 mg/l	–	–	–	119533
Formaldehyde Cell Test, 114500	HCHO	–	5.00 mg/l*	± 0.50 mg/l	–	–	–	see prep. instr.
Hardness, see Total Hardness Cell Test								
Iron Cell Test, 114549	Fe	CombiCheck 30, 114677	1.00 mg/l	± 0.15 mg/l	–	–	–	119781
Iron Cell Test, 114896	Fe	–	25.0 mg/l*	± 2.5 mg/l	–	–	–	119781
Lead Cell Test, 114833	Pb	CombiCheck 40, 114692	2.00 mg/l	± 0.20 mg/l	–	–	–	119776
Magnesium Cell Test, 100815	Mg	–	40.0 mg/l*	± 4.0 mg/l	–	–	–	see prep. instr.
Manganese Cell Test, 100816	Mn	CombiCheck 30, 114677	1.00 mg/l	± 0.15 mg/l	–	–	–	119789
Nickel Cell Test, 114554	Ni	CombiCheck 40, 114692	2.00 mg/l	± 0.20 mg/l	–	–	–	109989

* Self prepared, recommended concentration

Spectroquant® CombiCheck and Standard Solutions

Test kit, Cat. No. or method	Evalu- ation as	CombiCheck, Cat. No.	Confidence interval		Diluted and ready-to-use standard solutions, CRM			Ready-to-use standard solution, Cat. No.
			Spec. value for the standard	max. working tolerance	Cat. No.	concen- tration	expanded measurement uncertainty	
Nitrate Cell Test, 114542	NO ₃ -N	CombiCheck 20, 114675	9.0 mg/l	± 0.9 mg/l	125037	2.50 mg/l	± 0.06 mg/l	
					125038	15.0 mg/l	± 0.4 mg/l	119811
Nitrate Cell Test, 114563	NO ₃ -N	CombiCheck 20, 114675	9.0 mg/l	± 0.9 mg/l	125037	2.50 mg/l	± 0.06 mg/l	
					125038	15.0 mg/l	± 0.4 mg/l	119811
Nitrate Cell Test, 114764	NO ₃ -N	CombiCheck 80, 114738	25.0 mg/l	± 2.5 mg/l	125037	2.50 mg/l	± 0.06 mg/l	
					125038	15.0 mg/l	± 0.4 mg/l	
					125039	40.0 mg/l	± 1.0 mg/l	119811
Nitrat Cell Test, 100614	NO ₃ -N	–	100 mg/l*	± 10 mg/l	125039	40.0 mg/l	± 1.0 mg/l	
					125040	200 mg/l	± 5 mg/l	119811
Nitrite Cell Test, 114547	NO ₂ -N	–	0.300 mg/l*	± 0.030 mg/l	125041	0.200 mg/l	± 0.009 mg/l	119899
Nitrite Cell Test, 100609	NO ₂ -N	–	45.0 mg/l*	± 5 mg/l	125042	40.0 mg/l	± 1.3 mg/l	119899
Nitrogen (total) Cell Test, 114537 N		CombiCheck 50, 114695	5.0 mg/l	± 0.7 mg/l	125043	2.50 mg/l	± 0.06 mg/l	
					125044	12.0 mg/l	± 0.3 mg/l	see prep. instr.
Nitrogen (total) Cell Test, 100613 N		CombiCheck 50, 114695	5.0 mg/l	± 0.7 mg/l	125043	2.50 mg/l	± 0.06 mg/l	
					125044	12.0 mg/l	± 1.0 mg/l	see prep. instr.
Nitrogen (total) Cell Test, 114763 N		CombiCheck 70, 114689	50 mg/l	± 7 mg/l	125044	12.0 mg/l	± 0.3 mg/l	
					125045	100 mg/l	± 3 mg/l	see prep. instr.
Oxygen Cell Test, 114694	O ₂	–	–	± 0.6 mg/l	–	–	–	see the website
pH Cell Test, 101744	pH	–	7.0	± 0.2	–	–	–	109407
Phosphate Cell Test, 100474	PO ₄ -P	CombiCheck 10, 114676	0.80 mg/l	± 0.08 mg/l	–	–	–	119898
Phosphate Cell Test, 114543	PO ₄ -P	CombiCheck 10, 114676	0.80 mg/l	± 0.08 mg/l	125046	0.400 mg/l P	± 0.016 mg/l	
					125047	4.00 mg/l P	± 0.08 mg/l	119898
Phosphate Cell Test, 100475	PO ₄ -P	CombiCheck 80, 114738	15.0 mg/l	± 1.0 mg/l	–	–	–	
		CombiCheck 20, 114675	8.0 mg/l	± 0.7 mg/l	–	–	–	119898
Phosphate Cell Test, 114729	PO ₄ -P	CombiCheck 80, 114738	15.0 mg/l	± 1.0 mg/l	125047	4.00 mg/l P	± 0.08 mg/l	
		CombiCheck 20, 114675	8.0 mg/l	± 0.7 mg/l	125048	15.0 mg/l P	± 0.4 mg/l	119898
Phosphat Cell Test, 100616	PO ₄ -P	–	50.0 mg/l*	± 5.0 mg/l	–	–	–	119898
Phosphat Cell Test, 100673	PO ₄ -P	–	50.0 mg/l*	± 5.0 mg/l	125047	4.00 mg/l P	± 0.08 mg/l	
					125048	15.0 mg/l P	± 0.4 mg/l	
					125049	75.0 mg/l P	± 1.6 mg/l	119898
Phosphate Cell Test, 114546	PO ₄ -P	–	15.0 mg/l*	± 1.0 mg/l	–	–	–	119898
Potassium Cell Test, 114562	K	–	25.0 mg/l*	± 4.0 mg/l	–	–	–	170230
Potassium Cell Test, 100615	K	–	150 mg/l*	± 15 mg/l	–	–	–	170230
Residual Hardness Cell Test, 114683	Ca	–	2.50 mg/l*	± 0.30 mg/l	–	–	–	119778
Sodium Cell Test, 100885	Na	–	100 mg/l*	± 10 mg/l	–	–	–	see prep. instr.
Sulfate Cell Test, 102532	SO ₄	–	25.0 mg/l*	± 3.0 mg/l	–	–	–	119813
Sulfate Cell Test, 114548	SO ₄	CombiCheck 10, 114676	100 mg/l	± 15 mg/l	125050	40 mg/l	± 6 mg/l	
					125051	125 mg/l	± 6 mg/l	119813
Sulfat Cell Test, 100617	SO ₄	CombiCheck 10, 114676	100 mg/l	± 15 mg/l	125051	125 mg/l	± 6 mg/l	
					125052	400 mg/l	± 20 mg/l	119813
Sulfate Cell Test, 114564	SO ₄	CombiCheck 20, 114675	500 mg/l	± 75 mg/l	125051	125 mg/l	± 6 mg/l	
					125052	400 mg/l	± 20 mg/l	
					125053	800 mg/l	± 27 mg/l	119813
Surfactants (nonionic) Cell Test, 101787	n-Ten	–	4.00 mg/l*	± 0.40 mg/l	–	–	–	see prep. instr.
TOC Cell Test, 114878	TOC	–	40.0 mg/l*	± 3.0 mg/l	–	–	–	109017
TOC Cell Test, 114879	TOC	–	400 mg/l*	± 30 mg/l	–	–	–	109017
Total Hardness Cell Test, 100961	Ca	–	75 mg/l*	± 7 mg/l	–	–	–	see prep. instr.
Volatile Organic Acids Cell Test, 101763	HOAc	–	1500 mg/l*	± 80 mg/l	–	–	–	see prep. instr.
Volatile Organic Acids Cell Test, 101749	C ₃ H ₇ COOH	–	1500 mg/l*	± 80 mg/l	–	–	–	see prep. instr.
Volatile Organic Acids Test, 101809	C ₃ H ₇ COOH	–	1500 mg/l*	± 80 mg/l	–	–	–	see prep. instr.
Zinc Cell Test, 100861	Zn	–	0.500 mg/l*	± 0.050 mg/l	–	–	–	119806
Zinc Cell Test, 114566	Zn	CombiCheck 40, 114692	2.00 mg/l	± 0.40 mg/l	–	–	–	119806

* Self prepared, recommended concentration

Instructions for the Preparation of Standard Solutions

Standard solution of acid capacity

Preparation of a standard solution:

A sodium hydroxide solution of 0.1 mol/l (corresponds to 100 mmol/l) is used.

Further investigational concentrations may be prepared from this standard solution by diluting accordingly with distilled water.

Stability:

When stored in a cool place (refrigerator), the diluted investigational solutions remain stable for one week.

Reagents required:

1.09141.1000	Sodium hydroxide solution 0.1 mol/l Titripur®
1.16754.9010	Water for analysis EMSURE®

Standard solution of calcium

Preparation of a standard solution:

Dissolve 2.946 g of calcium nitrate tetrahydrate with distilled water in a calibrated or conformity-checked 500-ml volumetric flask and make up to the mark with distilled water.

The standard solution prepared according to this procedure has a concentration of 1000 mg/l calcium.

Further investigational concentrations may be prepared from this standard solution by diluting accordingly with distilled water.

Stability:

The standard solution of 1000 mg/l remains stable for one week. The diluted standard solutions (investigational concentrations) remain stable for one day.

Reagents required:

1.02121.0500	Calcium nitrate tetrahydrate for analysis EMSURE®
1.16754.9010	Water for analysis EMSURE®

Instructions for the Preparation of Standard Solutions

Standard solutions of free chlorine

All standard solutions described here for free chlorine yield equivalent results and are identically suited for the determination of chlorine.

Standard solution of free chlorine

Preparation of a standard solution:

Dissolve 1.85 g of dichloroisocyanuric acid sodium salt dihydrate GR with distilled water in a calibrated or conformity-checked 1000-ml volumetric flask and make up to the mark with distilled water.

The standard solution prepared according to this procedure has a concentration of 1000 mg/l free chlorine.

Further investigational concentrations may be prepared from this standard solution by diluting accordingly with distilled water.

Stability:

When stored in a cool place (refrigerator), the standard solution of 1000 mg/l and the diluted standard solutions (investigational concentrations) remain stable for one day.

Note:

This is a standard solution that can be prepared particularly rapidly and easily.

Reagents required:

1.10888.0250	Dichloroisocyanuric acid sodium salt dihydrate GR for analysis
1.16754.9010	Water for analysis EMSURE®

Instructions for the Preparation of Standard Solutions

Standard solution of free chlorine analogous to DIN EN ISO 7393

Preparation of a KIO₃ stock solution:

Dissolve 1.006 g of KIO₃ in 250 ml of distilled water in a calibrated or conformity-checked 1000-ml volumetric flask. Subsequently make up to the mark with distilled water.

Preparation of a KIO₃/KI standard solution:

Transfer 15.00 ml (5.00 ml) of the KIO₃ stock solution to a calibrated or conformity-checked 1000-ml volumetric flask, add approx. 1 g of KI and make up to the mark with distilled water.

1 ml of this solution is equivalent to 0.015 mg (0.005 mg) of free chlorine.

Preparation of the chlorine standard solution:

Pipette 20.0 ml (10.0 ml) (full pipette) KIO₃/KI standard solution into a calibrated or conformity-checked 100-ml volumetric flask, add 2.0 ml of H₂SO₄ 0.5 mol/l, leave to stand for 1 min, and then add NaOH 2 mol/l dropwise (approx. 1 ml) until the solution just loses its color. Subsequently make up the solution to the mark with distilled water.

The concentration of the solution is 3.00 mg/l (0.500 mg/l) free chlorine.

Stability:

The KIO₃ stock solution remains stable for 4 weeks when stored in a cool place (refrigerator). The KIO₃/KI standard solution can be used for 5 hours when stored in a cool place (refrigerator). The diluted chlorine standard solution is not stable and must be used immediately.

Note:

This procedure involves the preparation according to a standardized method.

Reagents required:

1.02404.0100	Potassium iodate, volum. standard
1.05043.0250	Potassium iodide for analysis EMSURE®
1.09072.1000	Sulfuric acid 0.5 mol/l Titripur®
1.09136.1000	Sodium hydroxide solution 2 mol/l Titripur®
1.16754.9010	Water for analysis EMSURE®

Instructions for the Preparation of Standard Solutions

Standard solution of free chlorine

Preparation of a stock solution:

First prepare a 1:10 dilution using a sodium hypochlorite solution containing approx. 13% of active chlorine. For this pipette 10 ml of sodium hypochlorite solution into a calibrated or conformity-checked 100-ml volumetric flask and then make up to the mark with distilled water.

Precise assay of the stock solution:

Pipette 10.0 ml of the stock solution into a 250-ml ground-glass-stoppered conical flask containing 60 ml of distilled water. Subsequently add to this solution 5 ml of hydrochloric acid 25% and 3 g of potassium iodide. Close the conical flask with the ground-glass stopper, mix thoroughly, and leave to stand for 1 min.

Titrate the eliminated iodine with sodium thiosulfate solution 0.1 mol/l until a weakly yellow color emerges. Add 2 ml of zinc iodide-starch solution and titrate from blue to colorless.

Calculation and preparation of a standard solution:

Consumption of sodium thiosulfate solution 0.1 mol/l (ml) x 355 = content of free chlorine, in mg/l

Further investigational concentrations may be prepared from the stock solution prepared according to the procedure described above by diluting accordingly with distilled water.

Stability:

When stored in a cool place (refrigerator), a standard solution remains stable for approx. one week. The diluted standard solutions (investigational concentrations) are stable for approx. 2 hours.

Note:

This is a standard solution that is absolutely necessary for the preparation of the monochloramine standard.

Standard solution of total chlorine

Preparation of a standard solution:

Dissolve 4.00 g of chloramine T GR with distilled water in a calibrated or conformity-checked 1000-ml volumetric flask and make up to the mark with distilled water.

The standard solution prepared according to this procedure has a concentration of 1000 mg/l total chlorine.

Further investigational concentrations may be prepared from this standard solution by diluting accordingly with distilled water.

Stability:

When stored in a cool place (refrigerator), the standard solution of 1000 mg/l and the diluted standard solutions (investigational concentrations) remain stable for one day.

Reagents required:

1.00316.1000	Hydrochloric acid 25 % for analysis EMSURE®
1.05614.9025	Sodium hypochlorite solution techn. approx. 13% active chlorine
1.09147.1000	Sodium thiosulfate solution 0.1 mol/l Titripur®
1.05043.0250	Potassium iodide GR for analysis
1.05445.0500	Zinc iodide-starch solution GR for analysis
1.16754.9010	Water for analysis EMSURE®

Reagents required:

1.02426.0250	Chloramine T trihydrate GR for analysis
1.16754.9010	Water for analysis EMSURE®

Instructions for the Preparation of Standard Solutions

Standard solution of COD

Preparation of a standard solution:

Dissolve 0.850 g of potassium hydrogen phthalate GR with distilled water in a calibrated or conformity-checked 1000-ml volumetric flask and make up to the mark with distilled water.

The standard solution prepared according to this procedure has a concentration of 1000 mg/l COD.

Further investigational concentrations may be prepared from this stock solution by diluting accordingly with distilled water.

Stability:

When stored in a cool place (refrigerator), the standard solution remains stable for one month. When stored under appropriate cool conditions (refrigerator), the diluted standard solutions (investigational concentrations) remain stable – depending on the respective concentration – for approx. one week to one month.

Reagents required:

1.02400.0080	Potassium hydrogen phthalate GR for analysis, volum. standard
1.16754.9010	Water for analysis EMSURE®

Standard solution of COD/chloride

Preparation of a chloride dilution solution:

Dissolve 32.9 g of sodium chloride GR with distilled water in a calibrated or conformity-checked 1000-ml volumetric flask and make up to the mark with distilled water.

The dilution solution prepared according to this procedure has a concentration of 20 g/l Cl⁻.

Preparation of a COD/Cl⁻ standard solution:

Dissolve 0.850 g of potassium hydrogen phthalate GR with **dilution solution** in a calibrated or conformity-checked 100-ml volumetric flask and make up to the mark with **dilution solution**.

The standard solution prepared according to this procedure has a concentration of 1000 mg/l and 20 g/l Cl⁻.

Further investigational concentrations may be prepared from this stock solution by diluting accordingly with **dilution solution**.

Stability:

When stored in a cool place (refrigerator), the dilution solution of 20 g/l Cl⁻ and the standard solution of 10000 mg/l COD / 20 g/l Cl⁻ remain stable for one month. When stored under appropriate cool conditions (refrigerator), the diluted standard solutions (investigational concentrations) remain stable - depending on the respective concentration - for approximately one week to one month.

Reagents required:

1.02400.0080	Potassium hydrogen phthalate GR for analysis, volum. standard
1.06404.0500	Sodium chloride for analysis EMSURE®
1.16754.9010	Water for analysis EMSURE®

Instructions for the Preparation of Standard Solutions

Standard solution of formaldehyde

Preparation of a stock solution:

In a calibrated or conformity-checked 1000-ml volumetric flask make up 2.50 ml of formaldehyde solution min. 37% GR to the mark with distilled water.

The stock solution prepared according to this procedure has a concentration of approx. 1000 mg/l formaldehyde.

Precise assay of the stock solution:

Pipette 40.0 ml (full pipette) of the formaldehyde stock solution into a 300-ml ground-glass conical flask and add 50.0 ml (buret) of iodine solution 0.05 mol/l and 20 ml of sodium hydroxide solution 1 mol/l.

Leave to stand for 15 minutes and subsequently add 8 ml of sulfuric acid 25%. Subsequently titrate with sodium thiosulfate solution 0.1 mol/l until the yellow iodine color has disappeared, add 1 ml of zinc iodide-starch solution, and continue to titrate until a milky, pure white color emerge.

Calculation and preparation of a standard solution:

$C1 =$ consumption of sodium thiosulfate solution 0.1 mol/l (ml)

$C2 =$ quantity of iodine solution 0.05 mol/l (50,0 ml)

$$\text{mg/l formaldehyde} = (C2 - C1) \times 37.525$$

Further investigational concentrations may be prepared from the stock solution exactly determined according to the procedure described above by diluting accordingly with distilled water.

Stability:

When stored in a cool place (refrigerator), the stock solution of approx. 1000 mg/l remains stable for one week. After this time, the stock solution must be determined anew. The diluted standard solutions (investigational concentrations) must be used immediately.

Reagents required:

1.04003.1000	Formaldehyde solution min. 37% GR for analysis
1.09099.1000	Iodine solution 0.05 mol/l Titripur®
1.09147.1000	Sodium thio-sulfate solution 0.1 mol/l Titripur®
1.09137.1000	Sodium hydroxide solution 1 mol/l Titripur®
1.00716.1000	Sulfuric acid 25% for analysis EMSURE®
1.05445.0500	Zinc iodide-starch solution GR for analysis
1.16754.9010	Water for analysis EMSURE®

Instructions for the Preparation of Standard Solutions

Standard solution of magnesium

Preparation of a standard solution:

Dissolve 1.055 g of magnesium nitrate hexahydrate with distilled water in a calibrated or conformity-checked 100-ml volumetric flask and make up to the mark with distilled water.

The standard solution prepared according to this procedure has a concentration of 1000 mg/l magnesium.

Further investigational concentrations may be prepared from this standard solution by diluting accordingly with distilled water.

Stability:

The standard solution of 1000 mg/l remains stable for one week. The diluted standard solutions (investigational concentrations) remain stable for one day.

Standard solution of nitrogen (total)

Preparation of a standard solution:

Dissolve 5.36 g of glycine GR with distilled water in a calibrated or conformity-checked 1000-ml volumetric flask and make up to the mark with distilled water.

The standard solution prepared according to this procedure has a concentration of 1000 mg/l total nitrogen.

Further investigational concentrations may be prepared from this standard solution by diluting accordingly with distilled water.

Stability:

When stored in a cool place (refrigerator), the standard solution of 1000 mg/l remains stable for one week. The diluted standard solutions (investigational concentrations) must be used immediately.

Standard solution of sodium

Preparation of a standard solution:

A chloride standard solution of 1000 mg/l is used.
1000 mg/l chloride corresponds to 649 mg/l sodium.

Further investigational concentrations may be prepared by diluting accordingly with distilled water.

Stability:

When stored in a cool place (refrigerator), the diluted standard solutions (investigational concentrations) remain stable for one month.

Reagents required:

1.05853.0500	Magnesium nitrate hexahydrate for analysis EMSURE®
1.16754.9010	Water for analysis EMSURE®

Reagents required:

1.04201.0100	Glycine GR for analysis
1.16754.9010	Water for analysis EMSURE®

Reagents required:

1.19897.0500	Chloride standard solution Certipur®
1.16754.9010	Water for analysis EMSURE®

Instructions for the Preparation of Standard Solutions

Standard solution of surfactants (nonionic)

Preparation of a standard solution:

Dissolve 1.00 g of Triton® X-100 with distilled water in a calibrated or conformity-checked 1000-ml volumetric flask and make up to the mark with distilled water.

The standard solution prepared according to this procedure has a concentration of 1000 mg/l non-ionic surfactants.

Further investigational concentrations may be prepared from this standard solution by diluting accordingly with distilled water.

Stability:

When stored in a cool place (refrigerator), the standard solution of 1000 mg/l remains stable for one week. The diluted standard solutions (investigational concentrations) must be used immediately.

Reagents required:

1.12298.0101	Triton® X-100
1.16754.9010	Water for analysis EMSURE®

Standard solution of total hardness

Preparation of a standard solution:

Dissolve 2.946 g of calcium nitrate tetrahydrate with distilled water in a calibrated or conformity-checked 500-ml volumetric flask and make up to the mark with distilled water.

The standard solution prepared according to this procedure has a concentration of 1000 mg/l calcium (corresponds to 175 °e).

Further investigational concentrations may be prepared from this standard solution by diluting accordingly with distilled water.

Stability:

The standard solution of 1000 mg/l remains stable for one week. The diluted standard solutions (investigational concentrations) remain stable for one day.

Reagents required:

1.02121.0500	Calcium nitrate tetrahydrate for analysis EMSURE®
1.16754.9010	Water GR for analysis

Standard solution of volatile organic acids

Preparation of a standard solution:

Dissolve 2,05 g of sodium acetate anhydrous with distilled water in a calibrated or conformity-checked 1000-ml volumetric flask and make up to the mark with distilled water.

The standard solution prepared according to this procedure has a concentration of 1500 mg/l acetic acid.

Stability:

When stored in a cool place (refrigerator), the standard solution remains stable for one week.

Reagents required:

1.06268.0250	Sodium acetate anhydrous for analysis EMSURE®
1.16754.9010	Water GR for analysis

Merck KGaA, 64271 Darmstadt, Germany, Tel. +49(0)6151 72-2440
www.analytical-test-kits.com

EMD Millipore Corporation, 290 Concord Road, Billerica,
MA 01821, USA, Tel. +1-978-715-4321
