B-150D Series

INSTRUCTION MANUAL

| Model |
|-------------|
| B-150D-MRPL |
| B-150D-BRPL |

Ver. 2.0 2019

Summary

| | |
|--|-----------|
| 1. Warning | 3 |
| 2. Symbols and conventions | 3 |
| 3. Safety Information | 3 |
| 4. Intended use | 3 |
| 5. Overview | 4 |
| 5.1 B-150D-MRPL | 4 |
| 5.2 B-150D-BRPL | 5 |
| 6. Unpacking | 6 |
| 7. Assembling | 6 |
| 7.1 B-150D-MRPL | 6 |
| 7.2 B-150D-BRPL | 7 |
| 7.3 Assembling the microscope | 8 |
| 7.4 Polarizing set (optional) | 9 |
| 8. Use of the microscope | 10 |
| 8.1 Light intensity adjustment | 10 |
| 8.2 Coarse focus tension adjustment | 10 |
| 8.3 Stage | 10 |
| 8.4 Adjust the interpupillary distance | 10 |
| 8.5 Diopter adjustment | 11 |
| 8.6 Use of oil immersion objective | 11 |
| 8.7 Aperture diaphragm | 12 |
| 8.8 Use with rechargeable batteries | 12 |
| 8.9 Use of the polarizer (optional) | 12 |
| 9. Maintenance | 13 |
| 10. Troubleshooting | 14 |
| Equipment disposal | 15 |

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

2. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

3. Safety Information

Avoiding Electrical Shock


Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

4. Intended use


For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

5. Overview

5.1 B-150D-MRPL

5.2 B-150D-BRPL

6. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

 Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

7. Assembling

Once opened the box, the microscope parts are the following:

7.1 B-150D-MRPL

- | | |
|-------------------------------|------------------|
| ① Frame | ⑥ Dust cover |
| ② Monocular observation head | ⑦ Green filter |
| ③ Eyepiece | ⑧ Frosted filter |
| ④ Tension adjustment tool | ⑨ Power supply |
| ⑤ Objectives (4X / 10X / 40X) | ⑩ USB cable |

7.2 B-150D-BRPL

- ① Frame
- ② Binocular observation head
- ③ Eyepieces
- ④ Tension adjustment tool
- ⑤ Objectives (4X / 10X / 40X / 100X)
- ⑥ Dust cover
- ⑦ Green filter
- ⑧ Frosted filter
- ⑨ Power supply
- ⑩ USB cable
- ⑪ Immersion oil

7.3 Assembling the microscope

1. Remove the dust cap from the microscope frame and from the bottom of the observation head.
 2. Insert the optical head above the stand and tighten the screw ①. (Fig. 1)
- **Hold the head with one hand during the locking in order to avoid that the head falls.**

3. Insert both eyepieces into the tubes of the optical head. (Fig. 2)
- **For model B-150D-MRPL: only one eyepiece is used as this model is monocular.**

4. Insert the power supply jack in the socket placed at the rear side of the microscope. (Fig. 3)


7.4 Polarizing set (optional)


1. Place the polarizer on the light exit ① at the base of the microscope. (Fig. 4)

2. Loosen the head fixing knob ② and remove the head from the microscope frame. (Fig. 5)

3. Insert the analyzer into the hole inside the frame ③. (Fig. 6)
4. Put back the head into its original position and lock the fixing knob.

8. Use of the microscope

8.1 Light intensity adjustment

Operate on the light intensity adjustment dial to turn ON / OFF the microscope and to increase / decrease the illumination voltage ①. (Fig. 7)

8.2 Coarse focus tension adjustment

- Adjust the tension using the provided tool.

The coarse knob tension is pre-set in the factory. To modify the tension according to personal's needs, rotate the ring ② using the provided tool (Fig. 8). Clockwise rotation increases the tension. If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

8.3 Stage

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm with coverslide 0,17mm. (Fig. 9)

1. Open the spring arm of the slide holder ③ and place the slide from the front on the stage.
 2. Gently release the spring arm of the slide holder.
- **A sudden release of the spring arm could cause the falling of the slide.**

8.4 Adjust the interpupillary distance

Only B-150D-BRPL

Hold the right and left parts of the observation head using both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. (Fig. 10)

- The graduation on the interpupillary distance indicator ④, pointed by the spot "." on the eyepiece holder, shows the distance between the operator's eyes.

The range of the interpupillary distance is 48- 75 mm.

8.5 Diopter adjustment

Only B-150D-BRPL

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
 2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the dioptic adjustment ring ① to compensate. (Fig. 11)
- **The adjustment range is ± 5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator's dioptic correction.**

Fig. 11

8.6 Use of oil immersion objective

Only B-150D-BRPL

1. Focus the specimen with a low power objective.
 2. Lower the stage.
 3. Put a drop of oil (provided) on the area of the specimen to be observed. (Fig. 12)
- **Make sure that there are no oil bubbles. Air bubbles in the oil damage the image quality.**
 - To check for bubbles: remove an eyepiece, fully open the aperture diaphragm and observe the objective exit pupil. (The pupil must be circular and bright).
 - To remove the bubbles, gently move the nose-piece to the right and left to move the immersion objective a few times and allow the air bubbles to move.
4. Insert immersion objective.
 5. Return the stage to the upper focusing point and obtain an optimal focus using the fine focus knob.
 6. After use, gently remove the oil with a soft paper towel or a lightly moistened optic paper with a mixture of ethyl ether (70%) and absolute ethyl alcohol (30%).
- **The immersion oil, if not immediately cleaned, could crystallize creating a glass-like layer. In this situation the observation of the specimen would be difficult (even not impossible) due to the presence of an additional thickness on the objective.**

Fig. 12

8.7 Aperture diaphragm

The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image. Move the diaphragm lever ① (Fig. 13) toward left or right to decrease or increase the N.A. value.

- With low contrast specimens set the numerical aperture to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's diaphragm in order to obtain an image like the one in Fig. 14


Fig. 13


Fig. 14

8.8 Use with rechargeable batteries

When the microscope is plugged with the power supply, the LED indicator for the battery recharge ② is lit. (Fig. 15)

- LED red: battery under charge
- LED green: battery fully charged.

When the microscope is unplugged, the LED is off.

During the normal use with batteries, LED is always OFF.

Fig. 15

8.9 Use of the polarizer (optional)

1. Remove the specimen from the stage.
2. Looking inside the eyepieces, rotate the polarizer until the darkest position is achieved.
3. Once the dark is achieved ("extinction" or "Crossed Nicol" position) it is possible to begin the observation.

9. Maintenance

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

10. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

| PROBLEM | CAUSE | SOLUTION |
|--|---|---|
| I. Optical Section: | | |
| LED operates, but field of view remains dark. | Power supply is unplugged. | Connect |
| | Brightness is too low | Set brightness to a proper level |
| | Batteries are uncharged | Fully charge the batteries |
| Dirt or dust is visible in the field of view. | Dirt/dust on the specimen | Clean the specimen |
| | Dirt/dust on the eyepieces | Clean the eyepieces |
| Image looks double | Aperture diaphragm is stopped down too far | Open aperture diaphragm |
| Visibility is poor. <ul style="list-style-type: none"> • Image is not good. • Contrast is poor. • Details are indistinct. • Image glares | Revolving nosepiece is in an incorrect position | Move the nosepiece to a click stop |
| | Aperture diaphragm is too closed or too open | Adjust aperture diaphragm |
| | Dust or dirt on lenses (condenser, objectives, eyepieces and slide) | Clean thoroughly |
| One side of the image is out of focus. | The nosepiece is not in the center of the light path | Turn the nosepiece to a click stop |
| | The specimen is out of place (tilted) | Place the specimen flat on the stage. |
| | The optical performance of the sample cover glass is poor | Use a cover glass of better quality |
| II. Mechanical Section: | | |
| The coarse focus knob is hard to turn. | The tension adjustment collar is too tight | Loosen the tension adjustment collar |
| The focus is unstable. | The tension adjustment collar is too loose | Tighten the tension adjustment collar |
| III. Electric section | | |
| The LED doesn't turn on. | No power supply | Check the power cord connection |
| | Batteries are uncharged | Fully charge the batteries |
| The brightness is not enough | The brightness adjustment is low | Adjust the brightness |
| The light blinks | The power cord is poorly connected | Check the power cord |
| IV. Observation tube | | |
| Field of view of one eye does not match that of the other. | Interpupillary distance is incorrect. | Adjust interpupillary distance. |
| | Incorrect diopter adjustment. | Adjust diopter. |
| | Your view is not accustomed to microscope observation. | Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope. |

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. “According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal.”

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.