EnzChek[®] Phospholipase A₂ Assay Kit

Catalog nos. E10217, E10218

Table 1. Contents and storage information.

| Material | E10217 (2-plates) | E10218 (10-plates) | Concentration | Storage | Stability |
|---|----------------------|-----------------------|---|---|---|
| Phospholipase A ₂ substrate (Red/Green BODIPY [®] PC-A2) (Component A) | 2 vials | 10 vials | | ≤-20°C Desiccate Protect from light | |
| Dioleoylphosphatidylcholine (DOPC) (Component B); MW = 785.59 | 800 µg | 4 mg | | | |
| Dioleoylphosphatidylglycerol (DOPG) (Component C); MW = 797.04 | 800 µg | 4 mg | Not applicable | ≤-20°CDesiccate | When stored as directed the product |
| Phospholipase A ₂ from honey bee venom (Component D) | 40 Units | 200 Units | | | is stable for at least 1 year. |
| Dimethylsulfoxide (DMSO, Component E) | 200 µL | 1 mL | | Desiccate | |
| 5X Phospholipase A ₂ reaction buffer (Component F) | 10 mL | 50 mL | 250 mM Tris-HCl, 500 mM NaCl, 5 mM CaCl ₂ , pH 8.9 | ≤-20°C | |

Number of assays: For Cat. no. E10217, sufficient material is supplied for 200 reactions in 96-well microplates at a volume of 100 µL per well as described in the protocol below, or 800 reactions using low-volume 384-well microplates at a volume of \leq 25 µL per well. For Cat. no. E10218, sufficient material is supplied for 1,000 reactions in 96-well microplates at a volume of 100 µL per well as described in the protocol below or 4,000 reactions using low-volume 384-well microplates at a volume of \leq 25 µL per well as described in the protocol below or 4,000 reactions using low-volume 384-well microplates at a volume of \leq 25 µL per well.

Approximate fluorescence excitation/emission maxima: Excitation = 505 nm (typical plate reader setting $Ex \approx 460$ nm); Emission = 515 nm (intensity based); Emission = 515/575 nm (ratiometric).

Introduction

The EnzChek[®] Phospholipase A_2 Assay Kit provides a simple fluorometric method designed for continuous monitoring of phospholipase A_2 (PLA₂). PLA₂ represents a family of enzymes that hydrolyze the *sn*-2 ester linkage of phospholipids and play important roles in cardiovascular, inflammatory and nervous system disorders, and cancers. The EnzChek[®] Phospholipase A_2 substrate is selective for PLA₂, and provides sensitive and continuous rapid real-time monitoring of PLA₂ enzyme activities for rapid and direct analysis of PLA₂ using automated instruments. (Figure 1). The enzyme assay can be used in a intensity or ratiometric based detection mode providing a simple method with low background, high sensitivity, and high specificity for PLA₂. In intensity-based detection mode, the PLA₂ activity is monitored by the intensity increase of a single wavelength at

approximately 515 nm. In ratiometric analysis, which is based on the distinct fluorescence resonance energy transfer (FRET) emission of the substrate prior to and after cleavage, PLA_2 is detected by changes in the emission intensity ratio at 515/575 nm with excitation at ~460 nm.

The EnzChek[®] Phospholipase A_2 Assay Kit can detect bee venom PLA₂ at 0.05 Units/mL or lower (see Figure 2). The kit may be amenable for use with cells or cell lysates and for characterizing PLA₂ inhibition (Figure 3).




Figure 1. Fluorescence emission spectra of EnzChek[®] Phospholipase A_2 substrate incorporated in liposomes with addition of bee venom PLA₂ at ambient temperature.






Figure 2. Plot of fluorescence emission intensities versus concentration of PLA₂ per well at 10 minutes, run at ambient temperature with liposomes. Fluorescence was measured exciting at 460 nm on a Spectra Max M5 (Molecular Devices). **Panel A:** Ratiometric (FRET)-based (515/575 nm) EnzChek® Phospholipase A₂ substrate. **Panel B:** Intensity-based detection mode (515 channel)- EnzChek® Phospholipase A₂ substrate. Background fluorescence determined for the no-PLA₂ enzyme control reaction has been subtracted for each value.





Before You Begin

| Provided | Samples |
|--------------------|--|
| | Deionized water |
| | • Ethanol |
| | Plastic vials for reagent preparation |
| | Microplates, 96-well or 384-well |
| | • Magnetic stirrer, stir bar, and pipettor (an air displacement pipettor with 100 μ L capacit fitted with a narrow orifice gel-loading tip is suitable and required for the liposomes preparation at step 2.6) |
| General Guidelines | - The kit is useful for detecting PLA_2 activity in samples or for screening PLA_2 inhibitors. |
| | • The assay protocol below is optimized for use with 96-well microplates using a 100 μ L reaction volume per assay. For 384-well plates, adjust the reaction volumes accordingly 25 μ L per assay (recommended). |
| | • The assay protocol is designed for use with a fluorescence microplate reader. A SpectraMax M5 (Molecular Devices) was used throughout the development of this kit. |
| | spectrumant mo (morecular Demeco) was used infoughout the development of this kit. |
| | Allow the kit components to equilibrate to room temperature before use. |

1.1 1 mM EnzChek[®] Phospholipase A₂ Substrate: Allow one vial of PLA₂ substrate (Component A) and DMSO (Component E) to warm to room temperature. Dissolve the contents of one vial of PLA₂ substrate in 40 μL DMSO. One vial of PLA₂ substrate is sufficient for approximately 100 assays with a final reaction volume of 100 μL per assay.

Store this stock solution frozen at $\leq -20^{\circ}$ C, protected from light.

1.2 1X EnzChek° **PLA**₂ **Reaction Buffer:** Add 4 mL of 5X PLA₂ reaction buffer (Component F) to 16 mL deionized water. This 20 mL volume of 1X reaction buffer is sufficient for approximately 100 assays of 100 μ L each with 10 mL excess for making stock solutions. Store remaining solution at 4°C.

- **1.3 10 mM DOPC (Dioleoylphosphatidylcholine):** Dissolve the contents of DOPC vial (Component B) in 100 μL (Cat. no. E10217) or 500 μL (Cat no. E10218) ethanol. Store solution at ≤-20°C.
- **1.4 10 mM DOPG (Dioleoylphosphatidylglycerol):** Dissolve the contents of DOPG vial (Component C) in 100 μL (Cat. no. E10217) or 500 μL (Cat no. E10218) ethanol. Store solution at ≤-20°C.
- **1.5 500** Units/mL PLA₂ Stock Solution: Dissolve contents of the PLA₂ vial (Component D) in 80 μ L (Cat. no. E10217, 2 plates) or 400 μ L (Cat. no. E10218, 10 plates) of 1X PLA₂ reaction buffer prepared in step 1.2. Sufficient enzyme is supplied to prepare 80 (Cat. no. E10217) or 800 (Cat. no. E10218) positive control samples at 5 U/mL in an assay volume of 100 μ L. You may use other sources of PLA₂ as a positive control, but the sensitivity and dynamic range of the assay may be affected. Store stock solution at 4°C.

Experimental Protocol

Perform the standard assay protocol below using a total volume of 100 μ L per well. Mix the samples and controls with the substrate-liposome mix at a ratio of 1:1 (50 μ L sample/control + 50 μ L substrate-liposome mix), such that the concentration of each component is two-fold lower in the final reaction volume. You may also use other volumes may be used, provided that you maintain the ratio of samples/controls to substrate at 1:1.

The protocol was developed and tested using group III \mbox{PLA}_2 from bee venom and Group IB \mbox{PLA}_2 from bovine pancreas.

Assay Protocol

- **2.1** Prepare a PLA₂ standard curve by diluting the appropriate amount of 500 Units/mL PLA₂ stock solution to 10 Units/mL in 1X PLA₂ reaction buffer to produce PLA₂ concentrations of 0–10 Units/mL, each in a volume of 50 μ L. Final PLA₂ concentration is two-fold lower (0–5 Units/mL).
- **2.2** If no standard curve is to be used, prepare positive and negative controls. For a positive control, dilute the 500 Units/mL PLA₂ stock solution to 10 Units/mL in 1X PLA₂ reaction buffer. For a negative control, use 1X PLA₂ reaction buffer without PLA₂.
- **2.3** Dilute the PLA_2 -containing samples in PLA_2 reaction buffer. You will need 50 µL sample for each reaction. A variable dilution may be required depending on the total PLA_2 present in each sample.
- **2.4** Pipet 50 μL of the standard curve samples (step 2.1) or controls (step 2.2), and experimental samples (step 2.3) into individual wells of a microplate.
- **2.5** Prepare the Lipid Mix by mixing together 30 μL 10 mM DOPC (from step 1.3), 30 μL 10 mM DOPG (from step 1.4), and 30 μL 1 mM PLA₂ substrate (from step 1.1).
- **2.6** Add 5 mL 1X PLA₂ reaction buffer to a 20 mL beaker containing a small magnetic stir bar and place the beaker on a magnetic stirrer to form a vortex. To prepare 5 mL substrate-liposome for 100 assays, slowly and steadily (over about 1 minute) inject 50 μ L of Lipid Mix (from step 2.5) into the side of the vortex using a pipettor fitted with a narrow orifice gelloading tip.
- 2.7 Add 50 μ L of the substrate-liposome mix (from step 2.6) to each microplate well containing standards, controls, and samples to start the reaction.

| 2.8 | ³ Incubate at room temperature for 10 minutes, protected from light. Because the assay continuous (not terminated), you may measure fluorescence at multiple time points to the kinetics of the reactions. | |
|---------------------|--|--|
| 2.9 | Measure the fluorescence using a microplate reader equipped for excitation in the range of 450–490 nm and fluorescence emission at \sim 515 nm and \sim 575 nm. | |
| | For ratiometric (FRET) detection, use a ratio of the emission 515 nm/575 nm. For single wavelength intensity detection, use emission at \sim 515 nm (see Figure 2). | |
| 2.10 | For each point, subtract the value derived from the no-PLA $_2$ control to correct for background fluorescence. | |
| Inhibitor Screening | For screening PLA_2 inhibitors, select the concentration of enzyme which yields a signal increase at least 20% of maximal. | |
| | For IC ₅₀ determinations, prepare serial dilutions of the inhibitor(s) of interest and include the selected amount of PLA_2 in the substrate-liposome mix at 2X the final concentration. | |

Troubleshooting

| Problem | Cause | Solution |
|--|--|---|
| No response from the control enzyme | Low substrate concentration or substrate is contaminated | Substrate stock solution in DMSO appears purple. If no color is visible, the substrate is too dilute. If color is visible, but there is no response, the substrate is contaminated. Repeat the experiment with fresh substrate. |
| No response from samples | PLA_2 absent, inactivated, or is present in low quantities | Increase incubation time or enzyme amount. If no signal, repeat the experiment with a fresh vial of substrate. |
| Response not in the linear range | PLA_2 in the sample is highly active | Dilute sample until the response falls within linear range of the standard curve. |
| DOPG is not in solution | DOPG precipitates from ethanolic solution when stored at –20°C | Redissolve DOPG by warming the solution to room temperature. |

Reference

1. Biochem Pharmacol. 36, 733 (1987).

Product List Current prices may be obtained from our website or from our Customer Service Department.

| Cat. no. | Product Name | Unit Size |
|----------|---|-----------|
| E10217 | EnzChek [®] Phospholipase A ₂ Assay Kit *2 Plates* | 1 kit |
| E10218 | EnzChek [®] Phospholipase A ₂ Assay Kit *10 Plates* | 1 kit |
| A10072 | Red/Green BODIPY® PC-A2(1-O-(6-BODIPY®558/568-aminohexyl)-2-BODIPY®FL C5-Sn-glycero-3-phosphocholine) | 100 µg |

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road Eugene, OR 97402 Phone: (541) 465-8300 Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time) Phone: (541) 335-0338 Fax: (541) 335-0305 probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209 Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time) Phone: (541) 335-0353 Toll-Free (800) 438-2209 Fax: (541) 335-0238 probestech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd. 3 Fountain Drive Inchinnan Business Park Paisley PA4 9RF, UK Phone: +44 (0) 141 814 6100 Fax: +44 (0) 141 814 6200 Email: euroinfo@invitrogen.com Technical Services: eurotech@invitrogen.com

For country-specific contact information, visit www.invitrogen.com.

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation [®] are registered with the U.S. Patent and Trademark Office.

Copyright 2009, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.