**Instruction Manual**  
**Manuale di istruzioni**  
**Manuel d'instructions**  
**Manual de instrucciones**  
**Bedienungsanleitung**

**RC Heating Plate**

F20700431, F20710431

**RC2 Heating Plate, 2-Position**

F20700430, F20710430

**General Information / Informazioni Generali / Informations Générales / Información General /  
Allgemeine Hinweise**

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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**This unit must be used for laboratory applications only.**

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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**This unit has been designed and manufactured in compliance with the following standards:**

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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**IEC/EN 61010-1  
IEC/EN 61010-2-051**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Do not use with explosive and dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie and to stir combustible liquids that have a low combustion temperature. The minimum fire point of flammable solution is 750 °C. Only small amounts (< 50 ml) of flammable liquid can be used with the device.

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RC is fitted with two fuses (2xT5 A L 250 V (for 230V), 2xF8 A 250 V (for 115 V)), found in the socket on the back. To replace one or more disconnect the mains cable and, using a screwdriver, lift up the small cover on the fuse box.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it.

It is also responsibility of the user to use safety substances for cleaning or decontaminating, which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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## 1. Introduction

The heating plate type RC is designed for general use, for all those applications that require the heating of samples in specific containers positioned on the heating plate of the instrument.

The temperature of the heating plate is controlled by an analogical bulb thermoregulator fitted directly into the heating plate. This makes it possible to obtain a constant control of the heating plate temperature from 50 to 370 °C.

**NOTE:** The container of the product used must be compatible with the temperature used. Normally, in these cases, Pyrex glass containers are used. The use of the heating plate with high temperatures may cause variations of the superficial colour but this doesn't change the characteristics of thermal, mechanical and chemical resistance.

## 2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- RC Heating Plate
- Power supply cord
- Instruction manual

### 2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

### 2.2 Start-up

Rotate the temperature knob completely to the left. Then, set the temperature by turning the dedicated knob.

## 3. Operating controls

### TEMPERATURE KNOB

To select the temperature of the heating plate turn the corresponding knob (Heating). The temperature can be programmed between 50 and 370 °C.

### ON-OFF SWITCH

The on-off switch turns the unit on and off. If the switch is in the "OFF" position the unit is off; if the switch is in the "ON" position the unit is on.

Always turn the unit off after use.

### 3.1 Error message

The unit is fitted with safety devices which cut-off the power supply to the heating plate in the case of malfunctions.

#### Display

- Flashing Heating LED (once every second)
- Flashing Heating LED (2 times every second)
- Flashing Heating LED (once every 3 seconds)

#### Cause

- Overheated heating plate or thermocouple open
- Fault in the temperature reading circuit
- Thermocouple wires reversed or white wire of potentiometer board not connected

Should any of the above occur, please contact your nearest VELP Scientifica service centre.

## 4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual.

In compliance with the product guarantee law, replace the fuses and repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors and using Velp spare part only.

The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

### 4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

## 5. Technical data

Velp Code	<b>F20700431</b>	<b>F20710431</b>	<b>F20700430</b>	<b>F20710430</b>
Power supply	230V/50-60Hz +/-10%	115V/50-60Hz +/-10%	230V/50-60Hz +/-10%	115V/50-60Hz +/-10%
Dimensions WxHxD	165x115x280 mm	165x115x280 mm	340x107x246 mm	340x107x246 mm
Weight	1.9 Kg	1.9 Kg	3.1 Kg	3.1 Kg
Overall power	600W	600W	600W x 2 (1200W)	600W x 2 (1200W)
Construction material	Painted aluminum	Painted aluminum	Painted aluminum	Painted aluminum
Diameter of the heating plate	135mm	135mm	135mm x 2	135mm x 2
Programmable temp. range	50 – 370 °C	50 – 370 °C	50 – 370 °C	50 – 370 °C
Type of temperature control	Analog	Analog	Analog	Analog
Overtemperature protection	Yes	Yes	Yes	Yes
Room temperature range	+5...+40 °C	+5...+40 °C	+5...+40 °C	+5...+40 °C
Storage temperature range	-10...+60 °C	-10...+60 °C	-10...+60 °C	-10...+60 °C
Max humidity	80%	80%	80%	80%
Level of electrical protection	IP 42	IP 42	IP 42	IP 42
Pollution degree	2	2	2	2
Overvoltage category	II	II	II	II
Noise level	<< 80 dBa	<< 80 dBa	<< 80 dBa	<< 80 dBa
Max altitude	4000 m	4000 m	4000 m	4000 m

## 6. Accessories / Spare parts

A00001069	Support rod
A00000351	Handle for AluBlock removal
A00000382	Extension for support rod ARE – AREX
A00000375	Hemispheric bowl 5L flasks, plate 135
A00000374	Hemispheric bowl 3L flasks, plate 135
A00000373	Hemispheric bowl 25ml flasks, plate 135
A00000346	PTFE Safety cover for bowl 1000 ml
A00000345	PTFE Safety cover for bowl 500 ml
A00000344	PTFE Safety cover for bowl 250 ml
A00000343	PTFE Safety cover for bowl 100 ml
A00000342	PTFE Safety cover for bowl 50 ml
A00000341	MonoAluBlock, 40 pos. Ø12 x h 14 mm
A00000340	MonoAluBlock, 17 pos. Ø28 x h 43 mm
A00000339	MonoAluBlock, 17 pos. Ø28 x h 30 mm
A00000338	MonoAluBlock, 17 pos. Ø28 x h 24 mm
A00000337	MultiAluBlockTM, 11 pos. Ø12 x h 14 mm
A00000334	Hemispheric bowl 1000ml flasks, plate 135
A00000333	Hemispheric bowl 500ml flasks, plate 135
A00000332	Hemispheric bowl 250ml flasks, plate 135
A00000331	Hemispheric bowl 100ml flasks, plate 135
A00000330	Hemispheric bowl 50ml flasks, plate 135
A00000327	MultiAluBlockTM, 4 pos. Ø21 x h31 mm
A00000325	MultiAluBlockTM, 4 pos. Ø28 x h 30 mm
A00000324	MultiAluBlockTM, 4 pos. Ø28 x h 43 mm
10000239	Foot 13Dx5H embedded (for code <b>F20700431 – F20710431</b> )
10000232	Foot 14Dx8H for screw (for code <b>F20700430 – F20710430</b> )
10005213	Knob 24D blue

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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and satisfies the essential requirements of the following directives:

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- Machinery directive 2006/42/EC / Macchine 2006/42/CE / Machines 2006/42/CE / Máquinas 2006/42/CE / Maschinen 2006/42/EG

- Low voltage directive 2014/35/EU / Bassa tensione 2014/35/EU / Basse tension 2014/35/EU / Baja tensión 2014/35/EU / Niederspannung 2014/35/EU

- Electromagnetic compatibility directive 2014/30/EU / Compatibilità elettromagnetica 2014/30/EU / Compatibilité électromagnétique 2014/30/EU / Compatibilidad electromagnética 2014/30/EU / Elektromagnetische Verträglichkeit 2014/30/EU

- plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung.

## Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortex mixers  
Homogenizers  
COD Thermoreactors  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Overhead Shakers  
Turbidimeter  
Radiation Detector  
Open Circulating Baths  
Pumps

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