VWR RC-10 Basic Chiller

VWR Cat. No. 462-0137

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Instructions

Mode d'emploi

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Instructions

Read this document thoroughly before using the chiller, with special care concerning the safety regulations. Keep this document in the vicinity of the chiller to react appropriately in all cases.

Table of content

1. List of materials
2. Safety instruction
 - 2.1. Usage
 - 2.2. Disrespect of the safety rules
3. Installation and start-up
4. Usage
5. Maintenance and care
6. Problems
7. Service
8. Storage and transport
9. Technical data

VWR International BVBA
Researchpark Hassrode 2020
Geldenaaksebaan 464
B - 3001 Leuven
Belgium
www.vwr.com

Made in Germany

1. List of material

VWR RC10 Basic

Operating instructions GER/ENG/FRA/ESP/ITA

Power cable type Euro

Power cable type CH

Power cable type GB

Connection hose NW8 2 pieces

Screw nut M16x1 2 pieces

Hose PVC 1500 mm 2 pieces

Clamping piece for hose 4 pieces

2. Safety instructions

This device has been manufactured based on a state-of-the-art technology. It was tested extensively before leaving the production line: it was fully functional, passing every safety controls.

To ensure a safe and correct usage, the user has to read thoroughly the following document and take into

2.1. Usage

- The installation and start-up have to be carried out by qualified personnel
- The chiller can only be used by someone familiar with the present document
- Maintenance and care can only be performed by qualified personal
- Service operations have to be performed by a technician specified by the manufacturer

The operator of the chiller is responsible for the training of his operating personnel and has to take care that every user knows and understands the contents of this operating manual.

2.2. Disrespect of the safety rules

The chiller is conceived and built as a piece of laboratory equipment. It should be used in a coolant cycle connected to an external device (e.g. rotary evaporators). Other applications should be up-front agreed upon with the manufacturer.

Do not use the device in rooms which require explosion-proof apparatus.

3. Installation and start-up

Check whether the voltage of the socket outlet is the same as the voltage given on the device identification plate. The chiller can only be connected to a network with ground connection.

The chiller can only be operated for temperatures ranging between +15°C and +32°C.

1. Make sure that the Chiller RC10 is leveled, placed on a solid ground
2. Close the drain tap (Picture 1).
3. Connect the in- and out-valves of the RC 10 with the external device. Only use hoses which can withstand at least temperatures of -10°C and pressure of 1 bar and secure these with the hose clamps (enclosed).
4. Fill in the coolant using the opening on top of the Chiller (Picture 2). The coolant cycle should be filled up to the maximum level indication. The ball should be on the maximum mark (Picture 3).
Water must not be used as a cooling medium for temperatures below + 5°C. In case of operating temperatures below + 5°C please use a water-glycol mixture. We recommend a mixture of water and glycol (70-30), which does not freeze until -15°C.
5. Plug in the chiller. Use the enclosed power cable.
6. Switch on the chiller (switch on the front panel).
7. After the external circuit has been filled up it may be necessary to add further liquid. Please check the level indicator. Add some coolant if necessary as described above on Line 4.

Photo 1: Drain tap

Photo 2: Filling aperture


Photo 3: Liquid level indicator

4. Usage

The apparatus is not splash-proof.

The Chiller RC 10 Basic is adjusted to a factory-defined temperature of +7°C.

The temperature setpoint can be adjusted by using the keys on the command panel. Using the chiller as a heating system is not possible.

To start the device push the main switch on the front panel (Picture 4), a green lamp being led when the device is on.

The temperature is maintained in a range of $\pm 0.5^{\circ}\text{K}$ around the setpoint.

The Duo Chiller can cool two cooling cycles. Open the required loop taps (in and out) at the back of the chiller, depending on how many external devices that require cooling.

Under normal usage, the temperature fluctuates ± 2 Kelvin around the set point of +7°C.


Photo 4 : Main switch

5. Maintenance and care

Maintenance and care can only be performed by qualified personal

Turn off the chiller and unplug the power cable.

Clean up regularly (about twice a year) the air cooling apertures (in and out) at the front and back of the chiller. In case where the chiller is extremely dirty, remove the panels to clean up the air apertures.

Replace the coolant once a year.

Clean up the tool panels using only a damp cloth, never use a detergent solution.

After having completed the maintenance and care operations, the chiller can be plugged in and started up as described previously.

6. Service

Service operations have to be performed by a technician authorized by the manufacturer.

Parts of the chiller can only be exchanged if they are original parts. A more detailed instruction manual (about spare parts) is available from the manufacturer.

7. Problems

Problem description	Possible reason	Solution
The chiller is not cooling	The coolant is dirty	Maintenance and clean up of the chiller
The chiller is not cooling	Room temperature is too high	Change chiller location
The coolant pump is weak	The coolant is dirty	Remplace the coolant
The coolant pump is not functioning	The coolant pump is defective	Contact your service personnel

Before any servicing operation, please contact :

VWR International BVBA
 Researchpark Haasrode 2020
 Geldenaaksebaan 464
 B-3001 LEUVEN
 Belgium

www.vwr.com

8. Storage and transport

The packaging of the device was carried out with greatest care to make sure that damage by transport and storage is avoided. Upon receipt, both packaging and device have to be checked for damages.

Should such damages be noticed, please contact the manufacturer.

Should you need to move the chiller, drain the coolant before the transport. The chiller must be dry before it is moved.

Transport should be done only using the original packaging.

9. Technical data

Electrical network	230 V / 50 Hz
Power	max 600 W
Functioning temperature	+7°C (pre-defined)
Temperature variation	+/- 2 K
Cooling power at +15°C	500 W
Pump power	14 l/minutes *
Pump pressure	0.4bar *
Coolant solution	R134a
Size	280 x 520 x 360 mm
Weight	
Pannels	Powder coated Steel
Reservoir (material)	Stainless steel V2A
Reservoir (size)	3,5 Liters

* When use with a rotary evaporator the values change to 8 l/min and 0.6bar.

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WARRANTY

The manufacturer undertakes a 3 year warranty for the operational capability of the Chiller, when used under appropriate conditions and in accordance with appropriate operating instructions.

Except from the warranty are the following wearing parts :

- Fan motor
- Circulation pump
- Fuse

The manufacturer reserves the right to decide the validity of any warranty claim. In case of faults parts will be repaired or replaced free of charge, or a new replacement unit will be supplied.

Damages in transit are not case of warranty.

Any other compensation claims are excluded from this guarantee.

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Disposal

The directive 2002/96/EC of the European Parliament and of the Council requires electrical and electronic equipment marked with a crossed out trash can to be disposed of separately in an environmentally friendly manner. Contact an authorized waste management company in your country. Disposal with household waste or similar collections of municipal waste is not permitted!

CE Declaration of conformity

We, hereby declare under our sole responsibility that the products:

462-0137 RC-10 Basic Chiller
462-0138 RC-10 Digital Chiller
462-0139 RC-10 Duo Chiller

and to which this declaration relates, are in conformity with the following Council directives of the EU member states:

Council directive on harmonization of the laws of the member states
• relating to electrical equipment for use within certain limits 2006/95/EC
• relating to electromagnetic compatibility 2004/108/EC
• relating to pressure equipment 97/23

and meet the essential requirements of the mentioned directives.

This equipment has been tested accordingly to the following standards:

Title of the applied standards:	Standard reference number
Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements	EN 61010-1/2002-08
Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 2-010: Particular requirements for laboratory equipment for the heating of materials	EN 61010-2-010/A1
Electromagnetic Compatibility	EN 61316-1:2006
Refrigerating systems and heat pumps. Safety and environmental requirements. Operation, maintenance, repair and recovery	EN 378

The technical documentation required demonstrating that the products meet the requirements of these standards has been compiled and can be made available for inspection by the relevant enforcement authorities

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Ireland / Northern Ireland

VWR International Ltd / VWR
International (Northern Ireland)
Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail: sales@ie.vwr.com

Italy

VWR International s.r.l.
Via Stephenson 94
20157 Milano (MI)
Tel.: 02 332 03 11
Fax: 800 152 999
E-mail: info@it.vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland

Labart Sp. z o.o.
A VWR International Company
Limbowa 5
80-175 Gdansk
Tel.: 58 32 38 210
Fax: 58 32 38 205
E-mail: labart@labart.pl

Portugal

VWR International - Material de
Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: info@se.vwr.com

Switzerland

VWR International AG
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard
Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com

Asia Pacific

China

VWR International China Co.,
Ltd
Suite 1802 - 1803, Xing Ye
Bank Mansion, No 168,
168 Jiangning Road
Shanghai 200041, China
Tel.: +86- 21 521 388 22
Fax: +86- 21 521 33 933
E-mail: sales_china@vwr.com

India

VWR Lab Products Pte Ltd
2nd Floor, Front Wing, 135/12,
Brigade Towers
Brigade Road
Bangaluru 560025 India
Tel.: +91-2522-277876
(Mumbai)
Tel.: +91-80-41117124
(Bangalore)
Fax: +91-80-41117120
E-mail: vwr_india@vwr.com

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

**Go to www.vwr.com for the latest news, special offers and details
of your local VWR distribution partner.**