Lovibond® Water Testing

Tintometer® Group

T-CAL® • T-CAL_{plus}® • stabilized Formazin Stock Solution

Turbidity Standards

(EN) Instruction manual.....	3	(DE) Gebrauchsanleitung.....	12
(ES) Manual de instrucciones.....	22	(FR) Mode d'emploi.....	31
(PT) Manual de instruções.....	40	(IT) Manuale di istruzione.....	49
(NL) Gebruiksaanwijzing.....	58	(TR) Kullanım için talimatlar.....	67
(RU) Инструкции по эксплуатации.....	76	(ZH) 使用说明.....	85

Table of Contents

	Page
1 Introduction	4
1.1 General Information	4
1.1.1 Important Notice	4
1.1.3 Authorized use	4
1.1.5 Keep instructions	4
1.1.6 User qualification	4
1.1.7 Handling of hazardous chemicals	5
2 Incoming inspection and storage	6
2.1 Incoming inspection	6
2.2 Storage	6
3 Instructions	7
3.1 Handling	7
4 Troubleshooting	9
4.1 Technical support	9
5 Standards & Accessories	10
5.1 Turbidity Standards	10
5.2 List of Accessories	10
6 Abbreviations	11

1 Introduction

1.1 General Information

1.1.1 Important Notice

This manual provides important information on the safe operation of the product. Please read this manual carefully and familiarize yourself with the product before use.

These instructions describe the correct handling of T-CAL® and T-CAL $plus$ ® turbidity standards, as well as formazin stock solutions to achieve the best possible results. It should be noted that turbidity measurement is subject to many influences and requires careful work.

To carry out the verification or calibration of turbidity measuring instruments, follow the procedure described in the corresponding operating manual of the measuring instrument.

Turbidity standards are stabilised formazin suspensions available in bottles, bags and ready-to-use cuvettes.

1.1.3 Authorized use

The manufacturer's liability and warranty for damage is voided with improper use, failure to follow this manual, use by unqualified personnel, or unauthorized changes to the product.

The manufacturer is not liable for costs or damages that arise from the user or third parties due to the use of this product, especially in cases of improper use of the product or misuse or faults in the connection of the product.

The manufacturer assumes no liability for print errors.

Note the following points for safe use:

- The product may only be used according to the authorized use specified above.
- The product may only be used under the environmental conditions mentioned in this operating manual.

The product must not be used if:

- it is visibly damaged (e.g. after being transported)
- it was stored under adverse conditions for a lengthy period of time

1.1.5 Keep instructions

The manual must be kept in the vicinity of the product so you can always find the information you need.

1.1.6 User qualification

The operating personnel must be able to understand and correctly implement the safety labels and safety instructions on the packages and inserts of the products.

The user must be capable and able to read and understand this manual in order to familiarize themselves with the handling and to ensure safe use.

1.1.7 Handling of hazardous chemicals

Chemical and/or biological hazards may exist where this product is used. Abide by all governing laws, regulations and protocols when using this product.

For the development of products, Lovibond® pays close attention to safety. Some hazards from dangerous substances cannot be avoided. If self-produced tests or solutions are used, the responsibility concerning any risks caused by those tests or solutions lies with the user (personal responsibility).

Dispose of used or expired standards according to the information in the safety data sheets and local government regulations.

2 Incoming inspection and storage

2.1 Incoming inspection

Upon receipt of the turbidity standard, check the condition of the container for damage and leaks. If damaged, the standard should not be used for calibration or verification of measuring instruments.

Do not transfer the turbidity standard to other containers as this may cause contamination.

EN

2.2 Storage

Leave the standard in the original container. Do not decant it.

Store the standard for daily use in a dark, cool place protected from sunlight and at a temperature between 5°C and 20°C. Observe additional information on the standard. Short-term storage or transport up to 30°C is possible. Storage above 30°C will reduce the shelf life.

Storage above 40°C can lead to destruction of the standard.

If the standard is not to be used for a longer period of time, it is best to store it in the refrigerator. Freezing the standard will not destroy it.

The shelf life of the standards is 24 months after production. Note the expiry date of the suspension before using it.

3 Instructions

3.1 Handling

Before using the standard, wait until its temperature has adjusted to the measuring environment. The turbidity standards can generally be used in the temperature range from 15°C to 25°C. For details, please refer to the information on the standards.

Before using the standards, mix them sufficiently so that settled substances are evenly distributed again. Proceed as follows in detail:

T-CAL® turbidity standards and stock solutions in bottles.

Place the bottle on a flat surface and roll the bottles back and forth for 3 minutes. In this way you ensure that the sediment is evenly distributed in the solution. Never shake a bottle with turbidity standard, this creates air bubbles in the suspension. Air bubbles can lead to a disturbance of the measurement, and thus falsify the measurement results. If air bubbles have formed, let the standard sit for some time to allow them to dissolve and remix the standard.

T-CAL^{plus}® turbidity standards in bags

Knead the bag well for at least 2 min and make sure that it is evenly distributed within the suspension. If there are air bubbles in the bag, avoid incorporating them into the suspension. Small air bubbles in the bag do not enter the PTV meter when filling via the connection tube and do not interfere with the measurement. If in doubt, let the bag rest for some time and remix.

T-CAL® ready-to-use cuvettes

Gently swirl the cuvettes at least 30 times within 30 seconds and make sure that the suspension is evenly distributed (the < 0.1 NTU cuvette must never be swirled or shaken). Avoid shaking and the associated incorporation of air bubbles within the suspension. If necessary, let the vial rest for some time and remix.

If air bubbles have been incorporated into the standard, it may take several hours, especially with low standards, for the measured values to stabilise in the target range. After mixing, at least one minute should pass before the standards are used for verification or calibration (some instruments already take this into account with a built-in timer).

It is advisable to carry out several measurements directly one after the other and to pay attention to a stable measured value or, in the case of online measurements, to observe the course. Fluctuating or strongly falling measured values can be an indication of air bubbles in the solution.

After a short time, the turbidity substances begin to settle again. The higher the turbidity, the faster the settling will begin. Turbidity standards with < 2 NTU should be remixed after 2 hours at the latest, standards with 2 to 40 NTU after 20 minutes, 40 to 200 NTU standards after 5 minutes, and standards with higher turbidity after 3 minutes at the latest. For remixing, it is sufficient to swirl them carefully several times.

When using cuvettes, first rinse them with the standard. Fill cuvettes by pouring the suspension along the edge of the cuvette. This avoids the formation of air bubbles.

Cuvettes used for measurement should be treated on the outside with a silicone oil to compensate for scratches on the material that falsify the measurement. To do this, apply a small amount of silicone oil to the side of the cuvette and spread it evenly. Wipe off excess oil by turning the cuvette in a dust-free cloth. This is especially necessary for measuring low turbidity values (< 1 NTU). Only touch the cuvettes by the cap to avoid fingerprints.

Close the bottle or bag again immediately after withdrawing the suspension to avoid contamination.

Withdrawn standard is intended for single use. After verification or calibration, do not return the used suspension to the bottle, as contamination of the used suspension will render the standard unusable. For standard solutions in the bag, mark the bag after the suspension has been returned to avoid reuse.

It is recommended to clean cuvettes with diluted hydrochloric acid (1:1) and then rinse with ultrapure water (< 0.1 NTU) before using them.

After use, always store the standard according to the instructions in the chapter Storage until next use.

Preparation of low turbidity suspensions from stock solutions

The dilutions should be prepared fresh immediately before use; they have a maximum shelf life of one day.

Use only ultrapure water to prepare low turbidity suspensions from the formazin stock solutions. Measure the turbidity of the water before using it. The turbidity should be below < 0.1 NTU. If necessary, the water must be degassed to remove invisible air bubbles. If the turbidity value is still higher, the water must be filtered (e.g. using a membrane filter $0.1 \mu\text{m}$) and checked again.

The glassware used to prepare the dilutions should first be cleaned with diluted hydrochloric acid (1:1) and then rinsed several times, first with tap water and then with ultrapure water.

To prepare the dilution, fill the volumetric flask about halfway with ultrapure water. Then add the calculated amount of the prepared (see section Handling) stock solution. Mix the preparation and then fill the flask up to the mark and mix again. When mixing, avoid the introduction of air bubbles by shaking too vigorously. If air bubbles are incorporated, there may be a conflict with the shelf life because of the subsequent resting.

It is recommended to mix the suspension for at least 1 hour using a magnetic stirrer. To prevent excessive mechanical stress, the stirring speed should not exceed 300 rpm.

Do not prepare dilutions from other dilutions already prepared (multiple dilution).

The safety datasheets of the chemicals comprise all instructions on safe handling, occurring hazards, preventive actions and actions to take in hazardous situations.

4 Troubleshooting

4.1 Technical support

If you have further questions about the handling or use of the standards, please contact our technical support.

Lovibond contact: techsupport@lovibond.com

5 Standards & Accessories

5.1 Turbidity Standards

Part Number	Description	
48012050	T-CAL®-Standard	<0.1 NTU 500 ml
48011050	T-CAL®-Standard	0.3 NTU, 500 ml
48011150	T-CAL®-Standard	1.0 NTU, 500 ml
48012250	T-CAL®-Standard	5.0 NTU, 500 ml
48011250	T-CAL®-Standard	10 NTU, 500 ml
48012350	T-CAL®-Standard	20 NTU, 500 ml
48011550	T-CAL®-Standard	100 NTU 500 ml
48011650	T-CAL®-Standard	200 NTU 500 ml
48011750	T-CAL®-Standard	800 NTU 500 ml
48011850	T-CAL®-Standard	1000 NTU, 500 ml
48012950	T-CAL®-Standard	4000 NTU, 500 ml
48010035	T-CAL <i>plus</i> ®-Standard	0.3 NTU, 350 ml
48010135	T-CAL <i>plus</i> ®-Standard	1.0 NTU, 350 ml
48010235	T-CAL <i>plus</i> ®-Standard	5.0 NTU, 350 ml
48010335	T-CAL <i>plus</i> ®-Standard	20 NTU, 350 ml
194150	Set Turbidity Standards T-CAL	<0.1 / 20 / 200 / 800 NTU (in cuvettes)

5.2 List of Accessories

Part Number	Description
54010435	PTV Cleaning Solution
54011010	Detergent Solution
54013003	Scale Remover for PTV-Series
194295	Silicone Oil
197635	Cleaning Cloth
197655	Set of 12 sample cuvettes with black lid, Height 55 mm, ø 24 mm

6 Abbreviations

NTU	Nephelometric Turbidity Unit Measurement according to EPA 180.1 at 90° (scattered light)
FNU	Formazine Nephelometric Units Measurement according to ISO 7027 at 90° (scattered light)
FAU	Formazine Attenuation Units Measurement according to ISO 7027 at 0° (transmitted light)

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Trademarks of the Tintometer Group
Tintometer Inc.