EQUIPEMENT DCO 6 PLACES/COMPLET UNIT COD 6 <i>SAMPLES</i>	4000638
EQUIPEMENT DCO 12 PLACES/COMPLET UNIT 12 <i>SAMPLES</i>	4000639
EQUIPEMENT DCO 20 PLACES/COMPLET UNIT 20 <i>SAMPLES</i>	4000640
BLOC-DIGEST MACRO 6 PLACES/6 <i>SAMPLES</i>	4000629
BLOC-DIGEST MACRO 12 PLACES/ <i>SAMPLES</i>	4000630
BLOC-DIGEST MACRO 20 PLACES/ <i>SAMPLES</i>	4000631
BLOC-DIGEST MICRO 12 PLACES/ <i>SAMPLES</i>	4001047
BLOC-DIGEST MICRO 24 PLACES/ <i>SAMPLES</i>	4001048
BLOC-DIGEST MICRO 40 PLACES/ <i>SAMPLES</i>	4001049
BLOC MACRO/ <i>BLOCK</i> 6 PLACES/ <i>SAMPLES</i>	4000507
BLOC MACRO/ <i>BLOCK</i> 12 PLACES/ <i>SAMPLES</i>	4000508
BLOC MACRO/ <i>BLOCK</i> 20 PLACES/ <i>SAMPLES</i>	4000509
BLOC MICRO/ <i>BLOCK</i> 12 PLACES/ <i>SAMPLES</i>	4001050
BLOC MICRO: <i>BLOCK</i> 24 PLAZAS/ <i>SAMPLES</i>	4001051
BLOC MICRO/ <i>BLOCK</i> 40 PLAZAS/ <i>SAMPLES</i>	4001052
RAT-2	4001538

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ACCESSORIES

Accessories that can be purchased separately:

* Metal heater blocks.

These metal heater blocks must work together with the control unit RAT 2. They never must be connected directly at mains.

Modèles Models	Code Code	N° de places Samples	Ø tube Ø tube	Dim. extérieures Overall dimensions			Température max. Maximum temperature	Consom. Consumption	Poids Weight (W)
				Haut. Height	Larg. Width	Prof. Depth			
MACRO	4000507	6	42	18	33	28	450	1500	18
	4000508	12			39	33		2000	25
	4000509	20			44	39		2500	31
MICRO	4001050	12	26		33	28		1500	16
	4001051	24			39	33		2000	22
	4001052	40			44	39		2500	27

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Deux sondes de température type K.

* Control unit R.A.T 2 Cod.: 4001538

Alphanumeric screen, programmable temperature between 50°C and 500°C. Programmable time periods between 1' and 10 h (600'). 20 four-stepped user programmes Two type-K temperature probes.

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* Tube rack with stand

In special chemically-treated duraluminium plate, with handles and side plates to prevent heat loss.

Modèles Models	Codes Code	Places Samples	Dim. extérieures Overall dimensions		
			Haut. Height	Larg. Width	Prof. Depth
MACRO	4005071	6	15	17,5	12,5
	4005081	12		23	18
	4005091	20		28,5	23,5
MICRO	4001053	12		17,5	12,5
	4001054	24		23	18
	4001055	40		28,5	23,5

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*** Fume exhaust assembly.**

*It comprises a multiple fume assembly and a stand.
Made of stainless steel, with borosilicate glass receivers.*

Modèles Models	Codes Code	Places Samples	Dim.extérieures Overall dimensions		
			Haut. Height	Larg. Width	Prof Depth
MACRO	005072	6	15	18	12
	005082	12		23	18
	005092	20		29	23
MICRO	001056	12		18	12
	001057	24		24	18
	001058	40		40	23

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*** Digestion and distilling tubes Cod.: 4042300**

*MACRO series of 250 ml of capacity.
Graduated up to 100 ml, Ø42 mm and 300 mm high.*

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*** Digestion and distilling tubes Cod.: 4001045**

*MICRO series of 100 ml of capacity.
Ø26mm and 300 mm high.*

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*** Water jet pump.****Cod.: 4000633**

In PVC, specially designed for the 20 and 40 unit fume exhaust assembly.

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*** Water jet pump.****Cod.: 7000923**

Metallic. Suitable for the 6, 12 and 24 unit fume exhaust assembly.

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EQUIPMENT DESCRIPTION

The BLOC-DIGEST digestion units are made, mainly, to do the digestion and determination of Kjeldahl nitrogen and measurement of C.O.D. Other possible applications are hydrolysis of samples in determination of fats, etc.

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The digestion units have a heater block, a control unit and a system to eliminate fumes.

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Heater block allocate inside digestion tubes in order to make an optimal heat transmission and reach a good temperature distribution among tubes, is necessary a suitable tubes diameter for the block that will be used.

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Heater block have two safety probes incorporated to avoid 550°C overtemperature.

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Heater block temperature is controlled by control unit which is apart of block in order to avoid its damaging by acid, vapours or splashing. This unit (RAT 2) allows work until 20 programmes with 4 steps (temperature/time) in each of them.

In digestion samples to Kjeldahl or hydrolysis stepped heating allows an optimal foam control because give the possibility to evaporate humidity of sample before make digestion at 400°C.

To remove the fumes produced during digestion you can use a fume collector together with a water jet pump, an extraction vent or a scrubber.

Union between fume collector and sample tubes is not hermetic. This allows that small explosions do not break glassware. Fume extraction conveys fume to the waste pipe because of the vacuum that produced by the water jet pump. So it is essential good operation of it wich have to do a vacuum level of -0,6 bar.

On C.O.D. technique fumes extraction is not used, but a special glassware condensate the steam produced and liquid returns to sample tube.

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TECHNICAL FEATURES

Voltage supply 115-230 V 50/60 Hz as indicated on the characteristics plate.

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INSTALLATION

Place metallic block and control unit on an horizontal, flat and out levelled surface.

With the connection cable supplied connect the RAT unit and the corresponding block. Connect the RAT unit at mains.

Be in mind that if you have to use a water jet pump you need install equip near a water tap an a waste pipe. In other case it is advisable install metallic block under vent air extractor

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Be sure that the voltage supply is the same as the one indicated on the characteristics plate.

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Do not use the apparatus if it is not earthed.

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If you have to change the plug cap be in mind the following:

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
OPERATION RAT-2 UNIT


CONTROL PANEL:

1. Main switch (POWER) with green indicator lamp.
2. Alphanumeric screen
3. "Menu" key
4. "Validate" key
5. "Increase value" key
6. "Reduce value" key
7. "Start/Stop" key
8. Overtemperature pilot lamp

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
PROGRAMMING R.A.T. 2 TEMPERATURE CONTROLLER.

1. Information on the display:

Note: Pay attention where the "cursor" is.

2. Status:


There are 4 status:

- STP: "Stop". This is the default state. Controller is waiting from user orders.
- RUN: On this state the controller is running a temperature / time program.
- PRG: Programming: State for program creating or editing.
- CFG: Configuration basic options as language, or calibration.

Note: On states: PRG and CFG, the controller returns to STP state if any key is pressed after 10 seconds.

Note: Once a program is running can be stopped by pressing R/S key.

3. Main menu:

4. Programming menu:

```
S °C MIN PR STA
1 125 30 1 STP
```

Select the **program** to edit. Use keys ▲▼
(Default: program 1)

 + 5 seg.

```
S °C MIN PR STA
1 125 30 1 PRG
```

Select the **step** to edit. Use keys ▲▼
(Default: step 1)

Validate ◀▶

```
S °C MIN PR STA
1 125 30 1 PRG
```

Set temperature for the step. Use keys ▲▼
(Temperature must be always higher than the previous step)
(Default: previous step temperature)

Validate ◀▶

```
S °C MIN PR STA
1 125 30 1 PRG
```

Set time (minutes) for the step. Use keys ▲▼
(Default: OFF)

Validate ◀▶

```
S °C MIN PR STA
2 125 20 1 PRG
```

Automatically jump to the next step.

⋮

```
S °C MIN PR STA
3 145 OFF 1 PRG
```

Program finish when there is an step with the time set at value OFF.

Validate ◀▶

Once, a step with the time set at value OFF, is validated, it left programming status and returns to STOP status.

```
S °C MIN PR STA
1 145 31 1 STP
```

It jumps to STOP state. Ready to RUN the program.

5. Checking program parameters:

The temperature / time parameters can be checked from the STP state.

```
S °C MIN PR STA
1 125 30 1 STP
```

Select the **program** to check. Use keys ▲▼

Use keys ◀▶ to move at **step** position on the display

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WORKING OF THE SAFETY THERMOSTAT

If a bad function of the regulation produce an abnormal increase of the temperature, the electronic safety thermostat with its second temperature probe, avoid that temperature reach more than 550°C. If this temperature exceed 550°C for 10 seconds, an acoustic alarm sounds and the indicator lamp (8) of the control panel lights up. If temperature low under 550°C pilot (8) lights off but alarm condition remain on display.

In case that the failure persists, call up the authorized technical service.

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ADVANTAGES OF A STEPPED TEMPERATURE CONTROL OVER A CONTINUOUS ONE

Continuous process has the problem that when sample contain humidity this cause foam. This may bring about losses of sample and damages on metallic block because of spilt or splashes of acid.

Processus Process		Pas 1 Step 1	Pas 2 Step 2	Pas 3 Step 3
Kjeldahl	(viande, fromage) (meat, cheese)	125°C / 30'	270°C / 30'	400°C / 120'
	(céréales) (cereal)	125°C / 15'	300°C / 15'	400°C / 90'
	(eaux résiduelles) (waste water)	125°C / 60'	270°C / 120'	400°C / 120'
DQO COD		150°C / 30'	170°C / 60'	---
Graisses hydrolysées Fat hydrolysis		125°C / 30'	150°C 30'	---

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These parameters could differ depending on cain of sample and the water quantity that it containing.

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
DIGESTION PROCESS:

WARNING!

Digestion process needs strong chemical products. Pay attention to its security sheets about how use it appropriately.

Only staff suitable trained on chemical products operation should work on digestion process of samples.

- 1.- *Connect the control unit and the digestion one.*
- 2.- *Place sample tubes into the tube rack.*


- 3.- *Introduce into the tubes: samples the catalyst and the digestion acid.*
- 4.- *Move the set at digestion block, using the handles.*

- 5.- *Insert the fume exhaust set and the water jet pump or C.O.D. coolers.*
- 6.- *Start up a stored programme in RAT 2 being in mind page 12 table values .*

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MAINTENANCE

Before removing the cover disconnect the apparatus from the mains.

The manipulation by unauthorized people can cause irreparable damage. Take it to one of the J.P.SELECTA authorized technical services.

CLEANING:

For cleaning of apparatus use the following products:

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Cleaning of stainless steel: Alcohol

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Cleaning of plastic: Alcohol with cotton duster.

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It is advisable to clean the glasses of the fume collector after each digestion.

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Be careful when manipulating the fume collector. After digestion it can contain acid from condensation. Use gloves.

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For the cleaning, remove the glasses and rinse just with water. If it is possible, use an ultrasonic cleaner for ten minutes.

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
To guarantee the safety of the equipment, all spare parts must be purchased from J.P.SELECTA, s.a.

Description / Description	Code/Code
Résistance/Heater element 1600W (Macro 6)	39070
Résistance/Heater element 1600W (Micro 12)	
Résistance/Heater element 2100W (Macro 12)	39069
Résistance/Heater element 2100W (Micro 24)	
Résistance/Heater element 2500 W (Macro 20)	39068
Résistance/Heater element 2500 W (Micro 40)	
Sonde type K / K thermocouple probe	43026
Joint viton blanc (Macro) /White viton gasket (Macro)	21157
Joint torique viton vert (Macro)/Green viton o-ring gasket (Macro)	21158
Joint viton noir (Micro)/Black viton gasket	21167
Tube viton / Viton tube	46084
Rampe verre (Macro 6)/Ramp glassware (Macro 6)	47003
Collecteur verre(Macro 6)/Colector glassware (Macro 6)	47007
Rampe verre (Macro 12)/Ramp glassware (Macro 12)	47002
Collecteur verre (Macro 12)/Colector glassware (Macro 12)	47006
Rampe verre (Macro 20)/Ramp glassware (Macro 20)	47001
Collecteur verre (Macro 12)/Colector glassware (Macro 20)	47005
Rampe verre (Micro 12)/Ramp glassware (Micro 12)	47013
Collecteur verre (Micro 12)/Colector glassware (Micro 12)	47014
Rampe verre (Micro 24)/Ramp glassware (Micro 24)	47015
Collecteur verre(Micro 24)/Colector glassware (Micro 24)	47016
Rampe verre 40)/Rampa glassware (Micro 40)	47011
Collecteur verre (Micro 40)/Colector glassware (Micro 40)	47012

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
Color funda: (Cover color)	(+)	(-)
Amarillo (Yellow)	Amarillo (Yellow)	Rojo (Red)
Verde (Green)	Rojo (Red)	Verde (Green)
Verde (Green)	Verde (Green)	Blanco (White)
Rojo (Red)	Rojo (Red)	Azul (Blue)
Magenta	Amarillo (Yellow)	Magenta

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KJ-02 KJELDAHL NITROGENUM DETERMINATION METHOD. (Método DIRECT METHOD).

1. Overview:

This methods consists on mineralizing the sample with concentrated sulphuric acid and alkanizing with NaOH. The ammonium liberated is carried by distillation and recovered on boric acid. The subsequently titration with Cl H allows the calculation of the amount of the ammonium on the sample.

2. Reagents:

- Sulphuric acid 96% (d=1.84).
- NaOH, dilution 35% (p/v).
- mixed of indicator, special for ammonium titration.
- Kjeldahl Catalyst
- Boric Acid 4% (p/v).
- HCl 0.25N.
- Distilled water
- Pumice stone, granules

Note: It is very important that all reagent will be free of nitrogen.

3. Material:

- Balance of resolution 0.1 mg.
- Digestion Unit (Bloc-Digest).
- Fume Removal.
- Distillation Unit Pro-Nitro I or II
- Titration Burette.

4. Digestion:

Weigh near 1 g. of sample perfectly milled and homogenized into a paper exempt of nitrogen and introduced in a digestion tube.

Add to sample tube 10 g. of Kjeldhal catalyst, 25 ml of sulphuric acid at 96% (d=1.84), and some granules purnice stone.

Put the digestion tubes with the sample into the Bloc-digest with the fume removal working.

between 350° and 420°C and in a period of time variable between 1 and 2 hours (see page 12 table).

At the end of the process, the obtained liquid is of green colour or transparent blue depending on the used catalyst.

Let cool down the sample at ambient temperature.

It is important to avoid the precipitation shaking from time to time.

Dosse slowly 50 ml of distilled water in each tube of sample. (Pay attention with the violence of the reaction).

Let cool down the sample at ambient temperature.

If there is a precipitation, shake or hot slightly.

5. Distillation.

Dosse 50 ml of Boric acid in an Erlenmeyer flask and some drops of mixed of indicator. Put the Erlenmeyer below the refrigerant paying attention to let the output tube be immersed under the Boric acid.

The distillation has to be extended the enough time in order to be distilled a minimum of 150 ml, approximately between 5 and 10 minutes.

6. Blank:

After the distillation of a sample, make a blank test, applying the described method, but using 5 ml of distilled water instead of sample.

7. Titration:

Titrate with Hydrochloric acid 0.25N the distilled obtained until the solu-

tion changes from green to violet colour.

Calculate the quantity of nitrogen detected by means of formulae:

$$\% \text{ Nitrogen} = \frac{1.4 \times (V_1 - V_0) \times N}{P}$$

$$\% \text{ Protein} = \% \text{ Nitrogen} \times F$$

Where:

P = Weigh g of sample.

V₁ = HCl consumption on tritration. (ml)

N = Normality del HCl

V₀ = HCl Consumption on blank tritration. (ml)

F = Conversion factor to pass from content in nitrogen to content in proteines. For the proteine it is commonly used one value of 6.25. For more accuracy, in order to distinguish the quality of the proteine according to the nature of the sample, there can be used other factors of conversion.

Do t

**CET ESSAI DEMANDE LA MANIPULATION DE SUBSTANCES DANGEREUSES
MANIPULATION OF CORROSIVE PRODUCTS IS REQUIRED**

DG-04 PREVIOUS HYDROLISIS FOR THE FAT DETERMINATION: Cheese, egg and meat.

Material needed (for 6 samples):

- 6 glass funnels Ø100mm
- 6 filter papers Ø100 normal, (Not slow).
- 6 digestion tubes (SELECTA Code 4042300).
- Digestion Bloc + RAT + Fume removal. (SELECTA Codes 4000507 + 4000051 + 4005072)
- 6 erlenmeyer 500ml.
- Bath in order to hot the water at 60°C

Reagents:

- 1L. of HCl 4N (350ml de HCl 35% and complete until 1L. with distillated water).
- Granules of pumice stone.

Procedure:

- Number the digestion tubes.
- Weigh around 3gr. of sample and to put it into the digestion tube.
- Add 150 ml of HCl 4N in each tube.
- Add a spoonful of granules of pumice stone.
- Put this into the bloc and to digest at 50°C - 200°C during 1 hour (see page 12 table). It has to boil slowly. See the below photo.
- Retire the digestion Bloc and let cool down.
- Prepare the Erlenmeyer with the funnels and the filter paper. See the below photo.
- Number the erlenmeyer.
- Empty the content of the sample tube into the funnel with filter paper.
- Clean the tube with hot water (±60°C) y and empty into the funnel. To repeat it until water that goes out of the tube is transparent.
- Put hot water (60°C) on the funnel.
- Empty the erlenmeyers.
- Repeat the operation until the water of erlenmeyer be transparent (at minimum 10 times).
- Let dry the paper filter.
- You will observe that keeps the fat of the sample.
- Proceed to the extraction with the det-gras equipment. Let dry 2-3h, and to extract in two phases: Boiling, rinsing, let dry 30 minutes extraer en dos fases: Boiling, rinsing, dry 30 min

**CET ESSAI DEMANDE LA MANIPULATION DE SUBSTANCES DANGEREUSES
MANIPULATION OF CORROSIVE PRODUCTS IS REQUIRED**

DQO-01 DETERMINATION OF CHEMICAL REQUEST OF OXYGEN (C.O.D.)**Introduction**

The most common methods in order to estimate the organic contamination are based on the global determination of oxidizable materia existing in the water and that are liable to oxide from the oxygen of the medium, since the methods that determine each of organic substances are complicated. The methods of global test of organic materia contaminated used currently are as follows:

- The chemical request of oxygen (C.O.D.). Measure of the oxidizable materia by the potassic dychromatic in an acid medium.

- The biological request of oxygen (DBO). Measure of the organic materia biodegradable.

Some times, the limit of contamination is determinated through the expression: $(\text{DBO} + \text{C.O.D.}) / 2$.

Overview:

The chemical request of oxygen, C.O.D., can be considered like a measure approximative of the theoretical request of oxygen, that is the quantity of oxygen consumed by the whole oxidation of the organical constituents to inorganic products. The C.O.D. is defined like the quantity of oxygen that is equivalent to the dychromat consumed by the dissolved materia and in suspension, when the sample is treated by the above oxidant in the defined conditions.

Target and field of application:

Determination of the C.O.D. for residual waters with value higher than 30 mg/l and less than 700mg/l (without dissolution)

This method can be applied to waters that after being dissolved contain more than 2000 mg/l of chlorure ions, like the water of sea and brackish waters.

Description of the method:

Most of the organical materia is oxidizable by the effect of a hot mixture of dychromat and sulfure.

Take the mixture and put back flow together with the acid solution with an overrun of the chrome salt .

After making the digestion, tritate the overrun with an iron amonium sulphate (FAS), that allows to know the oxidizable organical materia, calculated in terms of oxygen equivalent (mg O /l).

This relation is of 1.5 mol of oxygen by each one of dychromat.

Reagents:

Sulphuric acid 4 mol/l: to add to 500 ml. of water, 220 of sulphuric acid . Let cool down and dilute until 1000 ml.

Sulphuric acid - silver acid: Add 10 g of silver sulphate to 35 ml of water. Add 965 ml. of sulphuric acid. Let lie 1 or 2 days in order to be dissolved. You can help the dissolution shaking it.

Iron sulphate (II) and of amonium (FAS), titrated solution, $(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$ 0.12 mol/l Dissolve 47 g of FAS in water. Add 20 ml. of sulphuric acid ($r = 1.84$ g/ml), Let cool down and dissolve with water until 1000 ml.

This solution has to be contrasted each day as follows:

- Dilue 10 ml of patron solution of potassic dichromat until 100 ml. with sulfuric acid 4 mol/l.

The solution of iron sulphate (II) and amonium is titred, at the presence of 2 or 3 drops of ferroine as indicator.

The concentration c, expressed in moles by liter of the iron sulphate solution (II) and of amonium, is obtained by the follow formule:

$$C_{\text{FAS}} = \frac{0.04 \text{ M} \times 6 \times 10 \text{ ml}}{V_{\text{FAS}} \times 1} = \frac{2.4}{V_{\text{FAS}}}$$

4. Potasic dichromat, standard solution $\text{K}_2\text{Cr}_2\text{O}_7$ 0.040 mol/l that contains a mercury salt (II): to dissolve 80 g. of mercury sulphate (II) (HgSO_4) into 800 ml of water. Carefully add 100 ml. of sulphuric acid ($r = 1.84$ g/ml). Let cool down and dissolve into the solution 11.768 g of potasic dichromat, previously let dry at 105°C during 2 hours. Transfer the solution to a gauged balloon flask and dilute until 1000 ml.

The solution is stable during at least 1 month.

Potassium Hydrogen Phthalate, standard patron. Dilute 0.4251 g previously dried at 105°C, in the water and fill until 1000 ml with water.

The solution has a theoretical value of COD of 500 mg/l. This solution remains stable at least one week if it is conserved at 4°C.

Ferroine, indicated solution: 0.7 g of iron sulphate are dissolved (II) heptahydrated $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$ on water. Add 1.5 g. of fenantroine - 1,10 monohydrated and it is shaken until its dissolution. Then dilute until 100 ml.

Balloons regulators of ebullition Ø 3 mm.

**CET ESSAI DEMANDE LA MANIPULATION DE SUBSTANCES DANGEREUSES
MANIPULATION OF CORROSIVE PRODUCTS IS REQUIRED**

Material:

C.O.D. equipment composed by:

- Metallic bloc heater.
- Temperature controller
- Tubes for digestion C.O.D.
- Stand for tubes holder.
- Refrigerants for C.O.D.
- Support for refrigerated tubes.
- Burette of precision for titration.

Sampling:

1- The samples for the laboratory have to be collected preferably into glass or polyethylene flasks. The samples are analysed as soon as possible, never more than 5 days of their collection. If it is necessary to conserve the samples before making the tests, add 10 ml. of sulfuric acid by 1 liter of sample and to conserve between 0°C and 5°C. The flasks are shaken and ensure that the content is well homogenized before taking some part of the test.

2. To introduce in the digestion tube of D.C.O. 10,0 ml of the sample to be analysed. If the value of the D.C.O. of the sample is higher to 700 mg/l, this is made a dilution of the original sample in order to obtain a value of D.C.O. between 350 mg/l and 700 mg/l.

Blank test:

3. A blank test is done at the same time that the determination, following the same operational procedure, but instead using the sample, you use 10 ml. of distilled water.

The theoretic demand of oxygen of this solution is of 500 mg/l; the experimental procedure is satisfactory if it is obtained at least 96 % of this value.

Determination:

5. Add 5.0 ml of potassic dichromat and some glass ballons to regulate ebullition to the sample and it is homogenized carefully.

6. Add slowly 15 ml. of silver sulphate - sulfuric acid, shaking carefully the tube with a circular movement, to cool down with cool water of the tap or an ice bath, in order to avoid any loss of organical volatile substances.

7. Put the tubes dully prepared into the stand tubes holder on their support. With the help of a stand, to introduce into the heating bloc.

8. Accoply on each tubes one refrigerant in order to joint correctly the polished gaskets.

9. Programm the temperature of working (170-200°C) and the time of lenght (around 2 hours) of the back flowing into the temperature regulator and RAT time. Put the selector of ramp (on the back panel of the unit RAT) into the position C.O.D.

10. Push on «START» and at this moment the metallic bloc will hot where the tubes are put. This will be by a progressive way and when it reaches the preselectionated temperature, it will start to count the time of back flowing.

11. After been finished the time of back flowing, the heating in the metallic bloc stops automatically and you hear an acustic indicator. At this moment, you should push on «STOP» in order to finish the process.

12- To let cool down the tubes for a while and clean the projections that have been produced in the inner wall of the refrigerants into the recipient of sample with a little volume of distilled water.

13. To separate the refrigerants from the tubes, and the tubes go outside with the help of the stand, that is put on the support. The sample is dilued with 75 ml. of distilled water and let cool down the ambien temperature.

14. The excess of dichromat is titred with iron sulphate (II) and of amonium at the presence of 1 or 2 drops of ferroine. (It has to change from green-blue colour to red colour).

NOTES:

When the white is titred, you should obtain a demand around 10 ml of FAS.

When you titre the core test you should obtain a demand of around 4.8ml of FAS.

**CET ESSAI DEMANDE LA MANIPULATION DE SUBSTANCES DANGEREUSES
MANIPULATION OF CORROSIVE PRODUCTS IS REQUIRED**

DQO-01 DETERMINATION OF THE CHEMICAL DEMAND OF OXYGEN (C.O.D.)

Calculation of the results:

15- The chemical demand of oxygen, C.O.D., is expressed in mg of oxygen per litter, it is calculated by the following formulae:

$$\text{C.O.D.} = \frac{8000 \times c \times (V_1 - V_2)}{V_0} \text{ (mgO}_2\text{/l)}$$

where:

c is the concentration in mol/l of the solution of the iron sulphate (II) and of amonium already calculated (3).

V is the demanded volume, in ml, of the sample used for the determination before the dilution (if it exists).

V1 is the demanded volume, in ml, of the solution FAS (3) in the titration to the white.

V2 is the demanded volume, in ml, of the solution FAS (3) in the titration of the samples.

Volume of reagents for other ranges of C.O.D. determination

Note: for some of these ranges it could be necessary another shapes of the tubes and metallic bloc. These are specified here for your information.

	Macro/50ml	Macro/10ml	Micro/10ml	Micro/low COD
Range:	50-800 mg/l COD	50-800 mg/l COD	50-800 mg/l COD	50-800 mg/l COD
Sample:	50ml	10ml	2ml	2ml
H2SO4:	75ml	15ml	2.5ml	2.5ml
K2Cr2O7:	0.3065g *	0.0613g *	0.00245g	0.00245g
AgSO4:	0.75 g	0.15 g	0.03 g	0.03 g
HgSO4:	1.0 g	0.2 g	0.03 g	0.03 g

These values can be reduced 1/10 for the test with samples which COD is on the range 0-80.

**CET ESSAI DEMANDE LA MANIPULATION DE SUBSTANCES DANGEREUSES
MANIPULATION OF CORROSIVE PRODUCTS IS REQUIRED**

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GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damage caused by incorrect use or causes beyond the control of J.P. SELECTA, S.A.

Any manipulation of the apparatus by personnel not authorized by J.P. SELECTA, S.A. automatically cancels the guarantee.

DECLARATION DE CONFORMITE "CE" "EC" CONFORMITY DECLARATION

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