**Symbolen**

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CE 0120 mark of conformity for the following models, which are supplied in individual, sterile packages:

16534-GUK
16534-K
16555-GUK
16555-K
16596-HYK
17597-K
17598-K

CE mark of conformity for the following models available in bulk in non-sterile, PE bags:

16534-Q
16555-Q
16596-HYQ
17597-Q
17598-Q

Applications

- Sterile filtration and ultracleaning of liquids: 0,2 µm pore size
- Particulate removal from liquids; microbial reduction filtration and particulate removal from dry air and gases: 0,45 µm pore size

Minisarts® can also be used for sterile filtration and particulate filtration of liquid media for medical purposes.

Important Note

To ensure sterility, do not use a Minisart® unit if the individual packaging is damaged. Never reuse a Minisart® filter unit. They are intended for single use only. Minisarts® are designed to filter in either direction (bidirectional). However, once you start filtering in one direction, do not reverse it.

During handling, make sure to avoid contamination of the filtrate outlet and the attached accessories (e.g., needle).

Take care when using syringes with a volume less than 10 ml since they can generate a pressure greater than the maximum operating pressure of the Minisart® unit.

Minisart® has been tested for its suitability for in vivo use in conformance with the applicable regulations and guidelines. It does not contain any toxins.

Specifications

| | |
|--|---|
| Housing material | MBS-copolymer |
| Membrane filter material | Wetting-agent free cellulose acetate or PTFE if HY-coded |
| Filter diameter | 28 mm |
| Filtration area | 6,2 cm ² |
| Inlet connection | Female luer lock |
| Outlet connection | Male luer lock; 17597 and 17598 have a luer slip |
| Minimum effective burst pressure for housing | > 6 bar |
| Hold-up volume | 0,1 ml |
| Sterilization | Presterilized with ETO; GU-coded units presterilized by gamma irradiation |

| | |
|-------------------------|---|
| Pyrogenicity | Passed the USP test |
| Bacteria challenge test | with: 0,2 µm pore size passed the 10 ⁷ Brevundimonas diminuta/cm ² test |
| | No inhibition with L-929 or MRC-5 cells |

Directions for Use

- Remove a Minisart® from the box containing individual, sterile blister units.
- For liquid filtration, draw liquid into a syringe that has a luer lock connector.
- Starting at one corner, peel the protective backing off the blister part.
- Hold the blister part in one hand. Connect the syringe to the luer inlet of the Minisart® filter unit.
- If necessary, attach a sterile needle to the outlet of the Minisart® and remove the protective cap.
- Hold the syringe and Minisart® assembly with the needle pointing upwards. Slowly press in the plunger until a bead of liquid exits from the outlet or needle. This assures proper wetting of the membrane filter and prevents air bubbles.
- Filter the contents of the syringe.
- Upon completion of filtration, discard the Minisart® (if you are working with radioactive or infectious material, observe the regulations for disposal).
- Minisart®, non-sterile; can be sterilized with ethylene oxide gas (acc. to EN 550) or by gamma irradiation (acc. to EN 552); however, Minisart® with a PTFE membrane (16596-HYQ) may not be gamma irradiated.

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