CLEAN TEST

ITALIANO

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CLEAN TEST

Rapid test for the hygiene control of work surfaces

DESCRIPTION

CLEAN TEST is a swab for the detection of protein, fat and sugar residues suitable for the check of work surfaces cleaning level.

KIT CONTENT

The product consists of plastic tubes containing a swab, a vial (reagent A) and a disc (reagent B).

Each kit contains 50 CLEAN TEST swabs (REF 86001) or 3 CLEAN TEST swabs (REF 86001/3) and it is accompanied by 1 instruction sheet.

REQUIRED AIDS NOT INCLUDED

PHYSIOLOGICAL SOLUTION (REF 20142).

PRINCIPLE OF THE METHOD

The test highlights the presence of proteinic residues and of other reducing substances on the analyzed surface. The method, based on the reaction of bicinchoninic acid (reagent A) with cupric sulfate (reagent B) in alkaline conditions, produces the complexation of Cu ions with the peptidic bonds of proteins. Such complex takes on a purple colouring directly proportional to the concentration of proteins, fats and sugars which are present on the checked surface.

PROCEDURE

1. Take the swab contained in the tube and wet the tip with 3-4 drops of physiological solution (or sterile distilled water). Streak the swab on a small plot of surface measuring 10 x 10 cm horizontally and vertically.
2. Insert the swab into the tube, push the tube in order to lower the swab tip close to the disc.
3. Break the glass vial contained in the tube near the swab tip by pushing the tube gently, subsequently crush the vial.
4. Shake strongly the tube maintaining it in a vertical position with the cap turned upwards for 15-20 seconds.
5. Wait 8-10 minutes and compare the final color of the solution with the colourimetric pattern reported on the label of the tube.

INTERPRETATION OF RESULTS

A pale green, colourless or clear light grey colouring assures the proper cleaning of the analyzed surface. In case of a dark grey solution it is suggested to rinse the surface and/or repeat the test. However, if the solution turns to light or dark purple, a new cleaning cycle of the surface with suitable detergents is needed. False positive results may occur in presence of peroxides, ascorbic acid and of other strongly reducing agents. Otherwise, false negative ones may occur in presence of high concentrations of alkaline detergents.

QUALITY CONTROL FOR THE USER

1. Control of general characteristics, label and print
2. Chemical control

Product efficiency control through the following solutions :

COLORIMETRIC REACTION	STANDARD PROTEIN SAMPLES (BSA*)
Green	15 µg / 100 µL
	30 µg / 100 µL
Dark grey	60 µg / 100 µL
	80 µg / 100 µL
Light purple	120 µg / 100 µL
	300 µg / 100 µL
Dark purple	500 µg / 100 µL
	1000 µg / 100 µL

*BSA, Bovine Serum Albumin

WARNING AND PRECAUTIONS

The product not is classified as hazardous by current legislation; however for its use it is suggested to consult the safety data sheet.

The sampling can be performed on surfaces located in environments at temperature not lower than 0°C. The test is recommendable to be performed (break of vial, reaction and interpretation of results) at 10 to 25°C temperature range.

STORAGE

10-30°C, from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

PRESENTATION

Product	REF	
CLEAN TEST	86001	50 tests
CLEAN TEST	86001/3	3 tests

TABLE OF SYMBOLS

LOT Batch code	 Do not reuse	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult accompanying documents	

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