

FatExtractor E-500 Technical data sheet

The declaration of the total fat or crude fat is required for most foods and feeds. The FatExtractor E-500 is designed for solid-liquid extraction for quick and compliant fat determination. The design of the glass assemblies and the high-speed heaters combined with sophisticated process control allow for the fastest and most reproducible extraction processes with full compliance.

By simply changing the glass assembly, the FatExtractor E-500 complies with standard methods such as Soxhlet, Hot Extraction (HE) or Twisselmann (ECE). The glass assembly is interchangeable for a later conversion to another extraction method. Suitable for extractions with petroleum ether, diethyl <u>ether</u>, hexane or chloroform.



Description of function

The FatExtractor E-500 is designed to carry out the following solid-liquid extraction methods:

- Soxhlet Extraction (see Chapter "Soxhlet Extraction")
- Hot Extraction (see Chapter "Hot Extraction")
- Economic Continuous Extraction (see Chapter "Economic Continuous Extraction")

Soxhlet Extraction

Step 1 extraction

- The sample is located in the extraction chamber.
- The beaker contains the solvent.
- The solvent is heated, vapor rises up to the condenser, condenses and drops into the extraction chamber with the sample.
- The magnetic value is closed, the solvent is collected up to the optical sensor and extracts the analyte.
- When the optical sensor is reached, the magnetic valve opens and the solvent containing the analyte flows back into the beaker.

Step 2 rinsing

- The solvent is heated, vapor rises up to the condenser, condenses and drops into the extraction chamber with the sample.
- The magnetic valve is open, the solvent flows back into beaker, the solvent is not collected.
- This flushes traces of fat into the beaker.

Step 3 drying

- The solvent is heated, vapor rises up to the condenser, condenses and flows into tank.
- The analyte remains in the beaker.

Hot Extraction

Step 1 extraction

- The sample is located in the beaker.
- The beaker contains the solvent.
- The solvent is heated, vapor rises up to the condenser, condenses and drops into the beaker with the sample.

Step 2 rinsing

- The solvent in the beaker is heated up and evaporated.
- The vapor rises up to the condenser.
- The condensed solvent flows into the beaker with the sample.
- The tank bottle valve opens periodically and condensed solvent flows in the tank bottle.
- The solvent level reduce.

Step 3 drying

- The solvent in the beaker is heated up and evaporated.
- The vapor rises up to the condenser.
- The tank bottle valve is open and condensed solvent flows in the tank bottle.
- The analyte remains in the beaker.





Economic Continuous Extraction

Step 1 extraction

- The sample is located in the extraction chamber.
- The beaker contains the solvent.
- The solvent is heated, vapor rises up around the sample to the condenser, condenses and drops back into extraction chamber through the sample into the beaker.

Step 2 drying

- The solvent is heated, vapor rises up around the sample to the condenser, condenses and flows into tank.
- The analyte remains in the beaker.



Order code FatExtractor E-500

The FatExtractor E-500 can be ordered in the following configurations: Soxhlet, Hot Extraction (e.g. Randall) or Economic Continuous Extraction (e.g. Twisselmann). For large sample volumes, use the Soxhlet LSV version, suitable for extraction thimbles with a diameter of 43 mm.

The FatExtractor E-500 can be upgraded with the Analyte protection sensor for increased user safety and to prevent any loss of analyte during the extraction process. The Pro Interface consists of a 7 inch color touch display for even more intuitive use and overview on the interface.

Choose the configuration according to your needs:

1 1 E 5 0 0	0	0	0
Glass assembly			
S Soxhlet			
L Soxhlet LSV			
H Hot Extraction			
E Cont. Extraction ECE			
Analyte protection			
0 Without			
1 Analyte protection			
Interface			
r Pro			

S Standard

Order code FatExtractor E-500 System

The FatExtractor E-500 can be bundled with the HydrolEx H-506 for total fat determination and a Recirculating Chiller F-305/F-308.

HydrolEx H-506 performs convenient and safe acid hydrolysis as a step prior to extraction. It is an essential work step of the total fat determination where matrix structures enclosing the fat fraction of food and feed samples are broken up.

The Recirculating Chiller F-305/F-308 is an environmental friendly and water saving alternative to using tap water for cooling. The FatExtactor E-500 has to be operated either with a recirculating chiller or a tap water valve.

Choose the configuration according to your needs:

	1 1 S E 5 0 0		0	0	
			_	_	
Glass	assembly				
S So>	xhlet				
L So>	xhlet LSV				
H Hot	t Extraction				
E Cor	nt. Extraction ECE				
Analyt	te protection				
0 Wit	hout				
1 Ana	alyte protection				
Interfa	ace				
P Pro)				
S Sta	andard				
Coolir	ng				
0 Wit	hout				
1 Tap	o water valve				
2 Chi	iller F-305				
3 Chi	iller F-308				
Hydro					
0 Wit	hout				
1 H-5	506				
2 H-5	506 + waterjet pump				
Voltac					

1 220-240 V

2 100-120 V

Scope of delivery All configurations are supplied ready to use.

	Hot Extraction	Economic Continu- ous Extraction	Soxhlet Extraction	Soxhlet Extraction LSV
FatExtractor E-500	1	1	1	1
Condenser E-500	6	6	6	6
Extraction glass chamber ECE	-	6	-	-
Extraction glass chamber Soxhlet	-	-	6	-
Extraction glass chamber Soxhlet LSV	-	-	-	6
Soxhlet assembly cpl.	-	-	6	6
Extraction beaker	-	6	6	6
Extraction beaker HE	6	-	-	-
Sealing E-500	6	12	12	12
Set of glass sample tube holder	6	6	6	-
Set of glass sample tube holder LSV	-	-	-	6
Set of holders for thimbles 25 mm	6	6	6	-
Set of holders for thimbles 33 mm	6	6	6	6
Set of holders for thimbles 43 mm	-	-	-	6
Extraction thimbles 25 x 100 mm	6	6	6	-
Extraction thimbles 33 x 94 mm	6	6	6	6
Extraction thimbles 43 x 118 mm	-	-	-	6
Solvent tank cpl.	1	1	1	1
Cooling water hose 3 m	2	2	2	2
Beaker tong	1	1	1	1
Extraction beaker carrier	-	1	1	1
Extraction beaker carrier HE	1	-	-	-
Funnel	1	1	1	1
Power cord	1	1	1	1
Pliers for glass sam- ple tube with frit	1	1	1	1
Operation manual	1	1	1	1

Technical data

FatExtractor E-500

Power consumption	1300 W
Connection voltage	100 - 240 ± 10 % VAC
Fuse	10 A
Frequency	50 / 60 Hz
Overvoltage category	II
Pollution degree	2
Protection Class	1
Dimensions (W x D x H) (with glassware Soxhlet Extraction)	638 x 595 x 742 mm
Dimensions (W x D x H) Dimensions with glassware HE	638 x 595 x 613 mm
Dimensions (W x D x H) (with glassware Economic Continuous Extraction)	638 x 595 x 622 mm
Weight (without glassware)	41.8 kg
Weight (with glassware)	49.6 kg
Total heating power (rated)	600 W
Total heating power (maximum)	1200 W
Hose connection	6 / 9 mm
Allowed water pressure (nominal value)	6 bar
Allowed water pressure (maximum)	8 bar
Number of extraction positions	6
Solvent tank volume	2 L

Ambient conditions

For indoor use only.	
Max. altitude above sea level	2000 m
Ambient temperature	5–40 °C
Maximum relative humidity	80% for temperatures up to 31 °C decreasing linearly to 50 % relative humidity at 40 °C
Storage temperature	max. 45 °C

Conversion and upgrade kits

The FatExtractor E-500 can be converted into another configuration with an easy change of glass assembly. Choose the needed conversion kit according to the matrix.

from to	Soxhlet Extraction	Hot Extraction	Economic Continuous Extraction	Soxhlet Extraction LSV
Soxhlet Extraction		11068489	11068492	
Hot Extraction	11068487		11068487	11068487
Economic Continuous Extraction	11068488	11068491		11068488
Soxhlet Extraction LSV		11068490	11068493	

	Order no.
Conversion kit from Soxhlet / ECE to HE	11068487
incl.: 6 beakers HE (3 x 11067475) , beaker rack with beaker brackets, set of draining tubes for HE (11067480)	
Conversion kit from Soxhlet / Universal to ECE	11068488
Includes 6 Extraction glass chamber ECE (11062499), set of draining tubes for ECE (11067479)	
Conversion kit from HE to Soxhlet	11068489
Includes 6 Extraction glass chamber Soxhlet (11062496), 6 beaker (3 x 11067474), 6 Soxhlet assemblies cpl. (11067065) and a chamber rack (11067077), set of draining tubes for Soxhlet (11067478), two safety shields top, cpl. (2 x 11067832)	
Conversion kit from HE to ECE	11068491
Includes 6 Extraction glass chamber ECE (11062499), 6 beaker (3 x 11067474), set of draining tubes for ECE (11067479) and a chamber rack (11067077)	
Conversion kit from ECE to Soxhlet	11068492
Includes 6 Extraction glass chamber Soxhlet (11062496), 6 Soxhlet assemblies cpl. (11067065) and set of draining tubes for Soxhlet (11067478)	
Conversion kit HE to Soxhlet LSV	11068490
Includes 6 Extraction glass chamber Soxhlet LSV (11062497), 6 beaker (3 x 11067474), 6 Soxhlet assemblies cpl. (11067065) and a chamber rack (11067077), set of draining tubes for Soxhlet (11067478)	
Conversion kit ECE to Soxhlet LSV	11068493
Includes 6 Extraction glass chamber Soxhlet LSV (11062497), 6 Soxhlet assemblies cpl. (11067065) and a set of draining tubes for Soxhlet (11067478)	

	Order no.
Upgrade kit Analyte protection	11068524
For later upgrade of a FatExtractor E-500 with analyte protection sensor	
Upgrade kit Interface Pro	11068525
For later upgrade of a FatExtractor E-500 with a Pro interface (7 " touch screen)	

Chiller

No. Units	Ambient Temperature	Chiller
1	< 30 °C	Recirculating Chiller F-305
1	< 40 °C	Recirculating Chiller F-308
2	< 30 °C	Recirculating Chiller F-308
2	< 40 °C	Recirculating Chiller F-314

	Order no.
Recirculating Chiller F-305	11F30501
550 W at 15 °C, Display, 230 V	
Recirculating Chiller F-305	11F30502
550 W at 15 °C, Display, 115 V	
Recirculating Chiller F-308	11F30801
900 W at 15 °C, Display, 230 V	
Cooling capacity 900 W at 15 °C, for temperatures from -10 to 25 °C	
Recirculating Chiller F-308	11F30802
900 W at 15 °C, Display, 115 V	
Cooling capacity 900 W at 15 °C, for temperatures from -10 to 25 °C	
Recirculating Chiller F-314	11F31401
1400 W at 15 °C, Display, 230 V	
Cooling capacity 1400 W at 15 °C, for temperatures from -10 to 25 °C	
Recirculating Chiller F-314	11F31402
1400 W at 15 °C, Display, 115 V	
Cooling capacity 1400 W at 15 °C, for temperatures from -10 to 25 °C	

Spare parts

	Order no.	Image
Extraction glass chamber Soxhlet	11062496	
Extraction glass chamber Soxhlet LSV	11062497	

	Order no.	Image
Extraction glass chamber ECE	11062499	
Set of beakers, 2 pcs.	11067474	
Set of beakers HE (2pcs)	11067475	
Condenser E-500 cpl.	11067063	
Condenser tank bottle	11065966	
Tank adapter	11064590	
Tank bottle 2 L	11065983	
Ball joint clamp. For BJ35/20 To fasten receiving flask on condenser/secondary condenser.	003275	
Set of gliding elements including magnets, 10 pcs.	11067827	

	Order no.	Image
Soxhlet assembly cpl.	11067065	
One part constiting of magnetic valve and level sensor for extraction glass chamber Soxhlet		
Protection shield top, cpl.	11067832	
Protection shield bottom, cpl.	11067831	
Set of seals E-500, EKM, 6 pcs.	11069012	
Membrane with anchor for magnetic valve unit	037534	
Silicone hose D6/9 L=3 m	048355	
Set of draining tubes SOX, FEP	11067478	
The draining tubes connect the receiving funnel in the condensers with the tank valve to drain solvent into the tank.		
Set of draining tubes HE, FEP	11067480	
The draining tubes connect the receiving funnel in the condensers with the tank valve to drain solvent into the tank.		
Set of draining tubes ECE, FEP	11067479	

Accessories

	Order no.	Image
Holder for glass sample tubes, stainless steel	11067219	
Holder for glass sample tubes, PTFE	11067220	

	Order no.	Image
Holder for extraction thimbles (diameter 25 - 43 mm)	11068443	
Extraction beaker carrier	11067042	ð
Allows to carry 6 beakers (11067474)		
Extraction beaker HE carrier	11067493	
Allows to carry 6 beakers HE (11067475)		
Set condenser insulations E-500, 6 pcs.	11069078	
The insulation of the condensers prevent condensing water and is recom- mended in high humidity environment		
Set insulation cooling water hoses	11069079	
The insulation of the water hoses prevent condensing water and is recom- mended in high humidity environment.		
Support solvent supply	11068306	\frown
Allows to fix the tubes of solvent dispensers to the condensers for convenient solvent addition.		
Cooling water valve. 24VAC	031356	
Valve opens cooling water feed during distillation. Meant to be used with a vacuum controller/interface.		
Turning disk	11067985	
Allows for turning the instrument for easier access.		

Consumables

	Order no.
Quartz sand 0.3 - 0.9 mm, 2.5 kg	037689
Celite® 545, 1 kg	11068920

Glass sample tubes and extraction thimbles

	Order no.	Image
Glass sample tubes with frit, 6 pcs.	11067497	
Glass sample tubes LSV with frit, 6 pcs.	11067814	
Extraction thimbles 25 x 100 mm, 25 pcs.	018105	
Extraction thimbles 33 x 94 mm, 25 pcs.	11058983	
Extraction thimbles, Set. 25 pcs, 43 x 118 mm, cellulose For Soxhlet extraction unit.	018106	

Holder for extraction thimbles

	Order no.
Set of holders for glass sample tubes with frit, PTFE, 3 pcs.	11067485
Set of holders for LSV glass sample tubes, PTFE, 3 pcs.	11067486
Holders for thimbles d25, PTFE, 3 pcs.	11067488
Holders for thimbles d33, PTFE, 3 pcs.	11067490
Holders for thimbles d43, PTFE, 3 pcs.	11067491