
52 33

pH

CRISON

52 33

Puntos clave. *Punti chiave.* Points clé. *Key points.*

1 Tapón protector.
Tappo protettore.
Bouchon protecteur.
Protective cup.

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pH electrode 52 33

Introduction

The **52 33** is a puncture electrode with plastic body (POM). Suitable for measurements in semisolid samples.

Electrode preparation

Remove, with care, the storage protector containing electrolyte by unscrewing it (6).

Calibration

pH electrodes must be calibrated with buffer solutions before measuring.

⚠ See pH-meter instruction manual.

Measurements. Recommendations

- Do not introduce the electrode directly in thick samples or samples containing solid materials, bones, etc. It is necessary to open a path using a stainless steel perforator, introduce the electrode and take the measurement.
- During measurements, hold the measured product still.
- The sample must cover the electrode's diaphragms (4).
- After measurement, clean the electrode with cotton soaked in alcohol.
- Between measurements the electrode's membrane should be immersed in aqueous solution. The most appropriate place is the storage protector with some electrolyte.
- This electrode allows measure with the membrane in any position.

⚠ WARNING

- Do not leave the electrode immersed in distilled water.
- Do not bend the electrode. The internal body is made from glass.

Maintenance

- Maintain the membrane's surface (5) clean. Do not rub it. Avoid knocks and scratches.
- After measurement, clean the electrode with cotton soaked in alcohol.
- At the end, clean with lukewarm soapy water.

⚠ Never use cleaning solutions with $\text{pH} < 2$. It provokes alteration in the reference electrolyte and damages irreversibly the electrode.

Troubleshooting

... the pH-meter always measures the same pH value.

Short circuit in the electrode or connector.

Substitute the electrode.

... unstable reading.

Broken cable.

Substitute the cable.

... impossible to perform calibration.

... very slow response.

Dirt on the diaphragm or the membrane.

De-hydrated membrane.

Electrode ageing or wear.

Clean, regenerate or substitute the electrode.

Storage

The connector should be protected from humidity by the protective cap ①.

Place the electrode in the storage protector ⑥ with appropriate electrolyte (Crisolyt G, KCl + glycerine).

Life of an electrode

The average life expectancy of a pH electrode is approximately one year depending on the working conditions.

An electrode used for more than six months that does not respond appropriately after regeneration, must be replaced by a new one.

 *If an electrode has too short life span, less than three months, ask us about it.*

Specifications

Measuring range	2 ... 14 pH
Operating temp.	0 ... 60 °C
Reference	encapsulated Ag/AgCl
Diaphragms	1 ceramic + 1 open
Electrolyte	solid
Body material	glass / plastic
Min. puncture depth	12 mm

Applications:

- Cheeses, meats, fish, fruits, bread doughs, etc. It can be used for measurements in aqueous solutions too.

Limits

pH values < 2.

Temperatures > 60 °C.

Warranty

CRISON guarantees this electrode against manufacturing defects.

Validity: 6 months from shipment date.

Limitations:

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

QUALITY CERTIFICATE

The attached electrode 52 33 has passed the quality test and follows the next specifications:

- Asymmetry potential < ± 15 mV.
- Sensitivity, pH 4...7 (at 25 °C) > 98%.
- Response time, pH 4...7 < 20 s.

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