## Garantie

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## Warranty


Laboratoires Humeau guarantees this electrode against manufacturing defects.

**Validity:** 6 months from shipment date.

**Limitations:**

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

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## pH Electrode

### Introduction

The 2 Pore F is a penetration electrode with plastic body. It can be used in any position.

### Description

- 1-Plastic body  
Do not buckle it, inside is glass made.
- 2-Diaphragms, 2 open.
- 3- Sensitive membrane. Protect from friction or knocking
- 4-Storage protector with some drops of electrolyte (KCL+glycerine)

### Specifications

pH measuring range	0...14
Operating temp.	0...60°C
Reference	Encapsulated Ag/AgCl
Diaphragms	2 open
Electrolyte	Solid
Body material	Glass / Plastic
Min. puncture depth	15 mm
Connection	S7

### Applications

Cheeses, meats, fish, fruits, bread doughs, etc.

### Limits

pH values < 0 - Temperatures > 60 °C

### Warning

- **Do not introduce the electrode directly** in the thick samples or samples containing solid materials, bones etc. It is necessary to open a path using a stainless steel perforator.
- **Cleaning:** Between samples clean the electrode with lukewarm water. At the end, clean with lukewarm soapy water. Never use cleaning solutions pH < 2
- **Store** the electrode in the protector with KCL + Glycerine.

### Certificate of quality

The attached electrode 2 Pore F has passed the quality test and follows the next specifications:

- Asymmetry potential < ± 20mV
- Sensitivity pH 4...7 (at 25°C) > 98%
- Response time, pH 4...7 < 20 s