## Oxidase Test Stick

Rapid test for detection of cytochrome oxidase enzymatic activity.

### DESCRIPTION

Oxidase Test Stick is a diagnostic test used for differentiation and microbial identification, particularly of Gram-negative bacteria, on the basis of the presence of enzyme cytochrome oxidase.

The product matches with recommendations of EN ISO 16266 and ISO 9308-1 for detection of *Pseudomonas aeruginosa* and for confirmation of *Escherichia coli* and coliform bacteria, respectively.

### CONTENTS OF THE PACKAGES

Each package contains 50 Oxidase Test Stick.

### METHOD PRINCIPLE

Oxidase-positive bacteria produces the enzyme cytochrome oxidase (indophenol oxidase) that catalyzes the transport of electrons from donor compounds (NADH) to electron acceptors (usually oxygen).

Tetramethyl-p-phenylenediamine dihydrochloride contained in Oxidase Test Stick acts as an artificial electron donor and is oxidized by oxidase-positive bacteria forming the coloured compound indophenol blue.

### COMPOSITION

Oxidase Test Stick is made of a special paper with a zone impregnated with a solution of N,N,N',N'-tetramethyl-p-phenylenediamine dihydrochloride.

### TEST PROCEDURE

1. Allow container to come to room temperature before opening, for minimizing condensation on the stick.
2. Pick up one or more than one well isolated colony and wipe off by the zone of the stick indicated by arrows. Alternatively, transfer a drop of the suspension of the test organism to the reaction zone of the stick or dip directly the point of the stick into the microbial suspension.
3. Examine for an immediate color change (within 60 seconds) at the position of the inoculated area (NB. The usage of very dilute microbial suspensions may result in longer reactions time)

### INTERPRETING RESULTS

The development of a blue-purple color indicates a positive reaction. No color change corresponds to a negative test, i.e. the organism under investigation does not produce the enzyme cytochrome oxidase.

### LIMITATIONS

The most suitable cultures for the oxidase test are those from culture media without dyes, indicators or inhibitors. Bacterial colonies taken from media with pH values below 5.5 (e.g. after the metabolism of carbohydrates with subsequent acidification of the culture medium) can give a false negative oxidase reaction. Colonies taken from media containing nitrate may give unreliable results. Do not use steel, nichrome or iron containing loops to pick the colony. A platinum or plastic loop, or wooden applicator stick is recommended.

### STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

### SHELF LIFE

2 years.

### QUALITY CONTROL

Control strains are indicated in the QC table.

#### QC Table.

Microorganism	WDCM	Oxidase reaction
<i>Escherichia coli</i>	WDCM 00013	Negative, no color change
<i>Pseudomonas aeruginosa</i>	WDCM 00025	Positive, deep blue-purple coloration

### WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use and must be used only by properly trained operators.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

### BIBLIOGRAPHY

- ISO 9308-1:2014. Water quality – Enumeration of *Escherichia coli* and coliform bacteria – Part 1: Membrane filtration method for waters with low bacterial background flora.
- EN ISO 16266:2008. Water quality – Detection and Enumeration of *Pseudomonas aeruginosa* – Method by membrane filtration (ISO 16266:2006).
- Steel K. J. (1962) J. Appl. Bact. 25:445-447.

### PRESENTATION

	Contents	Ref.
Oxidase Test Stick	50 sticks	88029

### TABLE OF SYMBOLS

<b>LOT</b> Batch code	<b>IVD</b> <i>In vitro</i> Diagnostic Medical Device	Manufacturer	Use by	Fragile, handle with care
<b>REF</b> Catalogue number	Temperature limitation	Contains sufficient for <n> tests	Caution, consult Instruction For Use	Do not reuse

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