UREMILKIT

Reagent strips - Urea in Whole Milk

I. Description

The dry chemical pathway utilized in the UREMILKIT Reagent Strips has been described by Ormerod, T.P. in 1968. This procedure is a quick, convenient way to semi-quantitatively approximate Milk Urea levels in whole milk.

Care must be taken to follow the directions in this insert to obtain accurate results. Any discrepancies between normal test results should be confirmed by more accurate laboratory methods. UREMILKIT Reagent Strips do not replace laboratory methods.

II. Reagent composition

0.10% w/w ammonium citrate, 0,15% w/w enhancer, 0.20% bromothymolblue, 0.75kU/ml urease.

III. Warnings and precautions

1. For *in vitro* diagnostic use only.
2. Do not ingest Urease and BTB Indicator.
3. Reagents may be irritating to the skin. Avoid contact.
4. Raw milk specimens should be considered infectious and handled appropriately.

IV. Storage and handling

NOTE: All warnings, precautions, and proper storage procedures must be followed to avoid deterioration of UREMILKIT Reagent strips.

1. Store bottle at temperatures between 15°-30°C and out of direct sunlight.
2. All unused strips must store in original bottle. Do not transfer to any other container.
3. Do not remove the white packet (desiccants) from bottle. Replace cap quickly and tightly after use.
4. Do not touch test areas of the strip.
5. Protect against exposure to light, moisture, and heat which will affect reagent activity.

V. Specimen collection

The milk required for the test can be taken at any time during milking. The sample should be tested immediately after collection; if not, it should be kept in a refrigerator and warmed to room temperature prior to testing. If repeated tests are conducted in a herd, sampling should always take place at the same milking.

VI. Directions

Remove from the bottle only enough strips for immediate use and replace cap tightly.

Apply a large drop of fresh raw milk to completely cover reagent pad. Make sure the milk does not touch the edge of the pad.

Wait exactly 90 seconds using a stop-watch to keep time.

After 90 seconds, tap the side of the strip firmly on a paper towel to remove excess milk from the reagent pad.

Immediately compare pad to color chart and record result. Care must be taken to read results before colors fade.

VII. Results

Results are given directly in urea (both mg/dL and g/L) on color chart. To convert to urea nitrogen divide urea value by 2.14. Normal values for urea vary between species and should be verified by laboratory determined normal range values.

However, it was found that urea results from dairies and laboratories are measured by the urease method of 340nm at 37°C.

It is therefore advisable to compare results, to put itself in the same conditions of temperature, temperature differences varying rates more or less important.

Consequently, you will find inside the box a specific colored chart that has been established from 191 cow specimen measured at 37°C.

It is therefore preferable to use this card as reference instead of the chart displayed on the test strip vial.

VIII. Limitations

A large drop of fresh raw milk is necessary to completely cover reagent area. It is not desirable to use the milk sample mixed with preservatives (sodium azide, boric acid, etc.).

Any fat adhering to the reagent part can interfere with the color reaction. If the sample being tested is colostrum, the reading of the color reaction is more difficult because of the milk's intrinsic color. The test should be carried out as soon as possible after removal of the test strip from the container.

The color reaction should be read under sufficient light. If different shades of color appear at the margin and in the center of the reagent part, the reading is based on the central color.

IX. Materials provided

1. One bottle containing 25 strips of UREMILKIT Reagent strips.
2. A card with the color chart.

X. Other required materials

1. Stop-watch for precise timing.
2. Paper towel to wipe off milk from strip.

XI. References

1. Ormerod, T.P.: A rapid bedside method for measuring blood urea; *Practitioner* 201: 921, 1968. Jonker, J.S., R.A. Kohn and J. High. 2002. Use of milk urea nitrogen to improve cow diets. *J. Dairy Sci.* 85:939-946.
2. Kauffman, A.J. and N.R. St-Pierre. 2001. The relationship of milk urea nitrogen to urine nitrogen excretion in Holstein and Jersey cows. *J. Dairy Sci.* 84:2284-2294.
3. Nousiainen, J., K.J. Shingfield and P. Huhtanen. 2004. Evaluation of milk urea nitrogen as a diagnostic of protein feeding. *J. Dairy Sci.* 87:386-398.

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Z.A de Gesvrine - 4 rue Képler - B.P.4125
44241 La Chapelle-sur-Erdre Cedex - France
t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00
commercial@humeau.com

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1. Ormerod, T.P.: A rapid bedside method for measuring blood urea; Practitioner 201: 921, 1968. Jonker, J.S., R.A. Kohn and J. High. 2002. Use of milk urea nitrogen to improve cow diets. J. Dairy Sci. 85:939-946.
2. Kauffman, A.J. and N.R. St-Pierre. 2001. The relationship of milk urea nitrogen to urine nitrogen excretion in Holstein and Jersey cows. J. Dairy Sci. 84:2284-2294.
3. Nousiainen, J., K.J. Shingfield and P. Huhtanen. 2004. Evaluation of milk urea nitrogen as a diagnostic of protein feeding. J. Dairy Sci. 87:386-398.