Sugars Fermentation

Test for the sugars fermentation

Instructions For Use

Contents	Page
Italiano	1
English	3
Français	5
Deutsch	7
Español	9
Português	11
Ελληνικά	13
Products List	15

Rev.9 / 04.05.2015

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BIBLIOGRAFIA

- Bayley and Scott's. *Diagnostic Microbiology*. 7th Edition Mosby Company (1986).
- E.H. Lennette. *Manual of Clinical Microbiology*. 4th Edition A.S.M. (1985)

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Test for the sugars fermentation.

DESCRIPTION

Sugars Fermentation are tests to detect the bacterial sugars fermentation. The tubes contain a desiccated protein mixture and the sugar substratum of the fermentation. Such reaction identifies the bacterial strains.

CONTENTS OF THE PACKAGES

Each package contains 30 tubes.

The package of the product Xylose Test, ref. 88220, contains 30 tubes and 30 Xylose Disc.

PRINCIPLE OF THE METHOD

The carbohydrate fermentations are very important tests for the identification of several bacterial strains of particular clinical importance, both Gram-positive and Gram-negative.

The sugar fermentation is indicated by the color change of the pH indicator in the medium from purple to yellow.

COMPOSITION

Every tube contains desiccated protein substrata, the pH indicator bromocresol purple and the sugar specific for each test.

The tube of the product Xylose Test contains desiccated protein substrata and the pH bromocresol purple.

GATHERING SAMPLES

The colonies to be subjected to the fermentation test must be taken from a culture medium used to isolate the bacteria from the sample under examination.

TEST PROCEDURE

1. Take the container from the refrigerator and leave it for a few minutes on the bench until it has reached room temperature.
2. Take one tube and add 0.3 mL of physiological solution for microbiological use.
In case of Xylose Test add one disc of Xylose Disc in the tube rehydrated with physiological solution.
3. Using a sterile loop, pick up a well isolated bacterial colony from the isolation medium and suspend it into the culture medium of the tube.
4. Take another tube and add 0.3 mL of physiological solution.
In case of Xylose Test add one disc of Xylose Disc in the second tube rehydrated with physiological solution.
Do not inoculate this tube and use it as negative control.
5. Cap the tubes and incubate at $35 \pm 2^\circ\text{C}$ for 18-24 hours. Incubate Glucose Test for 16-18 hours.

INTERPRETATION OF RESULTS

Observe the color of the medium and interpret the results as indicated:

Purple color		Negative	Yellow color		Positive
--------------	---	----------	--------------	---	----------

QUALITY CONTROL

Every batch of Sugars Fermentation is submitted to the quality control by using one bacterial culture as positive control and one as negative control, as indicated in the QC table.

QC Table.

Test	Positive control	Negative control
Adonitol Test	<i>Klebsiella pneumoniae</i> ATCC® 13883	<i>Escherichia coli</i> ATCC® 25922
Arabinose Test	<i>Escherichia coli</i> ATCC® 25922	<i>Enterococcus faecalis</i> ATCC® 19433
Arabitol Test	<i>Listeria monocytogenes</i> ATCC® 35152	<i>Proteus mirabilis</i> ATCC® 25933
Dulcitol Test	<i>Klebsiella pneumoniae</i> ATCC® 13883	<i>Escherichia coli</i> ATCC® 25922
Galactose Test	<i>Escherichia coli</i> ATCC® 25922	<i>Pseudomonas aeruginosa</i> ATCC® 27853
Glucose Test	<i>Escherichia coli</i> ATCC® 25922	<i>Plesiomonas shigelloides</i> ATCC® 14029
Inositol Test	<i>Klebsiella pneumoniae</i> ATCC® 13883	<i>Escherichia coli</i> ATCC® 25922
Inulin Test	<i>Streptococcus pneumoniae</i> ATCC® 6303	<i>Streptococcus pyogenes</i> ATCC® 19615
Lactose Test	<i>Escherichia coli</i> ATCC® 25922	<i>Salmonella typhimurium</i> ATCC® 14028
Levulose Test	<i>Streptococcus pneumoniae</i> ATCC® 6303	<i>Neisseria meningitidis</i> ATCC® 13090
Maltose Test	<i>Escherichia coli</i> ATCC® 25922	<i>Pseudomonas aeruginosa</i> ATCC® 27853
Mannitol Test	<i>Staphylococcus aureus</i> ATCC® 25923	<i>Pseudomonas aeruginosa</i> ATCC® 27853
Mannose Test	<i>Salmonella typhimurium</i> ATCC® 14028	<i>Proteus mirabilis</i> ATCC® 25933
Raffinose Test	<i>Klebsiella pneumoniae</i> ATCC® 13883	<i>Proteus mirabilis</i> ATCC® 25933
Rhamnose Test	<i>Salmonella typhimurium</i> ATCC® 14028	<i>Proteus mirabilis</i> ATCC® 25933
Salicin Test	<i>Klebsiella pneumoniae</i> ATCC® 13883	<i>Escherichia coli</i> ATCC® 25922
Sorbitol Test	<i>Escherichia coli</i> ATCC® 25922	<i>Streptococcus pyogenes</i> ATCC® 19615
Sucrose Test	<i>Klebsiella pneumoniae</i> ATCC® 13883	<i>Salmonella typhimurium</i> ATCC® 14028
Trehalose Test	<i>Escherichia coli</i> ATCC® 25922	<i>Pseudomonas aeruginosa</i> ATCC® 27853
Xylose Test	<i>Salmonella typhimurium</i> ATCC® 14028	<i>Staphylococcus aureus</i> ATCC® 25923

PRECAUTIONS

The product does not contain substances classified as hazardous under current legislation. Sugars Fermentation is for *in vitro* diagnostic use only. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8°C in its original package. In such conditions the product is valid until the expiry date shown on the label. Do not use them beyond that date. Dispose of them if they show signs of deterioration.

ELIMINATING USED MATERIAL

After use, the tubes and the material that has come into contact with the sample must be decontaminated and disposed of in accordance with the laboratory procedures for the decontamination and disposal of potentially infected material.

BIBLIOGRAPHY

- Bayley and scott's. *Diagnostic Microbiology*. 7th Edition Mosby Company (1986).
- E.H. Lennette. *Manual of Clinical Microbiology*. 4th Edition A.S.M. (1985)

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In Vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Consult instructions for use	 Do not reuse

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Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.
© Copyright LIOFILCHEM 2015