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Rev.3 /18.12.2012

## Rapid test to evidence arginine decarboxylation

### DESCRIPTION

**ARGININE DECARBOXYLASE TEST** is a rapid test to evidence arginine decarboxylation to be used for enterobacteria identification.

### CONTENT OF THE PACKAGES

Each package contains 30 tubes with desiccated medium and one instruction sheet.

### ITEMS NECESSARY BUT NOT INCLUDED IN THE PACKAGES

- PHYSIOLOGICAL SOLUTION (ref. 20095).
- VASELIN OIL (ref. 80278)
- Sundry microbiology laboratory materials.

### PRINCIPLE OF THE METHOD

The arginine decarboxylase is an enzyme that attacks the carboxylic group of the amino acid arginine, with formation of amine putrescine. The sugar contained in the medium is fermented by all the enterobacteria with consequent initial colour change of indicator system from purple to yellow. The acid medium supports the amino acid decarboxylation reactions with formation of putrescine. The amines production alkalize the medium and induces a new indicator colour change from yellow to purple. The negative reaction is shown by the appearance of a yellow colour in the tube, while the positive reaction is shown by a purple colour.

### COMPOSITION

Each tube contains desiccated proteic substrata, a sugar, the pH indicator bromocresol purple and the amino acid arginine.

### TEST PROCEDURE

- Take the package from the fridge, take one or more tubes and leave it for a few minutes on the bench until it has reached room temperature.
- Add 0.3 mL of PHYSIOLOGICAL SOLUTION (ref. 20095) to the tubes.
- With a sterile loop suspend in the culture medium of the tube a well isolated bacterial colony of 24 hours incubation chosen from the culture medium.
- Add 2-3 drops of VASELIN OIL (ref. 80278)
- Put the top to the tubes and incubate them at  $36 \pm 1$  °C for 4 hours till a maximum of 24 hours.

### INTERPRETATION OF RESULTS

Interpret the results according to the table 1:

Table 1

Colour of the medium	ARGININE DECARBOXYLASE TEST
Purple	Positive
Yellow	Negative

### QUALITY CONTROL FOR THE USER

1. Appearance: purple colour desiccated medium adhering to the tube base.
2. Microbiological control.

Every batch of **ARGININE DECARBOXYLASE TEST** is subjected to the quality control, using one bacterial culture of *Enterobacter cloacae* ATCC 13047 as positive control and one of *Proteus mirabilis* ATCC 25933 as negative control.

### PRECAUTIONS

The product **ARGININE DECARBOXYLASE TEST** is classified as dangerous, it is recommended to consult the safety data sheet for its correct use. **ARGININE DECARBOXYLASE TEST** is a disposable device to be used only for diagnostic use *in vitro*. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

### DISPOSAL OF USED MATERIAL

After the use, **ARGININE DECARBOXYLASE TEST** and the material that has come into contact with the sample must be decontaminated and disposed of in accordance with the laboratory procedures for the decontamination and disposal of potentially infected material.

### STORAGE

Store **ARGININE DECARBOXYLASE TEST** at 2-8 °C in its original package. In such conditions the product is valid until the expiry date shown on the label. Do not use them beyond that date. Dispose of them if they show signs of deterioration.

### REFERENCES

1. Edwin H.Lenette: *Manual of Clinical Microbiology* (1995).
2. BLAZEVIC, D.J., and EDERER, G.M.: *Principles of biochemical tests in diagnostic microbiology*. 63-67. New York, John Wiley & Sons, 1975.

### PRESENTATION

Product	REF	$\Sigma$
<b>ARGININE DECARBOXYLASE TEST</b>	88016	30 tubes

### TABLE OF SYMBOLS

<b>IVD</b> In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
<b>REF</b> Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	<b>LOT</b> Batch code

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**LIOFILCHEM® S.r.l.**