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BIBLIOGRAFIA

1. Shewan J.M. and Hodgkiss W. (1954) Nature 63:208-209.
2. Health Protection Agency (HPA), Microbiology Services Division, Standard Unit (2011) UK Standards for Microbiology Investigations, Identification of *Vibrio* species.

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Rev. 0 / 11.10.2013

For the differentiation of *Vibrio* species.

DESCRIPTION

O129 DISC is a test used for differentiation of *Vibrio* spp based on sensitivity to vibrio-static agent O129 (2,4-diamino-6,7-di-isopropyl-pteridine phosphate).

CONTENT OF THE PACKAGES

Each package contains 1 cartridge with 50 discs, packaged in a blister with a desiccant, and an instruction sheet.

PRINCIPLE OF THE METHOD

This test is useful in the differentiation of *Vibrio* from other gram-negative bacteria especially *Aeromonas*, which are characteristically resistant to O129. Even among the genus *Vibrio*, different species show different sensitivity to O129. Hence two different concentration discs can be simultaneously tested to determine the degree of sensitivity of the species.

COMPOSITION

O129 DISC is available at two concentrations: 10 µg and 150 µg.

TEST PROCEDURE

- Inoculate the plate* by streaking the sample on the agar surface to produce even growth.
- Aseptically place one 150 µg disc only onto the agar surface to differentiate between *Vibrio* spp and other gram-negative bacteria. Otherwise use both the 10 µg disc and the 150 µg disc to also differentiate within *Vibrio* species. Ensure suitable distance between the discs to avoid zones of inhibition overlapping.
- Incubate aerobically at 35°C for 24 hours.
- Observe for zones of inhibition.

* It is recommended to use non-selective salt containing media (0.5%) such as NUTRIENT AGAR, ref. 10044 or COLUMBIA AGAR (Horse Blood 5%), ref. 10025.

INTERPRETATION OF THE RESULTS

Most *Vibrio* species are sensitive with 150 µg but species differ with 10 µg discs (some strains of *V. cholerae* O1 and O139 may be resistant to both disc contents). Test organism could result sensitive (S), partially sensitive (PS) or resistant (R), to O129 as indicated in table 1.

Table 1.

Organism	Zone of inhibition around the disc		Result
	O129 10 µg	O129 150 µg	
<i>Vibrio cholerae</i>	Yes	Yes	S
<i>Vibrio vulnificus</i>	Yes	Yes	S
<i>Vibrio</i> spp	Yes	Yes	S
<i>Vibrio parahaemolyticus</i>	No	Yes	PS
<i>Vibrio alginolyticus</i>	No	Yes	PS
<i>Vibrio fluvialis</i>	No	Yes	PS
<i>Vibrio furnissii</i>	No	Yes	PS
<i>Vibrio cholerae</i> O1	No	No	R
<i>Vibrio cholerae</i> O139	No	No	R
<i>Aeromonas</i> spp	No	No	R

Sensitive (S): zone of inhibition around both O129 discs.

Partially Sensitive (PS): zone of inhibition around the 150 µg disc and no zone around 10 µg disc.

Resistant (R): no zone of inhibition around both O129 discs.

QUALITY CONTROL FOR THE USER

Appearance: white filter paper discs of 6 mm diameter.

Microbiological control: cultural response observed on inoculated NUTRIENT AGAR plates, with placed O129 DISC of 10 µg and 150 µg, after incubation at 35°C for 24 hours.

Microorganism

Expected Result

<i>Vibrio parahaemolyticus</i>	ATCC® 17802	Partially Sensitive
<i>Aeromonas hydrophyla</i>	ATCC® 7966	Resistant

PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. Each disc is a disposable device intended for *in vitro* diagnostic use only and must be used in a professional environment by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8°C in its original package. In such conditions the product is valid until the expiry date shown on the label. Do not use beyond that date. Eliminate if signs of deterioration are evident.

DISPOSAL OF USED MATERIAL

After the use, O129 DISC and the material that has come into contact with the sample must be decontaminated and disposed of in accordance with the laboratory procedures for the decontamination and disposal of potentially infected material.


REFERENCES

- Shewan J.M. and Hodgkiss W. (1954) Nature 63:208-209.
- Health Protection Agency (HPA), Microbiology Services Division, Standard Unit (2011) UK Standards for Microbiology Investigations, Identification of *Vibrio* species.

PRESENTATION

Product	Concentration (µg)	Ref.	Tests
O129 DISC	150	88043	50 discs
O129 DISC	10	88044	50 discs

TABLE OF SYMBOLS

LOT Batch code	IVD In Vitro Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult accompanying documents	 Keep away from heat sources

LIOFILCHEM® s.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

Rev. 0 / 11.10.2013