

# Bacto™ Casitone • Trypticase™ Peptone Bacto™ Tryptone • BiTek™ Tryptone

## Intended Use

Bacto Casitone, Trypticase Peptone, Bacto Tryptone and BiTek Tryptone are used in preparing microbiological culture media.

Ingredients, where noted, meet United States Pharmacopoeia (USP) performance specifications.

## Summary and Explanation

The manufacturing process for an enzymatic digest of casein is not as destructive as an acid hydrolysis. Thus, the casein is not broken down as completely into its constituent components. In

of Analysis of AOAC International and meets specifications in the USP for pancreatic digest of casein.<sup>1,2</sup>

Bacto Tryptone was developed by Difco Laboratories while investigating a peptone particularly suitable for the elaboration of indole by bacteria. It is also notable for the absence of detectable levels of carbohydrates. Bacto Tryptone has been used in conjunction with caseinase acids in nutritional studies to determine amino acids vs. peptide utilization.<sup>3,4</sup> It is included in standard methods applications and is listed in the reagent section of the USP as meeting the specifications for pancreatic digest of casein, a component in many of the media listed.<sup>5,6,7,8</sup> The European Pharmacopoeia also lists pancreatic digest of casein as a component in many of the recommended media.<sup>9</sup> Bacto Tryptone also works well in fermentation

many cases this makes for a more nutritious hydrolysate, especially for those organisms that prefer peptides to amino acids.

Bacto Casitone can be used as a component in microbiological culture media or in fermentation applications. A recent publication has also reported that the stability of lyophilized influenza virus vaccine was augmented by the addition of 2% Casitone.<sup>7</sup>

Trypticase Peptone is the primary nitrogen source in Trypticase Soy Broth and Agar. This product is recommended for use in media formulations, where good growth of fungi and bacteria is required. Trypticase Peptone is referenced in Official Methods

applications. It has been used successfully with commonly used organisms, such as *Escherichia coli*,<sup>11</sup> as well as uncommon organisms, such as the diatom *Nitzschia laevis*.<sup>12</sup>

BiTek Tryptone is prepared similarly to Bacto Tryptone but the final product goes through fewer refinement steps during processing. This product provides some of the same benefits as Bacto Tryptone in instances where a less refined hydrolysate can be utilized.

## Principles of the Procedure

Bacto Casitone, Trypticase Peptone, Bacto Tryptone and BiTek Tryptone are pancreatic digests of casein. Casein is the main milk protein and a rich source of amino acid nitrogen.

## Typical Analysis

Refer to Product Tables in the Reference Guide section of this manual.

## Directions for Preparation from Dehydrated Product

Refer to the final concentration of Bacto Casitone, Trypticase Peptone, Bacto Tryptone and BiTek Tryptone in the formula of the medium being prepared. Add appropriate product as required.

## Procedure

See appropriate references for specific procedures using Bacto Casitone, Trypticase Peptone, Bacto Tryptone and BiTek Tryptone.

## Expected Results

Refer to appropriate references and procedures for results.

## References

1. Casitone, Trypticase and Tryptone. (2011). In: *Handbook of yeast*.
2. Hennessey, J.P. (2007). Official method of analysis of milk. *International milk code, volume 10: milk processing*, 262-263.
3. Hennessey, J.P. (2007). Official method of analysis of milk. *International milk code, volume 10: milk processing*, 262-263.
4. *Handbook of yeast*, 262-263.
5. *Handbook of yeast*, 262-263.
6. *Handbook of yeast*, 262-263.
7. *Handbook of yeast*, 262-263.
8. *Handbook of yeast*, 262-263.
9. *Handbook of yeast*, 262-263.
10. *Handbook of yeast*, 262-263.
11. *Handbook of yeast*, 262-263.
12. *Handbook of yeast*, 262-263.

## User Quality Control

NOTE: Differences in the identity specifications and Cultural Response testing for media affected with both Bacto™ and BiTek™ brands may reflect differences in the development and testing of media for industrial and clinical applications, per the referenced publications.

## Identity Specifications

### Bacto™ Casitone

Dehydrated Appearance: Fine, free-flowing granules.  
 Solution: 1.0%, 2.0% and 10.0% solutions, soluble in purified water. 1.0% solution is light amber, clear. 2.0% solution is light to medium amber, clear, may have a slight precipitate. 10.0% solution is medium to dark amber, clear to very slightly opalescent, may have a precipitate.

Reaction of 1.0% Solution at 25°C: pH 6.8-7.4

### Bacto™ Tryptone

Dehydrated Appearance: Light beige, free-flowing, homogeneous.  
 Solution: 1.0%, 2.0% and 10.0% solutions, soluble in purified water. 1.0% solution is very light to light amber, clear. 2.0% solution is light to medium amber, clear. 10.0% solution is medium to dark amber, clear to slightly opalescent, may have a slight precipitate.

Reaction of 2.0% Solution at 25°C: pH 6.5-7.5

### BiTek™ Tryptone

Dehydrated Appearance: Light beige, free-flowing, homogeneous.  
 Solution: 1.0%, 2.0% and 10.0% solutions, soluble in purified water. 1.0% solution is very light to light amber, clear. 2.0% solution is light to medium amber, clear. 10.0% solution is medium to dark amber, clear to slightly opalescent, may have a slight precipitate.

Reaction of 2.0% Solution at 25°C: pH 7.2 ± 0.2

### BiTek™ Trypticase™ Peptone

Dehydrated Appearance: Fine, homogeneous, free of extraneous material.  
 Solution: 2.0% solution, soluble in purified water. Solution is clear to slightly hazy.

Reaction of 2.0% Solution at 25°C: pH 6.5-7.5

Continued

## Cultural Response

### Biochemical Reactions

#### Bacto™ Casitone, Bacto™ Tryptone or BiTek™ Tryptone

Prepare a sterile solution as directed below. Adjust final pH to 7.2-7.4. Inoculate and incubate at 35 ± 2°C for 18-48 hours.

TEST	TEST SOLUTION	ORGANISM	ATCC <sup>®</sup>	INCUBATION (H)	RESULT
Fermentable Carbohydrates	1%	<i>Escherichia coli</i>	25922	-10 <sup>2</sup>	Negative
Indole Production	0.1%	<i>Escherichia coli</i>	25922	0.1 mL, unfiltered	Positive
Acetyl methylglutamate Production	0.1% with 0.5% dextrose	Enterobacter aerogenes	13688	0.1 mL, unfiltered	Positive
Hydrogen Sulfide Production	1%	<i>Salmonella enterica</i> subsp. <i>enterica</i> serotype <i>Schwarzinger</i>	14028	0.1 mL, unfiltered	Positive

### BiTek™ Trypticase™ Peptone

Prepare a sterile solution as directed below. Adjust final pH to 7.2-7.4. Inoculate and incubate at 35 ± 2°C for 18-48 hours.

TEST	TEST SOLUTION	ORGANISM	ATCC <sup>®</sup>	INCUBATION (H)	RESULT
Fermentable Carbohydrates	1%	<i>Escherichia coli</i>	25922	-10 <sup>2</sup>	Negative
Indole Production	0.1%	<i>Escherichia coli</i>	25922	0.1 mL, unfiltered	Positive
Acetyl methylglutamate Production	0.1% with 0.5% dextrose	Enterobacter aerogenes	13688	0.1 mL, unfiltered	Positive
Hydrogen Sulfide Production	1%	Clostridium freundii	8054	0.1 mL, unfiltered	Positive

## Growth Response

### Bacto™ Casitone, Bacto™ Tryptone or BiTek™ Tryptone

Prepare a sterile solution with 2.0% Bacto Casitone, Bacto Tryptone or BiTek Tryptone, 0.5% sodium chloride and 1.0% agar. Adjust final pH to 7.2-7.4. Inoculate and incubate plates at 35 ± 2°C for 18-48 hours.

ORGANISM	ATCC <sup>®</sup>	INCUBATION (H)	RECOVERY
<i>Escherichia coli</i>	25922	30-300	Good
<i>Staphylococcus aureus</i>	25923	30-300	Good

### BiTek™ Trypticase™ Peptone

1. Prepare a sterile solution of peptone agar without blend and with 1% sheep blood (DB) using 10 g Trypticase Peptone, 2.5 g sodium chloride and 6.5 g agar in 500 mL of purified water. Adjust final pH to 7.2-7.4. Inoculate and incubate plates at 35 ± 2°C for 3 days (incubate supplemented with CO<sub>2</sub>).

ORGANISM	ATCC <sup>®</sup>	INCUBATION (H)	RECOVERY PLAIN	RECOVERY WITH SB
Enterobacter aerogenes	13688	10 <sup>2</sup> -10 <sup>7</sup>	Good	N/A
Escherichia coli	25922	10 <sup>2</sup> -10 <sup>7</sup>	Good	N/A
Staphylococcus aureus	6539 <sup>†</sup>	10 <sup>2</sup> -10 <sup>7</sup>	Good	N/A
Staphylococcus epidermidis	12228	10 <sup>2</sup> -10 <sup>7</sup>	Good	N/A
Staphylococcus saprophyticus	12226	10 <sup>2</sup> -10 <sup>7</sup>	N/A	Good, beta hemolytic
Staphylococcus pneumoniae	6305	10 <sup>2</sup> -10 <sup>7</sup>	N/A	Good, alpha hemolytic
Staphylococcus pyogenes	49117	10 <sup>2</sup> -10 <sup>7</sup>	Good	Good, beta hemolytic

2. Prepare a sterile solution of chocolate peptone agar using Trypticase Peptone. Adjust final pH to 7.2-7.4. Inoculate and incubate plates at 35 ± 2°C for 3 days with CO<sub>2</sub>.

ORGANISM	ATCC <sup>®</sup>	INCUBATION (H)	RECOVERY
<i>Neisseria gonorrhoeae</i>	16624	10 <sup>2</sup> -10 <sup>7</sup>	Good

## Availability

### Bacto™ Casitone

100g 500g 1kg 2.5kg 5kg 10kg 25kg 50kg

Cat. No. 221030 Dehydrated - 100 g

221010 Dehydrated - 10 kg

### BiTek™ Trypticase™ Peptone

100g 500g 1kg 2.5kg 5kg 10kg 25kg 50kg

Cat. No. 211101 Dehydrated - 100 g

211102 Dehydrated - 2.5 (2.3 kg)

211103 Dehydrated - 25 (23.3 kg)

### Bacto™ Tryptone

100g 500g 1kg 2.5kg 5kg 10kg 25kg 50kg

Cat. No. 211105 Dehydrated - 500 g

211109 Dehydrated - 2 kg

211104 Dehydrated - 10 kg

### BiTek™ Tryptone

Cat. No. 201620 Dehydrated - 10 kg