**USER MANUAL
MANUAL DE USUARIO
MANUEL D'UTILISATION
BEDIENUNGSANLEITUNG**

**Microbiological Safety Cabinets
Cabinas de Seguridad Microbiológica
Cabines de sécurité microbiologique
Mikrobiologische Sicherheitswerkbänke**

BiOptima / BiOptima Cyto

USER'S MANUAL

Microbiological Safety Cabinet *BiOptima / BiOptima Cyto*

Failure to read, understand and follow the instructions in this manual may result in damage to the unit, injury to operating personnel and poor equipment performance.

SUMMARY

0.	SYMBOLS AND ABBREVIATIONS	4
1.	GENERAL DATA	5
2.	DESCRIPTION OF THE SAFETY CABINETS	6
2.1.	Working principle	6
2.2.	Applications	6
2.3.	Components description	7
2.4.	Standard equipment	7
2.5.	Certification	7
2.6.	Technical features	8
3.	INSTALLING INSTRUCTIONS	10
3.1.	Transport	10
3.2.	Placing on site	10
3.3.	Precautions for first connection	11
3.4.	Electrical connection and services	12
3.5.	Cabinet exhaust ducting	13
3.6.	Installation tests	14
4.	OPERATING INSTRUCTIONS	15
4.1.	Start-up	15
4.2.	Front panel positions	16
4.3.	General working recommendations	18
4.3.1.	Ergonomics	19
4.4.	Definition of the working area	22
4.5.	Control keyboard-display	23
4.6.	Disconnection	29
4.7.	Prolonged stoppage	29
4.8.	Optional accessories	30
5.	MAINTENANCE	31
5.1.	Maintenance table	31
5.2.	Absolute filters	32
5.3.	Spare parts	33
6.	CLEANING AND DISINFECTION	34
6.1.	Surface cleaning and disinfection	34
6.2.	Front glass cleaning and disinfection	34
6.3.	Disinfection with vaporized hydrogen peroxide	34
7.	MAINTENANCE TABLE	35
8.	ELECTRIC DIAGRAMS	35
9.	TESTS	35
9.1.	Downflow filters integrity test	35
9.2.	Exhaust filters integrity test	35
9.3.	Inflow filters integrity test (only Cyto model)	35
9.4.	Downflow air velocity test	35
9.4.1.	Set point alarm downflow air velocity test	35
9.5.	Inflow air velocity test	35
9.5.1.	Set point alarm inflow air velocity test	35
9.6.	Smoke test	35
10.	CERTIFICATES	36
10.1.	EU declaration of conformity	36
10.2.	Guarantee	36
11.	TROUBLESHOOTING	37

Translation of the original manual

0. SYMBOLS AND ABBREVIATIONS

List of stickers used in the biological cabinets.

"EARTH WIRE PROTECTED" STICKER.	
"FUNCTIONAL EARTH WIRE" STICKER.	
"BIOHAZARD" STICKER.	
"CLASS II CABINET CERTIFIED ACCORDING TO EN-12469-2000" STICKER.	
"CYTOTOXIC" STICKER	
"ATTENTION, ELECTRIC SHOCK RISK" STICKER.	
"SOCKET VOLTAGE" STICKER.	
"SOCKET MAXIMUM CURRENT" STICKER.	

"SAFETY WARNING PRIOR TO SERVICE OR REPAIR ACTIONS" STICKER

"DOP PORTS" STICKER

Abbreviations used in this manual:

V	Volts	s, h	Seconds, hours
A	Amperes	Hz	Hertz
W, kW, μ W	Watts, kilowatts, microwatts	dB	Decibel
mm, cm, m	Millimetres, centimetres, metres	$^{\circ}$ C	Celsius
kg	Kilogram	MSC	Microbiological Safety Cabinets
GS	"Tested safety" according German standards		

1. GENERAL DATA

Operating Instructions Manual for Microbiological Safety Cabinets Class II **BiOptima** according to EN 12469 standards, and version **BiOptima Cyto** for cytostatics, according to EN 12469 and DIN 12980 standards, which ensure (applying correct techniques in laboratory) high level protection for the operator, the environment and the sample.

WARNING

This Manual must be read and understood before its installation, maintenance and operation of the biological safety cabinet Class II. All the safety instructions given in the Manual must be followed.

Any use or application failing to comply with the recommendations expressed in the Manual may cause danger and loss of rights to manufacturing guarantee as well as freeing the manufacturer from all responsibility arising.


Any maintenance must be done by qualified staff according to the laboratory working correct conditions. In case these actions are not carefully carried out, as a result the safety cabinet can work in abnormal conditions.

All maintenance actions on this equipment must be performed by a qualified technician who is familiar with the proper maintenance procedure required for this equipment.

The European Union issued the Directive 2012/19/EU on Waste of Electrical and Electronic Equipment (WEEE).

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive. Mark symbol is shown at the right.

Azbil Telstar sells products through distributors throughout Europe. Contact your local distributor for recycling/disposal.

The name plate & CE mark sticker is located at one side of the cabinet.

MANUFACTURER: AZBIL TELSTAR TECHNOLOGIES, S.L.U

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E-mail: Telstar@telstar.com

<http://www.telstar.com>

2. DESCRIPTION OF THE SAFETY CABINETS

2.1. WORKING PRINCIPLE

CLASS II Biohazard Safety Cabinets are designed for microbiology work with the most advanced microprocessor control, ensuring a high level of protection for the operator, environment and sample.

Exterior of cabinet made of laminated, lacquered, oven-dried steel. Functional design with reduced external dimensions and front access for changing filters.

Spacious internal chamber with sliding-swinging front made from laminated safety glass. Segmented work surface and liquid collection tray made of polished stainless steel. Arm-rest is detachable.

Laminar down-flow and exhaust fans provide double security. In the case of a breakdown of one of the laminar down-flow fans, the other ones continues to provide an airflow to the laminar flow while the exhaust fan slightly increases its volume to ensure good negative pressure in the work area.

Principle diagrams in **Annex A1 & A2**.

As can be seen in **Annex A1**, the air propelled by the fans (1) is moved into the chamber or plenum (3) and, on the one hand via the Absolute HEPA filter (4), is filtered and flows according laminar regime in the working area (6). On the other hand, part of the air, also propelled by the exhaust fan (2) flows through an Absolute HEPA filter (5) before being expelled to the outside (8). The frontal air inlet (7) ensures the operator protection.

The unit is equipped with prefilters (9) to protect the interior of the unit against pollution by coarse materials and to collect spilled fluids. Furthermore the pre-filter will enhance the life span of the HEPA filter significantly.

In the BiOptima Cyto model, as can be seen in **Annex A2**, both, inflow and recirculated flow pass through a first filter step with several HEPA filters (10). Access to these filters is completely safe in order to replace them without any risk of contamination.

2.2. APPLICATIONS

The MSC models are appropriate for handling biologically pathogenic samples or products (up to micro organism risk BSL3), ensuring sterile air in the work zone, as well as a high level of protection for the operator and the environment. **BiOptima Cyto** models are also designed to work with chemically volatile samples or products (Cytostatics).

WARNING

Under no circumstances, is possible handling in corrosive, flammable or explosive gas ambient.

It is recommended that extracted air should be expelled to the exterior.

Right working conditions are obtained only when the display is green and without any alarm message on the display.

WARNING

BiOptima Cyto model is designed and manufactured to work with Cytostatic particles.
NEVER USE FOR GASES OR AEROSOLS

2.3. COMPONENTS DESCRIPTION

(See Annexes A1&A2)

2.4. STANDARD EQUIPMENT

- Two inlet sockets 2A - 460W each located on either side of the work zone. Insulation level IP44.
- UV kit

Easy-to-work UV germicidal lamp with a timer function. It can only be enabled in the *Disinfection* screen (**Chapter 4.5 point 2.4**) or by means of the **UV light Auto Power on** screen (**Chapter 4.5 point 3.5.1.2**). For safety reasons, this button is electrically linked with the light button, and it switches off if the front window is opened.

Average disinfection intensity over the work zone of 200 $\mu\text{W}/\text{cm}^2$ with 20% depreciation after a working life of 8000 hours.

WARNING

In the event of using this accessory, you should consider the normal precautions in the use of these types of lamps. Their use is only suitable for previous external sterilisation of materials introduced in the cabinet. So this lamp accessory should be used before the start of work in the cabinet and it is essential to keep it turned off when the operator is in front of the working area or working on it. The glasses provide protection from direct radiation but UV rays are reflected on polished metal parts and can seriously damage operator's eyes. That is why the glass front must be placed into the closed position to be able to switch on the U.V light, then the fans will automatically disconnect.

2.5. CERTIFICATION

Quality control and tests results certification carried out in line with the applicable standards EN 12469 & DIN 12980. See Certificate inside the cabinet.

2.6. TECHNICAL FEATURES

Features	Units.	BiOptima 4	BiOptima 6
External dimensions (LxWxH) (Width including armrest)	mm	1338 x 825 x 1450	1948 x 825 x 1450
Internal dimensions (LxWxH)	mm	1190 x 548 x 666	1800 x 548 x 666
Front opening height	mm	200	
Weight	Kg	220	285
Downflow velocity	m/s	0.35	0.35
Inlet velocity / Exhaust flow	m/s m ³ /h	0.65 552	0.6 773
Consumption*	W	189	240
Voltage	V	230	
Frequency	Hz	50/60	
Lighting	Lux	≥ 1000	
Noise	dB (A)	≤ 58 (Attached table)	
Vibration	mm RMS	< 0.005	
HEPA/ULPA H 14 Filters		Efficiency 99.995%, EN 1822 by MPPS method (99.999% D.O.P. test @ 0.3µm particles)	
Fans		High-efficiency fans. Flows measured by thermoanemometric probes.	

* Cabinet with factory settings and new filter.

DISSOCIATED DECLARED VALUES OF NOISE EMISSION Standard(s) to which conformity is declared ISO 4871		
	$V_{\text{downflow}} = 0.35\text{m/s}$ $V_{\text{inlet}} = 0.65\text{ m/s}$ (230 V, 50 Hz)	$V_{\text{downflow}} = 0.35\text{m/s}$ $V_{\text{inlet}} = 0.6\text{ m/s}$ (230 V, 50 Hz)
	BiOptima 4	BiOptima 6
L _{pA} (dB)	≤58	≤58
K _{pA} (dB)	2	2
Values measured according to acoustic test in standard EN 12469:2000 using standard ISO 11201		

Features	Units.	BiOptima Cyto 4	BiOptima Cyto 6
External dimensions (LxWxH) (Width including armrest)	mm	1338 x 832 x 2212	1948 x 832 x 2212
Internal dimensions (LxWxH)	mm	1190 x 548 x 666	1800 x 548 x 666
Front opening height	mm	200	
Weight	Kg	285	380
Downflow velocity	m/s	0.35	0.35
Inlet velocity / Exhaust flow	m/s m ³ /h	0.6 510	0.6 773
Consumption*	W	283	420
Voltage	V	230	
Frequency	Hz	50/60	
Lighting	Lux	≥ 1000	
Noise	dB (A)	≤ 60 (Attached table)	
Vibration	mm RMS	< 0.005	
HEPA/ULPA H 14 Filters		Efficiency 99.995%, EN 1822 by MPPS method filtration system (99.999% D.O.P. test @ 0.3µm particles)	
Fans		High-efficiency fans. Flows measured by thermoanemometric probes.	

* Cabinet with factory settings and new filter.

DISSOCIATED DECLARED VALUES OF NOISE EMISSION Standard(s) to which conformity is declared ISO 4871		
	$V_{\text{downflow}} = 0.35\text{m/s}$ $V_{\text{inlet}} = 0.6\text{ m/s}$ (230 V, 50 Hz)	
	BiOptima Cyto 4	BiOptima Cyto 6
L _{pA} (dB)	≤60	≤60
K _{pA} (dB)	2	2
Values measured according to acoustic test in standard EN 12469:2000 using standard ISO 11201		

3. INSTALLING INSTRUCTIONS

For a correct start-up procedure of the MSC, observe the following step-by-step instructions.

3.1. TRANSPORT

The cabinet is supplied in an individual package and in the suitable way to prevent it from being damaged.

This manual and the individual Quality Control certificate of each cabinet are included in the scope of supply.

3.2. PLACING ON SITE

- a) Special care should be taken when unpacking the front glass and the control display.
- b) Use the same packaging base for the transport to its final location and it is highly recommended doing it with a pallet jack.

WARNING 

Handle the pallet jack carefully to avoid turn over risks due to the high gravity center. Never open the frontal before the cabinet is fully fastened to its support stand.

WARNING 

We advise against the use of lift trucks for placing the cabinet directly by its base.

WARNING 

For any handling during assembly or dismantling of the cabinet, appropriate safety measures should be taken, such as wearing safety gloves or shoes.

- c) Cabinets should be placed away from entrances, corridors, ventilation supply points, etc., which might affect airflow in the work area.
- d) Ensure the cabinet is set perfectly on the table or workbench in order to prevent vibrations.
- e) To prevent the work area from overheating, we recommend that the room temperature is always within the comfort margins for working.

The optimal comfort margins are:

- Temperature: 18-28°C
- Relative humidity: 50% ± 20%

If any of these conditions is not fulfilled, optimum cabinet performance cannot be guaranteed.

- f) When necessary, air vents in the room should be modified or speed of air entry reduced so that air is not aimed directly at the cabinet, as this may produce turbulence. Side air vents are not appropriate for rooms where cabinets are installed.
- g) Due to the ceiling temperature affects air distribution, sites should be selected with enough space. We recommend 250 mm minimum free space between the exhaust and the ceiling.

- h) The location of the cabinet will directly influence the noise level. This level will rise considerably in small spaces with very reflective surfaces because of the proximity of the cabinets to the walls, and where fan speed has risen to compensate for air currents in the room or the resistance of clogged filters.
- i) Correct cabinet performance should be tested before putting it into operation if any modification has been made to the room or to the air intake system which may considerably change the airflow characteristics in the room.

3.3. PRECAUTIONS FOR FIRST CONNECTION


- Before connecting it to the mains supply, proceed to make a general cleaning to remove any dust particles that may be accumulated during transportation. Follow the instructions on cleaning and disinfection.
- Before using it for the first time, run the cabinet for around 6 hours to purge the filters. New filters give off a characteristic smell which tends to disappear in a short period of time.

WARNING 

The cabinet should work at the voltage shown on the name plate. It should be connected to a suitable socket base (minimum 16A) with an earthed plug. The cable section recommended for the base is 1.5 mm².

WARNING 

The plug connection must be located in such a place to prevent it from being unplugged intentionally or accidentally during operation, e.g., next to a wall or closed to another equipment.

3.4. ELECTRICAL CONNECTION AND SERVICES

The electrical characteristics of the cabinet are as follows:

- Electrical class: Class I.
- Insulation level: IP-20.


The cabinet is supplied with a plug cord according European standard (other plug cords available as optional). The on/off button is in the frontal control keyboard-display. (**Annex A3 pos. 1**).

The power socket in the wall should be located in a perimeter around the cabinet of:

- 1.25 m from the middle of the cabinet to the sides
- 2.5 m from the base of the cabinet upwards

On the top right-hand side of the cabinet are:

- Main fuse (1)
- Power supply (2)

WARNING

Unplug the cabinet from the electric current supply, in case of electrical maintenance.

The cabinet can only be activated from the  button in the frontal control keyboard.

In the event of a power failure the cabinet will stop, showing alarm message and providing enough energy to close the front by continuously pressing the additional keyboard key . Once the cabinet is closed, it can be disconnected.

The installation which the cabinet is connected to should meet the safety requirements in force in each country.

Exhaust air should be evacuated outside the lab if working with volatile substances not neutralised by absolute filters.

During normal work in these safety cabinets and because of a basic principle of recirculation and depression in the work zone, no perforated surface in the work zone should be covered with paper or sanitary fabric.

3.5. CABINET EXHAUST DUCTING

Ducting of the exhaust may be preferred to exhaust in the lab in order to, for example:

- Prevent from lab temperature rise.
- When volatiles and toxic vapours are involved during the process.

WARNING 

Ducting the exhaust to the outside of the lab creates a depression in the lab.


WARNING 

When the cabinet exhaust is connected to an assisted/air forced exhaust system, it is necessary to install a canopy between the cabinet and the exhaust duct. This will prevent from the loose of operator protection in case of any decrease of the exhaust airflow system (duct obstruction, blower control failure...).

If an additional fan is required to duct the air outside of the laboratory, it must be connected to the CPU Board in order to send an alarm if there is a breakdown. When installing and connecting the additional fan to the cabinet this option must be activated by the service technician. Velocity parameters must also be adjusted.

The following conditions must be fulfilled in order not to disturb the normal operation of the unit and its performance:

- Duct of at least 250 mm diameter.
- Outlet protected from back airflow, closure, entry of any element in the duct which may result in air flow restriction (animals...).

GS has not been tested for the canopy. It must be tested after installation.

3.6. INSTALLATION TESTS

According to standard DIN EN 12469 and in order to ensure correct cabinet operation, after the installation and connection of the cabinet, it is compulsory to perform the following tests (before the equipment is used for the first time)

- Impulsion filters integrity Test.
- Exhaust filters integrity Test.
- Inflow filters integrity test (Only model BiOptima Cyto)
- Downflow air velocity.
- Downflow set points alarm test.
- Inflow air velocity Test.
- Inflow set points alarm test.
- Smoke test.

TELSTAR also recommends making the following tests:

- Light Test.
- Noise Level Test

Telstar, its representatives and distributors put a qualified technical team at your disposition that uses control equipment with maximum health protection for carrying out this installation test under full safety guarantees.

4. OPERATING INSTRUCTIONS

4.1. START-UP

Once the above points have been controlled, start up the MSC by pressing the On/Off button (**Annex A3, pos. 1**). After entering the password correctly (by default the password is **0000**) or immediately if this function has been disabled, the fans will automatically start working.

Ensure that the aspiration area and air exhaust outlet are free (**Annexes A1 & A2, pos. 7 & 8**).

The electronic system of regulation by microprocessor via two thermoanemometer probe, laminar down flow and exhaust airflow, needs a few seconds to stabilise during connection in order to begin to work normally.


The cabinet can only be used by qualified staff who meets the correct work techniques in the laboratory. If these techniques are not fulfilled, faulty working of the equipment may happen.

WARNING 

The EN 12469 certification does not allow the deactivation of the password enquiry.

4.2. FRONT PANEL POSITIONS

By means of the central buttons, the front glass can be moved.

: It opens the front glass. Keep it pressed for some seconds to activate movement. Then release the button. When activated below working position, it will automatically stop just above it and go down till reaching the precise working position.

: It moves down the front glass. Due to security reasons, it is necessary to keep it pressed continuously to allow down direction of the glass. Depending on the initial glass position, it will automatically stop at the precise working position or in the closed position.

Note: Lights will blink during glass movement.

Closed position

Working position

Open position

WARNING 

The movement of the glass exerts a strong force. Check that no one has their hands, any part of body or obstacles below the glass when it is moving down, as it can cause injuries. Warn other people working near the cabinet.

In case it is needed to clean the glass by the inner part (Annex A5.2) or if it is necessary to place in the working area any big product or device, the front glass may be separated from the cabinet (after fully lifting it), providing full access to both the inner glass area and the working area. All needed precautions should be taken as the user protection decreases during that operation.

REMARK 

In case of power failure, user is able to close the glass thanks to the battery that is installed for that purpose.

WARNING 

The glass is required to be in full lifted position before closing the front glass towards the cabinet. In case this is not fulfilled, the glass may be damaged.

4.3. GENERAL WORKING RECOMMENDATIONS

- a) We recommend starting up the cabinet around 10 minutes before you begin to work. This will ensure particles are swept from the work zone and from material introduced.
- b) Wash your arms, hands and nails well, using an anti-bacterial soap, before and after work. Staff should avoid touching their mouths and eyes.
- c) We recommend wearing lab coats with adjusted cuffs and, in special works, protection gloves (or sleeve cuffs). Both the coat and cuffs should be made from a fabric that minimises the emission of fibres and particles.
- d) The work area inside the cabinet should not be used to store lab equipment. This can produce an accumulation of dust that is unnecessary and hazardous for working under sterile conditions. Only the fewest pieces of equipment should be placed in order to work correctly and with the minimum number of interruptions for the user (inserting and removing arms from the work zone).
- e) Although the work zone is under a clean and sterile laminar down flow, the area around the cabinet might be contaminated. The sources of contamination should be known in order to avoid them when inserting arms, material, etc.
- f) All the material needed for the work should be free of particles and carefully cleaned before it is placed on the work area.
- g) Materials such as paper, wood, cardboard, pencils, rubbers, etc. should not be placed in the work area as they shed a lot of particles.
- h) In the event of using pipettes, only mechanical aspiration type should be used. Blowing the pipette with the mouth should never be done because it is easy to inhale aerosols caused by suction.
- i) When platinum loops are used, it is advisable to use electric incinerators, and preferably single-use (disposable ones) loops.
- j) If the work requires the use of a gas flame from a Bunsen burner or similar, we recommend using a type where the flame is maintained by holding down a button. Please note that the use of a constant flame burner creates significant turbulence. Bear in mind that a flame too large could burn the absolute filters.
- k) Shadows and turbulence provoked by objects, equipment and materials placed in the work area should be studied before beginning to work in order to evaluate their possible effects. Bear in mind that the laminarity of the airflow is not recovered until a distance of 2.5 times the diameter of the object that causes the obstruction.
- l) When using bottles and tubes, they should preferably have a screw top instead of cotton stoppers which easily shed a lot of particles.
- m) When the work is finished, all disposable products (handles, Petri dishes, etc.) as well as means of cultivation, samples, tubes, bottles, etc. should be removed from the cabinet in impermeable bags suitable for sterilisation if necessary.
- n) When working in the cabinet and during cleaning operations, it is important to protect the HEPA filters from any damage by avoiding knocking them, projecting liquids or splashing, etc.
- o) Any equipment connected to the internal sockets inside the cabinet should meet all the electrical specifications.
- p) People traffic generates air movements that may cause unexpected loss of the containment barrier and, in the case of Class II cabinets, a reduced level of product

protection. Warning signs and other means should be used to avoid traffic during critical operations.

Cabinet locations can be divided to reduce the effects of air movement produced by other cabinets and people traffic.

- q) Suddenly doors opening and closing may disturb airflows in the cabinet and it should be avoided when cabinets are being used.

WARNING 

The UV light system can only be used when the front glass is closed.

Never connect equipment whose consumption exceeds the maximum one indicated on the internal sockets of the cabinets.

4.3.1. ERGONOMICS


Microbiological security cabinet BiOptima and its accessories have been designed following recommendations and specifications of the standard EN ISO 14738:2010.

Users should consider their position when working with the cabinet, for their own safety.


Adjusting properly the seat and adopting a correct position enables the user to work safely and to avoid any possible physical damage (sight, neck, back, etc.)

Below are examples of four positions **NOT** recommended and one that **IS** recommended.


In addition, in order to protect both the product and the user, it is necessary that the cabinet is placed in a location without air currents, which may cause turbulence in the laminar down flow.

4.4. DEFINITION OF THE WORKING AREA

To prevent from creating any turbulence in the laminar down flow and to protect the user and the product, it is necessary to establish a working area inside the cabinet.

- Do not work in any area with perforations and do not cover them up.
- Always work within the limits according the following image.

4.5. CONTROL KEYBOARD-DISPLAY

The following elements are found in the control keyboard:

A. On/Off Button  (**Annex 3 pos. 1**):

- Pressing when the cabinet is turned off:
It goes to the *Password Request* screen. When entering the correct password (**Chapter 4.5 point 1**), it goes to the *Main Operation* screen.
- Pressing when the cabinet is turned on:
It goes to the *Password Request* screen. When entering the correct password (**Chapter 4.5 point 1**), the cabinet switches off.

B. Navigation buttons ,  and  (**Annex 3 pos. 2**):

They move the cursor through all different function icons and select them. The icon pointed inverts its colours. If the selected item is a numeric field, the navigation buttons increase/decrease and validate the value.


C. Backlight graphic display (**Annex 3 pos. 3**):

It shows alarm messages, operating graphs, function icons and all the information required in cabinet-user communication. The background colour changes to:

- Green: Working mode. All the flows are working fine and there are no alarms activated.
- Red: Alarm mode. Any alarm has been activated or the cabinet is in unsafe mode. In case of alarm, an acoustic alarm is activated together with the red light.
- Yellow: Disinfection mode.
- White-blue: Information mode. It is used in Technical Service screens.

Below is an explanation of the different system screens:

1. **Password Request** screen (by default **0000**).

Increase the value of the selected field.

Decrease the value of the selected field.

Confirm the value of the selected field and move the selection to the next field. After pressing for the final time:

- Right Password: it accesses the *Main Operation* screen.
- Wrong Password: Cursor returns to the first field.

2. **Main Operation** screen.

It shows the direct access icons to the most common functions for the user.

Activate/Deactivate fans.

Activate/Deactivate lights.

* Press for 5 seconds the  button, to access the light intensity screen.

Activate/Deactivate sockets.

Activate/Deactivate gas solenoid valve (if assembled).*

Activate/Deactivate ECO mode: fans in low-velocity, lights switched off and acoustic alarm can be muted.

Go into Disinfection Menu.

Go to Configuration Menu.

***REMARK:** It is not possible to operate gas tap solenoid valve under unsafe conditions.

2.1. **Fan Deactivation Confirmation** screen.

Increase the value of the selected field.


Decrease the value of the selected field.

Confirm the value of the selected field and move the selection to the next field.

2.2. **Light intensity** screen.

Increase the value.


Decrease the value.

Confirm the value.

2.3. **ECO Mode Confirmation** screen.

Increase the value of the selected field.

Decrease the value of the selected field.

Confirm the value of the selected field and move the selection to the next field.

2.4.  **Disinfection Confirmation** screen.


Increase the value of the selected field.

Decrease the value of the selected field.

Confirm the value of the selected field and move the selection to the next field. After pin confirmation it goes to Disinfection Menu screen.

2.4.1. **Disinfection Menu** screen.

Go to *U.V. Menu* screen.

Go to *Vaporized hydrogen peroxide* screen.

Go to *Particles counter* screen (when applicable).

Back to *Main Operation* screen.

2.4.1.1. **U.V. Menu** screen.


Begin disinfection, if glass is closed, for the time determined by the timer.

Go to *U.V. Time Adjust* screen.

Back to *Disinfection Menu* screen.

2.4.1.1.1. **U.V. Lamp** screen.

Switch off the U.V. lamp and access to *Disinfection Menu* screen.

2.4.1.1.2. **U.V. Time adjust** screen.


Confirm the selected value and save it.*

Reset the numeric value.

Back to *U.V. Menu* screen.


* Value 000 means infinite time: start and stop will be made manually.

2.4.1.2. Vaporized hydrogen peroxide screen.

-  Start Vaporized hydrogen peroxide cycle.
-  Activate/Deactivate the sockets of the work zone.
-  Back to *Disinfection* Menu screen.

2.4.1.2.1. Vaporized hydrogen peroxide cycle screen.

-  Switch on manually the fans.
-  Switch off the fans.
-  Back to *Disinfection* menu screen (it requires confirmation).


3. Configuration Menu screen.

It shows the direct access icons to the configuration functions allowed to the user.

-  Go to Contrast screen.
-  Go to Maintenance screen.
-  Go to Language screen.
-  Go to PIN Configuration screen.
-  Go to Date & Hour screen.
-  Go to Sensors Information screen.
-  Go to Main Operation screen.


3.1. Contrast screen.

-  Increase the value.
-  Decrease the value.
-  Confirms the value and returns to *Configuration Menu* screen.

3.2. **Maintenance Menu** screen.

It shows information about the running hours of cabinet, fans, U.V. lamp and filters, together with the last Service visit. Besides it allows the following:

-  Record a new revision.
-  Get maintenance recommendations.
-  Show the historical of revisions and alarms.
-  Back to *Configuration Menu* screen.

3.3. **Language Selection** screen.

-  Go up.
-  Go down.
-  Confirms the selected language and goes back to *Configuration Menu* screen.

3.4. **PIN Configuration** screen.


-  Enable/disable the password request to Connect/Disconnect the cabinet.
-  Enable the user to change the user password.
-  Coded password. In case of forgotten password, by giving this code to Telstar the previous password will be recovered.
-  Back to *Configuration Menu* screen.

3.5. **Date-hour** screen.

-  Saves the value.
(format dd/mm/yy hh:mm).
-  Cancel and/or edit the entered value.
-  Access to *Auto Power-on* menu.
-  Back to *Configuration Menu* screen.

3.5.1. Auto Power on screen.

-  Cabinet Auto Power on.
-  UV light Auto Power on.
-  Back to *Date-hour* screen.

3.5.1.1. Cabinet Auto Power on screen.

-  Saves the value.
(format dd/mm/yy hh:mm).
-  Cancel and/or edit the entered value.
-  Back to *Auto Power on* screen.


3.5.1.2. UV light Auto Power on screen.

-  Saves the value.
(format dd/mm/yy hh:mm).
-  Cancel and/or edit the entered value.
-  Back to *Auto Power on* screen.

3.6. Sensors Information screen.

It shows information about the values of the sensors of the cabinet in real time.

-  Move window up.
-  Move window to the work-position.
-  Move window down.
-  Back to *Configuration Menu* screen.

D. Microprocessor:

The microprocessor keeps total control over the exhaust flow and impulsion flow via the thermoanemometer probes.

The microprocessor control system consists of:

- Backlight graphic display.
- Alarms shown on the display:

ALARM	SHOWN ON THE DISPLAY
Excessive laminar down flow	HIGH LAM. FLOW
Insufficient laminar down flow	LOW LAM. FLOW
Excessive inflow	HIGH INFLOW
Insufficient inflow	LOW INFLOW
Glass out of working position	WORK POSITION
Front opened	LIFTED FRONTAL
Redundant channel alarm	REDUNDANT FAIL.
Power supply failure	POWER FAILURE
Memory fail	MEMORY FAILURE
Exhaust duct obstructed	CLOGG.EXH.DUCT

- Messages shown on the display:

MESSAGE	SHOWN ON THE DISPLAY
Fans at reduced speed	ECO MODE
Uv light switched on	GERMICIDE
Battery needs replacement	LOW BATTERY
Glass in upper position	WINDOW OPEN
Window is moving	WINDOW MOVING

E. Auxiliary microprocessor or redundant circuit.

Should the main microprocessor fail, the auxiliary microprocessor will activate the optical alarm, and also the acoustic alarm. In such a case, move the glass to the closed position and cover the cabinet exhaust. Disconnect the equipment and contact the Telstar Technical Service.

4.6. DISCONNECTION

To stop the cabinet, it is necessary to push the On/Off button (**Annex A3 pos. 1**) to access the *Password Request* screen (**Chapter 4.5 point 1**). After entering the correct password, the cabinet will disconnect.

4.7. PROLONGED STOPPAGE

In the event of a prolonged stoppage, avoid environmental dust to go inside the cabinet by turning the cabinet off with the glass in the closed position. It is recommended to disconnect the battery to avoid its full discharge. When it is switched on again, proceed like the first time it was switched on. See **Chapters 3.3** and **4.1**.

4.8. OPTIONAL ACCESSORIES

- Support stand (two models for class II cabinets of 800 mm and 900 mm working area height. The cyto version cabinet includes model 900 mm working area height).
- Gas tap kit.
- Vacuum tap kit.
- Compressed air tap kit.
- Nitrogen tap kit.
- Oxygen tap kit.
- Water tap kit.

WARNING

Gas is a toxic and flammable substance which in the event of an accident can cause serious damage to the user and the cabinet.

The presence of a flame within the cabinet creates turbulence which affects its control and proper functioning, therefore we recommend you do **NOT** use burners unless strictly essential and if the user is prepared to assume this responsibility.

Gas connection should be done in line with the safety standards of each country, region and company.

For these reasons, the electrically controlled valve that enables gas entry can only be activated when the fans are working and it is automatically disconnected when the alarm indicating a lack of impulsion flow is activated.

5. MAINTENANCE

WARNING

All repair and maintenance work in the cabinet should be done by qualified technical staff with knowledge about the maintenance procedures required by the cabinet.

Before starting any maintenance work or dismantling the cabinet, the area should be wiped down with vaporized hydrogen peroxide (**Annex A5.4**). Once this has been done, disconnect the cabinet from the power socket.

In general, a regular control of the proper working of all cabinet components should be done once a year.

NOTICE

It is recommended to change the batteries every 2-3 years to guarantee the proper function. The disposal of exhausted batteries must be done in accordance with local requirements. The batteries constitute a chemical residue.

5.1. MAINTENANCE TABLE

The frequencies determined for each operation, and in particular with regard to changing filters, which appear on the Maintenance Table below are indicative values recommended for normal cabinet service. Very severe work conditions may reduce the time intervals shown on the Table.

TEST (See annex A9)	REGULARITY
<i>Downflow filters integrity Test</i>	<i>Minimum once a year or after filter changes</i>
<i>Exhaust filters integrity Test</i>	
<i>Aspiration filters integrity Test</i>	
<i>Downflow air velocity test</i>	
<i>Set point alarm downflow air velocity test</i>	
<i>Inflow air velocity Test</i>	
<i>Set point alarm inflow air velocity test</i>	
<i>Smoke test</i>	
<i>Light test</i>	<i>Recommended during the start up of the cabinet</i>
<i>Noise level Test</i>	

5.2. ABSOLUTE FILTERS

Under normal circumstances, the life span of a HEPA filter comprises several years. However, it will become clogged up gradually, until there is no longer a proper flow-through and the filter has to be replaced.

After this time, the propelled airflow falls because of the increased pressure drop of the filter as it is a mechanical filtering, although this does not mean that filtering efficacy is reduced. When it reaches the minimum flow measured by the thermoanemometric probes, the microprocessor activates both visual and acoustic alarms.

WARNING 

Do not replace absolute filters or directly handle fans in MSCs without having first wiped the area down with vaporized hydrogen peroxide. Opening the cabinet without taking the proper precautions can endanger your health. We recommend it is carried out by staff especially trained for this work.

To ensure correct installation and adjustment of the filter, as well as ensuring there are no leaks, filters should be replaced by specialised technical staff equipped with the corresponding monitoring instruments: anemometer, leak detector, aerosol generator, particle counter, etc.

The used clogged absolute filters from safety cabinets should be destroyed. This process is always the responsibility of the user. Our technical service does not have any responsibility for possibly contaminated parts.

WARNING 

In the cytotoxic cabinets, the replacement of the aspiration filters must be done with safety conditions keeping the cabinet working and following the instructions described in the Service Manual.

The pre-filter is assembled under the worktop.

According to EN12469:2000 Annex K the pre-filter needs to be replaced on every maintenance occasion. Replacement can be done easily following the instructions described in the Service Manual.

5.3. SPARE PARTS

Only use original spare parts from the manufacturer.

Please do not hesitate to contact the Customers Service pointing out the bench model, operating frequency (Hz), voltage (V) and serial number of the unit.

BiOptima all models

No.	Description
1	UV LAMP PHILIPS TUV 36T5 4P SE UNP
2	CPU CARD
3	BATTERY (2 UNITS)
4	SPEED PROBE 0.2-1.0 M/S (2 UNITS)

BiOptima 4 / Cyto 4

No.	Description
5	DOWNFLOW (2 UNITS) / EXHAUST FAN (1 UNIT)
6	PRE-FILTER BIOPTIMA 4 G3 (EN 779)
7	ABSOLUTE FILTER H14 1219X457X90
8	ABSOLUTE FILTER H14 610X457X90
9	LED LIGHTING BIOPTIMA 4

BiOptima 6 / Cyto 6

No.	Description
10	DOWNFLOW (3 UNITS) / EXHAUST FAN (1 UNIT)
11	PRE-FILTER BIOPTIMA 6 G3 (EN 779)
12	ABSOLUTE FILTER H14 1829X457X90
13	ABSOLUTE FILTER H14 1219X457X90
14	LED LIGHTING BIOPTIMA 6

BiOptima Cyto 4

No.	Description
15	V-SHAPE FILTER (9 UNITS)

BiOptima Cyto 6

No.	Description
16	V-SHAPE FILTER (14 UNITS)

6. CLEANING AND DISINFECTION

6.1. SURFACE CLEANING AND DISINFECTION

See **Annex A5.1**

6.2. FRONT GLASS CLEANING AND DISINFECTION

See **Annex A5.2**

6.3. DISINFECTION WITH VAPORIZED HYDROGEN PEROXIDE

It is recommended to do decontamination of the cabinet by using vaporized hydrogen peroxide emission, which ensures penetration in sealed and inaccessible areas of the Cabinet. It also guarantees the sterilisation of the working area, fans, plenum, HEPA filters and evacuation ducts.

Vaporized hydrogen peroxide systems provide rapid, safe, low-temperature decontamination methods for any enclosed area that may be contaminated with micro-organisms, including spore-forming bacteria.

See **Annex A5.4** for the disinfection with vaporized hydrogen peroxide procedure.

7. MAINTENANCE TABLE

See **Annex A6**

8. ELECTRIC DIAGRAMS

See **Annexes A7 (230V)**

9. TESTS

9.1. DOWNFLOW FILTERS INTEGRITY TEST

See **Annex A9.1**

9.2. EXHAUST FILTERS INTEGRITY TEST

See **Annex A9.2**

9.3. INFLOW FILTERS INTEGRITY TEST (ONLY CYTO MODEL)

See **Annex A9.3**

9.4. DOWNFLOW AIR VELOCITY TEST

See **Annex A9.4**

9.4.1. SET POINT ALARM DOWNFLOW AIR VELOCITY TEST

See **Annex A9.4.1**

9.5. INFLOW AIR VELOCITY TEST

See **Annex A9.5**

9.5.1. SET POINT ALARM INFLOW AIR VELOCITY TEST

See **Annex A9.5.1**

9.6. SMOKE TEST

See **Annex A9.6**

10. CERTIFICATES

10.1. EU DECLARATION OF CONFORMITY

See **Annex A10.1**

10.2. GUARANTEE

See **Annex A10.2**

11. TROUBLESHOOTING

No.	Type of error	Causes of error	Correction of error
1	Insufficient exhaust flow Alarm "LOW INFLOW" (see Chapter 4.5 point D)	Exhaust partially covered.	Visually check the exhaust to remove any type of obstruction.
		Front perforation in the working area totally/partially covered.	Visually check the perforation and remove any type of obstruction.
		Distance between cabinet exhaust and lab ceiling too small.	Increase space between cabinet exhaust and the ceiling. In some cases it is possible to modify or eliminate the false ceiling. Install a lower support table. This means the working area will be also lower. Duct the cabinet (Chapter 3.5).
		Clogged exhaust filter or aspiration filter (Model BiOptima Cyto).	Change the Prefilter (Annex A4.1). Check the date of the last filter change and the hours it has been under operation. (Chapter 4.5 point 3.2). Notify the technical service to change the filter.
		Ducted cabinet and ducting too long.	Fit an additional fan at the end of the duct (Chapter 3.5).
2	Insufficient laminar down flow Alarm "LOW LAM. FLOW" (see Chapter 4.5 point D)	Front and rear perforation in the work zone totally/partially covered.	Visually check the perforation and remove any type of obstruction.
		Clogged laminar down flow filter or aspiration filter (Model BiOptima Cyto).	Change the Prefilter (Annex A4.1). Check the date of the last filter change and the hours it has been under operation. (Chapter 4.5 point 3.2). Notify the technical service to change the filter.
3	Excessive laminar down flow Alarm "HIGH LAM. FLOW" (see Chapter 4.5 point D)	Air currents in front aperture.	Check and remove any airflow in front of the cabinet that distorts the reading of the sensor.

No.	Type of error	Causes of error	Correction of error
4	Laminar down flow not working properly Alarm "LOW LAM. FLOW" (see Chapter 4.5 point D)	Probe not working properly.	Check that the connection in the control card is correct (Annex A7).
		Laminar down flow fan stopped.	Check the fuse has not blown (Annex A7). Check that the connections in the control card are correct (Annex A7). Notify the Technical Service to change the fan.
5	Exhaust not working properly Alarm "LOW INFLOW" (see Chapter 4.5 point D)	Exhaust completely covered up.	Visually check the exhaust to remove any type of obstruction.
		Probe not working properly.	Check that the connection in the control card is correct (Annex A7).
		Exhaust fan stopped.	Check the fuse has not blown (Annex A7). Check that the connections in the control card are correct (Annex A7). Notify the Technical Service to change the fan.
6	Window does not move	Front is blocked.	Visually check the front movement to remove any type of obstruction. Check the fuse has not blown (Annex A7).
		Motor of the window does not work.	Notify the Technical Service to change the motor of the window.
		Inductive sensors unadjusted.	Notify the Technical Service to adjust the Inductive sensors position.
		Main control is damaged.	Call a specialised technician to replace the control.
7	Glass out of position Alarm "WORK POSITION" (see Chapter 4.5 pto. D)	Glass position.	Move the glass into the working position.
		Working position inductive sensors unadjusted.	Open the frontal and check that the working position inductive sensor is not loose. Notify the Technical Service to receive instructions to adjust the inductive sensor.
		Main control is damaged.	Call a specialised technician to replace the control.

No.	Type of error	Causes of error	Correction of error
8	Open frontal Alarm "LIFTED FRONTAL" (see Chapter 4.5 pto. D)	Lifted frontal inductive sensor unadjusted.	Open the frontal and check that the open frontal position inductive sensor is not loose. Notify the Technical Service to receive instructions to adjust the inductive sensor.
		Main control is damaged.	Call a specialised technician to replace the control.
9	Light does not turn on	Light not working properly.	Check that the connection in the control card is correct (Annex A7).
			Call a specialised technician to replace the Led bar.
10	UV germicidal lamp does not turn on	Closed frontal inductive sensor unadjusted.	Open the frontal and check that the closed frontal position inductive sensor is not loose. Notify the Technical Service to receive instructions to adjust the inductive sensor.
		Glass not in closed position.	Check the front slides properly.
		UV lamp does not work.	Check the lamp holders are properly connected. Check that the connection in the control card is correct (Annex A7). Change the UV lamp.
		Ballast is damaged.	Change the ballast.
11	Intermittent alarm activated	Inductive sensors are unadjusted.	Open the frontal and check that the inductive sensors are not loose. Notify the Technical Service to receive instructions to adjust the inductive sensors.
		Main control is damaged.	Call a specialised technician to replace the control.
12	Display does not turn on	Power cord cable disconnected.	Check that power cord cable is properly connected in socket on the wall.
		No power.	Check voltage on the socket of the wall.
		Main control is damaged.	Call a specialised technician to replace the control.

WARNING 

All repair and maintenance work in the cabinet should be done by qualified technical staff with knowledge about the maintenance procedures required by the cabinet.

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ADHESIVO "DOP PORTS"

DoP PORTS
DOWNFLOW LEFT DOWNFLOW RIGHT EXHAUST

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
EN Annex A3. CONTROL ALARM

1. On/Off button
2. Navigation buttons
3. Graphic screen

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
EN

This procedure must be done with a running fan. Then the negative pressure in the working area will stay. This is an extra safety aspect. Since this procedure can be done in a short time, the extra pollution of the HEPA filters will be low.




This procedure is described for a right-hand orientated person. The described left hand may be changed for a right hand and right for left hand as long as this is done consistently (this procedure needs consistence due to contamination danger).

Supplies needed:

- Plastic bag (minimum measurement 75x50 cm with closure).
- Two pair of gloves.
- Tape.
- Possible disposing requirements.
- New pre-filter.

Procedure:

- a) Switch on the unit.
- b) Make sure that the window is in the working position and that there are no alarms.
- c) Put on two pair of gloves.
- d) Lift up the segmented worktop to the left of the chamber, in order to reach the prefilter. (Fig. 1)
- e) Fold the edge of the plastic bag ± 10 cm and place the opening of the plastic bag inside the front aperture with your left hand. (Fig. 2)
- f) Position the plastic bag on the armrest.
- g) Pull the pre-filter from its frame with your right hand.
- h) Roll the filter up (still with your right hand). (Fig. 3)
- i) Put the filter into the plastic bag. (Fig. 4)
- j) Fold back the edge of the plastic bag, taking off your right glove and introducing it in the bag, not touching the inner side of the plastic bag. (Fig. 5)
- k) Press the air out of the bag (with its opening still in the working space).
- l) Seal the plastic bag airtight.
- m) Remove the plastic bag and place it in an appropriate container.
- n) Put the segmented worktop at the right and, using the right hand, introduce and place the new pre-filter in the frame. (Fig. 6)
- o) Put the segmented worktop back in its position.
- p) Take off the gloves and dispose them away in the appropriate container.

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EN Annex A5.1. SURFACE CLEANING AND DISINFECTION

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EN

	<p>CYTOTOXIC DANGER Before starting work in cytotoxics production the hazard risks of cytotoxics handling need to be evaluated and documented (industrial safety act, hazardous substances regulations). Based on these findings the employees must be instructed. For a number of cytotoxics alternatives cleaning procedures could be applied.</p>
	<p>BIOLOGICAL DANGER Never clean a used unit without safety precautions</p>

The unit has to be cleaned regularly. Preferably, each time after the unit has been used.

- Put on gloves and sleeves.
- Use glasses and mouth protection.

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EN

For worktable and sides, use sterile or non-use disposable cloths, which do not shed particles or fibres.

For disinfection, these cloths will be used, slightly dampened with a non-corrosive disinfectant solution (70% alcohol solutions), that do not damage or affect the stainless steel and the glass. Cleaning can also be carried out with water and soap prior to disinfection.

To disinfect correctly, the appropriate disinfectant solution should be used for the type of work that is being done. Many of these solutions require certain application times depending on the type of micro-biological agent used in the cabinet; we recommend that you consult the instructions for application and precautions for use given by the manufacturer of the disinfectant.

WARNING

It is important keep the glass in working position and the cabinet in working mode when performing the cleaning and disinfection of the working area of the cabinet. In case that it is needed to open the frontal in order to clean areas with difficult access, all needed cautions must be taken because the user protection decreases when the frontal is open.

It is recommended NOT to use disinfectant products with chlorinated derivatives or containing halogens or their compounds, due to they can damage the stainless steel and painted surfaces.

Should some bleach is poured, the surface must be always thoroughly rinsed and completely dried with a spatula (glass cleaner type).


It is required to dry completely all surfaces always after cleaning or disinfecting in order to preserve the original qualities of the surfaces.

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EN

To accede to corners or difficult access areas, another type of tools will be used as shown in the image.

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EN

The perforated stainless steel worktable is divided into several parts to facilitate cleaning of the interior container and provide for sterilisation of these parts into a vessel.

Remember that every segment or shelf has to be cleaned and/or disinfected both in its lower and upper sides. The parts will be cleaned inside the cabinet while this is running.

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EN

Also be sure of carrying out the cleaning/disinfection in the lower part of the working area (surface and corners).

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EN

Inside, safety cabinet cleaning is recommended on the following occasions:

- a) *Before starting any work in the safety cabinet.*
 - b) *After working in the safety cabinet.*
 - c) *Whenever there is a change of work program.*
 - d) *In event of liquid spilling on the work table.*
 - e) *Before carrying out mechanical or biological control tests in the working area.*
-

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EN Annex A5.2. CLEANING AND
DISINFECTION OF FRONT GLASS

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EN

For the cleaning of front glass you must use sterile textures or disposable ones without particle or fibres loosing.

For disinfection, use the previous textiles slightly dampened with a non-corrosive disinfectant solution (70% alcohol solutions), that do not damage or alter the glass. Cleaning can also be carried out with water and soap prior to disinfection.

To disinfect correctly, the appropriate disinfectant solution should be used for the type of work that is being done. Many of these solutions require certain application times depending on the type of micro-biological agent used in the cabin; we recommend that you consult the instructions for application and precautions for use given by the manufacturer of the disinfectant.

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EN

For the glass exterior part cleaning, it is recommended to move the glass into closed position. In this position the major part of the glass is accessible.

For cleaning the most accessible interior parts of the glass, it is recommended to leave the glass in its working position in order to assure the user's protection.

WARNING

It's important to clean and disinfect the cabinet's working area, with the glass in working position and the cabinet in performance. If it's necessary to separate the frontal-glass mechanism to clean the areas of difficult access, bring the glass to its fully lifted position to set free the frontal mechanism. By separating the frontal-glass, the inner surface of the glass is totally accessible. There must be taken the necessary precautions because the user's protection is lower with the glass separated from cabinet.

WARNING

The movement of the glass exerts a strong force. Check that no one has their hands, any part of the body or obstacles below the glass, as it can cause injuries. Warn other people working near the cabinet.

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EN

Gas decontamination of a cabinet that has handled biologically pathogenic products is necessary in the following cases:

- *Before undertaking any safety cabinet maintenance work.*
- *Before removal of HEPA filters for replacement.*
- *Before carrying out mechanical or biological control tests in the safety cabinet.*
- *Before moving the safety cabinet from one work area to another.*
- *Before changing a work program.*
- *Whenever liquid is spilled on the work table.*

The gas used in this procedure is the vaporized hydrogen peroxide.

The necessary material to perform decontamination by gas is as follows:

- *1 vaporized hydrogen peroxide decontamination kit accessory.*
- *1 Plastic bag 1800x600 mm.*
- *1 Roll of packaging tape.*
- *1 Roll of adhesive tape.*
- *1 Allen key.*
- *Vaporization hydrogen peroxide system.*
- *Vaporized hydrogen peroxide sensor.*

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

EN

Steps to follow:

- a) *Fumigation must be performed by a qualified person with a good knowledge of the procedure and the required precautions. At any case, a visible warning should indicate the cabinet is being disinfected.*
- b) *Install the vaporized hydrogen peroxide accessory kit in side glass (input port), and the exhaust cover on top (output port). Seal the exhaust cover with tape.*

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
EN

- c) *Separate the frontal of the cabinet.*
- d) *To guarantee the tightness of the equipment, seal the cabinet front opening by means of one plastic bag and packaging tape, covering any gap that could leak hydrogen peroxide to the ambient.*
- e) *Close the front window and place the glass in working position.*
- f) *Switch on the cabinet and set up in ECO mode (low velocity mode). Lights will be automatically switched off, screen will be red and alarm will sound. Mute the alarm.*

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
EN

- g) *Connect the vaporization hydrogen peroxide system hoses to the proper port. Check that input of vapor is fixed in side glass, and output is fixed on top.*
- h) *Cabinet is ready to start with decontamination cycle.*

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EN

Decontamination process.

It must be performed by specialized personnel of the vaporization hydrogen peroxide system.

After cycle is finished, high concentration may remain inside the cabinet. Once the plastic sealing the front and the exhaust cover are removed, concentration will be diluted with room air and be reduced rapidly to accepted levels.

ATTENTION 

During process the hydrogen peroxide concentration in the room must be monitored for safety reasons, so it will be necessary to place a sensor for that purpose. Ambient limit is 1 ppm, and at that concentration 8 h of continuous work is allowed. With ambient concentration of 5 ppm* 2h of continuous work is allowed.*

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EN

- a) *Even after the treatment, the HEPA filters can only be considered to be 'safe to handle when using the appropriate protective clothing' instead of no sterile. They should be treated in an autoclave or incinerated on removal.*
- b) *Suitable protective respiratory equipment designed for being used with hydrogen peroxide should be available in case of an emergency.*
- c) *Once the process has finished, switch off the equipment through the control panel.*

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
EQUIPMENT: BIOPTIMA CABINET 230V		PROJECT REFERENCE: 977083		CUSTOMER: AZBIL TELSTAR TECHNOLOGIES, S.L.	
DRAWING TYPE: ELECTRIC		DRAWINGS: 11		MATERIAL O.F.:	
PROJECT LEADER: Victor Lazaro		AUTOMATION LEADER: X. Cañellas		DATE: 23/04/2015	
APPROVED: X. Cañellas		DATE: 23/04/2015			
REVIEWED: X. Cañellas		DATE: 23/04/2015			
AS BUILT:		DATE:			
REVISIONS:					
ID:U	DATE: 28/05/2019	DRAWN BY: X. Cañellas			
ID:V	DATE: 21/06/2019	DRAWN BY: X. Cañellas			
ID:R	DATE: 20/11/2018	DRAWN BY: X. Cañellas			
ID:S	DATE: 21/01/2019	DRAWN BY: X. Cañellas			
ID:T	DATE: 20/05/2019	DRAWN BY: X. Cañellas			

Drawing	Description	Creation	Modification
01	SUMMARY	23/04/2015	13/02/2018
02	POWER	23/04/2015	21/06/2019
03	CONTROL BOARD	23/04/2015	13/02/2018
04	BLOWERS	23/04/2015	19/02/2018
05	LIGHTING AND UV KIT	23/04/2015	30/05/2018
06	WINDOW MOTOR, SOCKETS & GAS VALVE	23/04/2015	20/11/2018
07	WINDOW DETECTORS	23/04/2015	21/01/2019
08	AIR FLOW SENSORS	23/04/2015	20/02/2018
09	DISPLAY AND BUTTONS	23/04/2015	13/02/2018
10	ACCESSORIES BOARD (OPTIONAL PARTS)	23/04/2015	20/05/2019
11	BILL OF MATERIALS	23/04/2015	21/06/2019

PROPERTY OF TELSTAR REPRODUCTION FORBIDDEN		REPLACES:	REPLACED BY:	PROJECT NAME: 1014747
---	--	------------------	---------------------	--

Date 23/04/2015 Report of Drawings 1014747 Página 1

AZBIL TELSTAR
TECHNOLOGIES, S.L.


ELECTRICAL CHARACTERISTICS

MODEL	BIOPTIMA 6	BIOPTIMA 5	BIOPTIMA 4	BIOPTIMA 3
VOLTAGE	230VAC	230VAC	230VAC	230VAC
POWER	2.54 KVA	2.2 KVA	2.2 KVA	1.89 KVA
CURRENT	11 A	9.6 A	9.6 A	8.2 A
FREQUENCY	50.60 HZ	50.60 HZ	50.60 HZ	50.60 HZ

Revision	Date	Draw by	Date	Project description:	Situation:	Zone:	Drawings:	Previous:	Next:
V	21/06/2019	M.Diaz	23/04/2015	BIOPTIMA CABINET	AZBIL TELSTAR	POWER	2 / 11	1	3
M	13/02/2018	M.Diaz	Draw by: X. Canellas		TECHNOLOGIES, S.L.			1014747	
U	25/05/2019	M.Diaz	Approved: X. Canellas					977083	


CONNECTORS	
J1 - DISPLAY	J10 - SOCKETS 230VAC 4A
J2 - WINDOW MOTOR	J11 - POWER SUPPLY EXHAUST FANS
J3 - ELECTROVALVE	J12 - EXHAUST FAN CONTROL SIGNAL
J4 - LED LIGHTING	J13 - DOWNFLOW FAN CONTROL SIGNAL
J5 - EXHAUST AIR FLOW SENSOR	J14 - POWER SUPPLY DOWNFLOW FANS
J6 - DOWN AIR FLOW SENSOR	J15 - MAIN POWER SUPPLY 230/110VAC
J7 - ACCESSORY BOARD CONNECTION	J16 - BATTERY
J8 - UV LIGHT	J17 - 28.5VDC POWER SUPPLY
J9 - not used	J18 - WINDOW PUSHBUTTONS


FUSES	
F1 - 250VAC 8A/5x20R DC POWER	F5 - 250VAC 1.6A/5x20R LIGHT (not used)
F2 - 250VAC 10A/5x20R IMP FAN	F6 - 250VAC 1A/5x20R UV LIGHT
F3 - 250VAC 2.5A/5x20R EXT FAN	
F4 - 250VAC 5A/5x20R SOCKETS	


Revision	Date	Draw by	Date	Project description:	Situation:	Zone:	Drawings:	Previous:	Next:
K	13/02/2018	M. Diaz	23/04/2015	BIOOPTIMA CABINET	AZBIL TELSTAR TECHNOLOGIES S.L.	CONTROL BOARD	3 / 11	2	4
E	23/05/2016	X. Cateñelas	X. Cateñelas						
J	23/01/2017	M. Diaz	M. Diaz						

Project reference:	977083
Drawing number:	1014747
Drawing number:	3


Revision		Date	Draw by	Date	Project description:	Situation:	Zone:	Drawings:	Previous:	Next:		
I		06/02/2018	M. Diaz	23/04/2015	BIOOPTIMA CABINET	AZBIL TELSTAR TECHNOLOGIES, S.L.	BLOWERS	4 / 11	3	5		
J		19/02/2018	M. Diaz					X. Canelias				
H		15/11/2016	X. Canelias					X. Canelias				
Project description:		BIOOPTIMA CABINET		Situation:		AZBIL TELSTAR TECHNOLOGIES, S.L.		Drawings:		4 / 11		
Zone:		BLOWERS		Drawings:		4 / 11		Previous:		3		
Drawing number:		1014747		Next:		5		Project reference:		977083		
Drawing number:		1014747		Drawing number:		4		Telstar				



Revision	Date	Draw by	Date	Project description:	Situation:	Zone:	Drawings:	Previous:	Next:	
R	20/11/2018	X. Canelias	23/04/2015	BIOPTIMA CABINET	AZBIL TELSTAR TECHNOLOGIES S.L.	WINDOW MOTOR SOCKETS GAS VALVE	6 / 11	5	7	
J	19/02/2018	M.Diaz	X. Canelias				Project reference:	977083	Drawing number:	1014747
O	12/04/2018	X. Canelias	X. Canelias				Project reference:	977083	Drawing number:	6


Revision	Date	Draw by	Date	Project description:	Situation:	Zone:	Drawings:	Previous:	Next:	
F	06/02/2018	M. Diaz	23/04/2015	BIOPTIMA CABINET	AZBIL TELSTAR TECHNOLOGIES, S.L.	WINDOW DETECTORS	7 / 11	6	8	
G	13/02/2018	M. Diaz	Draw by: X. Canelias				Project reference:	977083	Drawing number:	1014747
S	21/01/2019	X. Canelias	Approved: X. Canelias				Project reference:	977083	Drawing number:	7

Revision	Date	Draw by	Date	Project description:	Situation:	Zone:	Drawings:	Previous:	Next:
N	13/02/2018	M.Diaz	23/04/2015	BIOOPTIMA CABINET	AZBIL TELSTAR	DISPLAY AND BUTTONS	9 / 11	8	10
L	06/02/2018	M.Diaz	X. Canelias		TECHNOLOGIES, S.L.				
H	15/11/2016	X. Canelias							

Revision	Date	Drawn by	Date	Project description	Situation	Zone	Drawings	Previous	Next
T	20/05/2019	X. Canellas	23/04/2015	BIOPTIMA CABINET	AZBIL TELSTAR	ACCESSORIES BOARD	10 / 11	9	1014747
G	12/08/2016	X. Canellas	X. Canellas		TECHNOLOGIES S.L.	OPTIONAL PARTS	977083	10	10
H	06/02/2018	M. Diaz							

Symbol	Drawing	Column	Description	Refer. Numbr.	Catalogue	Reference
-F1	02	B	SLOW BLOW FUSE 16A SCHURTER	0001.2516.T	...	616733
-FL1	02	B	POWER LINE FILTER EMC 20A YANBIXINKEJI	YB12E4-20A-Q	...	634174
-I1	02	B	SINGLE PHASE LINE INDUCTOR POLYUX	R12.5	...	644186
-PS1	02	F	POWER SUPPLY 27VDC 200W MEANWELL	RSP-200-27	...	620899
-BT1	02	F	BATERIA 12V 1,3AH BAOTE	BT-12M1.3AT	...	636420
-BT2	02	G	BATERIA 12V 1,3AH BAOTE	BT-12M1.3AT	...	636420
-PCB1	03	C	ATLAS BSC CONTROL CARD	632679
-PCB2	03	F	ATLAS ACCESSORIES EXPANSION CARD	627589
-FAN1	04	C	EBM EC FAN	K3G225-RE07-03	...	611044
-FAN2	04	E	EBM EC FAN	K3G225-RE07-03	...	611044
-FAN3	04	F	EBM EC FAN	K3G225-RE07-03	...	611044
-FAN4	04	G	EBM EC FAN	K3G225-RE07-03	...	611044
-LED1	05	C	LED STRIP (BIOPTIMA 3)	1741339
			LED STRIP (BIOPTIMA 4)
			LED STRIP (BIOPTIMA 5)
			LED STRIP (BIOPTIMA 6)
-BL1	05	G	BALASTRO T5 OSRAM	QTP-OPTIMAL 1X18-40	...	1741347
-L1	05	G	UV TUBE PHILIPS TUV 25W (BIOPTIMA 3)	TUV 25T5 4P SE UNP	...	620816
			UV TUBE PHILIPS TUV 36W (BIOPTIMA 4, 5, 6)	TUV 36T5 4P SE UNP	...	38341
-M1	06	F	DC TUBULAR MOTOR 24VDC	DM445DS-20/15	...	32073
-SK1	06	C	ELECTRIC SOCKET	B-012H	...	630082
-SK2	06	D	ELECTRIC SOCKET	B-012H	...	510032
-SV1	06	F	ELECTROVALVE 24VDC 8W	286293 TYPE 6013	...	640621
-IS1	07	C	INDUCTIVE SENSOR OMBRON	TL-W5MB1	...	64287
-IS2	07	C	INDUCTIVE SENSOR OMBRON	TL-W5MB1	...	64287
-IS3	07	C	INDUCTIVE SENSOR OMBRON	TL-W5MB1	...	64287
-IS4	07	G	INDUCTIVE SENSOR OMBRON	TL-W5MB1	...	64287
-IS5	07	C	INDUCTIVE SENSOR OMBRON	TL-W5MB1	...	64287
-SP1	08	C	AIR VELOCITY SENSOR 0-1m/s 0-5V DEGREE	FS22500	...	614060
-SP2	08	E	AIR VELOCITY SENSOR 0-1m/s 0-5V DEGREE	FS22500	...	614060
-DIS1	09	D	DISPLAY PARA BSC ATLAS	632700
-WB1	09	F	WINDOW BUTTONS ATLAS	1384557

Revision	Date	Draw by	Date	Project description:	Situation:	Zone:	Drawings:	Previous:	Next:
S	21/01/2019	M.Diaz	24/04/2014	BIOPTIMA CABINET	AZBIL TELSTAR	BILL OF MATERIALS	11 / 11	10	
V	21/06/2019	M.Diaz	Draw by: X. Canelias		TECHNOLOGIES, S.L.				
O	12/04/2018	X. Canelias	Approved: X. Canelias						

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EN 1. PURPOSE

Check the integrity of the downflow filters.

2. SCOPE

- BiOptima / BiOptima Cyto

3. TOOLS

- Aerosol generator (cool or hot).
- ATI PAO-4, Shell Ondina fluid or Emery 3004 oils.
- Photometer (ATI 2H).

4. METHODOLOGY

Introduce fume by the generator in the working area (Figure 1a) or in the aspiration chamber (Figure 1b). Concentration in cabinet usually within 30-70µg/l.

- Turn on the photometer and adjust it for measurements in accordance with the manufacturer's instructions (Figure 2)

Fig. 1a (BiOptima)

Fig. 1b (BiOptima Cyto)

Fig. 2

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EN

- With the nozzle of the probe approximately 2.5 cm from the surface, scan the downstream side of the HEPA filters, including the perimeter of each filter pack, by passing the photometer probe in slightly overlapping strokes over the entire surface.

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EN

- Also scan the entire periphery of the filter and the junction between filter and filter mounting frame at a scanning rate that does not exceed 5 cm/s.

Fig. 3

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EN

5. ACCEPTANCE CRITERIA

- The photometer signal after the filter should be smaller or equal to a 0,01% of the specific concentration before the filter.

Fig. 4. Resultado correcto: no fuga
Correct result: filter not leaking

Fig. 5. Resultado incorrecto: fuga
Wrong result: filter leaking

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EN 1. PURPOSE

Prove the **not** existence of leaks in the HEPA exhaust filters.

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EN

- With the nozzle of the probe approximately 2.5 cm from the surface, scan the downstream side of the HEPA filters, by passing the photometer probe in slightly overlapping strokes over the entire surface.

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EN

- Also scan the entire periphery of the filter and the joint between filter and filter mounting frame at a scanning rate that does not exceed 5cm/s.

Fig. 3

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EN 5. ACCEPTANCE CRITERIA

- The photometer signal after the filter should be smaller or equal to a 0.01 % of the specific concentration before the filter.

Fig. 4. Resultado correcto: no fuga
Correct result: there isn't Leak

Fig. 5. Resultado incorrecto: fuga
Wrong result: there is a leak

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EN 1. PURPOSE

- Prove the **not** existence of leaks in the HEPA inflow filters.

2. SCOPE

- BiOptima Cyto

3. TOOLS

- Aerosol generator (cool or hot).
- ATI PAO-4, Shell Ondina fluid or Emery 3004 oils.
- Photometer (ATI 2H).

4. METHODOLOGY

- Connect the photometer sample port to the left side DOP port (Figure 1 & 2).

Fig. 1

Fig. 2

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EN

- Turn on the photometer and adjust it for measurements in accordance with the manufacturer's instructions (Figure 3)

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Fig. 4

EN

Introduce fume by the generator in the left side of the working area (Figure 4). Concentration in cabinet usually within 30-70µg/l.

- *Move slowly the smoke generator to the right until reaching the central part of the cabinet.*
- *Check the photometer display reading during the whole process (Figure 5).*

Fig. 5

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EN

- *In the same way, connect the photometer sample port to the right side DOP port (Figure 6), and introduce smoke into the cabinet from the right side.*
- *Move slowly the smoke generator to the left until reaching the central part of the cabinet.*
- *Check the photometer display reading during the whole process (Figure 5).*

Fig. 6

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EN 5. ACCEPTANCE CRITERIA

- *The photometer signal after the filter should be smaller or equal to a 0.02% of the specific concentration before the filter.*

Fig. 4. Resultado correcto: no fuga
Correct result: filter not leaking

Fig. 5. Resultado incorrecto: fuga
Wrong result: filter leaking

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EN 1. PURPOSE

- Check that the speed control in different points of the working area is suitable to achieve product protection. To determine the speed average as well as the uniformity range in the laminar flow zone.

2. SCOPE

- BiOptima / BiOptima Cyto

3. TOOLS

- The tool to use is an omnidirectional anemometer.

4. METHODOLOGY

- Inside the considered working area according to each cabinet model divide the working area in squares (see attached table). Then switch on the cabinet and wait until the velocity and the temperature becomes stable.


- Place the anemometer at 290mm high of the shelves, over each square (take all readings with the anemometer in the same direction and without covering any holes, whenever possible).

- Take note of the readings and the position of the anemometer in each point.


- For this calculation, we proceed to the following formula:

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


EN CYTO models


- In the following image it is shown the reading points located in a 10x10 cm grid for the downflow test.

Modelo / Model	Puntos / Points
BiOptima Cyto 4	60
BiOptima Cyto 6	90

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EN 5. ACCEPTANCE CRITERIA

- The average of the airflow velocity at the work surface should be ≥ 0.34 m/s and ≤ 0.36 m/s
- Every single value should not differ more than $\pm 20\%$ from the average of the airflow velocity.

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EN Annex A9.4.1. SET POINT ALARM DOWNFLOW AIR VELOCITY TEST

EN 1. PURPOSE

- Check correct work of the alarms when flow speed is lower or higher than the admissible.

2. SCOPE

- BiOptima / BiOptima Cyto

3. TOOLS

- Omnidirectional anemometer

4. METHODOLOGY

- Accede to Technical Service menu through Password Petition screen (User's Manual chapter. 4.5 Point 1), using for it random security password that you are asked for.
- Move to Inflow Test Set Points screen and wait some seconds until the stabilization of the cabinet.
- Choose work area's higher lecture obtained in Downflow velocity Test (Annex A9.4) (avoid near zones to the cabinet walls). Place the anemometer at that point and increase slowly downflow velocity through function button

until "Excessive Laminar Flow" alarm activated. Take note:

- Anemometer velocity.
- Display velocity.

- Repeat the process with lower lecture and decreasing slowly downflow velocity through

the function button until "Insufficient Laminar Flow" activated. Take note:

- Anemometer velocity.
- Display velocity.

EN 5. ACCEPTANCE CRITERIA

- Anemometer reading at point of lowest speed ≥ 0.28 m/s.
- Anemometer reading at point of highest speed ≤ 0.42 m/s.

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EN 1. PURPOSE

- The purpose of this test is to calculate the inflow velocity through the front aperture of the equipment to obtain the Class II to the operator protection.

2. SCOPE


- BiOptima / BiOptima Cyto

3. TOOLS

- Anemometer with a range till 7 m/s.
- Pressure gauge.


4. METHODOLOGY

- Switch on the cabinet and wait until the velocity and the temperature are stable.
- Use the anemometer to measure 9 points divided in 3 horizontal rows.

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
- Probe must be held perpendicular to the direction of the airflow. This can be checked by moving the probe horizontally and vertically. The highest reading shows that the probe is perpendicular. You should also be helped by smoke.
- First row (A) at 10 mm below the top of the working opening.

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

EN

- Second row at the centre of the working opening.

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
- Once you get the 9 values, calculate the average per row and overall average.
- Take note of the display voltage of the exhaust fan and the downflow fans.
- Measure the exhaust plenum pressure.
- Measure the downflow plenum pressure (left port + right port / 2).

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BiOptima Cyto 6: 85 ± 20 Pa

EN 5. ACCEPTANCE CRITERIA

- *The average inflow velocity must fulfill:*

BiOptima 4: 0.65 ± 0.02 m/s

BiOptima 6 / Cyto 4 / Cyto 6: 0.6 ± 0.02 m/s

- *The voltage of the exhaust fan must fulfill:*

BiOptima 4: $6 \pm 0,5$ V

BiOptima 6: 6.2 ± 0.5 V

BiOptima Cyto 4: $7.3 \pm 0,5$ V

BiOptima Cyto 6: 8.3 ± 0.5 V

- *The voltage of the downflow fan must fulfill:*

BiOptima 4: 5.5 ± 0.5 V

BiOptima 6: 5.5 ± 0.5 V

BiOptima Cyto 4: 7.2 ± 0.5 V

BiOptima Cyto 6: 7.1 ± 0.5 V

- *The exhaust plenum pressure must fulfill:*

BiOptima 4: 123 ± 20 Pa

BiOptima 6: 85 ± 20 Pa

BiOptima Cyto 4: 114 ± 20 Pa

BiOptima Cyto 6: 85 ± 20 Pa

- *The downflow plenum pressure must fulfill:*

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EN 1. PURPOSE

- Check the correct operation of inflow air velocity alarm when flow is lower or higher than the admissible.

2. SCOPE


- BiOptima / BiOptima Cyto

3. TOOLS

- Anemometer with a range till 7 m/s.


4. METHODOLOGY

- This test has to be performed after calibration of the cabinet.
- Accede to Technical Service menu through Password Petition screen (User's Manual chapter. 4.5 Point 1), using for it random security password that you are asked for.
- Move to Inflow Set Points Test screen and wait some seconds until the stabilization of the cabinet.
- Decrease slowly downflow velocity through the function button  until "Insufficient inflow" activated. Take note:
 - Anemometer velocity.
 - Display velocity.
- Use the anemometer to measure 9 points divided in 3 horizontal rows.

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

EN

- *Once you get the 9 values, calculate the average per row and overall average.*
- *The same inflow velocity average must be calculated for the alarm "Excessive inflow". Repeat the process increasing slowly the inflow velocity through the function button  until the alarm is activated.*

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EN 5. ACCEPTANCE CRITERIA

"Insufficient inflow" alarm:

- *The inflow average value should be the same as shown in the display ± 0.05 (measuring error).*

Display v.-0.05 < Real inflow < Display v.+0.05

- *Display entry velocity ≥ 0.55 m/s when alarm is activated.*

"Excessive inflow " alarm:

- *The inflow average value should be the same as shown in the display ± 0.05 (measuring error).*

Display v.-0.05 < Real inflow < Display v.+0.05

- *Display entry velocity ≤ 0.72 m/s when alarm is activated.*

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EN Annex A9.6. SMOKE TEST

EN 1. OBJECT

- *Verify the cabinet for personnel and product protection.*

2. SCOPE

- *BiOptima / BiOptima Cyto*

3. TOOLS

- *Smoke generator*

4. METHODOLOGY

- *Product protection: smoke is produced 5-10 cm from the protective glass (Fig. 1 and 2). The behaviour of the smoke is checked visually.*
- *Smoke is produced inside the cabinet, at about 5-10 cm over the lower edge of the glass (Fig. 3). The behaviour of the smoke is checked visually.*
- *User protection: smoke is produced inside the cabinet, as close as possible to the lower edge of the glass (Fig. 4). The behaviour of the smoke is checked visually.*

Fig. 1


Fig. 2

Fig. 3

Fig. 4

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EN 5. ACCEPTATION CRITERIA

Product protection:

- *The generated smoke in the outside of the front aperture will be dragged to the inflow area. No fume sign has to reach the working surface.*
- *Turbulences in the working area shall not be created. The smoke generated should be dragged by the laminar flow to the surface of the working area.*

User protection:

- *No smoke shall be refluxed out of the cabinet once drawn in, guaranteeing the protection of the operator and environment.*

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GUARANTEE CERTIFICATE

AZBIL TELSTAR TECHNOLOGIES, S.L.U. guarantees the substitution of all the defect pieces in the period of TWO YEARS, starting in the date of the expedition of the equipment.

This equipment has undergone satisfactorily all the tests and revisions. Therefore, any damage caused by transport, wrong installation or inappropriate use is excluded of this guarantee.

This guarantee does not include any of the fungible materials comprised in the equipment, such as filters, fluids, filaments, etc.

The inappropriate use of the equipment excludes AZBIL TELSTAR TECHNOLOGIES, S.L.U. of any responsibility for eventual direct or indirect damages of any nature.

Likewise, this guarantee will be void if pieces other than the ones supplied by AZBIL TELSTAR TECHNOLOGIES, S.L.U. are used for the repairing of the equipment.

NOTE. This guarantee is valid for all the equipment unless otherwise stated. In order to obtain any information about this guarantee, please contact with our Sales Department.

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© Model 2H - Digital Photometer