CHROMagar™ ECC

Chromogenic medium for the detection and enumeration of *Escherichia coli* and coliforms

article number EF320 : 1000 ml bottle	NT-EXT-016
article number EF322 : 5000 ml bottle	Version 5
article number EF323 : bulk size	19-Jun-12

STORAGE Store the powder at 15/30°C until the shelflife date indicated on the label.

COMPOSITION in g/L Agar 15; Chromogenic mix, 4,8; Peptone & yeast extract 8; NaCl 5. pH: 7,2 ± 0,2.

(Classical formula adjusted and/or supplemented as required to meet performance criteria).

PREPARATION According to quantities desired, weigh out powder and use in the proportion 32,8 g/L of purified water, or use full pre-weighed dose with corresponding volume of purified water.

Disperse powder slowly in water by rotating for swelling of the agar. Bring to a boil (100°C) by repeated heatings, swirling or stirring regularly. If using an autoclave, do so without pressure. DO NOT HEAT TO MORE THAN 100°C. Mixture may also be brought to a boil in a microwave oven. In this case, after initial boiling, remove from oven to stir gently, return to oven for short repeated heatings. Continue until complete fusion of agar grains (large bubbles replacing foam: about 2 minutes). Cool in a water bath to 48°C.

Note: in case of product samples containing a high load of Pseudomonas and/or Aeromonas, cefsulodin can be added at 7,5 mg/l.

INOCULATION

• **If using pouring technique procedure:** Prepare 90mm Ø sterile Petri dishes and add 1ml of inoculum in each. Then pour 10ml of melted medium. Mix and let solidify. Invert and incubate at 37°C* for 24 hours.

• **If using surface technique procedure:** Pour medium into Petri dishes. Store in the dark before use. Streak the sample or place the inoculated membranes on plate surface and incubate at 37°C* for 24 hours.

* If research is focused on **faecal coliform** bacteria, incubate at 44°C. If research is targeted to maximise **total coliform** detection, incubate at 30°C

Prepared media plates can be kept for one day at ambient temperature. Plates can be stored for up to two weeks under refrigeration (2/8°C) if properly prepared and protected from light and dehydration.

INTERPRETATION

Microorganism → **Typical colony appearance**

<i>E.coli</i>	→ blue
other (faecal*) coliform bacteria	→ mauve
other Gram negative bacteria	→ colourless

PERFORMANCE and LIMITATIONS Sensitivity for *E.coli* is 97% (Ogden *et al.* 1991).

DISPOSAL OF WASTE After interpretation all plates should be destroyed by autoclaving at 121°C for at least 20 minutes.

English. For laboratory use. Laboratory product to be used only by trained personnel.

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Z.A de Gesvrine - 4 rue Képler - B.P.4125
44241 La Chapelle-sur-Erdre Cedex - France
t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00
commercial@humeau.com

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BIBLIOGRAPHY

1991. Ogden *et al.* An evaluation of Fluorogenic and Chromogenic Assays for the Direct Enumeration of *Escherichia coli*. Letters in Applied Microbiology, **13**: 212-215.

Available from CHROMagar :

CHROMagar™ Candida

Differentiation of major pathogenic *Candida* species

CHROMagar™ Orientation

Differentiation of urinary tract pathogens

Rambach™ Agar

Detection of *Salmonella* spp

CHROMagar™ Salmonella

Detection of *Salmonella* including *S. Typhi*

CHROMagar™ Salmonella Plus

Detection of *Salmonella* according to the ISO 6579:2002 norm

CHROMagar™ O157

Detection of *E.coli* O157

CHROMagar™ E.coli

Detection and enumeration of *E.coli*

CHROMagar™ ECC

Detection and enumeration of *E.coli* and coliforms

CHROMagar™ Liquid ECC

Broth for pad technique for *E.coli*-coliforms

CHROMagar™ Staph aureus

Detection and enumeration of *Staphylococcus aureus*

CHROMagar™ MRSA

Detection of MRSA including low level MRSA

CHROMagar™ Listeria

Detection and enumeration of *Listeria monocytogenes*

CHROMagar™ Vibrio

Detection and enumeration of *Vibrio parahaemolyticus*, *Vibrio vulnificus* and *Vibrio cholerae*

CHROMagar™ VRE

Detection of *E.faecium* VRE & *E.faecalis* VRE

CHROMagar™ and Rambach™ are trademarks of Dr. A. Rambach

Visit CHROMagar on internet via <http://www.chromagar.com>

CHROMagar

The Chromogenic Media Pioneer

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NT-EXT-016
Version 5 19-Jun-12