
COMPACTDRY LS

Instructions

Sample Preparation

Viable count in water or liquid foodstuff

1 mL of sample (diluted if necessary) is deposited on the LS CompactDry plate.

Viable count in solid food

Add buffered peptone water to the food sample and homogenize by stomacher® (as per manufacturer instructions). 1 mL of sample (diluted if necessary) is dispensed in the centre on the LS CompactDry plate.

Viable count in swab test specimen

Use the swab to wipe the surface, put the swab back into the tube and screw tightly. Shake the tube to mix. 1 mL of sample (diluted if necessary) is dispensed in the centre on the LS CompactDry plate.

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
CompactDry Protocol

- 1 Open the cap
- 2 Dispense 1 mL of sample on the middle of the CompactDry plate.
- 3 Sample diffuses automatically and evenly across the sheet transforming the dry medium into a gel within seconds.
- 4 Replace the cap, label and turn the plate over.
- 5 Incubate for the required time and temperature.
- 6 Following incubation, count the number of coloured colonies.

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Incubation temperature 35 - 37 +/- 2 °C
Incubation time 24 ± 2 hours
plus a further 24 ± 2 hours

Interpretation

Listeria species form light blue/blue colonies.

General information

Listeria ivanovii tends to grow slowly and Listeria seeligeri could be inhibited by this plate.

A combination of selective agents in the medium inhibits growth of other bacteria.

For effective recovery of Listeria, it is recommended that samples should be inoculated after an 1 hour resuscitation step at 20 °C in peptone buffered water (BPW). This resuscitation step, according to ISO 11290-2, (horizontal method for detection and enumeration of Listeria monocytogenes), does not affect the bacterial number.

The growth area is 20 cm². The back of the plate has a 1 cm² grid to ease colony counting. If there are a large number of colonies, the total viable count can be obtained by multiplying the average number of colonies in several grids by 20.

High concentrations (> 300 CFU) will cause the entire growth area to become blue/green. In this case dilute the sample. After use please follow local disposal regulations.

Storage and shelf life

Store at room temperature (+1 to +30 °C). Shelf life 18 months after manufacturing.

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NISSUI PHARMACEUTICAL CO.,LTD. (Japan)