

Instructions

Storage up to 30°C

24 months

 **Sample Preparation**

Viable count in water or liquid foodstuff

1 mL of sample (diluted if necessary) is deposited on the YM/YMR CompactDry plate.

Viable count in solid food

Add 9 part of buffer solution (0.1 % peptone water or MRD) to 1 part of food sample and homogenize by stomacher® (as per manufacturer instructions). 1 mL of sample (diluted if necessary) is dispensed in the centre on the YM/YMR CompactDry plate.

Viable count in swab test specimen

Use the swab to wipe the surface, put the swab back into the tube and screw tightly. Shake the tube to mix. 1 mL of sample (diluted if necessary) is dispensed in the centre on the YM/YMR CompactDry plate.

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
CompactDry Protocol

1 Open the cap **2** Dispense 1 mL of sample on the centre of the CompactDry plate. **3** Sample diffuses automatically and evenly across the sheet transforming the dry sheet into a gel within seconds. **4** Replace the cap, label and turn the plate over and **5** Incubate for the appropriate time and at required temperature. **6** Following incubation, count the number of coloured colonies.

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Interpretation

Blue coloured colonies indicate yeasts. Mould form cottony colonies with characteristic colours.

General information

Certain yeasts do not form blue colonies.

Antibiotics in the medium inhibit bacterial growth.

The growth area is 20 cm². The back of the plate has a 1 cm² grid to ease colony counting. If there are a large number of colonies, the total viable count can be obtained by multiplying the average number of colonies in several grids by 20.

High concentrations (> 300 CFU) will cause the entire growth area to become green/blue. In this case dilute the sample.

After use please follow local disposal regulations.

Storage and shelf life

Store at room temperature (+1 - +30 °C).
Shelf life 24 months after manufacturing.

PACKAGES	PRODUCT NUMBER
40 plates	54054-0YM-0040
240 plates	54054-0YM-0240
1400 plates	54054-0YM-1400

PACKAGES	PRODUCT NUMBER
40 plates	54063-YMR-0040
240 plates	54063-YMR-0240
1400 plates	54063-YMR-1400

Inkubationstemperatur 25 +/- 1 °C
Inkubationszeit 48 - 72 Stunden (YMR)
3 - 7 days (YM)

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YM: 3-7 days
YMR: 48-72 hours

Interpretation

Yeast : Blue or white to cream colonies
Mold : Cottony colonies with mold colors

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
3D Growth of Molds

Photos Courtesy of NPO Center for Fungal Consultation Japan

NordVal

MicroVal

Distributed in Europe
by Nissui Pharma Solutions