Key Code TSMX5862D
www.thermofisher.com

Europe + 800 135 79 135 US 1 855 236 0910
CA 1 855 805 8539 ROW +31 20 794 7071

Oxoid™ Biochemical Identification System (O.B.I.S.) Salmonella

REF ID0570M..... ▽60

INTENDED USE

The Oxoid Biochemical Identification System (O.B.I.S.) Salmonella is a rapid, qualitative, colorimetric test for the determination of pyroglutamyl aminopeptidase (PYRase) and nitrophenylalanine deaminase (NPA) activity. The test is used for the presumptive identification of *Salmonella* species from plated media. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

PRINCIPLE OF THE TEST

The O.B.I.S. Salmonella Test offers a rapid screening method to distinguish *Salmonella* spp. from those organisms exhibiting similar colonial appearance on common selective Salmonella media. This reduces the need for a full biochemical identification of suspect colonies.

Identification of Salmonella relies on primary isolation of the organism on selective enteric media. However, there are several other genera among the Enterobacteriaceae also capable of growth on such media and which can have similar colony morphology to salmonellae.

The lack of PYRase and NPA activity in *Salmonella* spp. can be used to differentiate them from *Citrobacter* spp. which possess PYRase activity^(1,2,3) and *Proteus*, *Morganella* and *Providencia* spp. which have NPA activity⁽⁴⁾.

The PYRase area on the O.B.I.S. Salmonella Test Card is impregnated with L-pyroglutamic acid 7-amido-4-methylcoumarin⁽⁵⁾. Dimethylaminocinnamaldehyde is used as a colour development reagent. The enzymatic hydrolysis of the substrate produces a purple colour on addition of the O.B.I.S. PYR Developing Solution⁽⁹⁾.

The NPA area on the O.B.I.S. Salmonella Test Card is impregnated with nitrophenylalanine. Deamination of the reagent is shown by an orange-brown colour when the O.B.I.S. NPA Developing Solution (0.25 M sodium hydroxide) is added.

COMPONENTS OF THE O.B.I.S. SALMONELLA KIT (ID0570M)

Each O.B.I.S. Salmonella Kit contains the following reagents sufficient for 60 tests:

ID0571M O.B.I.S. Salmonella Test Cards:

2 resealable plastic pouches, each containing 10 Test Cards and a desiccant.

There are 3 tests on each Test Card, each consisting of 2 reaction areas.

Each reaction area designated as PYR is impregnated with 1% w/v L-pyroglutamic acid 7-amido-4-methylcoumarin.

Each reaction area designated as NPA is impregnated with 5% w/v nitrophenylalanine.

ID0200M O.B.I.S. PYR Developing Solution

1 dropper bottle containing 7 ml of 1 M hydrochloric acid and 0.5% w/v dimethylaminocinnamaldehyde.

ID0210M O.B.I.S. NPA Developing Solution

1 dropper bottle containing 7 ml of 0.25 M sodium hydroxide.

ID0100M O.B.I.S. Buffer Solution

1 dropper bottle containing 7 ml of Phosphate Buffered Saline (PBS).

Instructions for Use (IFU).

MATERIALS REQUIRED BUT NOT PROVIDED

Microbiological loop.

Positive and negative quality control organisms.

PRECAUTIONS

This product is for *in vitro* diagnostic use only.

Do not use O.B.I.S. Salmonella reagents beyond the stated expiry date.

Specimen material may contain pathogenic organisms, handle with appropriate precautions.

The O.B.I.S. PYR Developing Solution (ID0200M) contains an acid.

The O.B.I.S. NPA Developing Solution (ID0210M) contains a weak base. Avoid direct contact by wearing suitable protective equipment. If the reagents come into contact with the skin, mucous membranes or eyes immediately wash the area thoroughly with plenty of water.

Used O.B.I.S. Test Cards and microbiological loops should be disposed of as biohazardous waste and incinerated or autoclaved for 15 minutes at 121°C.

Please refer to the Safety Data Sheet (SDS) on company website and product labelling for information on potentially hazardous components.

Directions should be read and followed carefully.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

STORAGE AND OPENING

The O.B.I.S. Salmonella Kit must be stored at 2°–8°C. Allow the pouches to equilibrate to room temperature before opening to prevent the formation of condensation on the Test Cards.

Open the pouches by cutting at the notch between the end seal and the zip lock opening.

Once opened, remove the number of Test Cards required for immediate testing (testing within next 10 minutes) and reseal the pouch straight away.

Protect the Test Cards from light

If fewer tests are required, cut the Test Card and return the unused portions to the pouch. Do not return used portions to the pouch as they will be contaminated.

When stored as indicated, O.B.I.S. Salmonella reagents will retain their activity until the expiry date shown on the box.

QUALITY CONTROL PROCEDURE

Each day the Kit is used the following control procedures should be performed:

PYRase Test

- Positive control** – Use a known pyroglutamyl aminopeptidase positive strain such as *Citrobacter freundii* ATCC™8090 (Thermo Scientific™ Culti-Loops™ R4601800). Follow the method given in the test procedure.

Ensure that a purple colour forms within 20 seconds.

- Negative control** – Use a known pyroglutamyl aminopeptidase negative strain such as *Salmonella enterica* subsp. *enterica* serovar Enteritidis ATCC® 13076™ (Thermo Scientific™ Culti-Loops™ R4608200). Follow the method given in the test procedure.

Ensure that no purple colour forms within 20 seconds.

NPA Test

- Positive control** – Use a known phenylalanine deaminase positive strain such as *Proteus mirabilis* ATCC™43071 (Thermo Scientific™ Culti-Loops™ R4607073). Follow the method given in the test procedure.

Ensure that an orange-brown colour forms within 20 seconds.

- Negative control** – Use a known phenylalanine deaminase negative strain such as *Citrobacter freundii* ATCC™8090 (Thermo Scientific™ Culti-Loops™ R4601800) or *Salmonella* ser enteritidis ATCC™13076 (Thermo Scientific™ Culti-Loops™ R4608200). Follow the method given in the test procedure.

Ensure that no orange-brown colour forms within 20 seconds.

Do not use the reagents if reactions with control organisms are incorrect.

SPECIMENS

Colonies with a typical Salmonella morphology on the following media may be tested: Brilliant Green Agar (Modified) (CM0329), MLCB Agar (CM0783), XLD Medium (CM0469), SS Agar (CM0099), SS Agar Modified (CM0533), Hektoen Enteric Agar (CM0419)*, Desoxycholate Citrate Agar (CM0035) or Bismuth Sulphite Agar (CM0201).

Colonies tested should be Gram-negative and Oxidase negative.

Plates should be tested a maximum of one hour after being removed from the incubator.

*Modified test procedure required – see below.

TEST PROCEDURE AND INTERPRETATION OF RESULTS

Important Note: nitrophenylalanine impregnated on the NPA Test area is light sensitive. Exposure to light for more than 1 hour must be avoided.

- Select one suspect colony (0.5 mm or larger). Apply the colony to the PYR and NPA Test areas (sufficient colony should be applied to each area to make a visible smear).
- Moisten both Test areas with 1 drop of Buffer Solution (ID0100M).
- Incubate the inoculated Test Card at room temperature (15–30°C) for 5 minutes.
- Dispense 1 drop of the O.B.I.S. PYR Developing Solution (ID0200M) to the PYR Test area. Development of a vivid purple colour on and around the smeared colony within 20 seconds confirms PYRase activity.
- Dispense 1 drop of the O.B.I.S. NPA Developing Solution (ID0210M) to the NPA Test area. Development of an orange-brown colour on and around the smeared colony within 20 seconds confirms NPA activity.

Modified Test Procedure – Hektoen Enteric Agar

- Apply one suspect colony (0.5 mm or larger), and inoculate both PYR and NPA Test areas using a plastic loop. Sufficient colony should be applied to each area to make a visible smear.
- Moisten both Test areas with one drop of Buffer Solution (ID0100M).
- Place the inoculated card into a suitable plastic bag/sleeve (approximately 75mm X 100mm).
- Incubate at 37°C ± 2°C for 30 minutes.
- Dispense one drop of O.B.I.S. PYR Developing Solution (ID0200M) onto the PYR Test area. Development of a vivid purple colour on and around the smeared colony within 20 seconds denotes a POSITIVE reaction. (The purple colour may be slightly obscured by the black pigmentation of ferrous sulfide if a black colony is selected.)
- A NEGATIVE reaction is denoted by no change in colour of the inoculum smear.
- Development of a greenish colour within one minute is considered a NEGATIVE PYR reaction and is due to the presence of indole. These organisms are not *Salmonella* spp.
- Dispense one drop of NPA Developing Solution (ID0210M) onto the NPA Test area. Development of a vivid orange-brown colour on and around the smeared colony within 20 seconds denotes a POSITIVE reaction. (The orange-brown colour may be slightly obscured by the black pigmentation of ferrous sulphide if a black colony is selected.)
- No change in colour of the smear, or intensification of the green-blue colour due to acidification of the bromothymol blue in the colony, indicates a NEGATIVE NPA reaction.

Interpretation chart

Typical organism reactions:

Organism	PYRase	NPA
<i>Salmonella</i> spp.	–	–
<i>Citrobacter</i> spp.	+	–
<i>Proteus</i> , <i>Morganella</i> and <i>Providencia</i> spp.	–	+

Trial Results (Data on File)

Results from pure culture:

Organism	N ^o of Specimens Tested	Positive PYRase Reaction	Positive NPA Reaction
<i>Salmonella</i> spp.	150	0	0
<i>Citrobacter</i> spp.	56	55*	0
<i>Proteus</i> , <i>Morganella</i> and <i>Providencia</i> spp.	50	0	50

* One *Citrobacter* strain gave a negative PYRase reaction. This was confirmed as *Citrobacter youngae* thus showing a truly negative result.

In a trial of 470 faecal samples O.B.I.S. *Salmonella* gave the following results:

Sensitivity 100%

Specificity 81.5%–95.4% †

† Depending on the plating medium used.

LIMITATIONS OF THE TEST

O.B.I.S. *Salmonella* is intended for the detection of PYRase and NPA activity in Gram-negative, Oxidase negative microorganisms. It can be used as a screen to differentiate *Salmonella* spp. from *Citrobacter*, *Proteus*, *Providencia* and *Morganella* spp. from the following selective enteric media: Brilliant Green Agar (Modified) (CM0329), MLCB Agar (CM0783), XLD Medium (CM0469), Bismuth Sulphite Agar (CM0201), SS Agar (CM0099), SS Agar Modified (CM0533) and Hektoen Enteric Agar (CM0419) or Desoxycholate Citrate Agar (CM0035).

A few isolates of H₂S-positive *Escherchia coli* may appear as salmonellae on MLCB Agar. These organisms will show the same results as *Salmonella* spp. on the O.B.I.S. *Salmonella* Test as they also lack both PYRase and NPA activity.

The neutral red in Desoxycholate Citrate Agar (CM0035) may produce a pink colour on the O.B.I.S. Test Card. To avoid any confusion with the positive purple reaction in the PYR test the result should be compared to a positive control strain.

E. coli and indole-positive *Proteus* spp. obtained from media with a high tryptophan content may generate a blue-green colour development on the PYR Test area. This is a negative result.

Citrobacter youngae may produce a PYRase-negative result and show the same result as *Salmonella* spp.

O.B.I.S. salmonella provides a reliable presumptive identification of salmonellae, but does not replace the need for full biochemical testing.

REFERENCES

1. Chagla, A. H., Borczyk, A. A., Aldom, J. E., Dalla Rosa, S. and Cole, D. D. (1993). Evaluation of the L-pyrrolidonyl-beta-naphthylamide hydrolysis test for differentiation of members of the families Enterobacteriaceae and Vibrionaceae. *J. Clin. Microbiol.* 31. 1946–1948.
2. Mulczyk, M. and Szewczuk, A. (1970). Pyrrolidonyl Peptidase in Bacteria: A New Colorimetric Test for Differentiation of Enterobacteriaceae. *J. Gen. Microbiol.* 61. 9–13.
3. Inoue, K., Miki, K., Tamura, K. and Sakazaki, R. (1996). Evaluation of L-Pyrrolidonyl Peptidase Paper Strip Test for Dif-ferentiation of Members of the Family Enterobacteriaceae, Particularly *Salmonella* spp. *J. Clin. Microbiol.* 34. 1811–1812.
4. Giammanco, G., Pignato, S. and Agodi, A. (1985). A Simple Chromogenic Test For Rapid Screening of *Proteus* and *Providencia* Bacteria. *Microbiologica.* 8. 395–397.
5. Druggan, P., Roberts, P. and Swaine, D. (1999). A Rapid Chromogenic Method for the Differentiation of *Citrobacter* spp. and *Salmonella* spp. Directly from Enteric Media. abstr. C-444 p71 Abstr. 99th Annual Meet. Am. Soc. Microbiol. 1999

SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Batch Code (Lot Number)
	Temperature Limitations (Storage temp.)
	Use By (Expiration Date)
	Do Not Reuse
	Consult Instructions for Use (IFU)
	Do not use if package is damaged
	Contains sufficient for <n> tests
	Manufactured by
	European Authorised Representative
	UK Conformity Assessed
	CE Mark
	Protect from light

IFU X5862D 2022-09-05 Printed in the UK

Oxoid Ltd, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

©2022 Thermo Fisher Scientific Inc. All rights reserved.

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Schlüssel-Code TSMX5862D
www.thermofisher.com

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Oxoid Ltd, Wade Road, Basingstoke, Hampshire, RG24 8PW, Reino Unido
©2022 Thermo Fisher Scientific Inc. Reservados todos los derechos.
Todas las demás marcas comerciales son propiedad de Thermo Fisher Scientific Inc. y sus subsidiarias.