Sabouraud Dextrose Agar

[REF] CM0041B, CM0041K, CM0041M, CM0041R, CM0041T

EN

*This Instructions for Use (IFU) is intended to be read in conjunction with the IFU for Chloramphenicol Selective Supplements (SR0078E, SR0785H).

Intended Use

Sabouraud Dextrose Agar devices (CM0041B, CM0041K, CM0041M, CM0041R, CM0041T) devices are acidic pH media for the isolation and maintenance of dermatophytes, other fungi and yeasts from skin, hair, nails, genital, respiratory and urine samples from patients. Sabouraud Dextrose Agar (CM0041B, CM0041K, CM0041M, CM0041R, CM0041T) devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having microbial infections.

The devices are for professional use only, are not automated, nor are a companion diagnostic.

These devices can also be used with non clinical samples (environmental, food etc.)

Summary and Explanation

Sabouraud Dextrose Agar is suitable for cultivation and differentiation of fungi and may be used in the place of the American medium Hodges.¹ The fungi maintain their typical cultural appearance and thus may be readily identified according to the standard macroscopic characters described by Sabouraud.² The medium is often used with antibiotics for the isolation of pathogenic fungi and *Candida albicans* from material containing large numbers of other fungi and bacteria.³⁻⁸

Principle of Method

Sabouraud Dextrose Agar contains mycological peptone and glucose to support growth of fungi. It is slightly acidic which helps to inhibit bacteria whilst being favourable for fungi. The addition of antimicrobial agents can increase the selectivity of the medium. For example, chloramphenicol is a bacteriostatic antimicrobial with a broad spectrum of activity against both Gram-positive and Gram-negative bacteria. Agar is added as a solidifying agent.

Typical Formula*

CM0041B, CM0041K, CM0041M, CM0041R, CM0041T

	grams per litre
Mycological peptone	10.0
Glucose	40.0
Agar	15.0

*Adjusted as required to meet performance.

Materials Provided

Product name	Product code	
Sabouraud Dextrose Agar	CM0041B	500g
Sabouraud Dextrose Agar	CM0041R	2.5kg
Sabouraud Dextrose Agar	CM0041M	8x1.69Kg
Sabouraud Dextrose Agar	CM0041T	5kg
Sabouraud Dextrose Agar	CM0041K	25kg

Materials Required but Not Supplied

Inoculating loops.

- Swabs.
- Collection containers
- Incubators
- Quality control organisms

Storage

Store dehydrated media in its original packaging between 10°C and 30°C. Keep out of direct sunlight. Keep tightly closed.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical

advice/attention.

- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com). Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) Q5, B 11, , B 28, B 39, B 41, B 57 and ID 1.

Procedure

1. Suspend Sabouraud Dextrose Agar (65.0 grams / 1 litre) in distilled water.
2. Bring to a boil to dissolve completely.
3. Sterilize by autoclaving at 121°C for 15 minutes.
4. Mix well and pour into sterile Petri dishes.

Interpretation

Cream, domed colonies (2 – 6mm) indicate the presence of target organisms, including *Candida albicans* or *Saccharomyces carlsbergensis*

White mycelia, black spores (>10mm) indicate the presence of *Aspergillus brasiliensis*.

Quality Control

Incubation Conditions: 20 - 25°C for up to 120 hours Colony Count: ≥ 70% of the control medium count Inoculum Level: 10 – 100 cfu	
<i>Saccharomyces carlsbergensis</i> ATCC®2700™ <i>Candida albicans</i> ATCC®10231	2-6mm cream, domed colonies
<i>Aspergillus brasiliensis</i> ATCC®16404™	Greater than 10mm colonies, white mycelia, black spores
Testing performed in accordance with ISO11133:2014 Reactions after incubation at 25 ± 2°C for 120 hours Inoculum Level: 50 – 120 cfu A satisfactory result is represented by recovery of positive strains equal to or greater than 70% of the control medium.	
<i>Saccharomyces cerevisiae</i> ATCC® 9763™	2-6mm cream, domed colonies
<i>Aspergillus brasiliensis</i> ATCC® 16404™	Greater than 10mm colonies, white mycelia, black spores

Limitations

Some strains of dermatophytes, other fungi and yeasts may grow poorly or not at all. The pH of the medium is favourable for fungi, but many bacteria will be able to grow. The choice of culture medium, for example with the addition of antimicrobials, and incubation protocol will depend on the specimen type.

Some of the pathogenic fungi may produce infective spores which are easily dispersed into the laboratory. Such organisms should be examined only within a protective cabinet.

The combination of cycloheximide and chloramphenicol that can be added to Sabouraud Dextrose Agar (CM0041B, CM0041K, CM0041M, CM0041R, CM0041T) devices inhibit many pathogenic fungi. However, the mycelial phase of *Histoplasma capsulatum*, *Paracoccidioides brasiliensis*, *Sporothrix schoenckii* and *Blastomyces dermatitidis* is not inhibited by these antibiotics when incubated at 25-30°C.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct isolation and maintenance of dermatophytes, other fungi and yeasts from clinical samples is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Sabouraud Dextrose Agar (CM0041) devices was demonstrated by an overall pass rate of 80% obtained for the 10 most recently produced batches. An overall pass rate of 80% (based on data from the most recent 10 batches) provides sufficient data on accuracy considering the low risk of the device and it is not reflective of the product that is released for clinical use. With the normal variation in the manufacturing process, a small number of failures may occur, however, there is no risk to the end user as batches will not be released for sale if out of specification. This shows that the performance is reproducible, thus meeting the specifications set out in the Performance Evaluation Plan. These media devices have been tested in-house as part of the QC process since they were first manufactured. Sabouraud Dextrose Agar devices (CM0041B, CM0041R, CM0041T and CM0041K) have been on the market for at least 20 years. The user can grow microorganisms that meet the defined acceptance criteria by following the incubation conditions described in the IFU, product specification and inspection plan.

Bibliography

1. Hodges R. S. (1928) *Arch. Derm. Syph., New York*, 18. 852.
2. Sabouraud R/ (1910) '*Les Teignes*', Masson, Paris.
3. Georg Lucille K., Ajello L. and Papageorge Calomira (1954) *J. Lab. Clin. Med.* 44. 422-428.
4. Ajello Libero (1957) *J. Chron. Dis.* 5. 545-551.
5. Williams Smith H. and Jones J.E.T. (1963) *J. Path. Bact.* 86. 387-412.
6. Hantschke D. (1968) *Mykosen.* 11. 113-115.
7. Dolan C.T. (1971) *Appl. Microbiol.* 21. 195-197.
8. Pagnano J., Levin J.G. and Trejo W. (1957-58) *Antibiotics Annual* 1957-58, 137-143.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limitation (storage temp.)
	Use by (expiration date) YYYY-MM
	Keep away from sunlight
	Consult Instructions for Use
	Do not use if packaging damaged
	Manufacturer
	Authorized representative in the European Community
	European Conformity Assessment
	UK Conformity Assessment

	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2024-03-19

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien