# CAMPYLOBACTER Blaser Wang Supplement

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Rev.0 / 06.04.2005

# CAMPYLOBACTER Blaser Wang Supplement

Selective supplement for the isolation of *Campylobacter spp*

## DESCRIPTION

CAMPYLOBACTER Blaser Wang Supplement is a selective supplement for the isolation of *Campylobacter spp* comprising a freeze-dried mixture of Polymyxin B, Trimethoprim, Amphotericin B, Cephalothin and Vancomycin. CAMPYLOBACTER Blaser Wang Supplement is used for selective enrichment of the medium CAMPYLOBACTER AGAR BASE code 611007 or 621007.

## KIT CONTENTS

Each kit contains:

- 10 bottles of freeze-dried CAMPYLOBACTER Blaser Wang Supplement
- 1 Instruction sheet

## PRINCIPLE OF THE METHOD

CAMPYLOBACTER Blaser Wang Supplement is derived from Skirrow's formulation with the addition of Amphotericin B and Cephalothin. Polymyxin B is a broad-spectrum antibiotic, Trimethoprim interferes with the synthesis of nucleic acids and proteins and has a bacteriostatic action on the majority of both Gram-positives and Gram-negatives, Cephalothin inhibits both Gram-positive and Gram-negative bacteria, Amphotericin B is active against fungi, Vancomycin inhibits synthesis of the wall of Gram-positive bacteria.

## COMPOSITION

CAMPYLOBACTER Blaser Wang Supplement		
Contents / bottle		Contents / l of medium
Polymixin B	1250 IU	2500 IU
Trimethoprim	2.5 mg	5.0 mg
Amphotericin B	1.0 mg	2.0 mg
Cephalothin	7.5 mg	15.0 mg
Vancomycina	5.0 mg	10.0 mg

## PROCEDURE FOR USE

1. Aseptically reconstitute the contents of a bottle of CAMPYLOBACTER Blaser Wang Supplement with 4 ml of a solution of sterile distilled water and ethyl alcohol in the ratio 1: 1. Shake until completely dissolved, avoiding foam formation.
2. Aseptically add the entire contents of a bottle (4 ml) to 500 ml of Campylobacter Agar Base medium code 611007-621007 autoclaved and cooled to 45-50 °C with the addition of 10% of defibrinated ram's blood or with 5-7% of laked horse blood and with the contents of a bottle of Campylobacter growth Supplement.
3. Mix with care.
4. Distribute into Petri dishes.

## TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation for CAMPYLOBACTER AGAR BASE code 611007 or 621007.

## QUALITY CONTROL

1. Control of the appearance: freeze-dried product, colour orange.
2. Microbiological control.

Prepare the plates using as base CAMPYLOBACTER AGAR BASE code 611007 or 621007 supplemented with 10% of ram's blood, with CAMPYLOBACTER Growth Supplement (1 bottle in 500 ml of medium) and with CAMPYLOBACTER Blaser Wang Supplement (1 bottle in 500 ml of medium). The plates are seeded with the strains indicated in the microbiological control table.

Incubation conditions: 48 h at 37 °C, in a microaerophilic atmosphere for Campylobacter.

Microbiological control

Control strains		Growth
<i>Campylobacter jejuni</i>	ATCC 33291	Good
<i>Escherichia coli</i>	ATCC 25922	Inhibited

## PRECAUTIONS

The product CAMPYLOBACTER Blaser Wang Supplement is classifiable as hazardous under current legislation; it is recommended that the Safety Data Sheet be consulted on its use.

CAMPYLOBACTER Blaser Wang Supplement is a selective supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

## STORAGE

Store CAMPYLOBACTER Blaser Wang Supplement at 2-8°C in its original packaging. In such conditions CAMPYLOBACTER Blaser Wang Supplement will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

## REFERENCES

- Skirrow, M.D. 1977 Campylobacter enteritidis: A New Disease. Br. Med. J. 2:9-11.
- Blaser, M.J., V. Berkowitz, F.M. Laforce, 1979. Campylobacter enteritidis: Clinical and Epidemiologic Features. Ann. Intern. Med: 91:179-185.
- Association of Official Analytical Chemists. 1995. Bacteriological analytical manual, 8th Ed.

## PRESENTATION

product	REF	$\Sigma$
CAMPYLOBACTER Blaser Wang Supplement	81051	10 bottles

One bottle is sufficient to prepare 500 ml of medium

## TABLE OF SYMBOLS

<b>IVD</b> In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	$\Sigma$ Contains sufficient for <n> tests	 Temperature limitation
<b>REF</b> Catalogue number	 Fragile, maneggiare con cura	 Use by	 Caution, consult accompanying documents	<b>LOT</b> Batch code

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