

H₂O₂ Reagent

Reagent for the catalase test

DESCRIPTION

H₂O₂ Reagent is a reagent for microbiological use made of a H₂O₂ solution, for the performance of the catalase test.

CONTENT OF THE PACKAGES

Each package contains:

- 1 vial containing 5 mL of 7.5% H₂O₂ solution.
- 1 instructions sheet.

ITEMS NECESSARY BUT NOT INCLUDED IN THE PACKAGES

- CATALASI/OXY TEST (ref. 88023)
- PHYSIOLOGICAL SOLUTION (ref. 20095)
- Various material for laboratory of microbiology

PRINCIPLE OF THE METHOD

The catalase is an enzyme that resolves H₂O₂ in O₂ and H₂O. With the exception of Streptococci, most of aerobe and facultative anaerobe microorganisms have catalase activity. The production of catalase enzyme in the microorganism in examination is evidenced by the immediate bubbles formation after the addition of **H₂O₂ Reagent** to a suspension of the microorganism. The absence of bubbles formation indicates a microorganism unable to resolve H₂O₂, therefore lacking in catalase enzyme.

COMPOSITION

H ₂ O ₂ Reagent	(Content / vial):
H ₂ O ₂	0.375 mL
Distilled water	4.625 mL

TEST PROCEDURE

1. Reconstitute one tube of CATALASI/OXY TEST (ref. 88023) with 0.2 mL di PHYSIOLOGICAL SOLUTION (ref. 20095).
2. Dilute one colony well isolated from culture medium without blood by a sterile loop into the tube of CATALASI/OXY TEST (ref. 88023). Shake gently.
3. Add 4-5 drops of **H₂O₂ Reagent**.
4. Observe the almost immediate development of gas bubbles and the foam formation in the tube.

INTERPRETATION OF THE RESULTS

Interpret results according to the table n°1:

Table n°1

	Catalase test
Bubbles and foam formation	Positive
No bubbles and foam formation	Negative

QUALITY CONTROL FOR THE USER

Appearance control: limpid and colourless solution.

Microbiological control: every batch of **H₂O₂ Reagent** is subjected to quality control, using a bacterial culture of *Pseudomonas aeruginosa* ATCC 27853 as positive control and one of *Enterococcus faecalis* ATCC 29212 as negative control.

PRECAUTIONS

The product **H₂O₂ Reagent** is classifiable as hazardous under current legislation; it is recommended that the Safety Data Sheet be consulted on its use. **H₂O₂ Reagent** must be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store **H₂O₂ Reagent** at 2-8°C in its original package. In such conditions the product **H₂O₂ Reagent** is valid until the expiry date shown on the label. Do not use them beyond that date. Dispose of them if they show signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, used **H₂O₂ Reagent** and the material that has come into contact with the sample must be decontaminated and disposed of in accordance with the laboratory procedures for the decontamination and disposal of potentially infected material.

REFERENCES

1. Murray, Baron, Pfaller, Tenorev and Tenover: *Manual of Clinical Microbiology* (1995).
2. Bayley and Scott's: *Diagnostic Microbiology* (1986).
3. Edwin H.Lenette: *Manual of Clinical Microbiology* (1995).

PRESENTATION












Product	REF	
H ₂ O ₂ Reagent	80057	1 vial x 5 mL

TABLE OF SYMBOLS

 In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	 Batch code



LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY

Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@lioilchem.net



Rev.3 / 14.02.2013