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# CAMPYLOBACTER Preston Supplement

Selective supplement for the enrichment of CAMPYLOBACTER AGAR BASE medium for the isolation of *Campylobacter* spp.

## DESCRIPTION

CAMPYLOBACTER Preston Supplement is a selective supplement for the isolation of *Campylobacter* spp. made of a freeze-dried mixture of Polymyxin B, Trimethoprim, Rifampicin and Cycloheximide. CAMPYLOBACTER Preston Supplement is used for the selective enrichment of medium CAMPYLOBACTER AGAR BASE ref. 611007 or 621007.

## KIT CONTENTS

Each kit contains:

- 10 bottles of freeze-dried CAMPYLOBACTER Preston Supplement
- 1 Instruction sheet

## PRINCIPLE OF THE METHOD

CAMPYLOBACTER Preston Supplement is based on the formulation described by Bolton and Robertson, that permits an efficacious isolation of *Campylobacter* species from all types of samples (human, animal, and environmental). Polymyxin B and Rifampicin are broad-spectrum antibiotics, Trimethoprim interferes with the synthesis of nucleic acids and proteins and has a bacteriostatic action on the majority of both Gram-positives and Gram-negatives, la Cycloheximide is active against mycetes.

## COMPOSITION

CAMPYLOBACTER Preston Supplement		
	Contents / bottle	Contents / l of medium
Polymyxin B	2500 IU	5000 IU
Trimethoprim	5.0 mg	10.0 mg
Rifampicin	5.0 mg	10.0 mg
Cycloheximide	50.0 mg	100.0 mg

## PROCEDURE FOR USE

1. Aseptically reconstitute the contents of one bottle of CAMPYLOBACTER Preston Supplement with 2 ml of a solution of sterile distilled water and acetone in the ratio 1:1. Shake until completely dissolved, avoiding foam formation.
2. Aseptically add the entire contents of one bottle (2 ml) to 500 ml of *Campylobacter* Agar Base ref. 611007-621007 autoclaved, cooled at 45-50°C and with the addition of 5% of lysed horse blood and with the content of one bottle of *Campylobacter* growth Supplement.
3. Mix with care.
4. Distribute into Petri dishes.

## TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation for CAMPYLOBACTER AGAR BASE ref. 611007 or 621007.

## QUALITY CONTROL

1. Control of the appearance: freeze-dried product, colour orange.
2. Microbiological control.

Prepare the plates using as base CAMPYLOBACTER AGAR BASE ref. 611007 or 621007 supplemented with 5% of lysed horse blood, with CAMPYLOBACTER Growth Supplement (1 bottle in 500 ml of medium) and with CAMPYLOBACTER Preston Supplement (1 bottle in 500 ml of medium). The plates are inoculated with the strains indicated in the microbiological control table.

Incubation conditions: 48 h a 42°C, in atmosphere of microaerophilia for *Campylobacter*.

Microbiological control

Control strains		Growth
<i>Campylobacter jejuni</i>	ATCC 33291	Good
<i>Escherichia coli</i>	ATCC 25922	Inhibited

## PRECAUTIONS

The product CAMPYLOBACTER Preston Supplement is classifiable as hazardous under current legislation; it is recommended that the Safety Data Sheet be consulted on its use. CAMPYLOBACTER Preston Supplement is a selective supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

## STORAGE

Store CAMPYLOBACTER Preston Supplement at 2-8°C in its original packaging. In such conditions CAMPYLOBACTER Preston Supplement will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

## REFERENCES

- Bolton, F.J. and Robertson, L. (1982) J. Clin. Pathol. 35: 462-467.
- Bolton, F.J., Coates, D., Hinchliffe, P.M. and Robertson, L. (1983) J. Clin. Pathol. 36: 78-83.
- Bolton, F.J., Coates, D. and Hutchinson, D.N. (1984) J. Appl. Bact. 56: 151-157.

## PRESENTATION

Product	REF	$\Sigma$
CAMPYLOBACTER Preston Supplement	81004	10 bottles

One bottle is sufficient to prepare 500 ml of medium.

## TABLE OF SYMBOLS

<b>IVD</b> In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
<b>REF</b> Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	<b>LOT</b> Batch code

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