BD Difco™ Salmonella O Antisera
Difco™ Salmonella H Antisera
Difco™ Salmonella Antiserum Vi

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INTENDED USE

BD Difco™ Salmonella O Antisera are used in slide agglutination tests for the identification of *Salmonella* by somatic (O) antigens.

BD Difco Salmonella H Antisera are used in tube agglutination tests for the identification of *Salmonella* by flagellar (H) antigens.

BD Difco Salmonella Vi Antiserum is used in slide agglutination tests for the identification of *Salmonella* Vi.

SUMMARY AND EXPLANATION

Salmonella species cause a variety of human diseases called salmonellosis. The range of disease is from mild self-limiting gastroenteritis to more severe forms, possibly with bacteremia or typhoid fever, which can be life threatening. Severe disease and bacteremia are associated primarily with three serovars of *S. enterica* subsp. *enterica* (Choleraesuis, Paratyphi A and Typhi), while most of the other 2,300 or more strains are associated with gastroenteritis. The severity of the diarrheal disease depends upon the virulence of the strain and the condition of the human host.

Salmonella is found in nature and occurs in the intestinal tract of many animals, both wild and domestic. The microorganism can spread to man from contact with the environment or from eating meat or vegetable food products.

All *Salmonella* serovars belong to two species: *S. bongori*, which contains 18 serovars, and *S. enterica*, which contains the remaining 2,300 or more serovars divided among six subspecies.^{1,2}

The six subspecies of *S. enterica* are:

- S. enterica* subsp. *enterica* (I or 1)
- S. enterica* subsp. *salamae* (II or 2)
- S. enterica* subsp. *arizonae* (IIIa or 3a)
- S. enterica* subsp. *diarizonae* (IIIb or 3b)
- S. enterica* subsp. *houtenae* (IV or 4)
- S. enterica* subsp. *indica* (VI or 6)

Nomenclature and classification of these bacteria are constantly changing.³ *Salmonella* and the former *Arizona* should be considered a single genus, *Salmonella*.⁴ It is recommended that laboratories report the names of *Salmonella* serovars for the subspecies *enterica*. The serovar names are no longer italicized and the first letter is capitalized. For example, the strain that used to be identified as *Salmonella typhimurium* is now known as *Salmonella* Typhimurium.

Serovars of other subspecies of *S. enterica* (except some in the subspecies *salamae* and *houtenae*) and those of *S. bongori* are not named and are designated by their antigenic formula. For the most recent information on nomenclature, consult appropriate references.¹⁻¹⁰

PRINCIPLES OF THE PROCEDURE

Salmonella O antigens are somatic (O) heat-stable antigens and are identified first. The Vi antigen is a heat-labile envelope antigen that may surround a cell wall and mask somatic antigen activity. Microorganisms having the Vi Antigen will not agglutinate in O antisera. In order to determine the O antigen of these cultures, a suspension of the organism must be boiled to destroy the heat-labile envelope antigen and then tested with O antisera. The flagellar (H) antigens are heat labile and are usually associated with motility.

Complete serological characterization of *Salmonella* is not required for successful detection of the microorganism when it occurs as a pathogen. The use of adequate isolation procedures and differential biochemical tests is of primary importance. Because antigenic relationships exist between genera of the family *Enterobacteriaceae*, it is recommended that the isolate be biochemically identified as *Salmonella* prior to Serology testing. Possible *Salmonella* isolates can be presumptively identified with a minimum of serological identification. Isolates can be sent to laboratories that perform the level of testing necessary to completely identify the microorganism.

For a further discussion of the serological identification of *Salmonella*, consult appropriate references.^{1-3,9,11-14}

Identification of *Salmonella* species includes both biochemical and serological identification. Serological confirmation involves the procedure in which the microorganism (antigen) reacts with its corresponding antibody. This *in vitro* reaction produces macroscopic clumping called agglutination. The desired homologous reaction is rapid, does not dissociate (high avidity) and bonds strongly (high affinity).

Because a microorganism (antigen) may agglutinate with an antibody produced in response to another species, heterologous reactions are possible. Such unexpected and perhaps unpredictable reactions may lead to some confusion in serological identification. Therefore, a positive homologous agglutination reaction should support the morphological and biochemical identification of the microorganism.

Agglutination of the somatic antigen in the slide test appears as a firm granular clumping. Homologous reactions are rapid and strong (3+). Heterologous reactions are slow and weak.

Agglutination of flagellar antigens in the tube test appears as a loose flocculation that can easily be resuspended.

REAGENTS

BD Difco Salmonella O, H and Vi Antisera are lyophilized, polyclonal rabbit antisera containing approximately 0.2% sodium azide as a preservative.

BD Difco Salmonella O Poly Antisera are polyvalent antisera. Each antiserum is specific for certain serogroup antigens. When properly rehydrated and used as recommended, each vial of BD Difco Salmonella O or Vi Antisera contains sufficient reagent for 60 tests. BD Difco Salmonella O Poly A-I and Vi Antisera are prepared with representative strains of these serogroups and are not absorbed. They may cross-react because of shared common O antigens.

BD Difco Salmonella O Group Antisera are specific for the major factors present in the serogroup. BD Difco Salmonella O Factor Antisera are specific for the factors of the individual serogroups. When using BD Difco Salmonella O Group Antisera, cross-reactions are possible because serogroups may share non-major group antigens. BD Difco Salmonella O Factor Antisera are absorbed as necessary to render each antiserum as specific as practical without reducing the homologous reactions to an unsatisfactory level.

BD Difco Salmonella H Poly Antisera are polyvalent antisera specific for certain flagellar antigens. Each vial of BD Difco Salmonella H Antisera contains sufficient reagent to perform between 150-1500 tests, depending on the antiserum used. BD Difco Salmonella H Antisera are either absorbed or unabsorbed specifically for either phase 1 or phase 2 antigens. BD Difco Salmonella H Antisera Spicer-Edwards are pooled, polyvalent antisera and additional adjunctive antisera to identify the more commonly occurring H antigens.

Warnings and Precautions

For *in vitro* Diagnostic Use.

This product contains dry natural rubber.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, specimens, containers, slides, tubes and other contaminated material must be sterilized by autoclaving. Directions for use should be followed carefully.

WARNING: This product contains sodium azide. Sodium azide is toxic by inhalation, by skin contact, and if swallowed. Contact with acid liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Warning

H302 Harmful if swallowed. **P264** Wash thoroughly after handling. **H301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage: Store lyophilized and rehydrated BD Difco Salmonella O, H and Vi Antisera at 2-8 °C.

Expiration date applies to the product in its intact container when stored as directed.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

Discard any antiserum that is cloudy or has a precipitate after rehydration or storage unless it can be clarified by centrifugation or filtration and demonstrates proper reactivity with validated positive and negative controls. Gross turbidity indicates contamination and such antisera should be discarded.

SPECIMEN COLLECTION AND PREPARATION

Clinical specimens: *Salmonella* can be recovered from selective differential media such as Hektoen Enteric Agar or XLD agar. For specific recommendations, consult appropriate references.^{11,12} Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as a *Salmonella* species. After these criteria are met, serological identification can be performed.

Food samples: *Salmonella* can be recovered when samples are processed to recover injured microorganisms and prevent overgrowth of competing microorganisms. Consult appropriate references for recommended procedures for isolation of *Salmonella* from foods.^{13,14} Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as a *Salmonella* species. After these criteria have been met, serological identification can be performed.

The isolate for serological testing should be subcultured from selective media to a nonselective agar.

PROCEDURE

Materials Provided: BD Difco Salmonella O Antisera; BD Difco Salmonella H Antisera; or BD Difco Salmonella Antiserum Vi.

Materials Required But Not Provided: Slide Test: 0.85% NaCl solution, sterile; agglutination slides with 1 inch squares; applicator sticks; boiling water bath; centrifuge.

Tube Test: 0.85% NaCl solution, sterile; culture tubes, 12 x 75 mm, and rack; water bath, 50 ± 2 °C; serological pipettes, 1 mL; formaldehyde.

Reagent Preparation: Equilibrate all materials to room temperature prior to performing the tests. Ensure that all glassware and pipettes are clean and free of residues such as detergents.

To reconstitute, add 3 mL of sterile 0.85% NaCl solution and rotate gently to completely dissolve the contents. Reconstituted antisera are considered a 1:2 dilution.

Test Isolate for Autoagglutination

1. From the test culture on nonselective media, transfer a loopful of growth to a drop of sterile 0.85% saline on a clean slide and emulsify the organism.
2. Rotate the slide for 1 min, then observe for agglutination.
3. If agglutination (autoagglutination) occurs, the culture is rough and cannot be tested. Subculture to nonselective agar, incubate and test the organism again as described in steps 1 and 2.
4. If no agglutination occurs, proceed with testing the organism.

Choice of Antisera

Salmonella O Antisera: This serological identification Scheme I (See Table 1 below) begins with BD Difco Salmonella O Antisera Poly A through Poly G, which contain the following:

Salmonella Group Antisera	Somatic Groups Present
Salmonella O Antiserum Poly A	A,B,D,E ₁ (E ₂ ,E ₃)*, E ₄ ,L
Salmonella O Antiserum Poly B	C ₁ ,C ₂ ,F,G,H
Salmonella O Antiserum Poly C	I,J,K,M,N,O
Salmonella O Antiserum Poly D	P,Q,R,S,T,U
Salmonella O Antiserum Poly E	V,W,X,Y,Z
Salmonella O Antiserum Poly F	51-55
Salmonella O Antiserum Poly G	56-61

*Strains of groups E₂ and E₃ are lysogenized by phage 15, then by phage 34. These strains are now classified into group E₁.²

If agglutination occurs, use individual BD Difco Salmonella O Group Antisera to determine the specific serogroup to which the isolate belongs. For efficiency, test first with individual BD Difco Salmonella O Group Antisera B, D and C₁ (the most common serogroups).

If agglutination does not occur with Poly A or B, test the isolate with BD Difco Salmonella Antiserum Vi. If positive, heat and retest with BD Difco Salmonella Antiserum Vi. If agglutination occurs after boiling with BD Difco Salmonella Antiserum Vi, the isolate is not likely to be *Salmonella*. If agglutination does not occur with BD Difco Salmonella Antiserum Vi after boiling, test boiled culture with individual BD Difco Salmonella O Antisera. If test results are indeterminate, isolate should be sent to a reference laboratory for additional testing.

If agglutination does not occur with Poly C, D, E, F and G, the isolate is not likely to be *Salmonella*.

Table 1. – Scheme I for using Salmonella O Antisera Poly Groups A, B, C, D, E, F and G.

Test with	Salmonella O Antisera Poly Groups A, B, C, D, E, F and G				
Test Result	+	– with Poly A or B		– with Poly C, D, E, F and G	
Test with	Individual Salmonella O Antisera	Salmonella O Antiserum Vi			
Test Result	+ with one Salmonella O Antiserum (required)	+	–		
Test with		Heat and retest with Salmonella O Antiserum Vi			
Test Result		+	–		
Test Conclusion or Next Action	Determine the Salmonella H Antigen	Test isolate is not a <i>Salmonella</i>	Test boiled culture with individual Salmonella O Antisera	Test isolate is not a <i>Salmonella</i>	Test isolate is not a <i>Salmonella</i>

Salmonella O Antiserum Poly A-I and Vi Scheme II (See Table 2 below): This antiserum detects factors 1-16, 19, 22-25, 34 and Vi. This combination of factors represents the most frequently isolated Groups A-I and the Vi antigens and is used to screen possible *Salmonella* isolates.

A positive reaction indicates that further serological testing is needed to identify the isolate using **BD Difco** Salmonella O Group Antiserum. The most common serogroups are B, D and C₁. For efficiency, first use the **BD Difco** Salmonella O Group Antiserum for these serogroups.

If the isolate is positive with **BD Difco** Salmonella O Antiserum Poly A-I & Vi but negative with the specific somatic antiserum, test the isolate with **BD Difco** Salmonella Antiserum Vi. If positive with **BD Difco** Salmonella Antiserum Vi, heat and retest with **BD Difco** Salmonella Antiserum Vi. If positive with **BD Difco** Salmonella O Antiserum Vi after boiling, the isolate is not likely to be *Salmonella*. If agglutination does not occur with **BD Difco** Salmonella Antiserum Vi after boiling, test boiled culture with individual **BD Difco** Salmonella O Antiserum. If test results are indeterminate, isolate should be sent to a reference laboratory for additional testing.

A negative reaction with **BD Difco** Salmonella O Antiserum Poly A-I & Vi indicates the isolate is not in serogroups A-I. If the biochemical reactions are consistent with *Salmonella*, a serogroup other than A-I is possible. Further testing with antisera for other serogroup antigens is necessary.

Table 2. – Scheme II for using Salmonella O Antiserum Poly A-I & Vi.

Test with	Salmonella O Antiserum Poly A-I & Vi				
Test Result	+				–
Test with	Individual Salmonella O Antiserum				
Test Result	+	–			
Test with		Salmonella O Antiserum Vi			
Test Result		+	–		
Test with		Heat and retest with Salmonella O Antiserum Vi			
Test Result		+	–		
Test Conclusion or Next Action	Determine the Salmonella H Antigen	Test isolate is not a <i>Salmonella</i>	Test boiled culture with individual Salmonella O Antiserum	Test isolate is not a <i>Salmonella</i>	May be a <i>Salmonella</i> detectable by use of Salmonella O Antiserum Poly C, D, E, F or G

Salmonella O Group Factor and Single Factor Antiserum: Use selected **BD Difco** Salmonella O Group Antiserum for further identification of isolates. Cross-reactions may occur between serogroups that share O antigens. Consider this partial list of **BD Difco** Salmonella O Group Antiserum as an example:

Salmonella O Antiserum Group A Factors 1, 2, 12; Salmonella O Antiserum Group B Factors 1, 4, 5, 12; Salmonella O Antiserum Group B Factors 1, 4, 12, 27

Factors 1 and 12 occur in combination with other antigens and may cause cross-reactions. The strength of the reactions will help in interpretation. Rapidly forming 3+ or greater agglutination indicates a homologous reaction. Use selected **BD Difco** Salmonella O Factor Antiserum. Absorbed antiserum specific for an identifiable antigen in a given serogroup is used to identify the isolate further. In the example above, **BD Difco** Salmonella O Factor Antiserum could be used:

Salmonella O Antiserum Factor 2; Salmonella O Antiserum Factor 4; Salmonella O Antiserum Factors 4, 5; Salmonella O Antiserum Factor 5

Polyvalent Salmonella H Antiserum: Further identification of a *Salmonella* isolate includes characterization of the flagellar antigens. Serogrouping with the following **BD Difco** Polyvalent H Antiserum can be done:

Salmonella H Poly Group Antiserum	Flagellar Antigens Present
Salmonella H Antiserum Poly a-z	Groups EN,G,L,Z ₄ , 1 complexes and a-k,r-z,Z ₆ ,Z ₁₀ ,Z ₂₉
Salmonella H Antiserum Poly A	Groups a,b,c,d,i,Z ₁₀ ,Z ₂₉
Salmonella H Antiserum Poly B	Groups eh,en,ex,enz ₁₅ , G complex
Salmonella H Antiserum Poly C	Groups k,l,r,y,z,Z ₄
Salmonella H Antiserum Poly D	Groups Z ₃₅ ,Z ₃₆ ,Z ₃₇ ,Z ₃₈ ,Z ₃₉ ,Z ₄₁ ,Z ₄₂
Salmonella H Antiserum Poly E	1 complex, Z ₆

Absorbed H antiserum specific for single antigens or a complex of antigens can be used to identify the isolate further.

Unabsorbed and Absorbed Salmonella H Antiserum: Complete identification of a *Salmonella* isolate involves analysis of phase 1 and phase 2 antigens using H antiserum. For the complex pattern of analysis and procedures, consult appropriate references.⁹

Salmonella H Antiserum Spicer-Edwards: **BD Difco** Salmonella H Antiserum Spicer-Edwards is used for screening and identifying the most commonly encountered *Salmonella* using a combination of polyvalent and single complex antiserum.

Table 3. – Identification of Salmonella H using Salmonella H Antiserum Spicer-Edwards.

H Antigen(s)	Salmonella H Antiserum Spicer-Edwards				H Antigen(s)	Salmonella H Antiserum Spicer-Edwards			
	1	2	3	4		1	2	3	4
a	+	+	+	–	k	–	+	+	+
b	+	+	–	+	r	–	+	–	+
c	+	+	–	–	y	–	+	–	–
d	+	–	+	+	z	–	–	+	+
e, h	+	–	+	–	Z ₄ Complex**	–	–	+	–
G Complex*	+	–	–	+	Z ₁₀	–	–	–	+
i	+	–	–	–	Z ₂₉	–	+	+	–

*The G complex component of Salmonella H Antiserum Spicer-Edwards 1 and 4 reacts with antigens f,g; f,g,s; f,g,t; g,m; g,m,q; g,m,s; g,m,s,t; g,m,t; g,p; g,p,s; g,p,u; g,q; g,s,t; g,t; m,p,t,u and m,t.

**The Z₄ Complex component reacts with Z₄,Z₂₃; Z₄,Z₂₄ and Z₄,Z₃₂.

Note that no antigen is positive with all four Salmonella H Antiserum Spicer-Edwards. Any antigen that reacts with all four sera should be checked for smoothness.

Slide Test Procedure

Salmonella O and Vi Antiserum

Use this procedure to test the isolate with each selected antiserum.

- Dispense 1 drop (35 µL) of each antiserum to be tested on an agglutination slide.
- Negative control:** Dispense 1 drop of 0.85% sterile NaCl solution on an agglutination slide.
- From a solid agar medium, transfer a portion of a loopful of an isolated colony to each reaction area above and mix thoroughly.
- Positive control:** Dispense 1 drop of each **BD Difco** Salmonella O Antiserum to be tested on an agglutination slide. Add 1 drop of an appropriate **BD Difco** QC Antigen Salmonella or stock cultures of known serological identification.
- Rotate the slides for 1 min and read for agglutination. Results must be read within 1 min.

Tube Test Preparation

- 0.6% Formalized Saline:** Prepare by adding 6 mL formaldehyde per 1000 mL of sterile 0.85% NaCl solution.
- Test organism:** It is often necessary to increase the motility of the test organism. To accomplish this, make several consecutive transfers in Motility GI Medium.
 - Inoculate the tube slightly below the surface of the medium using the stab method.
 - Incubate at 35–37 °C for 18–20 h.
 - Transfer only those organisms that have migrated to the bottom of the tube.
 - When the organism successfully travels 50-60 mm through the medium in 18–20 h, it is ready for use.
 - An infusion broth such as Veal Infusion Broth is recommended for cultivating motile *Salmonella* prior to testing. It should be inoculated and incubated at 35 °C for 24 h. Brain Heart Infusion Broth may be used with incubation at 35 °C for 4–6 h. If Tryptic Soy Broth is used, incubate at 35 °C for 24 h.
 - Prepare the test organism suspension by using equal volumes of broth culture and 0.6% formalized saline. The final density of this test suspension should be that of a McFarland Turbidity Standard No. 3.
- Positive control:** Commercially prepared QC Salmonella H antigens are not available. The user must maintain stock cultures of known serological identification for use in quality control. Prepare the antigen by using known serotypes and following the procedure described above. (See **Test organism**, above.)
- Salmonella H Antiserum:** Reconstituted antisera are considered a 1:2 working dilution. Prepare dilutions as follows and use on the day prepared. Discard any unused portion.
 - Most Salmonella H Antiserum:** Prepare a 1:250 dilution by adding 0.1 mL reconstituted antiserum to 24.9 mL of 0.85% NaCl solution. After mixing equal amounts (0.5 mL) of diluted antiserum and test isolate, the final dilution is 1:1,000.
 - Salmonella H Antiserum x, z₁₅ and z₂₈:** Prepare a 1:125 dilution by adding 0.1 mL reconstituted antiserum to 12.4 mL of 0.85% NaCl solution. After mixing equal amounts (0.5 mL) of diluted antiserum and test isolate, the final dilution is 1:500.
 - Salmonella H Antiserum Poly a-z:** Prepare a 1:25 dilution by adding 0.1 mL reconstituted antiserum to 2.4 mL of 0.85% NaCl solution. After mixing equal amounts (0.5 mL) of diluted antiserum and test isolate, the final dilution is 1:100.

Tube Test Procedure

Salmonella H Antiserum

- Prepare a 12 x 75 mm culture tube for each organism to be tested.
- Diluted antiserum:** Dispense 0.5 mL in each tube.
- Test isolate:** Add 0.5 mL to the appropriate tube.
- Positive control:** Add 0.5 mL of antigen positive control to a tube containing 0.5 mL of antiserum.
- Negative control:** Add 0.5 mL of 0.85% NaCl solution to a tube containing 0.5 mL of test isolate.
- Incubate all tubes in a water bath at 50 ± 2 °C for 1 h.
- Read for flocculation (agglutination).
- Repeat the Tube Test using a phase-reversed test organism. (See the procedure for **Phase Reversal** below.)

Phase Reversal

- Prepare Motility GI Medium phase reversal medium according to directions.
- Prepare the antiserum opposite to the phase desired. For example, incubating *Salmonella* Typhimurium phase 1[1] in GI Motility Medium containing i antiserum allows growth and spread of S. Typhimurium phase 2 [1,2].
- Add 1 mL of a 1:10 dilution of antiserum to 25 mL of sterile GI Motility Medium and mix well. Pour into a sterile Petri dish and allow to solidify.
- Inoculate by punching the edge of the solidified medium.
- Incubate at 35–37 °C for 24 h.
- Transfer growth from the spreading edge opposite the inoculation site to a liquid medium for testing according to steps under **Tube Test Procedure – Salmonella H Antiserum**.
- If motility is not acceptable, pass through Motility GI Medium again.

Salmonella H Antiserum Spicer-Edwards

- Prepare the test organism and the 1:2 antiserum dilution as described above in **Tube Test Preparation**.
- Final 1:1,000 dilution of antiserum:** Prepare by adding 0.1 mL of reconstituted antiserum (1:2 working dilution) to 24.9 mL of 0.85% NaCl solution.
- Prepare 4 culture tubes (12 x 75 mm) for each test organism.
- Salmonella H Antiserum Spicer-Edwards 1-4:** Add 0.5 mL of the diluted antiserum to the culture tubes.
- Test organism:** Add 0.5 mL to each tube.
- Incubate tubes in water bath at 50 ± 2 °C for 1 h.
- Remove from the water bath. Avoid excessive shaking when the tubes are in the water bath or when removing them from the water bath prior to reading the reactions.
- Read for flocculation (agglutination).

User Quality Control: At the time of use, apply both homologous and heterologous controls to check performance of the antiserum, techniques and methodology. **BD Difco** QC Antigens Salmonella may be used as a homologous control. Refer to the **BD Difco** QC Antigens Salmonella package insert for further instructions.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Slide Test

- Read and record results as follows:

4+	100% agglutination, background is clear to slightly hazy.	1+	25% agglutination, background is cloudy.
3+	75% agglutination, background is slightly cloudy.	–	No agglutination.
2+	50% agglutination, background is moderately cloudy.		

- The positive control should show 3+ or greater agglutination.
- The negative control should show no agglutination.
- For test isolates, a 3+ or greater agglutination is a positive result.
- A partial (less than 3+) or delayed agglutination reaction should be considered negative.
- If an H antigen identification is required, proceed to the next section.

Tube Test

- Read and record results as follows:

4+	100% agglutination, background is clear to slightly hazy.	1+	25% agglutination, background is cloudy.
3+	75% agglutination, background is slightly cloudy.	–	No agglutination.
2+	50% agglutination, background is moderately cloudy.		

- The positive control should show 3+ or greater agglutination at the routine test dilution (RTD).
- The negative control should show no agglutination.

When using **BD Difco** Salmonella H Spicer-Edwards, compare results with the flocculation (agglutination) patterns for the Spicer-Edwards schema (See Table 3, above).

LIMITATIONS OF THE PROCEDURE

1. Complete O and H antigen characterization of a *Salmonella* isolate is required for final identification. Due to the complexity of the laboratory procedures, identification with polyvalent antisera may be sufficient for most laboratories.
2. Possible *Salmonella* isolates having inconsistencies in biochemical reactions and O and H antigen tests should be referred to a reference laboratory for further testing.
3. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent making a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False-positive reactions may occur.
4. Rough culture isolates do occur and will agglutinate spontaneously, causing agglutination of the negative control reaction (autoagglutination). Smooth colonies must be selected and tested in serological procedures.
5. In the slide agglutination procedure for O antigen testing, it is recommended that several colonies be tested and that unabsorbed polyvalent antisera be used followed by absorbed single factor antisera. For example, colonies of a 1,2,12 culture on an agar plate will have varying degrees of each antigen. A 1,2,12 antiserum absorbed of 1 and 12 antibodies will be highly specific but will show weak or no agglutination with colonies that have less of antigen 2 and more of antigens 1 and 12. Using unabsorbed **BD Difco** *Salmonella* O Antiserum Group A Factors 1,2,12 to test several suspicious colonies on a plate followed by testing with absorbed **BD Difco** *Salmonella* O Antiserum Factor 2 gives the needed balance of sensitivity and specificity.
6. Agglutination reactions of 3+ or greater are interpreted as positive reactions. Cross-reactions resulting in a 1+ or 2+ agglutination are likely since there are somatic antigens shared among different groups as non-major group antigens.
7. There may exist common antigens between various "O" serogroups of *Salmonella*. As an example, **BD Difco** *Salmonella* O Antiserum Poly A contains, among others, agglutinins for factor 1, since cultures possessing factor 1 were used in immunization. It may be expected that this polyvalent antiserum will react with cultures other than those contained in "O" serogroups A, B, D, E and L due to the common 1 antigen (those organisms in Group G₁, G₂, H, R, T, etc., which contain factor 1).
8. **BD Difco** *Salmonella* O Antiserum Poly A-I & Vi has been prepared with representative members of those somatic groups and has not been absorbed. It is obvious that this serum may and will react with higher O groups of *Salmonella*.
9. The tube agglutination technique is recommended for H antigen testing because cross-reactions with somatic antigens may occur at the dilutions used in the slide technique.
10. In the tube test, make certain that the proper dilution is prepared for a given antiserum. Various dilutions are used for various antisera. Refer to the **Tube Test Preparation** section for details.

PERFORMANCE CHARACTERISTICS

Salmonella O Antisera and Salmonella Antiserum Vi: Sensitivity of **BD Difco** *Salmonella* O Antisera and *Salmonella* Antiserum Vi are determined by demonstrating appropriate reactivity, as defined in the Results, Slide Test section, against a battery of homologous *Salmonella* cultures (See Table 4). Specificity is determined by demonstrating non-reactivity against non-related (heterologous) *Salmonella* groups.

Salmonella H Antisera

Sensitivity of **BD Difco** *Salmonella* H Antisera is determined by demonstrating appropriate reactivity, as defined in the Results, Tube Test section, against a battery of homologous *Salmonella* cultures (See Table 5). Specificity is determined by demonstrating non-reactivity against non-related (heterologous) *Salmonella* groups.

AVAILABILITY: See page 15.

REFERENCES

1. McWhorter-Murlin, A.C., and F.W. Hickman-Brenner. 1994. Identification and serotyping of *Salmonella* and an update of the Kauffmann-White Scheme. Centers for Disease Control and Prevention, Atlanta, Ga.
2. Popoff, M.Y., and L. LeMinor. 1997. Antigenic formulas of the *Salmonella* serovars. WHO Collaborating Centre for Reference and Research on *Salmonella*. Institut Pasteur, Paris, France.
3. Holt, J.G., N.R. Krieg, P.H. Sneath, J.T. Staley and S.T. Williams. 1994. Bergey's manual of determinative bacteriology, 9th ed. Williams & Wilkins, Baltimore, MD.
4. Old, D.C. 1992. Nomenclature of *Salmonella*. J. Med. Microbiol. 37:361-363.
5. Penner, J.L. 1988. International committee on systematic bacteriology taxonomic subcommittee on *Enterobacteriaceae*. Int. J. Syst. Bacteriol. 38:223-224.
6. LeMinor, L., and M.Y. Popoff. 1987. Request for an opinion. Designation of *Salmonella enterica* sp. nov., nom. rev., as the type and only species of the genus *Salmonella*. Int. J. Syst. Bacteriol. 37:465-468.
7. Wayne, L.G. 1991. Judicial Commission of the International Committee on Systematic Bacteriology. Int. J. Syst. Bacteriol. 41:185-187.
8. Wayne, L.G. 1994. Actions of the Judicial Commission of the International Committee on Systematic Bacteriology on requests for opinions published between January 1985 and July 1993. Int. J. Syst. Bacteriol. 44:177.
9. Ewing, W.H. 1986. Edwards and Ewing's identification of *Enterobacteriaceae*, 4th ed. Elsevier Science Publishing Co., Inc., New York, N.Y.
10. Farmer III, J.J., A.C. McWhorter, D.J. Brenner and G.D. Morris. 1984. The *Salmonella-Arizona* group of *Enterobacteriaceae*: nomenclature, classification and reporting. Clin. Microbiol. Newsl. 6:63-66.
11. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry and M.A. Pfaller. 2007. Manual of clinical microbiology, 9th ed American Society for Microbiology, Washington, D.C.
12. Isenberg and Garcia (ed.). 2004 (update, 2007). Clinical microbiology procedures handbook, 2nd ed., American Society for Microbiology, Washington, D.C.
13. U.S. Food and Drug Administration. 2001. Bacteriological analytical manual online. <<http://www.cfsan.fda.gov/~ebam/bam-mm.html>>.
14. Andrews, W.H., R.S. Flowers, J. Silliker and J.S. Bailey 2001. *Salmonella*. In F.P. Downes and K. Ito (ed.), Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington, D.C.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

BD Difco Salmonella O Antisera **Difco Salmonella H Antisera** **Difco Salmonella Antiserum Vi**

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Benex Limited