English: pages 1 – 4 Italiano: pagine 11 – 14
Français : pages 4 – 7 Español: páginas 14 – 17
Deutsch: Seiten 8 – 11

R_x Only

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőtől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Inštrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

BD Difco™ E. Coli O Antiserum O157 and **BD Difco E. Coli H Antiserum H7** are used in tube agglutination tests for identifying *Escherichia coli* O157:H7.

SUMMARY AND EXPLANATION

E. coli O157:H7 is a foodborne pathogen that can cause potentially fatal enteric-related disease in humans.¹⁻⁷ This disease is characterized by sudden onset of severe cramps and abdominal pain, followed by a watery stool that may become markedly bloody. A series of 107 outbreaks involving 387 persons was traced to imported Camembert cheese in the United States in 1971.⁸ *E. coli* O157:H7 was recognized as a cause of hemorrhagic colitis in 1982⁴ and hemolytic uremic syndrome in 1983.⁹ The 1982 outbreak was derived from ingested hamburgers.^{4,10}

The incidence of disease caused by this organism has increased significantly over the past decade.^{3,11} The largest outbreak of *E. coli* O157:H7 disease occurred during January 1993, in Washington State where more than 600 patients with hemorrhagic colitis were confirmed.¹ The source of the outbreak was identified as undercooked hamburger at multiple outlets of the same fast food restaurant chain.

E. coli O157:H7 is an enteric pathogen that requires only a low inoculum to cause disease. Transmission is usually via high-volume food items whose preparation is not always under stringent control and is served to a target audience (children and the elderly are most at risk for complications of the illness). The organism has been isolated from several foods, including undercooked hamburger, drinking water, new potatoes, turkey roll, raw milk and apple cider. Serotyping of the entero-hemorrhagic *E. coli* is useful in the epidemiological documentation of the spread of a particular strain in a foodborne outbreak.¹¹

PRINCIPLES OF THE PROCEDURE

BD Difco E. Coli O Antiserum O157 is used in the tube agglutination technique for O antigen titration. **BD Difco E. Coli H Antiserum H7** is used in the tube agglutination technique for detecting H antigens. These antisera are used to confirm the presence of *E. coli* O157:H7 after selective isolation.

The serological technique is based on the reaction of a specific antiserum with its homologous antigen. While the specificity of serological methods is not absolute, serotyping of *E. coli*, along with biochemical characteristics, can provide an accurate identification of the etiological agent.

REAGENTS

BD Difco E. Coli O Antiserum O157 and **BD Difco E. Coli H Antiserum H7** are lyophilized, polyclonal rabbit antisera containing approximately 0.2% sodium azide as a preservative.

Warnings and Precautions

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, specimens, containers, slides, tubes and other contaminated material must be sterilized by autoclaving. Directions for use should be followed carefully.

WARNING: This product contains sodium azide. Sodium azide is toxic by inhalation, by skin contact, and if swallowed. Contact with acid liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Warning

H302 Harmful if swallowed.

P264 Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P301+P312** IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell. **P330** Rinse mouth. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage: Store lyophilized and rehydrated **BD Difco** E. Coli Antisera at 2–8 °C. The expiration date applies to the product in its intact container when stored as directed.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

Discard any antiserum that is cloudy or has a precipitate after rehydration or storage unless it can be clarified by centrifugation or filtration and demonstrates proper reactivity with validated positive and negative controls. Gross turbidity indicates contamination and such antisera should be discarded.

SPECIMEN COLLECTION AND PREPARATION

For isolation of *E. coli* O157 from bloody stool samples, use MacConkey Sorbitol Agar. For isolation of this organism from food, use SMAC-CT (Sorbitol MacConkey Agar w/Cefixime and Tellurite).¹² However, high levels of contaminating coliform organisms will mask the O157 strains.¹³

Procedures for serological confirmation by **BD Difco** E. Coli Antisera require a pure culture of the test organism isolated on Veal Infusion Agar or other enriched solid medium.

PROCEDURES

Materials Provided: **BD Difco** E. Coli O Antiserum O157 and **BD Difco** E. Coli H Antiserum H7.

Materials Required But Not Provided: Veal Infusion Agar, Motility GI Medium, Water bath, 50 ± 2 °C, Formalin, McFarland Barium Sulfate Standard No. 3, Culture tubes 12 x 75 mm and rack, Serological pipettes, 1 mL, Sterile 0.85% NaCl solution.

Reagent Preparation

Equilibrate all materials to room temperature before performing the tests.

Ensure that all glassware and pipettes are clean and free of residues such as detergents.

BD Difco E. Coli antisera: To rehydrate, add 3 mL sterile 0.85% NaCl solution to each vial. Rotate gently to dissolve contents completely. The rehydrated antisera are considered a 1:2 working dilution.

USER QUALITY CONTROL

At the time of use, test both positive and negative control cultures to check performance of the antiserum, techniques and methodology.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

Test Procedure

Tube Test for O Antigen Titration

1. Prepare pure cultures of the test organism by plating onto Veal Infusion Agar and incubating at 35 ± 2 °C for 16–18 h.
2. Suspend the test isolate growth obtained from solid agar medium in 0.85% NaCl solution to give a homogeneous suspension.
3. Heat the organism suspension in a boiling water bath for 30–60 min. The culture should be homogeneous. Precipitation indicates a rough culture and the suspension should be discarded.
4. Allow the suspension to cool; dilute with 0.85% NaCl solution to a density approximating that of a McFarland Barium Sulfate Standard No. 3.
5. Add formalin to a final concentration of 0.5% by volume.
6. In a rack, prepare a row of 8 test tubes for each organism suspension to be tested.
7. Dispense 0.9 mL 0.85% NaCl solution to the first tube in each row and 0.5 mL to the remaining tubes.
8. Prepare serial dilutions using the rehydrated **BD Difco** E. Coli O Antiserum O157. Using a 1 mL serological pipette, add 0.1 mL of antiserum to tube 1 in each row and mix thoroughly. Transfer 0.5 mL from tube 1 to tube 2 and

mix thoroughly. In a like manner, continue transferring 0.5 mL through tube 7, discarding 0.5 mL from tube 7 after mixing. Tube 8 is a negative control tube containing only 0.85% NaCl solution.

9. Add 0.5 mL of the test organism suspension to each of the 8 tubes. Upon addition of the test organism suspension, final dilutions will be 1:40 through 1:2,560 for tubes 1 through 7, respectively.
10. Shake the rack to mix. Incubate in a 50 ± 2 °C water bath for 18–24 h.
11. Read for agglutination.

Tube Test for H Antigen Detection

1. Prepare an actively motile culture of the suspect *E. coli* culture by several successive transfers in Motility GI Medium. At least 2–3 passages through Motility GI Medium are necessary before attempting to establish the presence and identity of H antigens. Fresh isolates of *E. coli* generally have poorly developed flagella.
2. Inoculate a loopful of the Motility GI Medium culture into a tube of Veal Infusion Broth. Incubate 6–8 h at 35 ± 2 °C or overnight, if necessary.
3. Inactivate the culture by adding formalin to a final concentration of 0.3% (0.3 parts formaldehyde per 100 parts of the Veal Infusion Broth culture). If necessary, adjust the density of the suspension with formalized saline to approximate a McFarland Barium Sulfate Standard No. 3. This broth culture will be used as the test antigen in step 6.
4. Prepare a 1:500 dilution of **BD Difco** E. Coli H Antiserum H7 by adding 0.2 mL of the rehydrated 1:2 working solution to 49.8 mL of 0.85% NaCl solution.
5. Pipette 0.5 mL of the diluted antiserum into a test tube.
6. Add 0.5 mL of the test organism broth culture from step 3 and mix well. The resulting antiserum dilution will be 1:1,000.
7. Incubate in a 50 ± 2 °C water bath for 1 h.
8. Read for agglutination.

RESULTS

Observe test result with indirect lighting against a dark background. Record agglutination as follows:

- 4+ 100% agglutination of cells; supernatant fluid is clear to slightly hazy.
- 3+ 75% agglutination of cells; supernatant fluid is slightly cloudy.
- 2+ 50% agglutination of cells; supernatant fluid is moderately cloudy.
- 1+ 25% agglutination of cells; supernatant fluid cloudy.
- ± Less than 25% agglutination of cells.
- No agglutination

BD Difco E. Coli O Antiserum O157 – Tubes showing a 2+ or greater agglutination at 1:320 or greater are considered positive.

BD Difco E. Coli H Antiserum H7 – Tubes showing 2+ or greater agglutination are considered positive.

LIMITATIONS OF THE PROCEDURE

1. Serological techniques employing E. Coli Antisera serve as corroborative evidence for the identification of *E. coli* O157:H7. Final identification cannot be made without consideration of morphological, serological and biochemical characterization.
2. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent making a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False-positive reactions may occur.
3. Rough culture isolates occur and will agglutinate spontaneously, causing agglutination in a normal serum (autoagglutination). Smooth colonies must be selected and tested in serological procedures.
4. Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.
5. It is important in this test to use the recommended time and temperature of incubation. In addition, care should be taken to make certain that the water bath is in a location free of mechanical vibration.
6. The tube agglutination technique is recommended for testing because cross-reactions with non-related antigens may occur at the dilutions used in the slide technique.

PERFORMANCE CHARACTERISTICS¹⁴

The performance of the **BD Difco** E. Coli O Antiserum O157, manufactured by Difco Laboratories, was compared to a commercially available latex test in a study published by Chapman.¹⁴ Pure cultures of 91 bacterial isolates were tested.^{14,15} The table below lists the results obtained. All *E. coli* previously confirmed as O157 were positive by the **Difco** tube agglutination method. All other organisms tested were negative by the **Difco** tube agglutination method.

No. of Samples	Strain	BD Difco E. Coli Antiserum O157 Tube Agglutination Test Positives	Commercially Available Latex Test Positives
30	<i>E. coli</i> O157	30	30
30	<i>E. coli</i> non-O157	0	0
30	Non-sorbitol fermenting organisms	0	0
1	<i>Escherichia hermanii</i>	0	0

AVAILABILITY

Cat. No. Description

229701 **BD Difco™** E. Coli O Antiserum O157, 3 mL

221591 **BD Difco™** E. Coli H Antiserum H7, 3 mL

REFERENCES

- Centers for Disease Control. 1993. Emerging infectious diseases. *Morb. Mort. Wkly.* 42:257–260.
- Glass, K.A., J.M. Leffelholtz, J.P. Ford, and M.P. Doyle. 1992. Fate of *Escherichia coli* O157:H7 as affected by pH or sodium chloride and in fermented, dry sausage. *Appl. Environ. Microbiol.* 58:2513–2516.
- Padhye, N.V., and M.P. Doyle. 1992. *Escherichia coli* O157:H7 epidemiology, pathogenesis, and methods for detection in food. *J. Food Prot.* 55:555–565.
- Riley, L.W., R.S. Remis, and S.D. Helgerson, et al. 1983. Hemorrhagic colitis associated with a rare *Escherichia coli* serotype. *N. Engl. J. Med.* 308:681–685.
- Samadpour, M., L. Grimm, B. Desai, D. Alfi, J.E. Ongerth, and P.I. Tarr. 1993. Molecular epidemiology of *Escherichia coli* O157:H7 strains using bacteriophage A–restriction fragment length polymorphism analysis: application to a multistate foodborne outbreak and a day care center cluster. *J. Clin. Microbiol.* 31:3179–3183.
- Tarr, P.I. 1994. Review of 1993 *Escherichia coli* O157: H7 outbreak: western United States. *Dairy, Food, and Environmental Sanitation.* 14:372–373.
- Gray, L.D. 1995. *Escherichia*, *Salmonella*, *Shigella*, and *Yersinia*, p. 450–455. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C.
- Marrier, R., J.G. Wells, R.C. Swanson, W. Callahan, and I.J. Mehlman. 1973. An outbreak of enteropathogenic *E. coli* foodborne disease traced to imported French cheese. *Lancet* 2:1376–1378.
- Flowers, R.S., W. Andrews, C.W. Donnelly, and E. Koenig. 1993. Pathogens in milk and milk products, p. 168–176. In R. T. Marshall (ed.), *Standard methods for the examination of dairy products*, 16th ed. American Public Health Association. Washington, D.C.
- Morbidity and Mortality Weekly Reports. November, 1982.
- Lior, H. 1994. *Escherichia coli* O157: H7 and verotoxigenic *Escherichia coli* (VTEC). *Dairy, Food, and Environ. Sanit.* 14:378–382.
- Association of Official Analytical Chemists. 2001. FDA bacteriological analytical manual online. <<http://www.cfsan.fda.gov/~ebam/bam-mm.html>>.
- Meng, J.P. Feng and M.P. Doyle. 2001. Pathogenic *E. coli*, p. 331–341. In P. Downes and K. Ito (ed.), *Compendium of methods for the microbiological examination of foods*, 4th ed. American Public Health Association, Washington, D.C.
- Chapman, P.A. 1989. Evaluation of commercial latex slide test for identifying *Escherichia coli* O157. *J. Clin. Pathol.* 43:1109–1110.
- Chapman, P.A., D.J. Wright and P. Norman. 1989. Verotoxin-producing *Escherichia coli* infections in Sheffield: cattle as a possible source. *Epidem. Infect.* 102:439–45.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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For US: "For Investigational Use Only"

Lower limit of temperature / Долен лимит на температурата / Dolní hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperatuuripiir / Limite inférieure de température / Najniža dozvoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Температураның төменгі рұқсат шегі / 하한 온도 / Žemiausia laikymo temperatūra / Temperatūras zemākā robeža / Laagste temperatuurlimiet / Nedre temperaturgrænse / Dolna granica temperatury / Limite minimo de temperatura / Limitā minimā de temperaturā / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sicaklık alt sınırı / Мінімальна температура / 温度下限

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

Difco is a trademark of Difco Laboratories, Inc., a subsidiary of Becton, Dickinson and Company.

© 2016 BD. BD, and the BD Logo are trademarks of Becton, Dickinson and Company.