Legionella GVPC Supplement

Selective supplement for isolation of *Legionella* spp.

DESCRIPTION

Legionella GVPC Supplement is a lyophilized selective supplement used in combination with Legionella Growth Supplement (ref. 80056) for the preparation of Legionella BCYE Agar Base (ref. 610049, 620049). The complete medium is used for detection and enumeration of *Legionella* spp from clinical specimens and environmental samples, according to ISO 11731 (all parts).

KIT CONTENTS

Each kit contains:

- 10 vials of freeze-dried Legionella GVPC Supplement;
- 1 instruction sheet.

PRINCIPLE OF THE METHOD

Legionella GVPC Supplement is a selective mix inhibiting contaminant microflora.

COMPOSITION

	Content / vial	Content / liter of medium
Glycine	1.5 g	3.0 g
Vancomycin	0.5 mg	1.0 mg
Polymyxin B	40000 IU	80000 IU
Cycloheximide	40.0 mg	80.0 mg

PROCEDURE FOR USE

1. Aseptically reconstitute the content of one vial of Legionella GVPC Supplement with 5 ml sterile distilled water.
2. Mix to complete dissolution and add 1 ml Legionella GVPC Supplement to 90 ml Legionella BCYE Agar Base, autoclaved and cooled to 45-50°C. To prepare the complete medium, also add 10 ml Legionella Growth Supplement (ref. 80056), previously reconstituted.
3. Mix well and dispense in Petri dishes.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation of the medium being prepared.

QUALITY CONTROL

1. Control of the appearance: freeze-dried white in colour.
2. Microbiological control.
Prepare the medium per label directions. Inoculate the plates with the strains indicated below and incubate at $36 \pm 2^\circ\text{C}$ for 2-5 days in humidified atmosphere.

Control strains

Legionella pneumophila WDCM 00107
Escherichia coli ATCC® 25922

Growth

Good
Partially to totally inhibited

PRECAUTIONS

Legionella GVPC Supplement is classifiable as hazardous under current legislation; it is recommended that the Safety Data Sheet be consulted on its use. The product is a selective supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration. Once reconstituted, the product can be stored for a maximum duration of 30 days at -20°C, shielded from light.

REFERENCES

- ISO 11731-2:2004. Water quality – Detection and enumeration of Legionella. Part 2: Direct membrane filtration method for waters with low bacterial count.
- ISO 11731:1998. Water quality – Detection and enumeration of Legionella.
- Edelstein P.H. (1981) Improved semiselective medium for the isolation of *Legionella pneumoniae* from contaminated clinical and environmental specimens. J. Clin. Microbiol. 14(3):298.

PRESENTATION

Product	Ref.	Contents
Legionella GVPC Supplement	81008	10 vials

One vial is sufficient to prepare 500 ml of medium.

TABLE OF SYMBOLS

 Batch code	 In Vitro Diagnostic Medical Device	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	 Do not reuse

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BIBLIOGRAFIA

- ISO 11731-2:2004. Water quality – Detection and enumeration of *Legionella*. Part 2: Direct membrane filtration method for waters with low bacterial count.
- ISO 11731:1998. Water quality – Detection and enumeration of *Legionella*.
- Edelstein P.H. (1981) Improved semiselective medium for the isolation of *Legionella pneumoniae* from contaminated clinical and environmental specimens. J. Clin. Microbiol. 14(3):298.

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LIOFILCHEM® S.r.l.