V.P. TEST EP Reagent

Reagents for performing the Voges-Proskauer test.

DESCRIPTION

V.P. TEST EP Reagent contains a 6% (w/v) solution of alpha-naphthol and a 40% (w/v) solution of potassium hydroxide used to carry out the Voges-Proskauer test for the identification of Enterobacteriaceae.

KIT CONTENT

Each package contains:

- 5 vials of 10 ml of sodium alpha-naphthol;
- 5 vials of 10 ml of potassium hydroxide;
- 2 empty dropper bottles;
- instruction sheet.

METHOD PRINCIPLE

The Voges-Proskauer reaction is used as qualitative test to distinguish between members of Enterobacteriaceae on the basis of the capability to produce acetylmethylcarbinol (acetoin) as a fermentation product from glucose. The positive reaction to the Voges-Proskauer test is denoted by the development of a red color obtained from acetoin oxidation in the presence of alkali and oxygen.

COMPOSITION

Alpha-naphthol: 6% (w/v) solution of alpha naphthol in ethyl alcohol.

KOH 40%: 40% (w/v) solution of potassium hydroxide.

TEST PROCEDURE AND RESULTS INTERPRETATION

Use V.P. TEST EP on the EnteroPluti-Test (ref. 78618/78619) as indicated in the instructions for use that accompany the product. Briefly, lay EnteroPluti-Test with its flat surface pointing upward and, punching the plastic film, inject with a syringe 3 drops of alpha-naphthol and 2 drops of potassium hydroxide in the sector VP. The positive test is given by the development of a bright red colour within 20 minutes after addition of reagents.

QUALITY CONTROL FOR THE USER

1. Appearance control: the alpha naphthol solution appears limpid and beige; the potassium hydroxide solution appears clear limpid colourless.
2. Microbiological control.

Control strains		VP test	
<i>Enterobacter cloacae</i>	ATCC® 13047	Red colour	Positive
<i>Escherichia coli</i>	ATCC® 25922	No colour development	Negative

LIMITATIONS

One should allow at least 20 min for color to develop before considering the Voges-Proskauer test negative.

PRECAUTIONS

V.P. TEST EP Reagent, is classified as hazardous under current legislation; for its use it is recommend to consult the MSDS. V.P. TEST EP Reagent is a disposable device to be used only for in vitro diagnostic use, is intended for professional use and should be used in the laboratory by properly trained operator, with approved methods of aseptic and safety against pathogens.

STORAGE AND TRANSPORT CONDITIONS

Store V.P. TEST EP Reagent at 2-8°C away from light in its original package, until the expiry date shown on the label. However, our stability studies have shown that the storage or transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, do not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

ELIMINATION OF THE USED MATERIAL

After use, V.P. TEST EP Reagent and the material that comes in contact with the sample must be decontaminated and disposed in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infectious materials.

BIBLIOGRAPHY

1. Barritt, M.M. (1936) The intensification of the Voges-Proskauer reaction by the addition of α -naphthol. J. Pathol. Bacteriol. 42:441-454.
2. Edwin H. Lennette (1995) Manual of Clinical Microbiology.
3. Blazevic, D.J., and Ederer, G.M. (1975) Principles of biochemical tests in diagnostic microbiology. 63-67. New York, John Wiley & Sons.

PRESENTATION

Product	Ref.	Contents
V.P. TEST EP Reagent	80281	5 x 10 ml vials 5 x 10 ml vials

TABLE OF SYMBOLS

LOT Batch code	IVD In vitro Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains Sufficient for <n> tests	 Caution, consult accompanying documents	 Do not reuse

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Z.A de Gsvrine - 4 rue Képler - B.P.4125
44241 La Chapelle-sur-Erdre Cedex - France
t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00
commercial@humeau.com

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