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Z.A de Gesvrine – 4 rue Képler – B.P.4125
44241 La Chapelle-sur-Erdre Cedex – France
t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00
commercial@humeau.com

w w w . h u m e a u . c o m

FOOD SYSTEM

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FOOD SYSTEM

ENGLISH

**System for detection and presumptive identification
of pathogenic microorganisms from foodstuffs validated to ISO 16140 standard.**

DESCRIPTION

FOOD-SYSTEM is a 24-well system containing desiccated biochemical substrates and culture media for detection and presumptive identification of microorganisms from meat, milk and cheese, fish and other food products.

The system provides detection and presumptive identification of: *Salmonella* spp., *Proteus* spp. / *Providencia* spp., *Pseudomonas* spp., *Staphylococcus aureus*, *E. coli*, *Bacillus cereus*, *Listeria* spp., yeasts and moulds, and particularly it is validated to ISO 16140 standard for the detection of *Salmonella* spp. and *Listeria* spp.

The system is inoculated with a suspension of food sample and incubated at 36 ± 1 °C for 18-24 hours.

The tests for detection and presumptive identification of the microorganisms present in the sample are interpreted by assessing the colour change in the various wells.

KIT CONTENTS

The kit contains:

20 FOOD SYSTEM (40 tests)	1 <i>TEST RESULT FORM</i>
40 vials of <i>Physiological solution</i> (4.5 mL/vial)	20 Sealing films
1 Instruction sheet	

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

VASELIN OIL (ref. 80278)	H ₂ O ₂ Reagent (ref. 80057)
OXIDASE TEST STICK (ref. 88029)	BUFFERED PEPTONE WATER (ref. 24099, 412090, 611014)
KOVAC'S – REAGENT (ref. 80270)	RINGER'S <i>Solution</i> (ref. 81059)
Sundry materials used in microbiology laboratories; cover slides and slides; microscope.	

CONFIGURATION

The configuration of the system is shown in Table 1. Each system allows two samples to be examined or one sample in double test.

Table 1:

Well	PRESUMPTIVE IDENTIFICATION of <i>Salmonella</i> spp.
1-LDC	Decarboxylation of lysine
2-H₂S	Production of hydrogen sulphide
Well	PRESUMPTIVE IDENTIFICATION of <i>Proteus</i> spp.
3-UR	Hydrolysis of urea
4-PRO	<i>Proteus</i> spp. / <i>Providencia</i> spp.
Well	PRESUMPTIVE IDENTIFICATION of <i>Pseudomonas</i> spp.
5-PSE	<i>Pseudomonas</i> spp.
Well	PRESUMPTIVE IDENTIFICATION of <i>Staphylococcus aureus</i>
6-STA	<i>Staphylococcus aureus</i>
Well	PRESUMPTIVE IDENTIFICATION of <i>Escherichia coli</i> and enteropathogenic <i>E.coli</i> O157
7-ESC	<i>Escherichia coli</i>
8-IND	Indole test for confirmation of <i>Escherichia coli</i>
Well	PRESUMPTIVE IDENTIFICATION of <i>Bacillus cereus</i>
9-BCE	<i>Bacillus cereus</i>
Well	PRESUMPTIVE IDENTIFICATION of <i>Listeria</i> spp.
10-LIS	<i>Listeria</i> spp.
11-CAT	Catalase test for confirmation of <i>Listeria</i> spp.
Well	PRESUMPTIVE IDENTIFICATION of Yeasts and moulds
12-Y/M	Yeasts and moulds

PRINCIPLE OF THE METHOD

FOOD SYSTEM allows the detection and presumptive identification of the following microorganisms: *Salmonella* spp., *Proteus* spp. / *Providencia* spp., *Pseudomonas* spp., *Staphylococcus aureus*, *E. coli*, *E. coli* 0157, *Bacillus cereus*, *Listeria* spp., yeasts and moulds.

- The presence of ***Salmonella* spp.** is shown by a yellow to red-orange colour change of well **1-LDC** and by the yellow to black color change of well **2-H₂S**.

The confirmation of *Salmonella* spp. is performed by taking a drop of the culture broth from the well **2-H₂S** and plating it onto a selective agar medium (i.e. XLD agar, ref. 10056).

Proceed to identification of bacterial colonies according to ISO 6579:2004 ².

The yellow to black color change of well **2-H₂S** and no color change in the well **1-LDC** can indicate the presence of ***Citrobacter* spp.**

- The presence of ***Proteus* spp. / *Providencia* spp.** is shown by the yellow to red-fuchsia color change of the well **3-UR**, and by the yellow to brown-black color change of the well **4-PRO**.

- The presence of ***Pseudomonas* spp.** is shown by the yellow to turbid green colour change of the well **5-PSE**.

The confirmation of *Pseudomonas* spp. is performed by the oxidase test directly in the well (OXIDASE TEST STICK ref. 88029).

- The presence of ***Staphylococcus aureus*** is shown by the appearance of a black ring on the bottom of the well **6-STA**.

The confirmation of *Staphylococcus aureus* is performed by taking a drop of the culture broth from the well **6-STA** by means of a sterile loop, and plating it onto a selective agar medium (i.e. Baird Parker Agar + RPF, ref. 10521).

Proceed to identification of bacterial colonies according to ISO 6888-2:1999 ³.

- The presence of ***E. coli*** and ***E. coli* 0157** is shown by the red to blue color change of the well **7-ESC** and by the formation of a pink-red coloring following the addition of Kovac's reagent into the well **8-IND**.

- The presence of ***Bacillus cereus*** is shown by the yellow to turbid green color change of the well **9-BCE**.

The confirmation of *Bacillus cereus* is performed by taking a drop of the culture broth from the well **9-BCE** by means of a sterile loop, and plating it onto a selective agar medium (i.e. Bacillus Cereus Agar-MYP, ref. 10027).

Proceed to identification of bacterial colonies according to ISO 7932:2004 ⁴.

- The presence of ***Listeria* spp.** is shown by the yellow to black color change of the well **10-LIS** and by the development of bubbles following the addition of Reagent H₂O₂ into the well **11-CAT**.

The confirmation of *Listeria* spp. is performed by taking a drop of the culture broth from the well **10-LIS** by means of a sterile loop, and plating it onto a selective agar medium (i.e. O.A. Listeria Agar, ref. 10620).

Proceed to identification of bacterial colonies according to ISO 11290:2002 ⁵.

- The presence of **yeasts and moulds** is shown by the green to yellow color change of the well **12-Y/M** and by watching for chlamydospores and hyphae at the microscope.

COMPOSITION

Table 2:

Well	Contents
1-LDC	Culture medium with substrate for the decarboxylation of lysine
2-H ₂ S	Culture medium with substrate for the production of hydrogen sulphide
3-UR	Culture medium with substrate for hydrolysis of urea
4-PRO	Culture medium with substrate for the isolation of <i>Proteus/Providencia</i> spp.
5-PSE	Culture medium with substrate for the isolation of <i>Pseudomonas</i> spp.
6-STA	Culture medium with substrate for the isolation of <i>Staphylococcus aureus</i>
7-ESC	Culture medium with substrate for the isolation of <i>Escherichia coli</i>
8-IND	Culture medium to show the production of bubbles
9-BCE	Culture medium with substrate for the isolation of <i>Bacillus cereus</i>
10-LIS	Culture medium with substrate for the isolation of <i>Listeria</i> spp.
11-CAT	Culture medium to show the catalase reaction
12-Y/M	Culture medium with substrate for the isolation of yeasts and moulds

Physiological solution (g/l): Sodium chloride 9 g, Distilled water 1000.0 mL, pH 7.0 ± 0.2

COLLECTION AND STORAGE OF THE SAMPLES

The samples to be subjected to the test must be collected in accordance with the standard methods for samples to be subjected to microbiological examination and must be sent to the laboratory for seeding into the **FOOD SYSTEM** immediately after being taken.

TEST PROCEDURE**a) PREPARATION OF THE SAMPLE**

1. Homogenize a quantity of foodstuff (1, 10 or 25 g) in a proper volume (9, 90 or 225 mL) of BUFFERED PEPTONE WATER or RINGER'S Solution in accordance with the standard procedures for the food sample in question. The final dilution of the sample should be 1:10. Store the homogenized sample at 2-4 °C until test ending.
2. Take 10 mL of homogenized sample (dil. 1:10) and transfer into a suitable tube.
3. Cover the tube and incubate at 36 ± 1 °C for 12-18 hours.

b) INOCULATION OF THE SYSTEM

1. Take a system from its wrapper and allow it to come to room temperature.
2. Write down the sample data and the starting date of the examination.
3. Transfer 0.5 mL, from the pre-incubated tube containing sample suspension, into a vial of *Physiological Solution* contained in the kit.
4. Distribute 0.2 mL (4-5 drops) of suspension in each well of the system.
5. Cover the first 3 wells **1-LDC**, **2-H₂S** and **3-UR** with 2 drops each of Vaseline oil.
6. Cover the system with the lid provided and incubate at 36 ± 1 °C for 18-24 hours.

The system is configured for performance of two tests. In case of performance of one test, cover the unused half system with the sealing film provided in the kit.

INTERPRETATION OF THE RESULTS

At the end of the incubation:

- Add 2 drops of Kovac's Reagent to the well **8-IND** and watch for the appearance of a pink-red ring.
- Add 2 drops of Reagent H₂O₂ to the well **11-CAT** and watch for the formation of bubbles.
- Observe the color changes in of the wells and interpret the results according to the Table 3.
- Note the results on the *TEST RESULT FORM* (copy as many forms as necessary).
-

Table 3:

Well	Presumptive Identification of <i>Salmonella</i> spp.	Well color	
		Positive reaction	Negative reaction
1-LDC	Decarboxylation of lysine	red-orange	yellow
2-H₂S	Production of hydrogen sulphide	black	yellow
Well	Presumptive Identification of <i>Proteus</i> spp. / <i>Providencia</i> spp.	Well color	
		Positive reaction	Negative reaction
3-UR	Hydrolysis of urea	red-fuchsia	yellow
4-PRO	<i>Proteus/ Providencia</i> spp.	brown-black	yellow
Well	Presumptive Identification of <i>Pseudomonas</i> spp.	Well color	
		Positive reaction	Negative reaction
5-PSE	<i>Pseudomonas</i> spp.	turbid green	yellow
Well	Presumptive Identification of <i>Staphylococcus aureus</i>	Well color	
		Positive reaction	Negative reaction
6-StA	<i>Staphylococcus aureus</i>	black	yellow
Well	Presumptive Identification of <i>Escherichia coli</i> , <i>Escherichia coli</i> 0157	Well color	
		Positive reaction	Negative reaction
7-ESC	<i>Escherichia coli</i> , <i>Escherichia coli</i> 0157	blue	grey-red
8-IND	Indole test	pink-red	yellow
Well	Presumptive Identification of <i>Bacillus cereus</i>	Well color	
		Positive reaction	Negative reaction
9-BCE	<i>Bacillus cereus</i>	turbid green	yellow
Well	Presumptive Identification of <i>Listeria</i> spp.	Well color	
		Positive reaction	Negative reaction
10-LIS	<i>Listeria</i> spp.	black	yellow
11-CAT	Catalase test	presence of bubbles	absence of bubbles
Well	Presumptive Identification of Yeasts and Moulds	Well color	
		Positive reaction	Negative reaction
12-Y/M	Yeasts and Moulds	yellow	verde
	Microscopic observation	presence of chlamydospores and hyphae	absence of chlamydospores and hyphae

QUALITY CONTROL

Each batch of **FOOD SYSTEM** is subjected to quality control using the following reference strains of microorganisms:

<i>Salmonella typhimurium</i>	ATCC 14028	<i>Pseudomonas aeruginosa</i>	ATCC 27853	<i>Bacillus cereus</i>	ATCC 11778
<i>Citrobacter freundii</i>	ATCC 8090	<i>Staphylococcus aureus</i>	ATCC 25923	<i>Listeria monocytogenes</i>	ATCC 35152
<i>Proteus mirabilis</i>	ATCC 25933	<i>Escherichia coli</i> 0157	ATCC 35150	<i>Saccharomyces cerevisiae</i>	ATCC 9763

PERFORMANCE

The validation of **FOOD SYSTEM** has been represented by a comparative study with the procedure indicated in the ISO 11290 for the isolation of *Listeria monocytogenes*, and with the procedure indicated in the ISO 6579 for the isolation of *Salmonella* spp. Overall 712 samples of meat and dairy products have been tested in duplicate with **FOOD SYSTEM** and reference method. The obtained results are indicated in the following table:

	<i>Listeria</i> spp. detection		<i>Salmonella</i> spp. detection	
	Meat products	Dairy products	Meat products	Dairy products
Relative Accuracy	94,7 %	97,2 %	94,0 %	92,4%
Relative Sensitivity	87,5 %	96,8 %	87,5 %	85,7%
Relative Specificity	100,0 %	97,6 %	100,0 %	100 %
Selectivity	100%	100%	100%	100%
Detection Limits	5- 10 UFC/ g	50- 100 UFC/ g	5- 10 UFC/ g	50- 100 UFC/ g

FACTORS THAT MAY INVALIDATE THE RESULTS

- Inaccurate standardisation of the inoculum.
- Material examined unsuitable.
- Use of expired systems.
- Temperature and time of incubation other than those recommended.

LIMITS AND WARNINGS

The presence of *Salmonella* spp., *Listeria* spp., *Escherichia coli*, *Staphylococcus aureus* and *Bacillus cereus* must be confirmed using appropriate culture media and biochemical and serological tests. For definitive identification of these micro-organisms, it is necessary to perform supplementary confirmation tests.

PRECAUTIONS

The product **FOOD SYSTEM** cannot be classified as hazardous under current legislation, nor does it contain harmful substances in concentrations $\geq 1\%$; it therefore does not require a Safety Data Sheet to be available.

FOOD SYSTEM must be used in the laboratory by properly trained personnel using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes of temperature.

In such conditions the product will remain valid until the expiry date indicated on the label. Do not use beyond that date.

Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **FOOD SYSTEM** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Code	Package
FOOD SYSTEM	71680	40 test

TABLE OF SYMBOLS

LOT Batch code	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	

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1. Difco & BBL Manual, Manual of Microbiological Culture Media 2003.
2. **ISO 6579:2002**. Horizontal method for the detection of *Salmonella* spp.
3. **ISO 6888-2:1999**. Enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species). - Technique that uses Rabbit Plasma Fibrinogen (RPF) Agar Medium.
4. **ISO 7932:2004**. Horizontal method for the enumeration of presumptive *Bacillus cereus*.
5. **ISO 11290:2002**. Horizontal method for the detection and enumeration of *Listeria monocytogenes*.
6. **ISO 16140:2003**. Microbiology of food and animal feeding stuffs - Protocol for the validation of alternative methods.

Microbiology Products

LIOFILCHEM S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto D.A. (TE) - ITALY

Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net