# X FACTOR TEST V FACTOR TEST V+X FACTOR TEST

## DESCRIZIONE

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# X FACTOR TEST V FACTOR TEST V+X FACTOR TEST

## DESCRIPTION

X FACTOR TEST, V FACTOR TEST and V+X FACTOR TEST are constituted by paper discs with special features, containing the respective coagulation factors, used for differentiating *Haemophilus* spp.

## CONTENT OF THE PACKAGES

Each package contains:

- 2 cartridges with 50 discs each, packaged in a heat-sealed container
- 1 dryer
- 1 instruction sheet

## PRINCIPLE OF THE METHOD

Different strains of *Haemophilus* grow on a culture medium only in presence of the coagulation factor or factors (X, V, or both) which they need. These different requirements allow the differentiation and the identification of *Haemophilus* spp.

## COMPOSITION

- Each disc of X FACTOR TEST contains Hemin.
- Each disc of V FACTOR TEST contains NAD (Nicotinamide-Adenine-Dinucleotide).
- Each disc of V+X FACTOR TEST contains Hemin and NAD (Nicotinamide-Adenine-Dinucleotide).

## TEST PROCEDURE

1. Take the cartridges container from the refrigerator and leave it on the test bench until it reaches room temperature (about 30 minutes). This will prevent humidity being deposited on the discs when the package is opened, which could prejudice their long-term stability.
2. Using a sterile swab, evenly inoculate the surface of a plate of Tryptic Soy Agar (ref. 10037) or Mueller Hinton Agar (ref. 10031) with a pure suspension of the microorganism to test.
3. Using sterile tools, press one disc of X FACTOR TEST, V FACTOR TEST and V+X FACTOR TEST on the inoculated surface, at a distance of 120° from one to another and at 1-2 cm from the edge of the plate. Incubate at 36±1°C for 24-48 in 5-10% carbon dioxide atmosphere.

## INTERPRETATION OF THE RESULTS

If the organism requires X Factor alone, it will grow only at the edge of the X and the X+V FACTOR TEST discs; if it requires V Factor alone, it will grow only at the edge of the V and the X+V FACTOR TEST discs; if both X and V Factors are required, it will grow only at the vicinity of the X+V FACTOR TEST discs.

Some examples are indicated in the following table:

	Without Factors	Factor X	Factor V	Factors V+X
<i>Haemophilus influenzae</i>	-	-	-	+
<i>Haemophilus aegyptius</i>	-	-	-	+
<i>Haemophilus parainfluenzae</i>	-	-	+	+
<i>Haemophilus ducreyi</i>	-	+	-	+
<i>Bordetella pertussis</i>	+	+	+	+

## QUALITY CONTROL

Each batch of X FACTOR TEST, V FACTOR TEST and V+X FACTOR TEST is subjected to microbial control, inoculating pure suspensions of *Haemophilus influenzae* ATCC 19418 and *Haemophilus parainfluenzae* ATCC 7901 on plates of Tryptic Soy Agar.

## PRECAUTIONS

X FACTOR TEST, V FACTOR TEST and V+X FACTOR TEST cannot be classified as being hazardous according to the current legislation. X FACTOR TEST, V FACTOR TEST and V+X FACTOR TEST is a disposable device to be used only for diagnostic use *in vitro*. It must be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

## STORAGE

Store X FACTOR TEST, V FACTOR TEST and V+X FACTOR TEST at 2-8°C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, X FACTOR TEST, V FACTOR TEST and V+X FACTOR TEST will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

## DISPOSAL OF USED MATERIAL

After use, X FACTOR TEST, V FACTOR TEST and V+X FACTOR TEST and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

## BIBLIOGRAPHY

- Kilian M. (1980) *Haemophilus*. in Manual of Clinical Microbiology. Eds. Lennette et al. Amer. Soc. for Microbiol. 3<sup>rd</sup> edn. Washington.

## PRESENTATION

Product	REF	µg	Σ
X FACTOR TEST	9503	5	100
V FACTOR TEST	9504	4	100
V+X FACTOR TEST	9505	4+5	100

## TABLE OF SYMBOLS

 IVD	In Vitro Diagnostic Medical Device		Do not reuse
 REF	Catalogue number		Fragile, handle with care
	Manufacturer		Contains sufficient for <n> tests
	Use by		Caution, consult accompanying documents
	Temperature limitation		Batch code

LIOFILCHEM® S.r.l.

Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

Tel. +390858930745

Fax +390858930330

Website: [www.liofilchem.net](http://www.liofilchem.net)

E-Mail: [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)

F00005

Rev. 1/ 28.03.2013