

3 surface sockets Unswitched

Ref. : 4668

Description :

- 3 single phase sockets 2P+E 220V 10/16A for fixing to skirting, compartment for services or fascia strip
- Supply of 3m 1000 Ro2V x 2.5² core cable
- Supply of fixing screws

4 surface sockets Switched, contactor with indicator

Ref. : 4669

Description :

- 4 single phase sockets 2P+E 220V 10/16A for fixing to skirting, compartment for services or fascia strip
- Supply of 3m 1000 Ro2V x 2.5² core cable
- Supply of fixing screws

3 surface sockets Switched, Shuko pattern (German standard)




Ref. : 4641

Description :

- 3 single phase sockets 2P+E 220V 10/16A for fixing to skirting, compartment for services or fascia strip
- Supply of 3m 1000 Ro2V x 2.5² core cable
- Supply of fixing screws

The accessories electricity and

14



3 surface sockets Unswitched, Shuko pattern (German standard)

Ref. : 4640

Description :

- 3 single phase sockets 2P+E 220V 10/16A for fixing to skirting, compartment for services or fascia strip
- Supply of 3m 1000 Ro2V x 2,5² core cable
- Supply of fixing screws

Surface socket Single, unswitched

Ref. : 4635

Description :

- 1 single phase sockets 2P+E 220V 10/16A for fixing to skirting, compartment for services or fascia strip
- Supply of fixing screws


Electrical bracket 2 unswitched sockets

Ref. : 4631

Description :

- 1 single phase sockets 2P+E 220V 10/16A for fixing to skirting, compartment for service, fascia strip or on a weighing table
- With cable

fitting



Taps on table ❶

Swan neck mixer CW/HW at 250	Ref. : 4774
Swan neck CW at 150	Ref. : 4770
Swan neck CW at 250	Ref. : 4778
Swan neck DW at 150	Ref. : 4777
Swan neck DW at 300	Ref. : 4779
CA	Ref. : 4771
Vacuum	Ref. : 4772
Gas single tap	Ref. : 4773
Gas double tap at 90°	Ref. : 4766

Taps on skirting-shelf ❷

CW outlet	Ref. : 4760
Swan neck mixer CW/HW at 250	Ref. : 4759
Swan neck CW at 150	Ref. : 4764
Swan neck CW at 250	Ref. : 4768
Swan neck DW at 160	Ref. : 4767
Swan neck DW at 250	Ref. : 4769
CA	Ref. : 4761
Vacuum	Ref. : 4762
Gas single tap	Ref. : 4763
Gas double tap at 90°	Ref. : 4765

Taps under central compartment ❸

CW double tap (1/2")	Ref. : 4780
CA double tap (3/8")	Ref. : 4781
Vacuum double tap (3/8")	Ref. : 4782
Gas double tap (3/8")	Ref. : 4783

Remote controlled taps ❹

Swan neck CW at 150	Ref. : 4790
Swan neck DW at 150	Ref. : 4797
CW outlet	Ref. : 4795
CA outlet	Ref. : 4791
Vacuum outlet	Ref. : 4792
Gas outlet	Ref. : 4793
Liquid jet vacuum pump	Ref. : 4775

Key
 CA : Compressed air
 CW : Cold water
 HW : Hot water
 DW : Distilled water