

TITOLCHIMICA S.P.A.
Via Dell'Artigianato, 2
45030 Villamarzana (RO)
Italia/Italy
C.F. e P.IVA: IT 00748970290
R.E.A. RO n.92331
Reg. Imprese C.C.I.A.A. RO n.4362

ACIDO SOLFORICO purissimo 96%

Code	TC43480
Batch Number	22602512
Production date	26/02/2025
Expiry date	march 2027
Appearance	Clear liquid

Feature	u.m.	Result	Specification
Assay	% p/p	95,8	95,5 - 96,5
Density (20° C)	g/ml	1,836	1,834 - 1,842

This product has been prepared and controlled in accordance to the Quality Management System UNI EN ISO 9001 of Titolchimica S.p.A. certified by Certiquality. The internal control method used is the MC1315.

Each “Result” is the result of a measurement, the value attributed to the measurand obtained by measurement. (Second Edition of the Guide EURACHEM/CITAC CG4 - Appendix B - Definitions).

Titolchimica S.p.A., in order to maintain appropriate homogeneity of the values found for each production lot, has established specification limits ("Specification") within which these values must lie at the time of production.

The scales used for weighing raw materials are calibrated by the ACCREDIA LAT N. 044 calibration center and are daily adjusted by using mass standards calibrated by the aforementioned calibration center, in reference to mass standards calibrated at the National Institute of Research and Metrology (INRIM).

Date of certificate
26/02/2025

QUALITY CONTROL LABORATORY

Ferris Vittoria

THE CERTIFICATE REFERS TO THE PRODUCT CODE. THE LETTERS PRESENT ON THE LABEL AFTER THE CODE ARE OMITTED HEREIN AS THEY EXCLUSIVELY IDENTIFY THE PACKAGING SIZE.