

ASTOR 20 COLONY COUNTER

The innovative counter based on pressure impulses and LEDs



The new Astor 20 Colony Counter is very easy to use. Thanks to its Petri dish holder, back-lit and sensitive to pressure, the enumeration of the microbial grown colonies on the agar surface can be performed by a simple and common marker pen. The Petri holder is back-lit by LEDs: this new technology grants a better light diffusion, a higher reliability and a minimum power consumption (2 W, as a total). Astor 20 does not require any particular marker, connected by any cable.

Technical characteristics:

- ☺ Easy-to-read display with automatic counter and Reset key
- ☺ Acoustic signal to confirm every colony count
- ☺ Transparent and back-lit Petri dish holder, sensitive to pressure
- ☺ Every kind of commercial marker can be used: no need of a special one
- ☺ Petri dishes with diameters up to 110 mm can be used
- ☺ 90-mm Petri dish adapter included
- ☺ Minimum power consumption: 2 W
- ☺ It can be used either with 110 or 220 V

Optional accessory:

- ☺ Benchtop lens with lamp and adjustable arm

Distribué par :

Z.A de Gesvrine - 4 rue Képler - B.P.4125
44241 La Chapelle-sur-Erdre Cedex - France
t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00
commercial@humeau.com

